Histostain[®]-Plus Kits

Invitrogen[®] 2nd Generation

LAB-SA Detection System

<u>CAT. NO.</u>	<u>FOR PRIMARY ANTIBODY</u>	<u>CHROMOGEN</u>	<u>GOOD FOR</u>
85-9643	Broad Spectrum*	DAB	150 slides

*Will react with mouse, rabbit, guinea pig and rat primary antibodies.

INTENDED USE

For In Vitro Diagnostic Use. CAUTION: Not for human or animal therapeutic use. Uses other than the labeled intended use may be a violation of local law.

Interpretation must be made within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Histostain[®]-Plus kits are designed to reveal antigens that have reacted with a user-supplied primary antibody on human tissue or cell samples.

Frozen or paraffin-embedded tissues, freshly prepared lymphocytes, and fixed cultured cells can be used.

BACKGROUND

In 1987 Invitrogen Laboratories, Inc. introduced the Histostain[®]-SP detection kit—the first commercial kit to use the labeled-(strept)avidin-biotin (LAB-SA) method. Since then, LAB-SA (also known as Streptavidin-Biotin Amplification) has been acknowledged as the premier detection method for IHC⁽¹⁻⁴⁾. LAB-SA methodology is now widely used in basic research and routine testing. Histostain[®]-Plus kits use Invitrogen's second generation LAB-SA technology. This enhanced LAB-SA system provides enhanced sensitivity and performance in immunohistochemistry (IHC) applications. DAB stained slides can be dehydrated and permanently mounted using Histomount[™] mounting medium. Histostain[®]-Plus kits are designed to reveal the presence of antigens in human tissue or cell preparations. These kits can be used with all common sample preparation methods, including: frozen tissues, paraffin-embedded sections, and cell smears.

After formalin-fixed, paraffin-embedded tissue sections have been deparaffinized in xylene and dehydrated in a graded series of ethanol, an endogenous peroxidase quenching step is performed. This endogenous peroxidase activity can be quenched by either using 3% hydrogen peroxide in methanol or Invitrogen's Peroxo-Block[™] (Cat. No. 00-2015). Non-immune serum is used to eliminate non-specific background. After a user-supplied primary antibody has been incubated with a sample, Invitrogen's LAB-SA method uses an enhanced, human-absorbed, biotinylated, affinity-purified secondary antibody to react with the primary antibody. Enhanced horseradish peroxidase conjugated streptavidin (HRP-SA) is then bound to the biotinylated secondary antibody. The DAB then is used to create an intense brown deposit around the antigen/antibody/enzyme complex in the sample. Histostain[®]-Plus kits can also be used with automated stainers. Call Invitrogen's technical service for details.

Invitrogen offers prediluted 2ndGen primary antibodies that are ready-to-use with Histostain[®]-Plus kits. For a complete list of our primary antibodies, call for a free copy of Invitrogen Pathology Catalog or call Invitrogen's technical service for details.

STORAGE AND SHELF LIFE

Store kit at 2-8°C. All performance claims are void after kit expiration date.

REAGENTS SUPPLIED

Reagent A. One dropper bottle (15 mL) of ready-to-use serum blocking solution, 10% non-immune serum (goat)

Reagent B. One dropper bottle (15 mL) of ready-to-use biotinylated secondary antibody

Reagent C. One dropper bottle (15 mL) of ready-to-use streptavidin-peroxidase conjugate

Reagent D1. One dropper bottle (3 mL) of concentrated substrate buffer (20x)

Reagent D2. One dropper bottle (3 mL) of concentrated chromogen solution, DAB (20x)

Reagent D3. One dropper bottle (3 mL) of 0.6% hydrogen peroxide (20x)

REAGENTS AND MATERIALS REQUIRED BUT NOT PROVIDED

– Xylene, ethanol, and absolute methanol

– Distilled or deionized water

– 30% hydrogen peroxide

– 10 mM phosphate-buffered saline (PBS), pH 7.4 or 50 mM Tris-buffer saline (TBS), pH 7.8. Optional: 0.05% Tween 20.

– Mini PAP Pen (Invitrogen Cat. No 00-8877)

– Primary antibody

– Hematoxylin (Invitrogen Cat. No. 00-8001)

– Mounting Media: GVA (Invitrogen Cat. No. 00-8000), Clearmount[™] (Invitrogen Cat. No. 00-8110), or Histomount[™] (Invitrogen Cat. No. 00-8030)

SUGGESTED STAINING PROTOCOL

A. PRELIMINARY PREPARATION OF SLIDES

<u>1. SPECIMEN PREPARATION</u>	<p>-Appropriate tissue and antigen fixation is required to obtain reproducible performance and reliable interpretations. Suitable fixatives include: 10% neutral buffered formalin, B5, Bouin's, Zinc formalin or alcohol-base fixatives.</p> <p>-10% neutral buffered formalin, B5, Bouin's, Zinc formalin or alcohol-base fixatives are regarded as suitable fixatives for most antigens of clinical significance.</p> <p>- Formalin-fixed tissues post-fixed in B5 before paraffin embedding may show improved stain. Cell smears prepared from body fluids should be made to assure a monolayer of cells. Multilayers of cells can trap staining reagents and interfere with the interpretation of the results. Smears should be fixed immediately after preparation. Depending on the properties of the antigen, cell smears are usually stable for one to two weeks when stored at 4°C.</p> <p>- Fixation of cytospin or frozen sections can be done with acetone (100%) at 4°C for a period of 10 minutes.</p>
<u>2. SLIDE PREPARATION</u>	<p>-Precoat slides with HistoGrip™ (Invitrogen Cat. No. 00-8050). As an alternative, precoat with 0.1% poly-L-lysine in water, then air dry.</p>
<u>3. DEPARAFFINIZATION AND REHYDRATION</u>	<p>-Paraffin sections are deparaffinized with xylene and rehydrated in a graded series of ethanol. Wash cell smears or tissue in a PBS bath for 10 minutes before starting the staining procedure. Note: Do not allow specimen or tissue sections to dry from this point on.</p>
<u>4. CONTROL SLIDES</u>	<p>Three control slides are necessary for the interpretation of results.</p> <ul style="list-style-type: none"> • Positive Tissue Control: A specimen, processed in the same way as the unknown, contains the antigen to be stained. • Reagent Control/Negative Control-1: An additional slide that will be treated with a non-immune serum instead of same concentration of primary antibody. Any staining observed on the specimen is probably due to non-specific protein binding or non-specific binding of other reagents. • Negative Control-2: A specimen, processed in the same way as the unknown, does not contain the antigen to be stained [optional].

B. STAINING PROTOCOL

Reagent	Specimen Preparation	Incubation Time (Min.)
<u>1. PEROXIDASE QUENCHING SOLUTION</u> (Not provided)	<p>This step is optional. Do step 1 if elimination of endogenous peroxidase activity is necessary.</p> <ol style="list-style-type: none"> a. For paraffin-embedded tissues, add 1 part of 30% hydrogen peroxide to 9 parts of absolute methanol. Mix well. b. Submerge slides in PEROXIDASE QUENCHING SOLUTION for 10 minutes. c. Wash with PBS 2 min., 3 times. Proceed to Step 2. <p>For frozen tissue, treat with Peroxo-Block™ (cat. no. 00-2015) for 45 sec. Wash immediately.</p>	10 (Optional)
<u>2. SERUM BLOCKING SOLUTION</u> Ready-to-use. Reagent A	<p>This step can be omitted if primary antibody contains a protein blocker (e.g. 10% non-immune serum).</p> <ol style="list-style-type: none"> a. Add 2 drops (100 µL) or enough to completely cover tissue, of SERUM BLOCKING SOLUTION to each section. Incubate for 10 min. b. Drain or blot off solution. DO NOT RINSE. 	10
<u>3. PRIMARY ANTIBODY</u> (Not provided).	<p>Optimal dilution and incubation times should be determined by the investigator. 2ndGen Primary antibodies are optimized for use with Histostain-Plus detection kits. (Dilution and incubation time are dependent on sample preparation, antibody affinity, amt. of antigen present, and antigen accessibility).</p> <ol style="list-style-type: none"> a. Apply 2 drops (100 µL) or enough to completely cover tissue, of PRIMARY ANTIBODY to each section. b. Incubate in moist chamber for 30-60 min. c. Rinse with PBS for 2 min., 3 times. 	30-60
<u>4. BIOTINYLATED SECOND ANTIBODY</u> Ready-to-use. Reagent B	<ol style="list-style-type: none"> a. Apply 2 drops (100 µL) or enough to completely cover tissue, of BIOTINYLATED SECOND ANTIBODY to each section. Incubate for 10 min. (For rat primary antibodies, incubate for 20 min.) b. Rinse with PBS for 2 min., 3 times. 	10
<u>5. ENZYME CONJUGATE</u> Ready-to-use. Reagent C	<ol style="list-style-type: none"> a. Apply 2 drops (100 µL) or enough to completely cover tissue, of ENZYME CONJUGATE to each section. Incubate for 10 min. b. Rinse with PBS for 2 min., 3 times. 	10
<u>6. DAB CHROMOGEN</u> Reagents D1, D2 and D3 20x concentrated components.	<ol style="list-style-type: none"> a. Add 1 drop of Reagent D1, 2 drops of Reagent D2 and 1 drop of Reagent D3 to 1 mL distilled or deionized water. Mix well. Protect from light and use within one hour. b. Apply 2 drops (100 µL) or enough to completely cover tissue, of DAB CHROMOGEN to each section. Incubate for 3-10 min. c. Rinse well with distilled water. 	3-10

From this point on, the investigator should proceed with his or her own established laboratory protocols for counterstaining and mounting. Invitrogen recommends using Hematoxylin (Cat. No. 00-8001) as a counterstain and either Histomount™ (Cat. No. 00-8030) for DAB or Clearmount™ (Cat. No. 00-8010) as a mounting solution for AEC or DAB.

www.invitrogen.com

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288 • E-mail: techsupport@invitrogen.com

PI 85-9643 (PIN 31296)

(Rev 09/09) DCC-09-1008

© 2007 Invitrogen Corporation. All right reserved. These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses.

INTERPRETATION OF RESULTS

The following table interprets results which might be achieved using various staining controls.

Case No.	+ control	reagent control	(-) control	Sample	Analysis
1	-	-	-	-	Procedure performed incorrectly.
2	+	+	+	+	Non-specific staining due to protein binding or endogenous peroxidase activity.
3	+	-	+	+/-	Negative control contains the antigen.
4	-	-	-	+	Positive control does not contain the antigen.
5	+	-	-	-	Sample does not contain the antigen.
6	+	-	-	+	Sample contains the antigen.

TROUBLESHOOTING

Possible causes for negative staining on positive slides:

1. Steps in the staining protocol were performed in incorrect sequence.
2. Primary or secondary antibody incubation steps were omitted.
3. Labile antigens were destroyed.
4. Specimen was improperly fixed and/or processed.
5. Specimen dehydrated during staining.

Possible causes for weak staining on all slides:

1. Specimen retained excess liquid after rinsing steps.
2. Incubation times were insufficient.
3. Substrate prepared improperly.
4. Deparaffinization was incomplete (staining may be accompanied by high background).

Possible causes for high background staining:

1. Endogenous peroxidase activity was incompletely blocked.
2. Deparaffinization was incomplete.
3. Excessive application of tissue adhesive.
4. Inadequate rinsing of slides.
5. Over-development of substrate.
6. Dehydration of specimen during staining.

REFERENCES

1. Elias JM, et al. *AM J Clin Pathol* 92:62-67, 1989.
2. Hsu SM and Raine L. In: *Advances in Immunochem* R.A.DeLellis, ed., Mason Publishing USA, Inc., New York, 1984, pp. 31-42.
3. Warnke R and Levy R. *J Histochem Cytochem* 28: 771-776, 1980.
4. Shi Z, et al. *J. Histochem Cytochem* 36(3): 317-322, 1988.

TRADEMARKS

NBA™, CAS-Block™, Clearmount™, HistoGrip™, Histomount™, and Peroxo-Block™ are trademarks of Zymed Laboratories, Inc. Zymed® and Histostain® are registered trademarks of Zymed Laboratories, Inc.

For more detailed information please refer to www.invitrogen.com.

IF YOU HAVE ANY QUESTIONS ABOUT THIS PRODUCT PLEASE CONTACT YOUR LOCAL DISTRIBUTOR

Authorized Representative for IVDD 98/79/EC: Invitrogen Ltd.
Inchinnan Business Park
3 Fountain Drive
Paisley
PA4 9RF
UK

www.invitrogen.com

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288 • E-mail: techsupport@invitrogen.com

PI 85-9643 (PIN 31296)

(Rev 09/09) DCC-09-1008

© 2007 Invitrogen Corporation. All right reserved. These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses.

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Bevollmächtigter Repräsentant für IVDD 98/79/EC: Invitrogen Ltd.