

Lipofectamine™ LTX Optimization

Terms and Conditions of Use

Invitrogen Corporation ("Invitrogen") grants to you, its customer, a non-exclusive, non-transferable license to access the Invitrogen transfection protocol applicable to Invitrogen products, as set forth below (the "Protocol"). You acknowledge that the Protocol is protected by copyright. All rights not specifically granted to you are expressly reserved to Invitrogen. You may use the Protocol for personal, educational, or scientific research or professional use, but in no case for a fee. You may not remove, obscure or modify any copyright or proprietary notices, author attribution or any disclaimer contained in the Protocol.

INVITROGEN PROVIDES THE PROTOCOL "AS IS", WITHOUT WARRANTIES OF ANY KIND, EITHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, WARRANTIES OF TITLE, OR IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR NON-INFRINGEMENT; YOUR USE OF THE PROTOCOL IS AT YOUR OWN RISK; NEITHER INVITROGEN NOR ITS AFFILIATES, EMPLOYEES, OFFICERS OR DIRECTORS SHALL BE LIABLE TO YOU OR ANY THIRD PARTY FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, CONSEQUENTIAL OR PUNITIVE DAMAGES ARISING OUT OF THIS LICENSE OR ANY USE OF THE PROTOCOL. Some jurisdictions do not allow the exclusion or limitation of liability for certain damages, so the above limitation may not apply to you to the extent prohibited by such local laws; if so, then Invitrogen's liability for damages hereunder shall not exceed an amount equal to the amounts paid by you hereunder, or one hundred dollars (\$100.00), whichever is less.

Confidentiality

You acknowledge the Protocol is Invitrogen's confidential information and constitutes Invitrogen's proprietary and valuable trade secrets. You agree to take reasonable care to protect the Protocol from disclosure to third parties, except as may be expressly authorized by Invitrogen.

Materials Needed

- Lipofectamine™ LTX Reagent (Cat. nos. 15338-100, 15338-500)
- PLUS™ Reagent (3 µg/µL) (Cat. no. 11514-001; also supplied with Cat. no. 15338-100)
- Opti-MEM® Reduced Serum Medium (Cat. nos. 319850, 519850, or 22600)
- 20 µg of test plasmid DNA (GFP plasmid or other)
- 24-well plate containing cells of interest at optimal confluency in 500 µL of growth medium, and ready for transfection
- Round or V-bottom 96-well plate for reagent preparation (× 1)
- 1.5 mL polypropylene microcentrifuge tubes (× 5)

Protocol

1. In four microcentrifuge tubes labeled A–D, dilute DNA and PLUS™ Reagent in Opti-MEM® media to obtain the final concentrations and volumes indicated below.
 - A) 398 µL of 5 ng/µL DNA and 2 µL PLUS™ Reagent
 - B) 396 µL of 10 ng/µL DNA and 4 µL PLUS™ Reagent
 - C) 394 µL of 15 ng/µL DNA and 6 µL PLUS™ Reagent
 - D) 392 µL of 20 ng/µL DNA and 8 µL PLUS™ Reagent

Mix well and incubate for 10 minutes at room temperature.

2. In a microcentrifuge tube labeled "LTX dilution," dilute the Lipofectamine™ LTX 1:10 in Opti-MEM® medium by adding 50 µL of Lipofectamine™ LTX to 450 µL of Opti-MEM® medium. Mix gently.
3. Following the table below, add Opti-MEM® media into the wells of the 96-well plate.

	1	2	3	4	5	6
A	46.25 µL	45 µL	43.75 µL	42.5 µL	40 µL	37.5 µL
B	42.5 µL	40 µL	37.5 µL	35 µL	30 µL	25 µL
C	38.75 µL	35 µL	31.3 µL	27.5 µL	20 µL	12.5 µL
D	35 µL	30 µL	25 µL	20 µL	10 µL	0 µL

4. Following the table below, add the diluted **Lipofectamine™ LTX** into the wells of the same 96-well plate.

	1	2	3	4	5	6
A	3.75 µL	5 µL	6.25 µL	7.5 µL	10 µL	12.5 µL
B	7.5 µL	10 µL	12.5 µL	15 µL	20 µL	25 µL
C	11.25 µL	15 µL	18.7 µL	22.5 µL	30 µL	37.5 µL
D	15 µL	20 µL	25 µL	30 µL	40 µL	50 µL

5. Add 50 µL of the DNA and PLUS™ Reagent dilution into each well across the row corresponding to the label on the tube. Mix gently by slowly pipetting up and down 3 times.


For example, add 50 µL from tube A into each well across the A row (i.e., A1, A2, A3, A4, A5, and A6). Repeat it for each of the remaining rows of B, C, and D, using the corresponding tubes so that each well contains 100 µL of the appropriate DNA-PLUS™-Lipofectamine™ LTX complex.

6. Incubate the plate containing the DNA-PLUS™-Lipofectamine™ LTX complexes for 30 minutes at room temperature.

7. Transfer the entire 100 µL of the DNA-PLUS™-Lipofectamine™ LTX complex from each well to the appropriate well of your 24-well plate containing your cells in 500 µL of growth medium (see **Example of a 24-well Optimization Plate** below). Mix gently by rocking the plate back and forth.

8. Return your cells to the incubator, and assay for gene expression in 18–48 hours.

Example of a 24-well Optimization Plate



Purchaser Notification

This product is covered by Limited Use Label License No. 5: Invitrogen Technology (see the Invitrogen catalog or our website, www.invitrogen.com). By the use of this product you accept the terms and conditions of the applicable Limited Use Label License.

©2009 Life Technologies Corporation. All rights reserved. For research use only. Not intended for any animal or human therapeutic or diagnostic use.

For technical support, email tech_support@invitrogen.com.
For country-specific contact information, visit www.invitrogen.com.