



- Add 40 µl 25 mM MES, pH 5 to final volume of 100 µl. Vortex to ensure good mixing.
- Incubate for at least 30 minutes at room temperature, or 2 hours at 4°C, with slow tilt rotation.
- After incubation, place the tube on the magnet for 4 minutes and remove the supernatant.
- Wash the coated Dynabeads as described below (see protocol 3.2).

### 3.2 Washing of Coated Beads

All immobilisation procedures require washing of the coated Dynabeads to remove excess ligand and to block un-reacted surface.

**NOTE:** In order to quench the non reacted activated carboxylic acid groups, incubate the Dynabeads coated with ligand with either 50 mM Tris pH 7.4 for 15 minutes or 50 mM ethanolamine in PBS pH 8.0 for 60 minutes, both at room temperature with slow tilt rotation.

- Wash the coated Dynabeads a total of four times with 100 µl PBS or 50 mM Tris. Blocking protein like BSA or skimmed milk powder may be added to a concentration of 0.1 - 0.5 % when it does not interfere with downstream applications of the Dynabeads. Also 0.1% Tween-20 or Triton X-100 can be added during washes to reduce nonspecific binding.
- Resuspend the coated Dynabeads to the desired concentration in PBS or a Tris storage buffer. The Dynabeads are now ready for use.

Store the coated Dynabeads at 2-8°C. Addition of 0.1 - 0.5 % protein (BSA) and/or 0.01 - 0.1 % Tween-20 or Triton X-100 is recommended to stabilize the immobilised ligand and increase the ease of handling. Coated Dynabeads can usually be stored for several months at this temperature, depending on the stability of the immobilised ligand. A final concentration of 0.02% (w/v) sodium azide (NaN<sub>3</sub>) may be added as a bacteriostatic agent. If the coated Dynabeads are stored for more than two weeks, they should be washed twice for five minutes with a buffer suitable for the application prior to use.

### 3.3 Isolation of Target Molecule

Efficient isolation of target molecules using Dynabeads is dependent on the bead-concentration, target molecule concentration, the ligand's affinity for the target molecule and time. Binding is performed from 10 minutes to 1 hour, at a recommended concentration of 1-10 x 10<sup>9</sup> beads/ml. Target-ligand equilibrium is reached after approximately 1 hour.

- Add sample containing target molecule to the coated Dynabeads (3 mg beads). For a 100 kD protein, use a volume containing approximately 25 µg target molecule to assure an excess of this molecule.
- Incubate the mixture with tilting and rotation for one hour to capture the target (incubation times as low as 10 minutes can be used with concentrated protein samples in volumes close to what was originally pipetted from the vial).
- Place the tube on the magnet for 4 minutes to collect the Dynabeads at the tube wall. For viscous samples, increase the time on the magnet. Pipette off the supernatant.
- Wash the Dynabeads 3 times using 1 ml PBS each time and exchanging buffers by the use of the magnet, according to protocol 3.2.

Efficient isolation of target molecules using Dynabeads is dependent on the bead-concentration, target molecule concentration, the ligand's affinity for the target molecule and the specific binding kinetics involved.

The concentration of required Dynabeads will depend on the size of your specific molecule. Also the salt-concentration and pH of the chosen binding, washing and elution buffers can be varied depending on the type of molecule to be immobilised. Similarly, the selected buffer used in the downstream application should be optimised for the specific application.

The size of the Dynabeads M-270 Carboxylic Acid presents a high surface area per mg beads and a corresponding high capacity for the target molecule. The effective binding capacity will depend on the size of the specific molecules to be immobilized.

As the Dynabeads M-270 Carboxylic Acid will not

inhibit enzymatic activity, bead-bound material can be used directly in downstream analysis. Alternatively, the target molecule can be eluted off the Dynabeads following conventional elution methods.

### 3.4 Target Protein Elution Procedure

Conventional elution methods can be applied for the elution of target protein from the Dynabeads. Low pH (2.8-3.5), change in ionic strength, affinity elution, electrophoresis, polarity reducing agents, deforming eluants can be applied, or even boiling the bead-target complex in SDS-PAGE application buffer for direct characterization of protein on SDS-PAGE. The method of choice depends on the affinity of the specific target molecule to the ligand coated onto the Dynabeads, the stability of the target molecule and the downstream application and detection method. Most proteins will be eluted off at pH 3.1 following the procedure described below, but some protein functionality might be lost under such harsh conditions. If maintaining functionality of the target molecule is important, try milder elution conditions first such as high salt (e.g. 2M NaI) or stepwise elution reducing pH from 6 down to 3. This is also recommended if the bead-bound ligand must remain functional to allow reuse of the Dynabeads.

- Add 30 µl 0.1 M citrate (pH 3.1) to the Dynabeads with immobilized target.
- Mix well by tilting and rotation for 2 minutes.
- Place the test tube on the magnet and transfer the supernatant, containing purified target, to a clean tube.
- Add additional 30 µl 0.1 M citrate (pH 3.1) to the Dynabeads to elute any remaining target.
- Mix well by tilting and rotation for 2 minutes.
- Place the test tube on the magnet, pipette off the eluate and pool the supernatants containing pure target molecule.

Total collected volume = 60 µl

To ensure reuse of the Dynabeads and functionality of the isolated target molecule, bring both the Dynabeads and the target molecules back to physiological pH (7.4) immediately after elution.

## 4. GENERAL INFORMATION

Invitrogen Dynal AS complies with the Quality System Standards ISO 9001:2000 and ISO 13485:2003.

### 4.1. Product Characteristics

Typical characteristics for any given batch of this product:

Diameter:	2.8 µm
Density:	1.6 g/cm <sup>3</sup>
Specific surface area:	2-5 m <sup>2</sup> /g beads
Active chemical functionality:	150 µmol/g beads
Concentration:	2 x 10 <sup>9</sup> beads/ml (approx. 30 mg/ml)

Certificate of Analysis (CoA) is available upon request. Material Safety Data Sheet (MSDS) is available at [www.invitrogen.com](http://www.invitrogen.com).

### 4.2. References

- Nakajima N and Ikade Y, "Mechanism of Amide Formation by Carbodiimide for Bioconjugation in Aqueous Media", *Bioconjugate Chem.* 1995, 6(1),123-130.
- Gilles MA, Hudson AQ and Borders CL Jr, "Stability of water-soluble carbodiimides in aqueous solution", *Anal Biochem.* 1990 Feb 1;184(2):244-248.
- Sehgal D and Vijay IK, "A method for the high efficiency of water-soluble carbodiimide-mediated amidation", *Anal Biochem.* 1994 Apr;218(1): 87-91.
- Szajani B et al, "Effects of carbodiimide structure on the immobilization of enzymes", *Appl Biochem Biotechnol.* 1991 Aug;30(2):225-231.

### 4.3. Additional Material Needed

- Magnetic device (Dynal MPC, Magnetic Particle Concentrator)
- Mixing/rotation device
- Test tubes, glassware and pipettes
- Ligands
- Buffers/solutions (see below)

### 4.4. Recommended Buffers/Solutions

Listed below are some recommended buffers for use with Dynabeads M-270 Carboxylic Acid.

0.01 M NaOH: 0.4 g NaOH (MW 40.0) dissolved in 1,000 ml distilled water.

100 mM MES pH 5: 2.13 g MES (2-[N-morpholino]ethane sulfonic acid, MW 213.25). Dissolve in 90 ml distilled water, adjust to pH 5 and adjust to 100 ml.

0.05 M Tris pH 7.4: 0.79 g Tris HCl (MW 157.6) Dissolve in 90 ml distilled water, adjust to pH 7.4 and adjust to 100 ml.

0.1 M Citrate pH 3.1: 2.10 g citric acid (C<sub>6</sub>H<sub>8</sub>O<sub>7</sub> x H<sub>2</sub>O, MW 210.14). Dissolve in 90 ml distilled water, adjust to pH 3.1 and adjust to 100 ml.

2 M NaI: 3 g NaI (MW 149.9) to 10 ml distilled water.

PBS pH 7.4 (phosphate buffered saline): 0.26 g NaH<sub>2</sub>PO<sub>4</sub> x H<sub>2</sub>O (MW 137.99) 1.44 g Na<sub>2</sub>HPO<sub>4</sub> x 2H<sub>2</sub>O (MW 177.99) 8.78 g NaCl (MW 58.5). Dissolve in 900 ml distilled water, adjust pH if necessary and adjust to 1,000 ml. PBS with 0.1% (w/v) BSA/HSA/skimmed milk: Include 0.1% (w/v) BSA/HSA/skimmed milk (0.1 g) in 100 ml PBS (above). PBS/Tween 20/Triton X: Include 0.5-1.0 % (w/v) Tween 20/Triton X (50-100 mg) in 100 ml PBS (above).

If a preservative is needed for storage of coated Dynabeads, a final concentration of 0.02% (w/v) sodium azide (NaN<sub>3</sub>) may be added to the storage buffer. This preservative is cytotoxic and must be carefully removed before use by washing. Required safety precautions must be followed when handling this material.

### 4.5. Storage & Stability

When stored in unopened vials at 2-8°C, Dynabeads M-270 Carboxylic Acid are stable until the expiration date printed on the label.

Dynabeads M-270 Carboxylic Acid coated with antibody may be stored at 4°C for several months without loss of antigen binding capacity. Coated Dynabeads M-270 should be washed once before use.

Dynabeads M-270 Carboxylic Acid should not be autoclaved, but can be incubated with ethanol (70%, 1 hour) or gamma irradiated after freeze drying.

Dynabeads M-270 Carboxylic Acid should be washed once before use (see washing procedure above).

Precautions should be taken to prevent bacterial contamination of the antibody-coated Dynabeads. If cytotoxic preservatives are added, these must be carefully removed before use by washing.

Primary antibody coated Dynabeads stored for more than two weeks should be washed once for 5 min in PBS/BSA before use.

### 4.6. Warnings & Limitations

This product is for research use only. Not intended for any animal or human therapeutic or diagnostic use unless otherwise stated.

Sodium azide is toxic if ingested. Avoid pipetting by mouth.

Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. When disposing through plumbing drains, flush with large volumes of water to prevent azide buildup.

### Intellectual Property Disclaimer

Invitrogen Dynal will not be responsible for violations or patent infringements that may occur with the use of our products.

### 4.7. Trademarks

Dynal®, Dynabeads® and Dynal® MPC™ are either registered trademarks or trademarks of Invitrogen Dynal AS, Oslo, Norway. Any registration or trademark symbols used herein denote the registration status of trademarks in the United States. Trademarks may or may not be registered in other countries.

Invitrogen Dynal will not be responsible for violations or patent infringements which may occur with the use of our products.

### 4.8 Limited Use Label License

No. 5: Invitrogen Technology – The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of

the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Invitrogen Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact

Licensing Department,  
Invitrogen Corporation,  
1600 Faraday Avenue, Carlsbad,  
California 92008.  
Phone (760) 603-7200.  
Fax (760) 602-6500.  
Email: [outlicensing@invitrogen.com](mailto:outlicensing@invitrogen.com)

### 4.9 Warranty

The products are warranted to the original purchaser only to conform to the quantity and contents stated on the vial and outer labels for the duration of the stated shelf life. Invitrogen Dynal's obligation and the purchaser's exclusive remedy under this warranty is limited either to replacement, at Invitrogen Dynal's expense, of any products which shall be defective in manufacture, and which shall be returned to Invitrogen Dynal, transportation prepaid, or at Invitrogen Dynal's option, refund of the purchase price.

Claims for merchandise damaged in transit must be submitted to the carrier.

This warranty shall not apply to any products which shall have been altered outside Invitrogen Dynal, nor shall it apply to any products which have been subjected to misuse or mishandling. ALL OTHER WARRANTIES, EXPRESSED, IMPLIED OR STATUTORY, ARE HEREBY SPECIFICALLY EXCLUDED, INCLUDING BUT NOT LIMITED TO WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. Invitrogen Dynal's maximum liability is limited in all events to the price of the products sold by Invitrogen Dynal. IN NO EVENT SHALL INVITROGEN DYNAL BE LIABLE FOR ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES. Some states do not allow limits on warranties, or on remedies for breach in certain transactions. In such states, the limits set forth above may not apply.

**Invitrogen Dynal is a part of the Invitrogen Group.**

**Contact details for your local Invitrogen sales office/technical support can be found at <http://www.invitrogen.com/contact>**

© Copyright 2007 Invitrogen Dynal AS, Oslo, Norway.  
All rights reserved.  
SPEC-06043