

Fluorescent Deoxyribonuclease I Conjugates

Quick Facts

Storage upon receipt:

- $\leq -20^\circ$
- Protect from light

Ex/Em of Conjugates: See Table 1

Note: Do not vortex or shake solutions containing DNase I.

Introduction

Actin is one of the most ubiquitous and evolutionarily conserved structural proteins of muscle and nonmuscle cells and is generally present in two forms: G-actin (globular) and F-actin (fibrous). The transformation between these two forms relates to the need of the cell for structural rigidity or motion. By using selective fluorescent probes with very high affinity for G- or F-actin, one can obtain information on the localization of G- and F-actin in fixed cells.

DNase I has been used to detect and measure unpolymerized actin;^{1–5} it binds to G-actin with an affinity of about $5 \times 10^8 \text{ M}^{-1}$. Molecular Probes offers a broad selection of fluorescent conjugates of bovine pancreatic DNase I (Table 1) for selectively labeling monomeric G-actin. Phallotoxin conjugates, also available from Molecular Probes, are typically used for the staining of F-actin.

The influence of cytochalasins on actin structure in monocytes has been quantitated by flow cytometry using Texas Red® DNase I (D972) and BODIPY® FL phalloidin (B607) to stain the G-actin and F-actin pools, respectively.⁶ Fluorescent DNase I has also been used in a model system to study the interactions of nucleotides, cations, and cytochalasin D with monomeric actin.⁷

Table 1. Fluorescent DNase I conjugates and their spectral characteristics.

Cat #	Conjugate	Abs *	Em *
D12371	Alexa Fluor® 488	495	519
D12372	Alexa Fluor® 594	590	617
D972	Texas Red®	597	615

* Approximate fluorescence excitation (Ex) and emission (Em) maxima, in nm. Complete spectra for these dyes are available at our website, probes.invitrogen.com.

Materials

DNase I conjugates are supplied lyophilized in unit sizes of 5 mg. For long-term storage, reconstitute at 5 mg/mL (~161 μM) in one of the following buffers:

- phosphate-buffered saline, pH 7.4 (PBS) + 50% v/v glycerol
- 20 mM Tris, pH 7.6 + 50 mM NaCl + 1 mM DTT + 0.1 mg/mL BSA + 50% v/v glycerol

Glycerol at 50% in these solutions will maintain the liquid phase at $\leq -20^\circ\text{C}$ without affecting the stability of the enzyme conjugate. To avoid nonspecific binding, use only freshly prepared solutions or solutions stored at $\leq -20^\circ\text{C}$ in glycerol. Pancreatic DNase I is *unusually sensitive* to physical denaturation by shaking. Mixing should be done by gentle inversion. The lyophilized enzyme is stable for at least one year when stored as recommended, protected from light. Reconstituted enzyme can also be stored at 2–6°C, but activity may decrease when stored longer than one day at this temperature.

Properties

The molecular weight of the DNase conjugates is about 31,000 daltons. Peak excitation and emission wavelengths are shown in Table 1.

Staining G-Actin with Fluorescent DNase

1. Grow adherent cells on the surface of a coverslip.
2. Wash cells with PBS to remove excess medium.
3. Fix cells in 3.7% formaldehyde in PBS at room temperature for 10–15 minutes. Methanol disrupts the cytoskeletal structure or dye binding, which results in the absence of filament staining.
4. Wash three times with PBS.
5. Permeabilize cells in cold 100% acetone at -20°C for 5 minutes. Permeabilization can also be achieved by incubation with 0.1% Triton X-100 in PBS for 5 minutes.
6. Air-dry the samples or immediately rehydrate in PBS for 5–10 minutes.
7. Wash three times with PBS.

8. To the coverslip, add 200 μ L of a 9 μ g/mL (or 0.3 μ M) solution of fluorescent DNase I in buffer. *Optional:* To simultaneously label F-actin, also add 1 unit (200 μ L of a 0.165 μ M solution) of a fluorescent phallotoxin. Stain for 15–20 minutes.

9. Wash three times with PBS.

10. Air-dry the samples. *Optional:* Add a drop of *SlowFade*[®] Gold or ProLong[®] Gold antifade reagent to reduce photobleaching. Mount and seal.

References

1. J Cell Sci 66, 39 (1984); 2. Anal Biochem 135, 22 (1983); 3. Exp Cell Res 147, 240 (1983); 4. Eur J Biochem 104, 367 (1980); 5. J Biol Chem 255, 5668 (1980); 6. J Biol Chem 269, 3159 (1994); 7. Eur J Biochem 182, 267 (1989).

Product List *Current prices may be obtained from our Web site or from our Customer Service Department.*

Cat #	Product Name	Unit Size
D12371	deoxyribonuclease I, Alexa Fluor [®] 488 conjugate	5 mg
D12372	deoxyribonuclease I, Alexa Fluor [®] 594 conjugate	5 mg
D972	deoxyribonuclease I, Texas Red [®] conjugate.....	5 mg
P36930	ProLong [®] Gold antifade reagent.....	10 mL
S36930	<i>SlowFade</i> [®] Gold antifade reagent.....	10 mL

Contact Information

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Please visit our website — probes.invitrogen.com — for the most up-to-date information.

Molecular Probes, Inc.

29851 Willow Creek Road, Eugene, OR 97402
Phone: (541) 465-8300 • Fax: (541) 335-0504

Customer Service: 6:00 am to 4:30 pm (Pacific Time)

Phone: (541) 335-0338 • Fax: (541) 335-0305 • probesorder@invitrogen.com

Toll-Free Ordering for USA:

Order Phone: (800) 438-2209 • Order Fax: (800) 438-0228

Technical Service: 8:00 am to 4:00 pm (Pacific Time)

Phone: (541) 335-0353 • Toll-Free (800) 438-2209

Fax: (541) 335-0238 • probestech@invitrogen.com

Invitrogen European Headquarters

Invitrogen, Ltd.
3 Fountain Drive
Inchinnan Business Park
Paisley PA4 9RF, UK
Phone: +44 (0) 141 814 6100 • Fax: +44 (0) 141 814 6260
Email: euroinfo@invitrogen.com
Technical Services: eurotech@invitrogen.com

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Limited Use Label License

For research use only. Not intended for any animal or human therapeutic or diagnostic use. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402. Tel: (541) 465-8300. Fax: (541) 335-0504.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2005, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.