

Dynabeads® Biotin Binder

For research use only

1. PRODUCT DESCRIPTION

- 1.1 Intended Use
- 1.2 Principle of Isolation
- 1.3 Description of Materials

2. PROTOCOLS

- 2.1 Dynabeads Washing Procedure
- 2.2 Sample Preparation
- 2.3 Critical Steps for Cell Isolation
- 2.4 Cell Isolation - Indirect Technique
- 2.5 Cell Isolation - Direct Technique

3. TECHNICAL ADVICE

- 3.1 Indirect versus Direct Technique
- 3.2 Antibody Selection
- 3.3 Labelling Cells with Biotinylated Antibodies
- 3.4 Isolation and Depletion of Target Cells

4. GENERAL INFORMATION

5. REFERENCES

1. PRODUCT DESCRIPTION

1.1 Intended Use

Dynabeads Biotin Binder in combination with biotinylated antibodies are ideal for depletion or positive isolation of cells from any species depending




Fig. 1: Cell isolation using indirect or direct technique. The primary antibody is biotinylated prior to use.

on the specificity of the antibody. Other biotinylated molecules (e.g. peptides/proteins, lectins or nucleic acids) may also be used depending on the target. Cells can be directly isolated from any sample, such as whole blood, bone marrow, mononuclear suspensions or tissue digests.

Note:

If you want to perform positive cell isolation with biotinylated antibodies followed by bead detachment, please use CELLection™ Biotin Binder (Cat. no. 161.02).

If you want to use biotinylated ligands for molecular applications such as protein purifications, DNA/RNA binding protein isolations, preparing single-stranded templates etc., please see www.invitrogen.com for other streptavidin coated beads.

1.2 Principle of Isolation

Either add the biotinylated antibody to the cell sample (indirect technique) or pre-coat onto the beads (direct technique) prior to cell isolation. Mix the Dynabeads with the cell sample in a tube. The Dynabeads will bind to the target cells during a short incubation, and then separate the bead-bound cells by a magnet (see fig.1).

- Positive Isolation - discard the supernatant and use the bead-bound cells for downstream applications.
- Depletion - discard the bead-bound cells and use the remaining, untouched cells for any application.

1.3 Description of Materials

Dynabeads Biotin Binder are uniform, superparamagnetic polystyrene beads (2.8 µm diameter) coated with recombinant streptavidin. The streptavidin coated onto Dynabeads will bind most biotinylated ligands. Unwanted binding of cells to streptavidin via lectin-like receptors or other adhesive receptors is avoided since the recombinant streptavidin contains neither sugar nor the RYD-sequence.

Materials Supplied

5 ml Dynabeads Biotin Binder

4 x 10⁸ beads/ml in phosphate buffered saline (PBS), pH 7.4, containing 0.1% BSA and 0.02% sodium azide (NaN₃).

This product will process up to 2 x 10⁹ cells

Additional Materials Required

- Biotinylated ligand (e.g. antibodies, peptides/proteins, lectins or nucleic acids).
- Magnet (Dynal MPC™): See www.invitrogen.com/magnets-selection for magnet recommendations.
- Mixer allowing both tilting and rotation.
- Buffer 1: PBS w/0.1% BSA and 2 mM EDTA, pH 7.4.

Important Notes:

Serum or serum albumin fractions may contain free biotin and should not be used. Dynal recommends Bovine Serum Albumin fraction V from Sigma (cat. no. A4503 or equivalent), which does not contain free biotin.

Do not use culture media as it may contain free biotin.

BSA can be replaced by other biotin-free blocking proteins or 1 mg/ml pluronic (F68).

EDTA can be replaced by sodium citrate.

PBS containing Ca²⁺ or Mg²⁺ is not recommended.

2. PROTOCOLS

2.1 Dynabeads Washing Procedure

Dynabeads should be washed before use.

1. Resuspend the Dynabeads in the vial.
2. Transfer the desired volume of Dynabeads to a tube.
3. Add the same volume of Buffer 1, or at least 1 ml, and mix.
4. Place the tube in a magnet for 2 min and discard the supernatant.
5. Remove the tube from the magnet and resuspend the washed Dynabeads in the same volume of Buffer 1 as the initial volume of Dynabeads (step 2).

2.2 Sample Preparation

2.2.1 Whole Blood and Buffy Coat

Whole blood and buffy coat may be washed before use to remove interfering factors.

1. Dilute the whole blood or buffy coat in Buffer 1 (1+2).
2. Centrifuge at 600 x g for 10 min at RT (18-25°C).
3. Discard the plasma fraction/upper layer.
4. Resuspend to the original volume in Buffer 1 (2-8°C).

2.2.2 MNC Preparation from Whole Blood, Buffy Coat or Bone Marrow

1. Resuspend at 1 x 10⁷ cells per ml in Buffer 1 (2-8°C).

Please contact technical centre for a list of recommended human and mouse sample preparation procedures.

2.3 Critical Steps for Cell Isolation

- Use a mixer that provides tilting and rotation of the tubes to ensure Dynabeads do not settle at the bottom of the tube.
- When incubating Dynabeads and cells, the incubation temperature must be 2-8°C to reduce phagocytic activity and other metabolic processes.
- Never use less than 25 µl (1 x 10⁷) Dynabeads per ml cell sample and at least 4 Dynabeads per target cell.

Table 1: Volume of Dynabeads added per ml of cell sample. The volumes can be scaled up as required.

	Positive isolation	Depletion
Sample volume (1 x 10 ⁷ cells/ml*)	1 ml Max 2.5 x 10 ⁶ target cells	1 ml Max 2.5 x 10 ⁶ target cells
Volume of Dynabeads	25 µl	50 µl

* If the concentration of cells is increased or the target cell concentration exceeds 2.5 x 10⁶, the Dynabeads volume must be increased accordingly. Cell concentration can be up to 1 x 10⁸ cells per ml.

2.4 Cell Isolation - Indirect Technique

2.4.1 Labelling Cells with Biotinylated Antibodies or Ligands

Label cells with biotinylated antibody as recommended by the manufacturer. If necessary, optimize incubation time, antibody titer and cell concentration for best signal to noise ratio.

Use approximately 1 µg biotinylated antibody per 10⁶ target cells (if no other recommendation is stated) and ≥ 1 x 10⁷ cells/ml.

1. Add primary antibody to the cell suspension and mix.
2. Incubate for 10 min at 2-8°C.
3. Wash the cells by adding 2 ml Buffer 1 per 1 x 10⁷ cells and centrifuge at 300 x g for 8 min. Discard the supernatant.
4. Resuspend the cells in Buffer 1 at 1 x 10⁷ cells per ml.
5. Proceed to Isolation or Depletion of Cells (2.4.2).

2.4.2. Isolation or Depletion of Cells

1. Add Dynabeads to the prepared sample according to table 1.
2. Incubate for 20 min (positive isolation) or 30 min (depletion) at 2-8°C with gentle tilting and rotation.
3. Optional: double the volume with Buffer 1 to limit trapping of unbound cells.
4. Place the tube in a magnet for 3 min.
5. Depletion: Transfer the supernatant containing the unbound cells to a fresh tube for further experiments.
6. Positive isolation: Discard the supernatant and gently wash the bead-bound cells 4 times, using the following procedure:
 - i) Add 1 ml Buffer 1 per 1 x 10⁷ Dynabeads.
 - ii) Place the tube in the magnet for 2 min and discard the supernatant.
7. Resuspend the cells in buffer/medium for downstream application.

2.5 Cell Isolation - Direct Technique

2.5.1 Pre-coating Dynabeads

- Use 0.5-1.5 µg biotinylated antibody per 25 µl (1 x 10⁷) Dynabeads. It is recommended to titrate the amount of antibody.
1. Transfer washed Dynabeads to a tube.
 2. Add biotinylated antibodies.
 3. Incubate for > 30 min with gentle tilting and rotation.
 4. Place the tube in a magnet for 2 min and discard the supernatant.
 5. Wash the beads twice using 2 ml of Buffer 1.
 6. Remove the tube from the magnet and resuspend the washed Dynabeads in the same volume of Buffer 1 as the initial volume of Dynabeads.

2.5.2 Isolation or Depletion of Target Cells

1. Add the pre-coated Dynabeads to the cells according to table 1.
2. Incubate for 20 min (positive isolation) or 30 min (depletion) at 2-8°C with gentle tilting and rotation.
3. Optional: double the volume with Buffer 1 to limit trapping of unbound cells.
4. Place the tube in a magnet for 3 min.
5. Depletion: Transfer the supernatant containing the unbound cells to a fresh tube for further experiments.
6. Positive isolation: Discard the supernatant and gently wash the bead-bound cells 4 times, using the following procedure:
 - i). Add 1 ml Buffer 1 per 1 x 10⁷ Dynabeads.
 - ii). Place the tube in the magnet for 2 min and discard the supernatant.
7. Resuspend the cells in buffer/medium for downstream application.

3. TECHNICAL ADVICE

3.1 Indirect versus Direct Technique

Use the indirect technique when:

- A cocktail of biotinylated antibodies is used.
- Very high depletion efficiency is needed.
- The affinities of biotinylated monoclonal antibodies are low.

- The cells express low number of target antigens.
- The direct technique gives unsatisfactory purity.

The direct technique may be used if:

- The affinity of the biotinylated antibody is high.
- The cells express a high number of target antigens.
- A stock preparation of pre-coated Dynabeads is desired.

3.2 Antibody Selection

The choice of biotinylated antibody is the most important factor for successful cell isolation. Note that some antibodies may show reduced antigen-binding efficiency when coated onto beads (direct technique), even though the antibody shows good results in other immunological assays.

Some antibodies may be delivered in buffers containing high levels of free biotin. To be able to pre-coat the Dynabeads (direct technique), the free biotin needs to be removed (e.g. dialysis or spin column with MW 10, 000 cut-off).

3.3 Labelling Cells with Biotinylated Antibodies

- To avoid non-specific binding of cells (e.g. monocytes, B cells), add aggregated IgG to block Fc receptors prior to adding the biotinylated antibody.
- Excess antibody must be removed by washing before cell isolation.

3.4 Isolation and Depletion of Target Cells

- Remove density gradient media (e.g. Ficoll): Wash cells prior to adding biotinylated antibodies or Dynabeads.
- Remove soluble factors in serum: Serum may contain soluble factors (e.g. antibodies or cell surface antigens), which can interfere with the cell isolation protocol. Washing the cells once may reduce this interference.

4. GENERAL INFORMATION

Invitrogen Dynal AS complies with the Quality System Standards ISO 9001:2000 and ISO 13485:2003.

Storage/Stability

This product is stable until the expiry date stated on the label when stored unopened at 2-8°C.

Store opened vials at 2-8°C and avoid bacterial contamination.

Keep Dynabeads in liquid suspension during storage and all handling steps, as drying will result in reduced performance. Resuspend well before use.

Technical Support

Please contact Invitrogen Dynal for further technical information (see contact details).

Warning and Limitations

This product is for research use only. Not intended for any animal or human therapeutic or diagnostic use unless otherwise stated.

Follow appropriate laboratory guidelines.

This product contains 0.02% sodium azide as a preservative, which is cytotoxic. Avoid pipetting by mouth! Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. When disposing through plumbing drains, flush with large volumes of water to prevent azide build up.

Certificate of Analysis (CoA) is available upon request.

Material Safety Data Sheet (MSDS) is available at 'www.invitrogen.com'.

Patents and Trademarks

Dynal[®], Dynabeads[®] and Dynal MPC[™] are either registered trademarks or trademarks of Invitrogen Dynal AS, Oslo, Norway. Any registration or trademark symbols used herein denote the registration status of trademarks in the United States. Trademarks may or may not be registered in other countries.

CELLlection[™] is a trademark of Invitrogen Dynal AS.

Intellectual Property Disclaimer

Invitrogen Dynal will not be responsible for violations or patent infringements that may occur with the use of our products.

Warranty

The products are warranted to the original purchaser only to conform to the quantity and contents stated on the vial and outer labels for the duration of the stated shelf life. Invitrogen Dynal's obligation and the purchaser's exclusive remedy under this warranty is limited either to replacement, at Invitrogen Dynal's expense, of any products which shall be defective in manufacture, and which shall be returned to Invitrogen Dynal, transportation pre-paid, or at Invitrogen Dynal's option, refund of the purchase price.

Claims for merchandise damaged in transit must be submitted to the carrier.

This warranty shall not apply to any products which shall have been altered outside Invitrogen Dynal, nor shall it apply to any products which have been subjected to misuse or mishandling. ALL OTHER WARRANTIES, EXPRESSED, IMPLIED OR STATUTORY, ARE HEREBY SPECIFICALLY EXCLUDED, INCLUDING BUT NOT LIMITED TO WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. Invitrogen Dynal's maximum liability is limited in all events to the price of the products sold by Invitrogen Dynal. IN NO EVENT SHALL INVITROGEN DYNAL BE LIABLE FOR ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES. Some states do not allow limits on warranties, or on remedies for breach in certain transactions. In such states, the limits set forth above may not apply.

Limited Use Label License

No. 5: Invitrogen Technology – The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Invitrogen Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its

components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact

Licensing Department,
Invitrogen Corporation,
1600 Faraday Avenue, Carlsbad,
California 92008.
Phone (760) 603-7200.
Fax (760) 602-6500.
Email: outlicensing@invitrogen.com

5. REFERENCES

1. Almeida ARM *et al* (2002) Homeostasis of peripheral CD4⁺ T cells: IL-2R α and IL-2 shape a population of regulatory cells that controls CD4⁺ T cell numbers. *J. Imm.* 169: 4850-4860.
2. Anjuere F *et al* (2003) Transcutaneous immunization with cholera toxin B subunit adjuvant suppresses IgE antibody responses via selective induction of Th1 immune responses. *J. Imm.* 170: 1586-1592.
3. Borrello I *et al* (2000) Sustaining the graft-versus-tumor effect through post-transplant immunization with granulocyte-macrophage colony-stimulating factor (GM-CSF)-producing tumor vaccines. *Blood* 95:3011-3019.
4. Callahan MK *et al* (2003) Phosphatidylserine on HIV envelope is a cofactor for infection of monocytic cells. *J. Imm.* 170: 4840-4845.
5. Cardona A *et al* (2002) γ/δ T cell-deficient mice exhibit reduced disease severity and decreased inflammatory response in the brain in murine neurocysticercosis. *J. Imm.* 169: 3163-3171.
6. Cardona A *et al* (2003) CC chemokines mediate leukocyte trafficking into the central nervous system during murine neurocysticercosis: Role of γ/δ T cells in amplification of the host immune response. *Inf. & Imm.* 71: 2634-2642.
7. Chao S *et al* (2004) Pre-clinical study: antioxidant levels and immunomodulatory effects of wolfberry juice and other juice mixtures in mice. *JANA* 7: 31-38.
8. Erlich S *et al* (1999) Fluorescence-based selection of gene-corrected hematopoietic stem and progenitor cells from acid sphingomyelinase-deficient mice: Implications for Niemann-Pick disease gene therapy and the development of improved stem cell gene transfer procedures. *Blood* 93: 80-86.
9. Faili A *et al* (2002) AID-dependent somatic hypermutation occurs as a DNA single-strand event in the BL2 cell line. *Nat. Imm.* 3:815-821.
10. Galkina E *et al* (2003) L-selectin shedding does not regulate constitutive T cell trafficking but controls the migration pathways of antigen-activated lymphocytes. *J. Exp. Med.* 198: 1323-1335.
11. Gaudin E *et al* (2004) Positive selection of B cells expressing low densities of self-reactive BCRs. *J. Exp. Med.* 199: 843-853.
12. Halin C *et al* (2005) The S1P-analog FTY720 differentially modulates T cell homing via HEV: T cell-expressed S1P1 amplifies integrin activation in peripheral lymph nodes but not in Peyer's patches. *Blood* 106: 1314-1322.
13. Hamad ARA *et al* (2001) Lack of coreceptor allows survival of chronically stimulated double-negative α/β T Cells: Implications for autoimmunity. *J. Exp. Med.* 193: 1113-1121.
14. Huang Y *et al* (2004) Matching at the MHC class I K locus is essential for long-term engraftment of purified hematopoietic stem cells: a role for host NK cells in regulating HSC engraftment. *Blood* 104: 873-880.
15. Kao KJ *et al* (2001) Characterization of immunologic tolerance induced by transfusion of UV-

B-irradiated allogeneic mononuclear leukocytes. *Blood* 98:1239-1245.

16. Krajewski S *et al* (1997) Immunolocalization of the ICE/Ced-3 – family protease, CPP32 (caspase-3), in non-Hodgkin's lymphomas, chronic lymphocytic leukemias, and reactive lymph nodes. *Blood* 89: 3817-3825.
17. Landais E *et al* (2004) Direct killing of Epstein-Barr virus (EBV)-infected B cells by CD4 T cells directed against the EBV lytic protein BHRF1. *Blood* 103:1408-1416.
18. Meissner U *et al* (2001) A soluble form of the murine common γ chain is present at high concentrations *in vivo* and suppresses cytokine signaling. *Blood* 97: 183-191.
19. Naper C *et al* (2005) Two structurally related rat Ly49 receptors with opposing functions (Ly49 stimulatory receptor 5 and Ly49 inhibitory receptor 5) recognize non-classical MHC class Ib-encoded target ligands. *J. Imm.* 174: 2702-2711.
20. Peichev M *et al* (2000) Expression of VEGFR-2 and AC133 by circulating human CD341 cells identifies a population of functional endothelial precursors. *Blood* 95:952-958.
21. Otsu M *et al* (2002) Reconstitution of lymphoid development and function in ZAP-70-deficient mice following gene transfer into bone marrow cells. *Blood* 100: 1248-1256.
22. Sevilla N *et al* (2000) Immunosuppression and resultant viral persistence by specific viral targeting of dendritic cells. *J. Exp. Med.* 192: 1249-1260.
23. Zhang TY *et al* (2005) The expression of 11 β -hydroxysteroid dehydrogenase type I by lymphocytes provides a novel means for intracrine regulation of glucocorticoid activities. *J. Imm.* 174: 879-889.

Invitrogen Dynal is a part of the Invitrogen Group.

Contact details for your local Invitrogen sales office/technical support can be found at <http://www.invitrogen.com/contact>

© Copyright 2007 Invitrogen Dynal AS, Oslo, Norway. All rights reserved.

SPEC-05724