

# One Shot® INV110 Competent Cells

Catalog no. C7171-03

Revision date: 20 August 2009

Manual part no. 25-0269

MAN0000103



**Important Kit Change!**

This kit now includes  
**Recovery Medium** instead of  
SOC Medium. See the note on  
page 4 for details.



# Table of Contents

Kit Contents and Storage .....4

Description of the Cells .....5

Transformation Protocol.....7

Technical Support .....9

# Kit Contents and Storage

---

## Shipping and Storage

One Shot® INV110 Competent Cells are shipped on dry ice. Upon receipt, remove the Recovery Medium and store at room temperature or at +4°C. Store the remainder of the kit at -80°C.

---

## Kit Contents

Each kit contains the following:

- 21 tubes of chemically competent INV110 *E. coli* in 50 µl aliquots (transformation efficiency =  $1 \times 10^6$  cfu/µg supercoiled DNA)
  - Supercoiled pUC19 plasmid DNA (100 pg/µl in 5 mM Tris-HCl, 0.5 mM EDTA, pH 8; 20 µl) for testing transformation efficiency
  - Recovery Medium (10 ml) for plating (see **Important Note** below)
- 



## Important

Effective 20 August 2009, this kit now includes **Recovery Medium** instead of SOC Medium. Recovery Medium is a slightly richer formulation than SOC for cultivating recombinant strains of *E. coli*, and provides better lot-to-lot consistency. This change has been validated and does not affect the product's protocol or specifications.

Please continue to perform your transformation procedure as before, simply substituting Recovery Medium for SOC Medium. **All volumes, protocols, and handling procedures are the same for the new medium.**

---

## Product Qualification

The Certificate of Analysis provides detailed quality control information for each product. Certificates of Analysis are available on our website at [www.invitrogen.com/support](http://www.invitrogen.com/support). Search for the Certificate of Analysis by product lot number, which is printed on the box.

---

## Information for European Customers

The INV110 *E. coli* strain is genetically modified and carries Tn10 (Tet<sup>R</sup>) to allow tetracycline selection of the host strain. As a condition of sale, this product must be in accordance with all applicable local legislation and guidelines including EC Directive 90/219/EEC on the contained use of genetically modified organisms.

---

# Description of the Cells

---

## Features of INV110

The INV110 *E. coli* strain contains the following features:

- *dam* and *dcm* deficiencies to allow restriction digestion with *dam*- and *dcm*-sensitive restriction enzymes (see below for more information)
  - The  $\Delta(mcrC-mrr)$  allele eliminates two restriction systems to allow more efficient transformation of DNA from highly methylated sources (i.e. eukaryotic DNA)
  - The *endA1* mutation to permit isolation of higher quality plasmid DNA
  - The *lacI<sup>q</sup>* allele for high expression of the Lac repressor
  - *Tn10* to permit selection of the host strain using tetracycline
- 

## Genotype

F' {*tra* $\Delta$ 36 *proAB lacI<sup>q</sup> lacZ* $\Delta$ M15} *rpsL* (Str<sup>R</sup>) *thr leu endA thi-1 lacY galK galT ara tonA tsx dam dcm supE44*  $\Delta(lac-proAB)$   $\Delta(mcrC-mrr)$ 102::*Tn10* (Tet<sup>R</sup>)

---

## *dam* and *dcm* Methylation

INV110 is a *dam* and *dcm* methylase-deficient *E. coli* strain. In *dam*<sup>+</sup> *dcm*<sup>+</sup> *E. coli* strains, the Dam and Dcm methylases methylate DNA at the following sites:

Methylase	Methylation Site
Dam	G <sup>m</sup> ATC
Dcm	C <sup>m</sup> C(A/T)GG

---

The *dam* and *dcm* deficiencies allow production of DNA that is unmethylated at these sites. Please refer to the catalog of your restriction enzyme supplier for information on methylation-sensitive restriction enzymes.

---

*Continued on the next page*

## Description of the Cells, Continued

---



The presence of the *dam* deficiency affects INV110 in a number of ways. Please be aware that INV110 cells have the following characteristics when compared to *dam*<sup>+</sup> *E. coli* strains:

- A higher mutation rate (greater than 250-fold)
- A higher recombination frequency and lower viability on plates
- A lower transformation efficiency

If you are isolating DNA from transformants, we recommend that you pick colonies and isolate DNA within one day after transformation. **INV110 should not be used to maintain your plasmid of interest for long-term storage.**

---

### General Handling

Be extremely gentle when working with competent cells. Competent cells are highly sensitive to changes in temperature or mechanical lysis caused by pipetting. Transformation should be started immediately following the thawing of the cells on ice. Mix reagents by swirling or tapping the tube gently, not by pipetting.

---

### Blue-White Screening

If blue-white screening is required to select for transformants, your LB agar plates must contain X-Gal and isopropyl -D-thiogalactoside (IPTG). INV110 One Shot<sup>®</sup> cells express the Lac repressor and require IPTG to induce expression from the *lac* promoter.

To perform blue-white screening, prepare the following stock solutions of X-Gal and IPTG:

- 40 mg/ml X-Gal in dimethylformamide (DMF)
- 100 mM IPTG in water (filter sterilize)

You may add X-Gal and IPTG to LB agar plates using one of the following methods:

1. Spread 40  $\mu$ l of 40 mg/ml X-Gal and 40  $\mu$ l of 100 mM IPTG on top of the agar. Let the X-Gal and IPTG diffuse into the agar for approximately 1 hour.
  2. Add X-Gal and IPTG directly to the autoclaved agar prior to pouring your plates. When the autoclaved agar has cooled to 55°C, add X-Gal to a final concentration of 20  $\mu$ g/ml and IPTG to a final concentration of 0.1 mM.
-

# Transformation Protocol

---

## Materials Supplied by the User

You will need the following items for transformation:

- 37°C shaking and non-shaking incubator
  - 10 cm diameter LB agar plates containing 10 µg/ml tetracycline, the appropriate antibiotic for selection of your plasmid of interest (i.e. 50–100 µg/ml ampicillin), and X-Gal and IPTG, if desired
  - Ice bucket with ice
  - 42°C water bath
  - A test tube rack to hold all transformation tubes so that they can all be put into a 42°C water bath simultaneously
- 

## Before Starting

- Equilibrate a water bath to 42°C
  - Warm the vial of Recovery Medium to room temperature. *Note:* If you are using an older kit that contains SOC Medium, simply follow the protocol in this manual, substituting the SOC Medium for the Recovery Medium specified here.
  - Spread X-Gal and IPTG onto LB agar plates containing tetracycline and the appropriate antibiotic, if desired
  - Warm the plates in a 37°C incubator for 1 hour (use one plate for each transformation)
- 

## Transformation Procedure

Follow the instructions provided below to transform your plasmid of interest into INV110.

1. Centrifuge the vial(s) containing the plasmid DNA briefly and place on ice.
2. Thaw, on ice, one 50 µl vial of One Shot<sup>®</sup> cells for each transformation.
3. Pipet 1–5 µl (10–100 ng) of each DNA sample directly into the competent cells and mix by tapping gently. **Do not mix cells by pipetting.** The remaining DNA sample(s) can be stored at -20°C.
4. Incubate the vial(s) for 30 minutes on ice.

*Procedure continued on next page*

---

*Continued on next page*

# Transformation Protocol, Continued

---

## Transformation Procedure, continued

*Procedure continued from previous page*

5. Transfer all vials at the same time (i.e., in a rack) to the 42°C water bath and incubate for exactly 30 seconds. Do not mix or shake.
  6. Remove vial(s) from the 42°C bath and place on ice.
  7. Add **250 µl** of pre-warmed Recovery Medium to each vial. **Important:** Use sterile technique when handling the medium to avoid contamination.
  8. Place the vial(s) in a microcentrifuge rack on its side and secure with tape to avoid loss of the vial(s). Shake the vial(s) at 37°C for 1 hour at 225 rpm in a shaking incubator.
  9. Spread 20 µl to 200 µl from each transformation vial on a separate, labeled LB agar plate containing X-Gal and IPTG (if desired), 10 µg/ml tetracycline, and the appropriate antibiotic for selection of the plasmid of interest.
  10. Invert the plate(s) and incubate at 37°C overnight.
  11. Select colonies and analyze by restriction digest, PCR, or sequencing.
- 



We recommend that you test the efficiency of the competent cells by using the supercoiled pUC19 plasmid supplied with the kit as described below.

- Prepare LB agar plates containing 10 µg/ml tetracycline and 50–100 µg/ml ampicillin.
  - Transform 1 µl (100 pg) into 50 µl of competent cells according to the transformation procedure above.
  - Incubate overnight at 37°C and count colonies. Calculate the transformation efficiency as transformants per 1 µg of plasmid using the formula below. The cells should have an efficiency of at least  $1 \times 10^6$  transformants/µg of supercoiled plasmid.
- 

## Calculation

Use the formula below to calculate transformation efficiency.

$$\frac{\text{\# of colonies}}{100 \text{ pg transformed DNA}} \times \frac{10^6 \text{ pg}}{\mu\text{g}} \times \frac{300 \mu\text{l total transformation volume}}{X \mu\text{l plated}} = \frac{\text{\# transformants}}{\mu\text{g plasmid DNA}}$$


---



# Technical Support

---

## World Wide Web



Visit the Invitrogen website at [www.invitrogen.com](http://www.invitrogen.com) for:

- Technical resources, including manuals, vector maps and sequences, application notes, MSDSs, FAQs, formulations, citations, handbooks, etc.
  - Complete technical support contact information
  - Access to the Invitrogen Online Catalog
  - Additional product information and special offers
- 

## Contact Us

For more information or technical assistance, call, write, fax, or email. Additional international offices are listed on our website ([www.invitrogen.com](http://www.invitrogen.com)).

### Corporate Headquarters:

5791 Van Allen Way  
Carlsbad, CA 92008 USA  
Tel: 1 760 603 7200  
Tel (Toll Free): 1 800 955 6288  
Fax: 1 760 602 6500  
E-mail: [tech\\_support@invitrogen.com](mailto:tech_support@invitrogen.com)

### European Headquarters:

Inchinnan Business Park  
3 Fountain Drive  
Paisley PA4 9RF, UK  
Tel: +44 (0) 141 814 6100  
Tech Fax: +44 (0) 141 814 6117  
E-mail: [eurotech@invitrogen.com](mailto:eurotech@invitrogen.com)

---

## MSDS

Material Safety Data Sheets (MSDSs) are available on our website at [www.invitrogen.com/msds](http://www.invitrogen.com/msds).

---

## Certificate of Analysis

The Certificate of Analysis provides detailed quality control and product qualification information for each product. Certificates of Analysis are available on our website. Go to [www.invitrogen.com/support](http://www.invitrogen.com/support) and search for the Certificate of Analysis by product lot number, which is printed on the box.

---

*Continued on next page*

# Technical Support, continued

---

## Limited Warranty

Invitrogen (a part of Life Technologies Corporation) is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, contact our Technical Support Representatives.

All Invitrogen products are warranted to perform according to specifications stated on the certificate of analysis. The Company will replace, free of charge, any product that does not meet those specifications. This warranty limits the Company's liability to only the price of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. The Company reserves the right to select the method(s) used to analyze a product unless the Company agrees to a specified method in writing prior to acceptance of the order.

Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore the Company makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, report it to our Technical Support Representatives.

**Life Technologies Corporation shall have no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.**

---

## Trademarks

Agilent® is a registered trademark of Agilent Technologies, Inc.

---

©2009 Life Technologies Corporation. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.





**Corporate Headquarters**

5791 Van Allen Way  
Carlsbad, CA 92008

T: 1 760 603 7200

F: 1 760 602 6500

E: [tech\\_support@invitrogen.com](mailto:tech_support@invitrogen.com)

For country-specific contact information, visit our web site at [www.invitrogen.com](http://www.invitrogen.com)