

AlignFlow™ and AlignFlow™ Plus Flow Cytometry Alignment Beads

Table 1. Contents and storage information.

Material	Amount	Concentration	Storage	Stability
AlignFlow™ * and AlignFlow™ Plus † beads	3 mL dropper vials	Suspensions in water containing 0.05% Tween 20 and 2 mM sodium azide	<ul style="list-style-type: none"> • 2–6°C • Protect from light • Sonicate before use • Do not freeze 	When stored as directed, reagents are stable for at least 1 year. ‡

* The AlignFlow™ suspensions contain $\sim 1.2 \times 10^8$ beads/mL (0.1% solids). † The AlignFlow™ Plus suspensions contain $\sim 1.7 \times 10^7$ beads/mL (0.2% solids). ‡ No leaching or degradation is expected.

Approximate fluorescence excitation and emission maxima: See Table 2.

Introduction

In order to ensure accurate and reproducible results, flow cytometers should be checked daily for proper performance. Reference standards are indispensable for conducting these tests. Molecular Probes AlignFlow™ and AlignFlow™ Plus flow cytometry alignment beads are reliable references for aligning, focusing, and calibrating flow cytometers. These fluorescently stained, polystyrene microspheres are highly uniform with respect to size and fluorescence intensity. The 2.5 μm -diameter AlignFlow™ beads and the 6.0 μm -diameter AlignFlow™ Plus beads are each available in four versions: for UV (350–370 nm) excitation, for 488 nm excitation, for 633 nm excitation, and for 630–660 nm excitation (Table 2). The fluorescent dyes have been carefully selected for optimal excitation by laser sources commonly used in flow cytometry. The 488 nm-excitables emit broadly from 515 nm to 660 nm. The UV-excitables emit from 400 nm to 470 nm; the 633 nm-excitables, from 645 nm to 680 nm; and the 630–660 nm-excitables, from 670 nm to 720 nm. Because the dyes are contained inside each microsphere, instead of merely on the bead surface, AlignFlow™ beads exhibit superior signal stability.

Using AlignFlow™ or AlignFlow™ Plus beads, an operator can reproducibly adjust parameters crucial to flow cytometry before an experiment is started. The beads approximate the size, emission wavelength, and intensity of many biological samples and permit the calibration of the flow cytometer's laser source, optics, and stream flow without wasting valuable and sensitive experimental material. Because of their excellent stability, these beads can also serve as a daily reference standard by which an instrument's performance and reliability are evaluated.

Table 1. Molecular Probes AlignFlow™ and AlignFlow™ Plus flow cytometry alignment beads.

Cat #	Bead Diameter	Ex *	Em range †
A7304	2.5 µm	350–370	400–470
A7305	6.0 µm		
A7302	2.5 µm	488	515–660
A7303	6.0 µm		
A7312	2.5 µm	633	645–680
A7313	6.0 µm		
A14835	2.5 µm	630–660	670–720
A14836	6.0 µm		

* Fluorescence excitation, in nm, for which the beads were designed. † Approximate useful emission range, in nm.

Guidelines for Use

Molecular Probes AlignFlow™ and AlignFlow™ Plus beads serve as a reference standard for calibrating flow cytometers. Experimental protocols depend somewhat on the flow cytometer and software used; please refer also to the reference materials applicable to your particular instrument. Before sampling, be sure that the polystyrene beads are uniformly suspended by vortex mixing and sonicating the suspension. Generally, one drop added to 1 mL of Haema-Line 2 sheath fluid or buffered saline solution provides an appropriate concentration for analysis; mix well before applying the sample.

During the alignment, set the approximate photomultiplier tube (PMT) voltage to place the sample in a convenient range for the flow cytometer, then maximize the signal and minimize the coefficient of variation (CV) by making the appropriate optical adjustments. Follow the instrument manufacturer’s manual carefully for optimal performance.

When using the AlignFlow™ or AlignFlow™ Plus beads as a daily reference source, maximize the peak heights and minimize the CVs, as outlined above. Once the settings are optimized, record all settings and print out the relevant data for your record. Each day pre-run the flow cytometer using these reference beads and compare the output to the reference record. Deviations from the baseline readings may indicate instrument malfunction.

References

1. *Flow Cytometry, A Practical Approach*, 2nd Edition, M.G. Ormerod, Ed. IRL Press (1994);
2. Givan A.L., *Flow Cytometry, First Principles*, Wiley-Liss (1992).

Product List Current prices may be obtained from our website or from our Customer Service Department.

Cat #	Product Name	Unit Size
A7304	AlignFlow™ flow cytometry alignment beads, 2.5 µm *for UV excitation*	3 mL
A7302	AlignFlow™ flow cytometry alignment beads, 2.5 µm *for 488 nm excitation*	3 mL
A7312	AlignFlow™ flow cytometry alignment beads, 2.5 µm *for 633 nm excitation*	3 mL
A14835	AlignFlow™ flow cytometry alignment beads, 2.5 µm *for 630-660 nm excitation*	3 mL
A7305	AlignFlow™ Plus flow cytometry alignment beads, 6.0 µm *for UV excitation*	3 mL
A7303	AlignFlow™ Plus flow cytometry alignment beads, 6.0 µm *for 488 nm excitation*	3 mL
A7313	AlignFlow™ Plus flow cytometry alignment beads, 6.0 µm *for 633 nm excitation*	3 mL
A14836	AlignFlow™ Plus flow cytometry alignment beads, 6.0 µm *for 630-660 nm excitation*	3 mL

Contact Information

Molecular Probes, Inc.
29851 Willow Creek Road
Eugene, OR 97402
Phone: (541) 465-8300
Fax: (541) 335-0504

Customer Service:
6:00 am to 4:30 pm (Pacific Time)
Phone: (541) 335-0338
Fax: (541) 335-0305
probesorder@invitrogen.com

Toll-Free Ordering for USA:
Order Phone: (800) 438-2209
Order Fax: (800) 438-0228

Technical Service:
8:00 am to 4:00 pm (Pacific Time)
Phone: (541) 335-0353
Toll-Free (800) 438-2209
Fax: (541) 335-0238
probestech@invitrogen.com

Invitrogen European Headquarters
Invitrogen, Ltd.
3 Fountain Drive
Inchinnan Business Park
Paisley PA4 9RF, UK
Phone: +44 (0) 141 814 6100
Fax: +44 (0) 141 814 6260
Email: euroinfo@invitrogen.com
Technical Services: eurotech@invitrogen.com

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Limited Use Label License No. 223: Labeling and Detection Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2006, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.