

## PCR SuperMix

Cat. No. 10572-014

Cat. No. 10572-063

Size: 100 reactions

Size: 5000 reactions

Store at **-20°C** (in a non-frost-free freezer)

### Description

PCR SuperMix provides qualified reagents for the amplification of nucleic acid templates by the polymerase chain reaction (PCR). PCR SuperMix contains  $Mg^{++}$ , dNTPs, and recombinant *Taq* DNA Polymerase at concentrations sufficient to allow amplification during PCR. PCR SuperMix is supplied at 1.1X concentration to allow approximately 10% of the final reaction volume to be used for the addition of primer and template solutions. Reagents sufficient for 100 and 5000 amplification reactions of 50  $\mu$ l each are provided.

PCR SuperMix may be stored at either -20°C or 4°C. Storage at 4°C avoids the necessity of thawing the mix before assembling the PCR. No detectable reduction of PCR performance or enzyme activity is observed after storage of PCR SuperMix for 8 months at 4°C. Repeated freeze-thaw cycles do not reduce performance or activity.

<u>Component</u>	<u>100-rxn size</u>	<u>5000-rxn size</u>
PCR SuperMix	4 × 1.125 ml	4 × 56.25 ml

### Components

22 mM Tris-HCl (pH 8.4), 55 mM KCl, 1.65 mM  $MgCl_2$ , 220  $\mu$ M dGTP, 220  $\mu$ M dATP, 220  $\mu$ M dTTP, 220  $\mu$ M dCTP, 22 U recombinant *Taq* DNA Polymerase/ml, stabilizers.

### **Guidelines and Recommendations**

Since PCR is a powerful technique capable of amplifying trace amounts of DNA, all appropriate precautions should be taken to avoid cross-contamination. Ideally, amplification reactions should be assembled in a DNA-free environment. Use of aerosol-resistant barrier tips is recommended. Take care to avoid contamination of PCR SuperMix with the primers or template DNA used in individual reactions. PCR products should be analyzed in an area separate from the reaction assembly area.

A standard 50- $\mu$ l reaction uses 45  $\mu$ l of PCR SuperMix and 5  $\mu$ l of primer and template solution. We have observed no decrease in product yield if the amount of template and primer solution added is between a fraction of a microliter and 20  $\mu$ l. Lower yields occur as the  $Mg^{++}$  concentration drops to a suboptimal level. If the final  $Mg^{++}$  concentration is adjusted to 1.5 mM, the volume of primer and template solution that can be added to 45  $\mu$ l of PCR SuperMix can exceed 50  $\mu$ l.

### **Quality Control**

The Certificate of Analysis provides detailed quality control and product qualification information for each product. Certificates of Analysis are available on our website. Go to [www.invitrogen.com/support](http://www.invitrogen.com/support) and search for the Certificate of Analysis by product lot number, which is printed on the box.

**Protocol**

The following protocol is suggested as a starting point and guideline when using PCR SuperMix. We recommend assembling reactions on ice from pre-chilled components. This protocol is for a reaction size of approximately 50  $\mu$ l. The reaction size may be adjusted as desired.

**Note:** For multiple reactions with common components, prepare a master mix of the components common to all reactions to reduce pipetting errors.

1. Set up reaction tubes/plates on ice.
2. Add the following components in any order to each reaction vessel.
  - 45- $\mu$ l PCR SuperMix
  - Primers (200 nM final concentration per primer is recommended)\*
  - Template DNA solution\*

\*Total volume of primer and template solution can be 0.5-20  $\mu$ l.

3. Mix contents and cover with mineral or silicone oil, if necessary.
4. Cap reaction vessels and load in thermal cycler at 80°C.
5. Run cycling program.

**Limited License No. 1: Thermostable Polymerases**

Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 5,789,224, 5,618,711, and 6,127,155. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim, no right to perform any patented method, and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

**Limited License No. 5: Invitrogen Technology**

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For products that are subject to multiple limited use label licenses, the terms of the most restrictive limited use label license shall control. Life Technologies Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Life Technologies Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Life Technologies is willing to accept return of the product with a full refund. For information about purchasing a license to use this product or the technology embedded in it for any use other than for research use please contact Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008 or [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com).

©2009 Life Technologies Corporation. All rights reserved.  
For research use only. Not intended for any animal or human therapeutic or diagnostic use.