

5-ethynyl uridine (EU)

Catalog no. E10345

Table 1. Contents and storage information.

Material	Amount	Storage	Stability
5-ethynyl uridine (EU)	5 mg	<ul style="list-style-type: none"> • $\leq -20^{\circ}\text{C}$ • Desiccate 	When stored as directed, the product is stable for at least 1 year

Introduction

5-Ethynyl uridine (EU) is a novel alternative for bromo uridine (BrU or BrUTP)¹ to directly image spatially and temporally nascent global RNA transcription both *in vitro* and *in vivo*. EU is a nucleoside analog of uracil, and is incorporated into RNA during active RNA synthesis. Detection of the incorporated EU is based on a powerful click chemical reaction that utilizes bio-orthogonal or biologically unique moieties, specifically azides and alkynes.²⁻⁵

EU contains an alkyne that can react with an azide-containing detection reagent in the Click-iT® Cell Reaction Buffer Kit (Cat. no. C10269) to form a stable triazole ring (Figure 1, Table 2). Because click chemistry-based label and detection tags react selectively and specifically with one another, you can apply click-labeled molecules to complex biological samples and subsequently detect with unprecedented sensitivity due to an extremely low background.




Figure 1. Click reaction between EU and azide-modified dye or hapten.

Table 2. Azide-modified fluorophores and haptens for detection.

Label	Ex/Em*	Cat. no.	Use	Detection platform
Alexa Fluor® 488	495/519	A10266	Fluorescent dye or hapten	Flow, HCS, Fluorescence microscopy
Alexa Fluor® 594	590/617	A10270	Fluorescent dye	HCS, Fluorescence microscopy
Alexa Fluor® 647	650/655	A10277	Fluorescent dye	Flow, HCS, Fluorescence microscopy
Biotin**	NA	B10184	Hapten	Flow, HCS, Fluorescence microscopy,
Oregon Green® 488	496/524	O10180	Fluorescent dye or hapten	Flow, HCS, Fluorescence microscopy
Tetramethylrhodamine (TAMRA)	555/580	T10182	Fluorescent dye or hapten	HCS, Fluorescence microscopy

*Approximate fluorescence excitation (Ex) and emission (Em) maxima, in nm. **Requires streptavidin.

Click chemistry fills the void when methods such as direct labeling with biotin or a fluorophore (MW ~300–1,000 daltons), or the use of antibodies (MW ~150,000 daltons) are not applicable or efficient, particularly with small molecules and complex samples (Figure 2). The click chemistry label is small enough that the tagged molecules (*e.g.*, nucleotides, sugars, and amino acids) are acceptable substrates for the enzymes that assemble these building blocks into biopolymers.

It has been demonstrated that EU can be efficiently incorporated into RNA by several polymerases, but it does not incorporate into DNA.⁶ The click detection reaction requires copper as the catalyst, and copper has been shown to induce some RNA degradation (scission). Because of this, we do not recommend RNA sequencing following the click reaction at this time. However, there are other novel applications for EU. For example, the small size of the detection reagent enables the simultaneous detection of RNA-interactive proteins with antibodies, which cannot be accomplished if antibodies are used for the detection of BrU.




Figure 2. Relative size of detection molecules commonly used in cellular analysis.

Before You Begin

- Preparing Stock Solution** EU is readily soluble in DMSO, alcohol, water, or aqueous buffers. For use with *in vitro* applications, we recommend preparing a 0.1–1 M EU stock solution in DMSO or an aqueous buffer. You may store any unused EU stock solution at $\leq -20^{\circ}\text{C}$ for up to one year.
- Handling and Disposal** EU is a nucleoside analog which can be incorporated into RNA. Handle and dispose of EU in compliance with all pertaining local regulations. When EU is dissolved in DMSO, which is known to facilitate the entry of organic molecules into tissue, use additional precautions appropriate for the hazards posed by such materials.

Experimental Protocols

- EU Labeling** In initial experiments, we recommend testing a range of EU concentrations to determine the optimal concentration. If currently using a BrU-based assay, a similar concentration and duration to BrU is a good starting concentration for EU. The optimal concentration may vary depending upon the duration of the pulse, with lower concentrations recommended for longer incubations. With cultured cells, acceptable EU incorporation has been observed with 0.5–5 mM EU for 0.5–24 hours.
- Note:** Do not exceed DMSO concentrations of 0.5% when labeling live cells.
- EU Detection** You may detect incorporated EU with an available dye- or hapten-containing azide (Table 2) together with the Click-iT[®] Cell Reaction Buffer Kit (Cat. no. C10269). If signal amplification is desired, in addition to the classic hapten, biotin, you may use three of the fluorescent dyes, Alexa Fluor[®] 488, Oregon Green[®] 488, and tetramethylrhodamine (TAMRA), with their corresponding anti-dye antibodies.
- Imaging and Analysis** See Table 2 for the appropriate fluorescence excitation/emission maxima for detection of EU following the click reaction.

References

1. RNA 6, 1750 (2000); 2. ChemBioChem 4, 1147 (2003); 3. J Am Chem Soc 125, 3192 (2003); 4. Angew Chem Int Ed Engl 41, 2596 (2002); 5. Angew Chem Int Ed Engl 40, 2004 (2001); 6. PNAS 105, 15779 (2008).

Product List Current prices may be obtained from our website or from our Customer Service Department.

Cat. no.	Product Name	Unit Size
E10345	5-ethynyl uridine (EU)	5 mg
Related Products		
A889	anti-fluorescein/Oregon Green [®] , rabbit IgG fraction *1 mg/mL*	0.5 mL
A982	anti-fluorescein/Oregon Green [®] , rabbit IgG fraction, biotin-XX conjugate *1 mg/mL*	0.5 mL
A6397	anti-tetramethylrhodamine, rabbit IgG fraction *1 mg/mL*	0.5 mL
A6413	anti-fluorescein/Oregon Green [®] , rabbit IgG Fab fragment *0.5 mg/mL*	0.5 mL
A6421	anti-fluorescein/Oregon Green [®] , mouse IgG _{2a} , monoclonal 4-4-20	0.5 mg
A10266	Alexa Fluor [®] 488 azide (Alexa Fluor [®] 488 5-carboxamido-(6-azidohexanyl), bis(triethylammonium salt) *5-isomer*	0.5 mg
A10270	Alexa Fluor [®] 594 azide (Alexa Fluor [®] 594 carboxamido-(6-azidohexanyl), triethylammonium salt) *mixed isomers*	0.5 mg
A10277	Alexa Fluor [®] 647 azide, triethylammonium salt	0.5 mg
A11090	anti-fluorescein/Oregon Green [®] , rabbit IgG fraction, Alexa Fluor [®] 488 conjugate *1 mg/mL*	0.5 mL
A11091	anti-fluorescein/Oregon Green [®] , rabbit IgG fraction, Alexa Fluor [®] 594 conjugate *1 mg/mL*	0.5 mL
A11094	anti-Alexa Fluor [®] 488, rabbit IgG fraction *1 mg/mL*	0.5 mL
A11095	anti-fluorescein/Oregon Green [®] , goat IgG fraction *1 mg/mL*	0.5 mL
A11096	anti-fluorescein/Oregon Green [®] , goat IgG fraction, Alexa Fluor [®] 488 conjugate *1 mg/mL*	0.5 mL
A21250	anti-fluorescein/Oregon Green [®] , rabbit IgG fraction, R-phycoerythrin conjugate *2 mg/mL*	250 µL
A21253	anti-fluorescein/Oregon Green [®] , rabbit IgG fraction, horseradish peroxidase conjugate	0.5 mg
B10184	biotin azide (PEG ₄ carboxamide-6-azidohexanyl biotin)	1 mg
C10269	Click-iT [®] Cell Reaction Buffer Kit.	1 kit
O10180	Oregon Green [®] 488 azide (Oregon Green [®] 488 6-carboxamido-(6-azidohexanyl), triethylammonium salt) *6-isomer*	0.5 mg
Q10137MP	Qdot [®] 565 goat anti-fluorescein conjugate *2 µM* *whole IgG*	500 µL
Q15421MP	Qdot [®] 655 goat anti-fluorescein conjugate *1 µM* *whole IgG*	200 µL
Q15431MP	Qdot [®] 565 goat anti-fluorescein conjugate *1 µM* *whole IgG*	200 µL
T10182	tetramethylrhodamine (TAMRA) azide (tetramethylrhodamine 5-carboxamido-(6-azidohexanyl)) *5-isomer*	0.5 mg

Contact Information

Molecular Probes, Inc.

29851 Willow Creek Road
Eugene, OR 97402
Phone: (541) 465-8300
Fax: (541) 335-0504

Customer Service:

6:00 am to 4:30 pm (Pacific Time)
Phone: (541) 335-0338
Fax: (541) 335-0305
probesorder@invitrogen.com

Toll-Free Ordering for USA:

Order Phone: (800) 438-2209
Order Fax: (800) 438-0228

Technical Service:

8:00 am to 4:00 pm (Pacific Time)
Phone: (541) 335-0353
Toll-Free (800) 438-2209
Fax: (541) 335-0238
probetech@invitrogen.com

Invitrogen European Headquarters

Invitrogen, Ltd.
3 Fountain Drive
Inchinnan Business Park
Paisley PA4 9RF, UK
Phone: +44 (0) 141 814 6100
Fax: +44 (0) 141 814 6260
Email: euroinfo@invitrogen.com
Technical Services: eurotech@invitrogen.com

For country-specific contact information,
visit www.invitrogen.com.

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Limited Use Label License No. 223: Labeling and Detection Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation [®] are registered with the U.S. Patent and Trademark Office.

Copyright 2009, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.