00-8120 100 mL
CAS-Block™
Lot No. See product label

CAS-Block™

INVITROGEN HISTOCHEMICAL REAGENT

FORM:
Ready-to-use

APPLICATION:
CAS-Block™ is a universal blocking agent for reducing nonspecific background staining in immunolabelling techniques. It has been shown to be superior to conventional blocking techniques, and eliminates the need for species-specific blocking serum. CAS-Block™ can be used for immunohistochemistry, ELISA, and western blotting procedures. CAS-Block™ also works well as a diluting reagent for primary antibodies.

- RECOMMENDED GUIDELINES FOR USE:**
1. Deparaffinize slides and hydrate to water.
 2. Block endogenous peroxidase activity (if needed).
 3. Wash tissue slides in several changes of PBS.
 4. Apply CAS-Block™ to tissue sections for 10 min. at room temp.
 5. Drain-off tissue section (do not wash) and apply primary antibody.

STORAGE:
Store at 2-8°C.

Authorized Representative for IVDD 98/79/EC: Invitrogen Ltd.
Inchinnan Business Park
3 Fountain Drive
Paisley
PA4 9RF
UK

Explanation of symbols			
Symbol	Description	Symbol	Description
	Catalogue Number		Batch code
	Research Use Only		<i>In vitro</i> diagnostic medical device
	Use by		Temperature limitation
	Manufacturer		European Community authorised representative
	Without, does not contain		With, contains
	Protect from light		Consult accompanying documents
	Directs the user to consult instructions for use (IFU), accompanying the product.		

General purpose reagent. For laboratory use

www.invitrogen.com
Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288
E-mail: techsupport@invitrogen.com

00-8120 100 mL
CAS-Block™
Lot No. See product label

CAS-Block™

INVITROGEN HISTOCHEMICAL REAGENT

FORM:
Ready-to-use

APPLICATION:
CAS-Block™ is a universal blocking agent for reducing nonspecific background staining in immunolabelling techniques. It has been shown to be superior to conventional blocking techniques, and eliminates the need for species-specific blocking serum. CAS-Block™ can be used for immunohistochemistry, ELISA, and western blotting procedures. CAS-Block™ also works well as a diluting reagent for primary antibodies.

- RECOMMENDED GUIDELINES FOR USE:**
1. Deparaffinize slides and hydrate to water.
 2. Block endogenous peroxidase activity (if needed).
 3. Wash tissue slides in several changes of PBS.
 4. Apply CAS-Block™ to tissue sections for 10 min. at room temp.
 5. Drain-off tissue section (do not wash) and apply primary antibody.

STORAGE:
Store at 2-8°C.

Authorized Representative for IVDD 98/79/EC: Invitrogen Ltd.
Inchinnan Business Park
3 Fountain Drive
Paisley
PA4 9RF
UK

Explanation of symbols			
Symbol	Description	Symbol	Description
	Catalogue Number		Batch code
	Research Use Only		<i>In vitro</i> diagnostic medical device
	Use by		Temperature limitation
	Manufacturer		European Community authorised representative
	Without, does not contain		With, contains
	Protect from light		Consult accompanying documents
	Directs the user to consult instructions for use (IFU), accompanying the product.		

General purpose reagent. For laboratory use

www.invitrogen.com
Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288
E-mail: techsupport@invitrogen.com

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PI008120

(Rev 02/10) DCC-10-0301

© 2007 Invitrogen Corporation. All right reserved. These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses.

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Bevollmächtigter Repräsentant für IVDD 98/79/EC: Invitrogen Ltd.