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Sterisart® Universal-Pumpe 16419

Contents

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	EC Declaration of Conformity	2
1.	Intended Use of the Sterisart® Universal Pump	3
2.	Important Notes on the Safe and Proper Usage of the Sterisart® Universal Pump	3
3.	Unpacking the Pump and Checking the Equipment Supplied	3
4.	Description of the Sterisart® Universal Pump	4
5.	Startup	5
6.	Application	6
6.1	Operating the Sterisart® Universal Pump 16419	6
7.	Operating the Sterisart® Pump 16419	7
8.	Maintenance	7
9.	Cleaning and Sanitization	8
10.	Installation in an Isolator	8
11.	Important Notes Concerning the Warranty	9
12.	Troubleshooting	9
13.	Technical Specifications	10
14.	Ordering Information for the Sterisart® Universal Pump	10
15.	Ordering Information for Consumables	11
16.	Literature	11
17.	Information and Instructions on Disposal	12
	Non-hazardous Declaration	14

EC Declaration of Conformity
 We hereby declare that the Sterisart® Universal pump 16419/16420 has been designed and manufactured according to the following directives:

73/23/ EEC – Low voltage
 98/37/EEC – Machinery
 89/336/ EEC – Electromagnetic compatibility
 2002/95/EEC – RoHS Directive

The Sterisart® Universal pump 16419 is labeled with the CE mark.

1. Intended Use of the Sterisart® Universal Pump

The new Sterisart® Universal pump is a peristaltic pump equipped with appropriate devices for attaching a closed sterility test system and a sample container.

The pump draws liquid from the sample container into the specially designed double tubing, distributing the liquid equally to both sample containers of the sterility test system. The sample is rinsed in a similar way. After the nutrient media is added, the sample containers are incubated for a prescribed time and evaluated.

2. Important Notes on the Safe and Proper Usage of the Sterisart® Universal Pump

- The Sterisart® Universal pump must be switched off before any electrical connections are made or broken (e.g. when plugging in or unplugging the power cable or foot switch).
- The Sterisart® Universal pump must only be used with the original current-carrying parts supplied with it (power cable with ferrite core, D-Sub connecting cable and foot switch). The D-Sub protective cap supplied must cover the port when the D-Sub connecting cable is not connected.
- The Sterisart® Universal pump is designed for use with Sterisart NF filtration systems. The suitability of any other tubing must be approved by Sartorius Stedim Biotech.
- Only authorized Sartorius Stedim Biotech service technicians are allowed to open the pump housing.
- Do not expose the Sterisart® Universal pump to extreme temperatures.
- Improper use or handling will invalidate the manufacturer's guarantee.
- Sterisart® Universal pump may only be used with a grounded, readily accessible and fuse-protected electrical outlet that complies with national standards.

3. Unpacking the Pump and Checking the Equipment Supplied

When unpacking the Sterisart® Universal pump 16419 make sure not to leave any parts in the packaging. The equipment supplied includes the following accessories (one of each):

- Sterisart® Universal pump 16419
- Power cable (depends on country rating)
- Foot switch
- Ampoule breaker
- Container for draining
- Pump cover for container for draining (Sterisart® units)
- Pump cover for container for draining (Millipore units)
- Drain tubing
- Fastening arm
- Holding ring for bottles
- Support rod
- Bow
- Clamping lever with sleeve
- Operating instructions German | English

4. Description of the Sterisart® Universal Pump

The Sterisart® Universal pump 16419 (see Fig. 1) consists of pump housing (1) with pump head (2), lever (3), bow (4) and handles (5). On the right handle, there is also a bolt with a tilting lock for the container for draining (11). The selecting knob (6) is at the front on the right. The stand (7) with the holding ring for bottles (8) and fastening arm (9) and the optional ampoule breaker (10) is located on the left side of the pump.

The container for draining (Fig. 1; 11) for the sterility test containers is located on the right side of the pump. The drain tubing can be attached to the appropriate connector (barb) on the container for draining (Fig. 1; 12).


Fig. 1

Located on the rear panel (see Figure 2) is the On | Off switch (1), AC power socket (2), socket for the foot switch (3) and the manufacturer's ID label (8). Sockets (4), (5), (6) and (7) are not active in the Sterisart® basic model 16419.

Fig. 2

5. Startup

Set up the Sterisart® universal pump ergonomically in a suitable area (cleanroom, laminar flow work bench, isolator).

Before initial start-up, the accessories supplied must be assembled on the Sterisart Universal pump. Please proceed in the following order:

- At the back left, push the stand (Fig. 1, 7) onto the pin of the upper part of the pump until it is firmly in place. The holes in the foot of the stand allow the positioning of the stand to be altered.
- Fasten the holding ring for bottles (Fig. 1, 8) and the fastening arm (Fig. 1, 9) to the support rod. If necessary, attach the ampoule breaker (10) beforehand.
- Insert the bow (Fig. 1, 4) and the clamping lever (Fig. 1, 3) on the corresponding bolts on the top of the pump. The lever is jacketed in a plastic sleeve, which can be removed for cleaning. If signs of wear appear on this sleeve after it has been in use for some time, it can be re-ordered as a spare part (see chapter "Ordering information for the Sterisart® Universal Pump").
- After connecting the drain tubing to the barb (Fig. 1, 12) of the container for draining (Fig. 1, 11); attach the container for draining to the appropriate bolts on the right handle (Fig. 1, 5). The container for draining and the cover plate are attached by moving the tilting lock.
- Fasten the end of the drain tubing to the appropriate drain.

Very Important Note

Always make sure that there is a sufficient liquid level differential (> 5 cm) between the container for draining and the drain. This makes safe drainage possible and prevents backflow into the container for draining. Make sure that the drain tubing is kept lower than the pump to prevent the drain pan from overflowing.

Treat the waste as a biohazard. Only dispose of it in accordance with current regulations.

- Plug the foot switch into the designated socket located on the rear panel. (see Fig 2, 3). Fasten the plug to the socket by turning it clockwise.
- Plug the power cable into the male power socket on the back of the Sterisart® Universal pump and into a grounded electrical outlet (Fig. 2, 2).
- Use the On|Off switch (Fig. 2, 1) at the rear side of the Sterisart® Universal pump to turn it on.

6. Application

Caution!

The pump must only be operated with the tubing inserted and only by one person at a time. No corrosive or explosive fluids may be pumped! Close-fitting clothing must be worn (no ties, etc.).

To be able to operate the pump, Sterisart sterility test systems must be installed beforehand. To do this, follow the steps given below:

Remove the sterility test systems from their respective packaging and place the sterility test sample containers in the container for draining on the right side of the pump.

- Now open the lever and the bow.
- Keep the tubing around the rotor of the pump (the arrow on the fixed part surrounding the pump head indicates the subsequent direction of rotation), so that no kinks occur and the tubing is not under tension (Fig. 3). The tubing is attached by simply inserting it (with it lying on the ground) and fastening the lever (Fig. 4). In this way, the tubing is automatically correctly positioned.

Important Note!

When threading the tubing, make sure that it is not kinked or twisted.

- To start the pumping process, the tubing must be fully inserted and the lever must be engaged. A safety function ensures that the pump does not start if the lever or the bow are open.

Fig. 3

- After threading the tubing and closing the lever, insert the removal spike/needle of the Sterisart® NF system as described in the instructions for the Sterisart® NF sterility test system. Then attach the test liquid container (e.g., bag or bottle) by hanging it on the support rod (bag) or by placing it in the holding ring on the support rod (bottle). The holding ring accepts bottles of various sizes. The holder must be turned so that the sterility test system tubing is not under tension after the spike | needle has been inserted and so that the spike | needle cannot slip out of the container.
- For the following pumping process, make sure that the Sterisart® NF sample containers are correctly positioned in the intended holders, and that the tubing has been inserted without tension or kinking.

Important Note!

When the lever is opened, the Sterisart® Universal pump is turned off immediately by a safety function. In this event, the Sterisart® Universal pump can be restarted after adequately closing all components.

For testing ampoules, an ampoule breaker is also attached to the support rod. Use the integrated knife to scratch the neck of the ampoule, then snap off the top. The broken-off top of the ampoule falls into the collection container (which can be emptied later).

7. Operating the Sterisart® Universal Pump 16419

8. Maintenance

Fig. 4

The Sterisart® Universal Pump 16419 can be started up either by pressing the selecting knob or by activating the foot switch. The speed can be smoothly adjusted using the selecting knob. The pump is switched off by again pressing the selecting knob or by releasing the foot switch.

The Sterisart® Universal pump requires no user maintenance.

9. Cleaning and Sanitization

Caution

Ports on the back panel are not tight.

The Sterisart® Universal pump is not suited for sanitization (gaseous or atomized) with peracetic acid, hydrogen peroxide, or other decontamination procedures commonly used in isolators, unless it is built into an isolator countertop. These types of disinfectant can be used when the pump has been professionally installed in the isolator.

- Before cleaning the pump rotor, turn off the main power switch and unplug the power cable.

The following parts of the pump have to be removed for cleaning:

- Bow, lever, support rod and container for draining
- Disconnect the pump head:
 - The rotor must be turned by hand until the marking on the upper side (black dot) is pointing backwards.
 - The marking indicates the position of the fastening screw that is located in the rotor mid-section.
 - Loosen this screw with a 4-mm Allen wrench. Do not turn the screw out completely.
 - Pull the rotor upwards.
- Pull the speed-selecting knob upwards.

General Instructions:

Comply with the IP protection rating! No liquid may enter the Sterisart® Universal pump. Disconnect the pump from line voltage before any cleaning, maintenance, or repair work.

Cleaning stainless steel surfaces:

The stainless steel housing can be cleaned at regular intervals. Clean stainless steel parts with a wet cloth. Customary cleaning agents suitable for stainless steel can be used without risk. For surface sanitization, use commercially available disinfectants that are suitable for stainless steel.

Do not use any cleaning agents on the stainless steel housing that contain sodium hydroxide or acetic, sulfuric or citric acid. The use of steel wool cleaning sponges is prohibited.

All stainless steel surfaces can be cleaned by wiping with a moist cloth and with the following products:

- Peracetic acid, H₂O₂ and acetic acid (2.5%)
- Quaternary ammonium
- Formaldehyde and glutaraldehyde
- Isopropyl alcohol
- Ethanol 70%

Autoclavable at 121°C for 20 min. are bow, lever, speed selecting knob, support rod, holding ring for bottles, counter-holding ring for bottles and container for draining (complete).

Never autoclave the pump rotor

After cleaning, all parts should be secured as described in Chapter "Startup" (5). The rotor is assembled as follows:

- The rotor shaft has a groove on the side turned towards the back.
- Push the rotor upwards onto the shaft.
- Turn the marking on the rotor towards the back, so that the groove on the rotor shaft aligns with the marking on the rotor
- Tighten the screw in the rotor to guarantee that the rotor is firmly seated.

10. Installation in an Isolator

The installation must be carried out by the isolator manufacturer.

To install the Sterisart® Universal pump in an isolator, it must be opened up first. Make sure to unplug it from AC current (mains supply) before opening it!

During assembly, please follow the installation instructions precisely (part of the installation kit for isolators, order number 1ZE---0041). The instructions are available separately from the manufacturer.

11. Important Notes Concerning the Warranty

The correct functioning and safe use of the Sterisart® Universal pump is only guaranteed when it is operated according to the information given in these instructions.

The equipment manufacturer is responsible only for effects on safety, dependability and performance of the instrument if assembly, extensions, new settings, changes or repairs are carried out by persons authorized by him to do so, if the electrical installation in the room is in conformance with IEC requirements and if the equipment is used in accordance with these operating instructions.

Any changes made to the Sterisart® Universal pump will invalidate the manufacturer's warranty unless they are previously agreed to by authorized Sartorius Service staff, or carried out by them.

12. Troubleshooting

This section lists a few problems that could arise in connection with the Sterisart® pump. Several solutions are given for each problem. The first solution given is usually the most likely one. If this does not solve the problem, go down the list of suggestions until the problem has been solved.

The pump must be switched off before electrical connections are unplugged or plugged in!

Should it not be possible to locate and solve the problem in this way, contact the appropriate service center (see the back cover of this manual).

Problem	Solution
No function; no display	<p>The pump is not switched on, or is not correctly connected to the electric supply.</p> <p>Plug the power cable into the electric outlet correctly and to the Sterisart® Universal pump. Switch on the on-off switch.</p> <p>Check the electric supply; switch on multiple outlet sockets; check the wiring.</p> <p>If these possibilities are not feasible, the primary fuse in the main line filter may be defective.</p> <p>Changing this, with the power cable unplugged, must only be performed by authorized personnel!</p>
The pump rotor does not rotate.	<p>The locking lever is open or the foot switch is not properly connected.</p> <p>The bow has been removed from the pump.</p> <p>Check the connection of the foot switch to the pump.</p>

13. Technical Specifications

Dimensions pump in mm:	(W × D × H) approx. 336 × 260 × 210 (with lever)
Dimensions pump with holding ring for bottles, container for draining in mm:	(W × D × H) 440 × 365 × 485
Weight: Basic model 16419:	Approx. 13.5 kg
Maximum operating temperature:	40°C
IP protection class:	41

Electrical connection

Voltage:	100–240 VAC (alternating voltage)
Frequency:	50–60 Hz
Power consumption:	100 VA

Pump

Pump flow rate:	70–650 ml/min (depending on the tubing)
-----------------	---

Tubing

Correct functioning of the Sterisart® Universal pump is only guaranteed when used with tubing that has been approved by Sartorius Stedim Biotech.

Materials

Base plate, lever, bow, rotor:	Stainless steel 316L (1.4435)
Housing, connection port, support rod, holding ring for bottles:	Stainless steel 316 (1.4401)
Connecting cable:	PVC
Container for draining, plastic sleeve for lever, hydrodynamic bearing for bow:	PET
Hydrodynamic bearing for lever:	Iglidur X
Seals, gaskets, drain tubing:	Silicone

14. Ordering Information for the Sterisart® Universal Pump

16419	Sterisart® Universal pump, basic model
16973	Ampoule tongs
1ZG---0009	Holding ring for bottles
1ZF---0004	Fastening arm
1ZA---0002	Drain tubing
1ZW---0002	Ampoule breaker
1ZE---0033	Foot switch
1ZF---0013	Stand
1ZG---0022	Container for draining
1ZB---0001	Plastic sleeve for lever
1ZF---0007	Wing nut
1ZS---0005	Bow
1ZS---0002	Lever
1ZL---0003	Rotor
1ZG---0023	Cover for container for draining (Sterisart® sterility test units)
1ZG---0024	Cover for container for draining (Millipore sterility test units)
1ZG---0014	Adaptor for Sterisart® container for draining for use with a Millipore Equinox pump
1ZG---0039	Transport trolley
1ZE---0041	Installation kit isolators

Additional accessories are available on request.

15. Ordering Information for Consumables

Recommended disposable sterility test units for use with pump

Order No.	Description
16466-ACD	Sterisart® NF alpha, dual-needle metal spike for closed containers (box of 10, individually sterilized with ETO; single-packed).
16467-ACD	Sterisart® NF alpha, 6-cm metal needle for open containers (box of 10, individually sterilized with ETO; single-packed)
16468-ACD	Sterisart® NF alpha, system for medical devices with Luer lock connectors (box of 10, individually sterilized with ETO; single-packed)
16466-GBD	Sterisart® NF gamma, dual-needle metal spike for closed containers (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16467-GBD	Sterisart® NF gamma, 6-cm metal needle for open containers (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16468-GBD	Sterisart® NF gamma, system for medical devices with Luer or Luer lock connectors (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16469-GBD	Sterisart® NF gamma, system with adapter for prefilled syringes (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16470-GBD	Sterisart® NF gamma, system for slightly soluble powders in unvented vials (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16475-GBD	Sterisart® NF gamma, system for lyophilized or freely soluble powders in unvented vials (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16476-GBD	Sterisart® NF gamma, system with short dual-needle metal spike for closed containers (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16477-GBD	Sterisart® NF gamma, system for liquids in plastic containers (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.
16478-GBD	Sterisart® NF gamma, system for medical containers with Luer lock male connectors (box of 10, individually gamma sterilized, double-packed). Optimal for use in isolators.

Additional Sterisart® units are available on request.

16. Literature

Data sheet Sterisart® Universal Pump Publ. No.: SLD2010-e	English Order No.: 85030-531-17
Brochure Sterisart® Universal Pump Publ. No.: SLD1003-e	English Order No.: 85030-531-14

17. Information and Instructions on Disposal

Attention must be paid to the applicable laws, rules and regulations, when disposing of the contaminated components.

The packaging is to be taken to a local waste disposal site if no longer required.

The packaging is made of environmentally friendly materials that can be used as secondary raw materials.

The equipment, including accessories and batteries, does not belong in your regular household waste. EU legislation requires its Member States to collect electrical and electronic equipment and disposed of it separately from other unsorted municipal waste with the aim of recycling it.

In Germany and several other countries, Sartorius Stedim Biotech GmbH itself assumes responsibility for the return and conformant disposal of its electronic and electrical products. These products may not be placed with household waste or brought to collection centers run by local public disposal operations – not even by small commercial operators.

For disposal in Germany and in the other member nations of the European Economic Area (EEA), please contact our local service technicians or our Service Center in Goettingen, Germany:

Sartorius Stedim Biotech GmbH
Servicezentrum Biotechnologie
August-Spindler-Strasse 11
37079 Goettingen, Germany

WEEE-Reg.-No. DE 89907997

In countries that are not members of the European Economic Area (EEA) or where no Sartorius Stedim Biotech subsidiaries or dealerships are located, please contact your local authorities or a commercial disposal operator.

Prior to disposal and/or scrapping of the equipment, any batteries should be removed and disposed of in local collection boxes.

Sartorius Stedim Biotech will not take back equipment contaminated with hazardous materials (ABC contamination) – either for repair or disposal. Please refer to the package leaflet | the included CD ROM or to our web site (www.sartorius-stedim.com) for more detailed information and addresses for repair service or disposal of your device.

For equipment disposal in the following countries of the European Economic Area (EEA), please contact your local Sartorius Stedim affiliate, subsidiary, dealer or distributor:

Sales and Service Contacts

For further contacts, visit www.sartorius-stedim.com

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Asia | Pacific

China

Sartorius Stedim Beijing
Representative Office
No. 33, Yu'an Road,
Airport Industrial Zone B, Shunyi District
Beijing 101300

Phone +86.10.80426516
Fax +86.10.80426580

Sartorius Stedim Shanghai
Representative Office
Room 618, Tower 1, German Centre,
Shanghai, PRC., 201203

Phone +86.21.28986393
Fax +86.21.28986392.11

Sartorius Stedim Guangzhou Office
Room 704, Broadway Plaza,
No. 233-234 Dong Feng West Road
Guangzhou 510180

Phone +86.20.8351.7921
Fax +86.20.8351.7931

India

Sartorius Stedim India Pvt. Ltd.
10, 6th Main, 3rd Phase Peenya
KIADB Industrial Area
Bangalore – 560 058

Phone +91.80.2839.1963|0461
Fax +91.80.2839.8262

Japan

Sartorius Stedim Japan K.K.
KY Building, 8-11
Kita Shinagawa 1-chome
Shinagawa-ku
Tokyo 140-0001

Phone +81.3.3740.5407
Fax +81.3.3740.5406

Malaysia

Sartorius Stedim Malaysia Sdn. Bhd.
Lot L3-E-3B, Enterprise 4
Technology Park Malaysia
Bukit Jalil
57000 Kuala Lumpur

Phone +60.3.8996.0622
Fax +60.3.8996.0755

Singapore

Sartorius Stedim Singapore Pte. Ltd.
10, Science Park Road, The Alpha
#02-25, Singapore Science Park 2
Singapore 117684

Phone +65.6872.3966
Fax +65.6778.2494

Australia

Sartorius Stedim Australia Pty. Ltd.
Unit 5, 7-11 Rodeo Drive
Dandenong South Vic 3175

Phone +61.3.8762.1800
Fax +61.3.8762.1828

Non-hazardous Declaration

Declaration of Safety and Acceptability as Non-Hazardous to Health

The safety and health of our employees, the legal regulations governing the handling of materials that pose hazards to health and the environment and the regulations governing occupational safety mandate and require that this Non-hazardous Declaration be filled out completely for all equipment returned to the respective Sartorius Stedim affiliate, subsidiary, dealer or distributor listed for your country. No repair or return is possible unless a completed Non-hazardous Declaration form has been submitted.

Fax or mail a copy of a completed Non-hazardous Declaration form in advance to the respective Sartorius Stedim affiliate, subsidiary, dealer or distributor listed for your country (for the fax number or address, please refer to our list of countries). We require this information before the equipment/component arrives. An additional copy must accompany the equipment | component. If necessary, the shipping company must be notified.

Incomplete entries or non-compliance with this procedure will automatically lead to considerable delays in processing.

In the event of return, we do not differentiate between

- still utilizable and | or repairable used electrical and electronic equipment and no longer utilizable and/or repairable used electrical and electronic equipment or
- between waste that is subject to monitoring by the supervisory authorities and waste that is not subject to any special monitoring by supervisory authorities (used electrical and electronic equipment).

Model: _____ Serial No.: _____

I | we guarantee that the following measures have been taken:

- The equipment | component has been freed of hazardous materials to ensure that persons involved with the handling | repair are not exposed to any hazard or danger whatsoever.
- The equipment has been packaged for safety and the complete labeling has been affixed to it.
- The shipping agent has been notified (as prescribed) about the hazardous nature of the shipment.

The person sending in the equipment | component is hereby notified that said person shall be liable to Sartorius Stedim Biotech GmbH , its affiliates, subsidiaries, dealers and distributors as well as to any third parties – in particular to any of employees of Sartorius Stedim Biotech GmbH, its affiliates, its subsidiaries, dealers or distributors involved with the handling/repair of the equipment | component – for any damage, caused by entries that are incomplete or incorrect on account of negligence, gross negligence or willful intent.

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Phone +49.551.308.0
Fax +49.551.308.3289

www.sartorius-stedim.com

Copyright by Sartorius Stedim Biotech GmbH, Goettingen, Germany. All rights reserved. No part of this publication may be reprinted or translated in any form or by any means without the prior written permission of Sartorius Stedim Biotech GmbH.

The status of the information, specifications and illustrations in this manual is indicated by the date given below. Sartorius Stedim Biotech GmbH reserves the right to make changes to the technology, features, specifications, and design of the equipment without notice.

Status:
November 2009,
Sartorius Stedim Biotech GmbH,
Goettingen, Germany