

Amplite™ Colorimetric Aldehyde Quantitation Kit

Blue Color

Ordering Information	Storage Conditions	Instrument Platform
Product Number: 10053 (200 assays)	Keep at -20 °C Avoid exposure to light	Absorbance microplate readers

Introduction

The formation, reactivity and toxicity of aldehydes resulted from the peroxidation of lipids of cellular membranes have received great attention in recent years. Rapid and accurate measurement of aldehydes is an important task for biological research, chemical research, food industry and environmental pollution surveillance. There are a few reagents or assay kits available for quantifying aldehydes. Most of the existing aldehyde test methods are based on separations either by the tedious and expensive HPLC-MS or GC-MS.

Our Amplite™ Colorimetric Aldehyde Quantitation Kit uses a proprietary sensor that generates a chromogenic product with an absorbance at 620 nm upon reacting with an aldehyde. This kit provides a sensitive mix-and-read method to detect as little as 0.3 nanomole of aldehyde in a 100 µL assay volume (3 µM). The assay can be performed in a convenient 96-well or 384-well microtiter-plate format and easily adapted to automation without a separation step and the signal can be read by an absorbance plate reader at the wavelength between 620 and 660 nm.

Kit Key Features

Broad Application:	Used for quantifying aldehydes in a variety of applications, such as enzyme reactions.
Sensitive:	Detect as little as 0.3 nanomoles of aldehyde in a 100 µL assay volume.
Continuous:	Easily adapted to automation without a separation step.

Kit Components

Components	Amount
Component A: AldeView™ Blue	2 bottles
Component B: Assay Buffer	1 bottle (25 mL)
Component C: AldeView™ Blue Enhancer	1 bottle (10 mL)
Component D: Aldehyde Standard	1 vial

Assay Protocol for One 96-Well Plate

Brief Summary

**Prepare Aldehyde standards and/or test samples (50 µL) → Add 2X AldeView™ Blue reaction mixture (50 µL)
→ Incubate at RT for 20 minutes → Add 50 µL of AldeView™ Blue Enhancer → Incubate
at RT for 20 minutes → Monitor absorbance increase at 620 nm**

Note: Thaw all the kit components to room temperature before starting the experiment.

1. Prepare 2X AldeView™ Blue reaction mixture:

Add 5 mL of Assay Buffer (Component B) into one bottle of AldeView™ Blue (Component A) to make AldeView™ Blue reaction mixture.

Note: 5 mL of 2 X AldeView™ Blue reaction mixture is enough for one plate. The reaction mixture is not stable, and best used within 2 hours.

2. Prepare serial dilutions of aldehyde standard (0 to 100 µM):

- 2.1 Add 1 mL of Assay Buffer (Component B) into the vial of Aldehyde Standard (Component D) to make a 10 mM aldehyde standard stock solution.

Note: The unused 10 mM aldehyde standard stock solution should be divided into single use aliquots and stored at -20 °C.

- 2.2 Take 100 µL of 10 mM aldehyde standard stock solution (from Step 2.1) to perform 1:100, and 1:2 serial dilutions to get 100, 50, 25, 12.5, 6.25, 3.125, 1.56 and 0 µM serial dilutions with Assay Buffer (Component B).
- 2.3 Add serially diluted aldehyde standards and aldehyde-containing test samples into a white 96-well microplate with clear bottom as described in Tables 1 and 2.

Table 1. Layout of Aldehyde Standards and test samples in a white 96-well microplate with clear bottom

BL	BL	TS	TS														
AS1	AS1														
AS2	AS2																		
AS3	AS3																		
AS4	AS4																		
AS5	AS5																		
AS6	AS6																		
AS7	AS7																		

Note: AS= Aldehyde Standards, BL=Blank Control, TS=Test Samples.

Table 2. Reagent composition for each well

Aldehyde Standards	Blank Control	Test Sample
Serial dilutions*: 50 µL	Assay Buffer: 50 µL	50 µL

**Note: Add the serially diluted aldehyde standards from 1.56 µM to 100 µM into wells from AS1 to AS7 in duplicate.*

3. Run aldehyde assay:

- 3.1 Add 50 µL of 2X AldeView™ Blue reaction mixture (from Step 1) into each well of aldehyde standard, blank control, and test samples (see Step 2.3) to make the total aldehyde assay volume of 100 µL/well.

Note: For a 384-well plate, add 12.5 µL of test sample and 12.5 µL of 2X AldeView™ Blue reaction mixture into each well.

- 3.2 Incubate the reaction mixture at room temperature for 20-30 minutes (protected from light).

- 3.3 Add 50 µL of AldeView™ Blue Enhancer (Component C) into each well.

Note: For a 384-well plate, add 25 µL of AldeView™ Blue Enhancer into each well.

- 3.4 Monitor the absorbance increase at around 620 to 660 nm (Max at 620 nm) using an absorbance plate reader.

Data Analysis

The absorbance in blank wells (with 0 μM Aldehyde Standard) is used as a control, and is subtracted from the values of those wells with the aldehyde reactions. An aldehyde standard curve is shown in Figure 1.

Note: The absorbance background increases with time, thus it is important to subtract the absorbance value of the blank wells for each data point.




Figure 1. Aldehyde dose response was measured in a white wall/clear bottom 96-well plate with Amplitude™ Colorimetric Aldehyde Quantitation Kit using a SpectraMax microplate reader (Molecular Devices). As low as $\sim 3 \mu\text{M}$ of aldehyde can be detected with 30 minutes incubation ($n=3$).

References

1. Trevor M. Kitson. (1985) High concentrations of aldehydes slow the reaction of cytoplasmic aldehyde dehydrogenase with thiol-group modifiers *Biochem. J.* 228, 765.
2. Crabb DW, Matsumoto M, Chang D, You M (2004). Overview of the role of alcohol dehydrogenase and aldehyde dehydrogenase and their variants in the genesis of alcohol-related pathology. *The Proceedings of the Nutrition Society* 63 (1): 49.
3. Steinmetz CG, Xie P, Weiner H, Hurley TD (1997). Structure of mitochondrial aldehyde dehydrogenase: the genetic component of ethanol aversion. *Structure* 5 (5): 701.
4. O'Donnell JM, Kudej RK, LaNoue KF, Vatner SF, Lewandowski ED. (2004) Limited transfer of cytosolic NADH into mitochondria at high cardiac workload. *Am J Physiol Heart Circ Physiol*, 286, H2237.
5. Zurek G, and karst U (2000). 2,4-Dinitro-3,5,6-trideuterophenylhydrazones for the quantitation of aldehydes and ketones in air samples by liquid chromatography-mass spectrometry. *J of chromatography A*, 869, 251.
6. Ou Z, Ogamo A, Guo L, Konda Y, Harigaya Y, and Nakagawa Y. (1995). Identification and quantitation of choline glycerophospholipids that contain aldehyde residues by fluometric high-performance liquid chromatography. *Analytical biochemistry* 227, 289.

Warning: This kit is only sold to end users. Neither resale nor transfer to a third party is allowed without written permission from AAT Bioquest®. Chemical analysis of the kit components is strictly prohibited. Please call us at 408-733-1055 or e-mail us at info@aatbio.com if you have any questions.