

## Amplite™ Fluorimetric Acetylcholine Assay Kit

### \*Red Fluorescence\*

| Ordering Information:              | Storage Conditions:                        | Instrument Platform:            |
|------------------------------------|--|---------------------------------|
| Product Number: 11403 (200 assays) | Keep in freezer<br>Avoid exposure to light | Fluorescence microplate readers |

### Introduction

Acetylcholine (ACh) and its metabolites are involved in three main physiological purposes: structural integrity and signaling roles for cell membranes, cholinergic neurotransmission (acetylcholine synthesis), and as a major source for methyl groups via its metabolite. Acetylcholine is a neurotransmitter in both the central and peripheral nervous systems. It is one of many neurotransmitters in the autonomic nervous system (ANS) and the only neurotransmitter used in the motor division of the somatic nervous system. It is involved in a number of biological events that are related to diabetic vasculopathy, hypertension, and Alzheimer's disease.

Our Amplite™ Fluorimetric Acetylcholine Assay Kit provides one of the most sensitive methods for quantifying acetylcholine. The kit uses Amplite™ Red to quantify acetylcholine through the choline oxidase-mediated enzyme coupling reactions. The fluorescence intensity of Amplite™ Red is proportional to acetylcholine formation. The kit is an optimized "mix and read" assay. It provides an ultrasensitive one-step fluorimetric assay to detect as little as 0.01 nanomoles ACh in a 100 µL assay volume (0.1 µM). Its signal can be easily read with a fluorescence microplate reader at Ex/Em = ~540/590 nm or an absorbance microplate reader at ~576 nm.

### Kit Key Features

|                           |  |
|---------------------------|--|
| <b>Broad Application:</b> | Can be used for quantifying acetylcholine in solutions and in cell extracts. |
| <b>Sensitive:</b>         | Detect as low as 0.01 nmoles of acetylcholine in solution.                   |
| <b>Continuous:</b>        | Easily adapted to automation without a separation step.                      |
| <b>Convenient:</b>        | Formulated to have minimal hands-on time.                                    |
| <b>Non-Radioactive:</b>   | No special requirements for waste treatment.                                 |

### Kit Components

| Components                          | Amount                         |
|-------------------------------------|--------------------------------|
| Component A: Amplite™ Red           | 1 vial                         |
| Component B: Acetylcholine Probe    | 2 bottles (lyophilized powder) |
| Component C: Acetylcholine Standard | 1 vial                         |
| Component D: Assay Buffer           | 1 bottle (25 mL)               |
| Component E: DMSO                   | 1 vial (100 µL)                |

### Assay Protocol for One 96-well Plate

#### Brief Summary

**Prepare ACh Assay mixture (50 µL) → Add ACh standards or ACh test samples (50 µL) → Incubate at room temperature for 10 - 30 minutes → Monitor fluorescence intensity at Ex/Em = 540/590 nm**

*Note: Thaw all the kit components at room temperature before starting the experiment.*

**1. Prepare stock solutions:**

- 1.1 **250X Amplite™ Red stock solution:** Add 40 µL of DMSO (Component E) into the vial of Amplite Red™ (Component A) to make a 250X Amplite™ Red stock solution.  
*Note: The unused Amplite Red™ stock solution should be divided into single use aliquots. Store at -20 °C and avoid exposure to light.*
- 1.2 **Acetylcholine stock solution:** Add 200µL of ddH<sub>2</sub>O into the vial of Acetylcholine Standard (Component C) to make 50 mM acetylcholine stock solution.  
*Note: The unused acetylcholine stock solution should be divided into single use aliquots and stored at -20 °C.*

**2. Prepare acetylcholine assay mixture:**

- 2.1 Add 5 mL of Assay Buffer (Component D) to the bottle of Acetylcholine Probe (Component B) and mix well.
- 2.2 Add 20 µL of 250X Amplite Red™ stock solution (from Step 1.1) into the bottle of Acetylcholine Probe solution (from Step 2.1) to make the acetylcholine assay mixture.  
*Note: The Assay mixture should be used promptly and kept from light. The assay background would increase with longer storage time.*

**3. Prepare serial dilutions of acetylcholine standard (0 to 100 µM):**

- 3.1 Add 20 µL of 50 mM acetylcholine standard stock solution (from Step 1.2) to 980 µL Assay Buffer (Component D) to generate 1000 µM acetylcholine standard solution.  
*Note: Diluted acetylcholine standard solution is unstable, and should be used within 4 hours.*
- 3.2 Take 200 µL of 1000 µM acetylcholine standard to perform 1:10 and 1:3 serial dilutions to get 100, 30, 10, 3, 1, 0.3, 0.1 and 0 µM serial dilutions of acetylcholine standard.
- 3.3 Add acetylcholine standards and acetylcholine containing test samples into a 96-well solid black microplate as described in Tables 1 and 2.  
*Note: Treat cells or tissue samples as desired.*

**Table 1** Layout of acetylcholine standards and test samples in a solid black 96-well microplate

|     |     |      |      |      |      |  |  |  |  |  |  |
|-----|-----|------|------|------|------|--|--|--|--|--|--|
| BL  | BL  | TS   | TS   | .... | .... |  |  |  |  |  |  |
| AS1 | AS1 | .... | .... | .... | .... |  |  |  |  |  |  |
| AS2 | AS2 |      |      |      |      |  |  |  |  |  |  |
| AS3 | AS3 |      |      |      |      |  |  |  |  |  |  |
| AS4 | AS4 |      |      |      |      |  |  |  |  |  |  |
| AS5 | AS5 |      |      |      |      |  |  |  |  |  |  |
| AS6 | AS6 |      |      |      |      |  |  |  |  |  |  |
| AS7 | AS7 |      |      |      |      |  |  |  |  |  |  |

*Note: AS= Acetylcholine Standards; BL=Blank Control; TS=Test Samples*

**Table 2** Reagent composition for each well

| Acetylcholine Standard   | Blank Control       | Test Sample |
|--------------------------|---------------------|-------------|
| Serial Dilutions*: 50 µL | Assay Buffer: 50 µL | 50 µL       |

*\*Note: Add the serial dilutions of acetylcholine standard from 0.01 to 100 µM into wells from AS1 to AS7 in duplicate.*


**4. Run acetylcholine assay:**

- 4.1 Add 50 µL of acetylcholine assay mixture (from Step 2.2) to each well of the acetylcholine standard, blank control, and test samples (see Step 3.3) to make the total acetylcholine assay volume of 100 µL/well.  
*Note: For a 384-well plate, add 25 µL of sample and 25 µL of acetylcholine assay mixture per well.*
- 4.2 Incubate the reaction at room temperature for 10 to 30 minutes, protected from light.
- 4.3 Monitor the fluorescence increase with a fluorescence microplate reader at Ex/Em = 540/590 nm.

## Data Analysis

The fluorescence in blank wells (with the assay buffer only) is used as a control, and is subtracted from the values for those wells with the acetylcholine reactions. An acetylcholine standard curve is shown in Figure 1.

*Note: The fluorescence background increases with time, thus it is important to subtract the fluorescence intensity value of the blank wells for each data point.*



**Figure 1** Acetylcholine dose response was measured in a 96-well solid black plate with Amplite™ Fluorimetric Acetylcholine Assay Kit (Cat. # 11403) using a Gemini fluorescence microplate reader (Molecular devices). As low as 0.01 nmoles/well (0.1 µM) of acetylcholine can be detected with 10 minutes incubation (n=3).

## References

1. Kovarik, Z et al. (2003). Acetylcholinesterase active centre and gorge conformations analysed by combinatorial mutations and enantiomeric phosphonates. *Biochem. J.* (2003) 373, 33–40.
2. Ordentlich, A. et al. (1996). The Architecture of Human Acetylcholinesterase Active Center Probed by Interactions with Selected Organophosphate Inhibitors. *J. Biol. Chem.* 271 (20):11953–11962.
3. Magnottl, RA. et al. (1987). Measurement of Acetylcholinesterase in Erythrocytes in the Field. *Clin. Chem.* 33/10, 1731-1 735.

**Warning: This kit is only sold to end users. Neither resale nor transfer to a third party is allowed without written permission from AAT Bioquest. Chemical analysis of kit components is strictly prohibited. Please call us at 408-733-1055 or e-mail us at [info@aatbio.com](mailto:info@aatbio.com) if you have any questions.**