

**Thermo Scientific 1.1X PCR Master Mix  
w/ 1.5 mM MgCl<sub>2</sub>**

**Description:** The PCR Master Mix is a ready-to-use master mix. It is a convenient way of amplifying DNA fragments without the need to thaw individual components, reducing the risk of contamination and pipetting errors. The ThermoPrime *Taq* DNA Polymerase, dNTPs, reaction buffer and magnesium chloride are all present in the mix.

**Enzyme Source:** *Thermus aquaticus*

**Associated Activities:** ThermoPrime has 5' to 3' polymerization and exonuclease activity but lacks 3' to 5' exonuclease activity (proofreading).

**Kit Contents**

Vial	Pack Size (cap color)	
	A	B
PCR Master Mix	2 x 1.8 ml (red)	20 x 1.8 ml (red)

Master Mix: Each vial contains 1.8ml of a 1.1X working concentration PCR Master Mix sufficient for 80 x 25µl reactions. The final reaction 1X contains:

0.625 units	ThermoPrime <i>Taq</i> DNA Polymerase
75 mM	Tris-HCl (pH 8.8 at 25°C)
20 mM	(NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub>
1.5 mM	MgCl <sub>2</sub>
0.01% (v/v)	Tween <sup>®</sup> 20
0.2 mM	each of dATP, dCTP, dGTP and dTTP

**Storage Conditions:** Store the master mix at -20°C until ready for use for up to 1 year. Avoid freeze thawing. The vial can be stored at 4°C for up to 1 month. Shipped on ice within the UK and on dry ice internationally and within the US.

**Example of Protocol:**

Mix gently and spin down the master mix prior to use

	Volume	Final Concentration 1X
PCR Master Mix	22.5 µl	1X
Primer forward (12.5 µM each)	1 µl*	0.5 µM*
Primer reverse (12.5 µM each)	1 µl*	0.5 µM*
Water (PCR Grade)	variable	
DNA Template	0.5 µl	0.5 – 125 ng
Total volume	25 µl	

\*Scale up or down the volume and concentration as appropriate

**Example of program:**

	Temp.	Time	Number of cycle
Initial Denaturation	95°C	2 min	1 cycle
Denaturation	95°C	20 sec	30 to 40 cycles
Annealing	50-65°C	30 sec	
Extension**	72°C	60 sec	
Final Extension	72°C	5 min	1 cycle

\*\*Increase length of time in proportion to size of amplicon, *Taq* DNA Polymerase extends at approximately 1000 bp/min.

**Ordering Information:**

AB-0575/A	PCR Master Mix	160 x 25µl rxns
AB-0575/B	PCR Master Mix	1,600 x 25µl rxns

**Troubleshooting**

For technical information or troubleshooting contact Thermo Scientific Genomics Tech Support:

	United States	Rest of World
<b>Email</b>	<a href="mailto:abgene.techservice@thermofisher.com">abgene.techservice@thermofisher.com</a>	<a href="mailto:abgene.techsupport@thermofisher.com">abgene.techsupport@thermofisher.com</a>
<b>Phone</b>	+1 (800) 235-9880	+44 (0)1372 840410

Use of this product is covered by US Patent No. 6,127,155. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim, no right to perform any patented method and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California

Revised April, 2011. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Literature Code: AB-0575-v2-0411