

# LIBRARY PREPARATION

## NEBNext<sup>®</sup> DNA Library Prep Master Mix Set for 454<sup>™</sup>


Instruction Manual

NEB #E6070S/L  
10/50 reactions

 NEW ENGLAND  
**BioLabs**<sup>®</sup> Inc.  
*enabling technologies in the life sciences*

AMPURE® is a registered trademark of Beckman Coulter, Inc.  
NEBNEXT® a registered trademark of New England Biolabs, Inc.  
454™ and GS FLX TITANIUM™ are trademarks of Roche.  
LITMUS™ is a trademark of New England Biolabs, Inc.





## Table of Contents:

---

Applications.....	2
Protocols.....	3–4
NEBNext End Repair Enzyme Mix.....	5
NEBNext End Repair Reaction Buffer.....	6
Quick T4 DNA Ligase.....	7
NEBNext Quick Ligation Reaction Buffer.....	8
<i>Bst</i> DNA Polymerase, Large Fragment.....	9
NEBNext Adapter Fill-in Reaction Buffer.....	10
Molecular Biology Grade Water.....	11
Hydrophilic Streptavidin Magnetic Beads.....	12
NEBNext Bead Binding Buffer.....	13
NEBNext Bead Wash Buffer.....	14

## The Reagent Set Includes:

---

*The volumes provided are sufficient for preparation of up to 10 reactions (NEB #E6070S) and 50 reactions (NEB #E6070L).*

### Box 1: Store at –20°C

- NEBNext End Repair Enzyme Mix **RR**
- NEBNext End Repair Reaction Buffer (10X)
- Quick T4 DNA Ligase **RR**
- NEBNext Quick Ligation Reaction Buffer (5X) **RR**
- Bst* DNA Polymerase, Large Fragment **RR**
- NEBNext Adaptor Fill-in Reaction Buffer (10X) **RR**
- Molecular Biology Grade Water **RR**

### Box 2: Store at 4°C

- Hydrophilic Streptavidin Magnetic Beads (4 mg/ml)
- NEBNext Bead Binding Buffer (2X)
- NEBNext Bead Wash Buffer (1X)

## Applications:

The NEBNext DNA Library Prep Master Mix Set for 454 contains enzymes and buffers in convenient master mix formulations that are ideally suited for sample preparation for next-generation sequencing (1), and for preparation of single stranded DNA for use in high density hybridization arrays (2) or for genomic subtraction hybridization methods (3). Each of these components must pass rigorous quality control standards and are lot controlled, both individually and as a set of reagents.

**Lot Control:** The lots provided in the NEBNext DNA Library Prep Master Mix Set for 454 are managed separately and are qualified by additional functional validation. Individual reagents undergo standard enzyme activity and quality control assays, and also meet stringent criteria in the additional quality controls listed on each individual component page.

**Functional Validation:** Each set of reagents is functionally validated together through construction and sequencing of a genomic DNA library by 454 GS FLX Titanium™ (Roche), and by generation of single stranded DNA.

For larger volume requirements, customized and bulk packaging is available by purchasing through the OEM/Bulks department at NEB. Please contact OEM@neb.com for further information.

## References:

1. Maricic, T. and S. Paabo (2009). "Optimization of 454 sequencing library preparation from small amounts of DNA permits sequence determination of both DNA strands." *Biotechniques*, 46, 51–52, 54–57.
2. Straus, D. and F.M. Ausubel (1990). "Genomic subtraction for cloning DNA corresponding to deletion mutations." *Proc. Natl. Acad. Sci. USA* 87, 1889–1893.
3. Zhou, X. and D.T. Wong (2007). "Single nucleotide polymorphism mapping array assay." *Methods Mol. Biol.* 396, 295–314.

## Protocols:

### NEBNext End Repair Module Protocol

Starting Material: 1–5 µg of DNA Fragmented to 100–1000 bp in ≤ 85 µl

1. Mix the following components in a sterile microfuge tube:

Fragmented DNA	1–85 µl
NEBNext End Repair Reaction Buffer (10X)	10 µl
NEBNext End Repair Enzyme Mix	5 µl
Sterile H <sub>2</sub> O for a final volume of 100 µl	variable
<hr/> total volume	<hr/> 100 µl

2. Incubate in a thermal cycler for 30 minutes at 20°C.
3. Purify DNA sample on one column and elute in 30 µl of sterile dH<sub>2</sub>O or elution buffer.

### NEBNext Quick Ligation Module Protocol

1. Mix the following components in a sterile microfuge tube:

End Repaired, Blunt or dA-Tailed DNA	30 µl
Quick Ligation Reaction Buffer (5X)	10 µl
DNA Adaptors (not provided please use adaptors appropriate to specific application)	5 µl
Quick T4 DNA Ligase	5 µl
<hr/> total volume	<hr/> 50 µl

2. Incubate in a thermal cycler for 15 minutes at 20°C.
3. Purify DNA sample on one column and elute in 25 µl of sterile dH<sub>2</sub>O or elution buffer.

### NEBNext Fill-in and ssDNA Isolation Module Protocol

*Recommended: Removal of small fragments using Agencourt AMPure® Beads (Beckman Coulter, Inc.) or gel size selection.*

1. Transfer 50 µl of Hydrophilic Streptavidin Magnetic Beads to a 1.5 ml tube.
2. Using a magnet, pellet the beads and remove the buffer
3. Wash beads twice with 100 µl of 2X Bead Binding Buffer, pelleting the beads with a magnet to remove the buffer after each wash.
4. Resuspend beads in 25 µl of 2X Bead Binding Buffer.
5. Add 25 µl adapter-ligated DNA fragments to the beads.
6. Vortex and place on a tube rotator at room temperature for 20 minutes.
7. Using a magnet, wash the beads twice with 100 µl (1X) Bead Wash Buffer, pelleting the beads with a magnet to remove the buffer after each wash.

8. Mix the following components in a separate sterile microfuge tube:
 

Molecular Biology Grade Water	42 $\mu$ l
Adapter Fill-in Reaction Buffer	5 $\mu$ l
<i>Bst</i> DNA Polymerase, Large Fragment	3 $\mu$ l
9. Transfer the 50  $\mu$ l Fill-in Reaction Mix to the beads.
10. Vortex lightly and incubate at 37°C for 20 minutes.
11. Wash beads twice with 100  $\mu$ l of 1X Bead Wash Buffer, pelleting the beads with a magnet to remove the buffer after each wash.
12. Prepare Melt Solution:
 

10 N NaOH	125 $\mu$ l
Water	9.875 ml
13. Prepare Neutralization Solution in a 1.5 ml tube.
 

3 M sodium acetate, pH 5.2	10 $\mu$ l
Column binding buffer with pH indicator	500 $\mu$ l
14. Add 50  $\mu$ l of Melt Solution to the beads.
15. Vortex well, pellet the beads with a magnet.
16. Carefully transfer the Melt Solution containing the ssDNA Fragment Library to 1.5 ml tube containing the Neutralization Solution.
17. Repeat steps 14–16, adding the second round of Melt Solution containing the ssDNA Fragment Library to the same 1.5 ml tube containing the Neutralization Solution and the first round of Melt Solution and ssDNA Fragments. Adjust pH if necessary by adding an additional 5  $\mu$ l of 3 M sodium acetate.
18. Purify DNA on one column without adding any additional column binding buffer. Wash column twice with column wash buffer to remove all residual salts. Elute in 25  $\mu$ l of sterile dH<sub>2</sub>O or elution buffer.

# NEBNext End Repair Enzyme Mix

#E6041A: 0.06 ml

#E6041AA: 0.3 ml



Store at -20°C

**Description:** NEBNext End Repair Enzyme Mix is optimized to convert 1 to 5 µg of fragmented to repaired DNA having 5'-phosphorylated, blunt ends.

## **NEBNext End Repair Enzyme Mix:**

10,000 units/ml T4 Polynucleotide Kinase

3,000 units/ml T4 DNA Polymerase

## **Storage Conditions:**

10 mM Tris-HCl

100 mM KCl

1 mM DTT

0.1 mM EDTA

50% Glycerol

0.1% Triton X-100

pH 7.4 @ 25°C

## Quality Control Assays

**SDS-PAGE Purity:** SDS-PAGE analyses of each individual enzyme indicates > 95% enzyme purity.

**Endonuclease Activity:** Incubation of a minimum of 10 µl of this enzyme mix with 1 µg of φX174 RF I DNA in assay buffer for 4 hours at 37°C in 50 µl reactions results in less than 10% conversion to RF II as determined by agarose gel electrophoresis.

**Phosphatase Activity:** Incubation of a minimum of 10 µl of this enzyme mix in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM MgCl<sub>2</sub>) containing 2.5 mM *p*-nitrophenyl phosphate at 37°C for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

**Functional Activity (Nucleotide Incorporation):** 0.2 µl of this enzyme mix incorporates 10 nmol of dNTP into acid-precipitable material in a total reaction volume of 50 µl in 30 minutes at 37°C in 1X T4 DNA Polymerase Reaction Buffer with 33 µM dNTPs including [<sup>3</sup>H]-dTTP, 70 µg/ml denatured herring sperm DNA and 50 µg/ml BSA.

**Functional Activity (Nucleotide Incorporation and Phosphorylation):** 5 µl of this enzyme mix repairs and phosphorylates the ends of > 95% of 10 µg of DNA fragments containing both 3' and 5' overhangs within 30 minutes at 20°C in 1X End Repair Buffer, as determined by capillary electrophoresis.

**Lot Controlled**

# NEBNext End Repair Reaction Buffer

#E6042A: 0.12 ml

Concentration: 10X

#E6042AA: 0.6 ml

Store at -20°C

## 1X NEBNext End Repair Reaction Buffer:

50 mM Tris-HCl

10 mM MgCl<sub>2</sub>

10 mM DTT

1 mM ATP

0.4 mM dATP

0.4 mM dCTP

0.4 mM dGTP

0.4 mM dTTP

pH 7.5 @ 25°C

## Quality Control Assays

**16-Hour Incubation:** 50 µl reactions containing this reaction buffer at 1X concentration and 1 µg of HaeIII digested φX174 RF I DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis. 50 µl reactions containing this reaction buffer at 1X concentration and 1 µg T3 DNA incubated for 16 hours at 37°C also results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis.

**Endonuclease Activity:** Incubation of this reaction buffer at a 1X concentration with 1 µg of φX174 RF I DNA for 4 hours at 37°C in 50 µl reactions results in less than 10% conversion to RF II as determined by agarose gel electrophoresis.

**RNase Activity:** Incubation of this reaction buffer at 1X concentration with 40 ng of a FAM-labeled RNA transcript for 16 hours at 37°C results in no detectable RNase activity as determined by polyacrylamide gel electrophoresis.

**Phosphatase Activity:** Incubation of this reaction buffer at a 1X concentration in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM MgCl<sub>2</sub>) containing 2.5 mM *p*-nitrophenyl phosphate at 37°C for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

**Lot Controlled**



# Quick T4 DNA Ligase

#E6047A: 0.06 ml

#E6047AA: 0.3 ml



Store at  $-20^{\circ}\text{C}$

Source: Purified from *E. coli* C600 pcl857 pPLc28 lig8 (2).

## Quality Control Assays

**SDS-PAGE Purity:** SDS-PAGE analysis of this enzyme indicates > 95% enzyme purity.

**16-Hour Incubation:** 50  $\mu\text{l}$  reactions containing a minimum of 2,000 units of this enzyme and 1  $\mu\text{g}$  of HaeIII digested  $\phi\text{X174}$  RF I DNA incubated for 16 hours at  $37^{\circ}\text{C}$  results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis. 50  $\mu\text{l}$  reactions containing a minimum of 2,000 units of this enzyme and 1  $\mu\text{g}$  T3 DNA incubated for 16 hours at  $37^{\circ}\text{C}$  also results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis.

**Endonuclease Activity:** Incubation of a minimum of 3,200 units of this enzyme with 1  $\mu\text{g}$  of  $\phi\text{X174}$  RF I DNA in assay buffer for 4 hours at  $37^{\circ}\text{C}$  in 50  $\mu\text{l}$  reactions results in less than 10% conversion to RF II as determined by agarose gel electrophoresis.

**Phosphatase Activity:** Incubation of a minimum of 20,000 units of this enzyme in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM  $\text{MgCl}_2$ ) containing 2.5 mM *p*-nitrophenyl phosphate at  $37^{\circ}\text{C}$  for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

**RNase Activity:** Incubation of a minimum of 2,000 units of this enzyme with 40 ng of a FAM-labeled RNA transcript for 16 hours at  $37^{\circ}\text{C}$  results in no detectable RNase activity as determined by polyacrylamide gel electrophoresis.

**Exonuclease Activity:** Incubation of a minimum of 3,200 units of this enzyme with 1  $\mu\text{g}$  sonicated [ $^3\text{H}$ ] DNA ( $10^5$  cpm/ $\mu\text{g}$ ) for 4 hours at  $37^{\circ}\text{C}$  in 50  $\mu\text{l}$  reaction buffer releases < 0.1% radioactivity.

**Functional Activity (Blunt End Ligation):** 50  $\mu\text{l}$  reactions containing a 0.5  $\mu\text{l}$  Quick T4 DNA Ligase, 18  $\mu\text{g}$  HaeIII digested  $\phi\text{X174}$  and 1X T4 DNA Ligase Buffer incubated at  $16^{\circ}\text{C}$  for 7.5 min results in > 95% of fragments ligated as determined by agarose gel electrophoresis.

**Functional Activity (Cohesive End Ligation):** 20  $\mu\text{l}$  reactions containing 0.5  $\mu\text{l}$  Quick T4 DNA Ligase, 12  $\mu\text{g}$  HindIII digested lambda DNA and 1X T4 DNA Ligase Buffer incubated at  $37^{\circ}\text{C}$  overnight results in > 95% of fragments ligated as determined by agarose gel electrophoresis. Redigestion of the ligated products, 50  $\mu\text{l}$  reactions containing 6  $\mu\text{g}$  of the ligated fragments, 40 units HindIII, and 1X NEBuffer 2 incubated at  $37^{\circ}\text{C}$  for 2 hours, results in no detectable undigested fragments as determined by agarose gel electrophoresis.

**Functional Activity (Adapter Ligation):** 50  $\mu\text{l}$  reactions containing 0.125  $\mu\text{l}$  Quick T4 DNA Ligase, 8 nmol 12 bp adapter, and 1X T4 DNA Ligase Buffer incubated at  $16^{\circ}\text{C}$  overnight results in no detectable unligated adapter as determined by agarose gel electrophoresis.

**Functional Activity (Transformation):** After a five-minute ligation of linearized, dephosphorylated LITMUS™ 28 (containing either blunt [EcoRV] or cohesive [HindIII] ends) and a mixture of compatible insert fragments, transformation into chemically competent *E. coli* DH-5 alpha cells yields a minimum of  $1 \times 10^6$  recombinant transformants per  $\mu\text{g}$  plasmid DNA.

## Lot Controlled

### References:

1. Engler, M. J. and Richardson, C. C. (1982). In P. D. Boyer (Ed.), *The Enzymes* Vol. 5, (p. 3). San Diego: Academic Press.
2. Remaut, E., Tsao, H. and Fiers, W. (1983) *Gene*, 22, 103–113.

# NEBNext Quick Ligation Reaction Buffer

#E6048A: 0.12 ml

Concentration: 5X

#E6048AA: 0.6 ml

Store at -20°C

## 1X NEBNext Quick Ligation Reaction Buffer:

66 mM Tris-HCl

10 mM MgCl<sub>2</sub>

1 mM dithiothreitol

1 mM ATP

6% Polyethylene glycol (PEG 6000)

pH 7.6 @ 25°C

## Quality Control Assays

**16-Hour Incubation:** 50 µl reactions containing this reaction buffer at 1X concentration and 1 µg of HaeIII digested φX174 RF I DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis. 50 µl reactions containing this reaction buffer at 1X concentration and 1 µg T3 DNA incubated for 16 hours at 37°C also results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis.

**Endonuclease Activity:** Incubation of this reaction buffer at a 1X concentration with 1 µg of φX174 RF I DNA for 4 hours at 37°C in 50 µl reactions results in less than 10% conversion to RF II as determined by agarose gel electrophoresis.

**RNase Activity:** Incubation of this reaction buffer at 1X concentration with 40 ng of a FAM-labeled RNA transcript for 16 hours at 37°C results in no detectable RNase activity as determined by polyacrylamide gel electrophoresis.

**Phosphatase Activity:** Incubation of this reaction buffer at a 1X concentration in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM MgCl<sub>2</sub>) containing 2.5 mM *p*-nitrophenyl phosphate at 37°C for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

**Lot Controlled**

# Bst DNA Polymerase, Large Fragment

#E6030A: 0.03 ml

#E6030AA: 0.15 ml



**Store at -20°C**

**Description:** *Bst* DNA Polymerase, Large Fragment is the portion of the *Bacillus stearothermophilus* DNA Polymerase protein that contains the 5' → 3' polymerase activity, but lacks 5' → 3' exonuclease activity.

**Source:** *Bst* Polymerase, Large Fragment is prepared from an *E. coli* strain containing a genetic fusion of the *Bacillus stearothermophilus* DNA Polymerase gene, lacking the 5' → 3' exonuclease domain, and the gene coding for *E. coli* maltose binding protein (MBP). The fusion protein is purified to near homogeneity and the MBP portion of the fusion is cleaved off *in vitro*. The remaining polymerase is purified free of MBP (1).

**Supplied in:** 50 mM KCl, 10 mM Tris-HCl (pH 7.5), 0.1 mM EDTA, 1 mM dithiothreitol, 0.1% Triton X-100 and 50% glycerol.

## Quality Control Assays

**SDS-PAGE Purity:** SDS-PAGE analysis of this enzyme indicates > 95% enzyme purity.

**16-Hour Incubation:** 50 µl reactions containing a minimum of 8 units of this enzyme and 1 µg of HaeIII digested φX174 RF I DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis. 50 µl reactions containing a minimum of 8 units of this enzyme and 1 µg T3 DNA incubated for 16 hours at 37°C also results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis.

**Endonuclease Activity:** Incubation of a minimum of 50 units of this enzyme with 1 µg of φX174 RF I DNA in assay buffer for 4 hours at 37°C in 50 µl reactions results in less than 10% conversion to RF II as determined by agarose gel electrophoresis.

**Phosphatase Activity:** Incubation of a minimum of 50 units of this enzyme in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM MgCl<sub>2</sub>) containing 2.5 mM *p*-nitrophenyl phosphate at 37°C for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

**RNase Activity:** Incubation of a minimum of 8 units of this enzyme with 40 ng of FAM-labeled RNA transcript for 16 hours at 37°C results in no detectable RNase activity as determined by polyacrylamide gel electrophoresis.

**Functional Activity (Nucleotide Incorporation):** One unit of this enzyme incorporates 10 nmol of dNTP into acid-precipitable material in a total reaction volume of 50 µl in 30 minutes at 37°C in 1X NEBuffer 2 with 33 µM dNTPs including [<sup>3</sup>H]-dTTP, 70 µg/ml denatured herring sperm DNA and 50 µg/ml BSA.

**Lot Controlled**

**References:**

1. Kong, H., Aliotta, J. and Pelletier, J.J., New England Biolabs, unpublished results.

U.S. Patent No. 5,814,506.

# NEBNext Adapter Fill-in Reaction Buffer

#E6035A: 0.05 ml

Concentration: 10X

#E6035AA: 0.25 ml

Store at  $-20^{\circ}\text{C}$

## 1X NEBNext Adapter Fill-in Reaction Buffer:

20 mM Tris-HCl

10 mM  $(\text{NH}_4)_2\text{SO}_4$

10 mM KCl

2.0 mM  $\text{MgSO}_4$

0.1% Triton X-100

0.4 mM dATP

0.4 mM dCTP

0.4 mM dGTP

0.4 mM dTTP

pH 8.8 @  $25^{\circ}\text{C}$

## Quality Control Assays

**16-Hour Incubation:** 50  $\mu\text{l}$  reactions containing this reaction buffer at 1X concentration and 1  $\mu\text{g}$  of HaeIII digested  $\phi\text{X174}$  RF I DNA incubated for 16 hours at  $37^{\circ}\text{C}$  results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis. 50  $\mu\text{l}$  reactions containing this reaction buffer at 1X concentration and 1  $\mu\text{g}$  T3 DNA incubated for 16 hours at  $37^{\circ}\text{C}$  also results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis.

**Endonuclease Activity:** Incubation of this reaction buffer at a 1X concentration with 1  $\mu\text{g}$  of  $\phi\text{X174}$  RF I DNA for 4 hours at  $37^{\circ}\text{C}$  in 50  $\mu\text{l}$  reactions results in less than 10% conversion to RF II as determined by agarose gel electrophoresis.

**Phosphatase Activity:** Incubation of this reaction buffer at a 1X concentration in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM  $\text{MgCl}_2$ ) containing 2.5 mM *p*-nitrophenyl phosphate at  $37^{\circ}\text{C}$  for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

**Lot Controlled**

# Molecular Biology Grade Water

#E6031A: 1 ml

#E6031AA: 5 ml

**Store at -20°C**

**Description:** Molecular Biology Grade Water is free of detectable DNA and RNA nucleases and phosphatases and suitable for use in molecular biology applications.

## Quality Control Assays

**16-Hour Incubation:** 50 µl reactions containing molecular biology grade water, 1X NEBuffer 2 (50mM NaCl, 10 mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 1 mM dithiothreitol, pH 7.9) and 1 µg of HaeIII digested φX174 RF I DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis.

50 µl reactions containing 43 µl of Molecular Biology Grade Water and 1 µg T3 DNA incubated for 16 hours at 37°C also results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis.

**Endonuclease Activity:** Incubation of molecular biology grade water, 1X NEBuffer 2 (50mM NaCl, 10 mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 1 mM dithiothreitol, pH 7.9) and 1 µg of φX174 RF I DNA for 4 hours at 37°C in 50 µl reactions results in less than 10% conversion to RF II as determined by agarose gel electrophoresis.

**RNase Activity:** Incubation of molecular biology grade water with 1X NEBuffer 2 (50mM NaCl, 10 mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 1 mM dithiothreitol, pH 7.9) and 40 ng of a FAM-labeled RNA transcript in a 10 µl reaction for 16 hours at 37°C results in no detectable RNase activity as determined by polyacrylamide gel electrophoresis.

**Phosphatase Activity:** Incubation of molecular biology grade water in protein phosphatase assay buffer (1M diethanolamine @ pH 9.8 and 0.5 mM MgCl<sub>2</sub>) containing 2.5 mM *p*-nitrophenyl phosphate in a 200 µl reaction at 37°C for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

**Lot Controlled**

## Hydrophilic Streptavidin Magnetic Beads

#E6032A: 0.5 ml

Concentration: 4 mg/ml

#E6032AA: 2.5 ml

**Store at 4°C**

**Description:** Hydrophilic Streptavidin Magnetic Beads are 2 µm supermagnetic particles covalently coupled to a highly pure form of streptavidin. The beads can be used to capture biotin labeled DNA.

**Supplied in:** 4 mg/ml suspension in phosphate buffer (PBS) (pH 7.4) containing 0.1% BSA and 0.02%  $\text{NaN}_3$ .

**Support Matrix:** 2 µm non-porous magnetic microparticle.

**Binding Capacity:** The beads will bind greater than 800 pmol of free biotin per mg and greater than 400 pmol of single-stranded 20 bp biotinylated oligonucleotide per mg.

**Lot Controlled**

# NEBNext Bead Binding Buffer

#E6034A: 2.25 ml

Concentration: 2X

#E6034AA: 11.25 ml

Store at 4°C

## 1X NEBNext Bead Binding Buffer:

5 mM Tris-HCl

0.5 mM EDTA

1 M NaCl

pH 7.5 @ 25°C

**Description:** NEBNext Bead Binding Buffer has been optimized for binding biotin-labeled DNA to Magnetic Streptavidin Beads.

## Quality Control Assays

**16-Hour Incubation:** 50 µl reactions containing this reaction buffer at 0.2X concentration and 1 µg of HaeIII digested φX174 RF I DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis. 50 µl reactions containing this reaction buffer at 1X concentration and 1 µg T3 DNA incubated for 16 hours at 37°C also results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis.

**Endonuclease Activity:** Incubation of this reaction buffer at a 0.2X concentration with 1 µg of φX174 RF I DNA for 4 hours at 37°C in 50 µl reactions results in less than 10% conversion to RF II as determined by agarose gel electrophoresis.

**RNase Activity:** Incubation of this reaction buffer at 0.2X concentration with 40 ng of a FAM-labeled RNA transcript for 16 hours at 37°C results in no detectable RNase activity as determined by polyacrylamide gel electrophoresis.

**Phosphatase Activity:** Incubation of this reaction buffer at a 1X concentration in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM MgCl<sub>2</sub>) containing 2.5 mM *p*-nitrophenyl phosphate at 37°C for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

**Lot Controlled**

# NEBNext Bead Wash Buffer

#E6033A: 4.0 ml

Concentration: 1X

#E6033AA: 20.0 ml

Store at 4°C

## 1X NEBNext Bead Wash Buffer:

5 mM Tris-HCl

0.5 mM EDTA

1 M NaCl

pH 7.5 @ 25°C

**Description:** NEBNext Bead Wash Buffer has been optimized to remove non-specific binding to Magnetic Streptavidin Beads.

## Quality Control Assays

**16-Hour Incubation:** 50  $\mu$ l reactions containing this reaction buffer at 0.1X concentration and 1  $\mu$ g of HaeIII digested  $\phi$ X174 RF I DNA incubated for 16 hours at 37°C results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis. 50  $\mu$ l reactions containing this reaction buffer at 0.1X concentration and 1  $\mu$ g T3 DNA incubated for 16 hours at 37°C also results in no detectable non-specific nuclease degradation as determined by agarose gel electrophoresis.

**Endonuclease Activity:** Incubation of this reaction buffer at a 0.1X concentration with 1  $\mu$ g of  $\phi$ X174 RF I DNA for 4 hours at 37°C in 50  $\mu$ l reactions results in less than 10% conversion to RF II as determined by agarose gel electrophoresis.

**RNase Activity:** Incubation of this reaction buffer at 0.1X concentration with 40 ng of a FAM-labeled RNA transcript for 16 hours at 37°C results in no detectable RNase activity as determined by polyacrylamide gel electrophoresis.

**Phosphatase Activity:** Incubation of this reaction buffer at a 0.5X concentration in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM  $\text{MgCl}_2$ ) containing 2.5 mM *p*-nitrophenyl phosphate at 37°C for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

**Lot Controlled**











## USA

New England Biolabs, Inc.  
240 County Road  
Ipswich, MA 01938-2723  
Telephone: (978) 927-5054  
Toll Free: (USA Orders) 1-800-632-5227  
Toll Free: (USA Tech) 1-800-632-7799  
Fax: (978) 921-1350  
e-mail: [info@neb.com](mailto:info@neb.com)  
[www.neb.com](http://www.neb.com)

## Canada

New England Biolabs, Ltd.  
Telephone: (905) 837-2234  
Toll Free: 1-800-387-1095  
Fax: (905) 837-2994  
Fax Toll Free: 1-800-563-3789  
e-mail: [info@ca.neb.com](mailto:info@ca.neb.com)

## China, People's Republic

New England Biolabs (Beijing), Ltd.  
Telephone: 010-82378265/82378266  
Fax: 010-82378262  
e-mail: [info@neb-china.com](mailto:info@neb-china.com)

## France

New England Biolabs France  
Free Call: 0800/100 632  
Free Fax: 0800/100 610  
e-mail: [info@fr.neb.com](mailto:info@fr.neb.com)

## Germany

New England Biolabs GmbH  
Telephone: +49/(0)69/305 23140  
Free Call: 0800/246 5227 (Germany)  
Fax +49/(0)69/305 23149  
Free Fax: 0800/246 5229 (Germany)  
e-mail: [info@de.neb.com](mailto:info@de.neb.com)

## Japan

New England Biolabs Japan, Inc.  
Telephone: +81 (0)3 5669 6191  
Fax +81 (0)3 5669 6192  
e-mail: [info@neb-japan.com](mailto:info@neb-japan.com)

## United Kingdom

New England Biolabs (UK), Ltd.  
Telephone: (01462) 420616  
Call Free: 0800 318486  
Fax: (01462) 421057  
Fax Free: 0800 435682  
e-mail: [info@uk.neb.com](mailto:info@uk.neb.com)