β-N-Acetylhexosaminidase,





1-800-632-7799 info@neb.com www.neb.com

P0721S

RX

500 units 5,000 U/ml Lot: 0011210 RECOMBINANT Store at -20°C Exp: 10/14

Description: β-N-Acetyl-hexosaminidase, is a recombinant protein fusion of β-N-Acetyl-hexosaminidase (1) and maltose binding protein. It has identical activity to β-N-Acetyl-hexosaminidase. β-N-Acetyl-hexosaminidase, catalyzes the hydrolysis of terminal β-D-N-acetyl-galactosamine and glucosamine residues from oligosaccharides.

*Note: Specificity Change

β-N-Acetylhexosaminidase,



1-800-632-7799 info@neb.com www.neb.com

P0721S

 \mathbf{R}^{*}

500 units 5,000 U/ml Lot: 0011210 RECOMBINANT Store at -20°C Exp: 10/14

Description: β-N-Acetyl-hexosaminidase, is a recombinant protein fusion of β-N-Acetyl-hexosaminidase (1) and maltose binding protein. It has identical activity to β-N-Acetyl-hexosaminidase. β-N-Acetyl-hexosaminidase, catalyzes the hydrolysis of terminal β-D-N-acetyl-galactosamine and glucosamine residues from oligosaccharides.

*Note: Specificity Change

*Specificity:

GICNAC $\begin{cases} \begin{matrix} \downarrow \\ \beta 1 = 3 \text{ R} \\ \beta 1 - 4 \text{ R} \\ \beta 1 - 6 \text{ R} \end{cases}$ GalNAc $\begin{cases} \downarrow \\ \beta 1 = 4 \text{ R} \end{cases}$

Source: Cloned from *Streptomyces plicatus* (1) and overexpressed in *E. coli* (2).

Supplied in: 50 mM NaCl, 20 mM Tris-HCl (pH 7.5 @ 25°C) and 5 mM Na₂EDTA.

Reagents Supplied with Enzyme:

10X G2 Reaction Buffer

Reaction Conditions:

1X G2 Reaction Buffer: 50 mM Sodium Citrate (pH 4.5 @ 25°C). Incubate at 37°C.

Optimal incubation times and enzyme concentrations must be determined empirically for a particular substrate.

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal β -D-N-acetyl-galactosamine from

*Specificity:

GICNAC $\begin{cases} \beta 1 \stackrel{\checkmark}{=} 3 R \\ \beta 1 - 4 R \\ \beta 1 - 6 R \end{cases}$ GaINAC $\begin{cases} \beta 1 \stackrel{\checkmark}{=} 4 R \\ \beta 1 \stackrel{\checkmark}{=} 6 R \end{cases}$

Source: Cloned from *Streptomyces plicatus* (1) and overexpressed in *E. coli* (2).

Supplied in: 50 mM NaCl, 20 mM Tris-HCl (pH 7.5 @ 25°C) and 5 mM Na₂EDTA.

Reagents Supplied with Enzyme:

10X G2 Reaction Buffer

Reaction Conditions:

1X G2 Reaction Buffer: 50 mM Sodium Citrate (pH 4.5 @ 25°C). Incubate at 37°C.

Optimal incubation times and enzyme concentrations must be determined empirically for a particular substrate.

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal β-p-N-acetyl-galactosamine from

1 nmol of GalNAcβ1-4Galβ1-4Glc-7-amino-4-methyl-coumarin (AMC), in 1 hour at 37°C in a total reaction volume of 10 ul.

Unit Definition Assay: Two fold dilutions of β -N-Acetyl-hexosaminidase, are incubated with 1 nmol AMC-labeled substrate in 1X G2 Reaction Buffer in a 10 μ l reaction. The reaction mix is incubated for 1 hour at 37°C. Separation of reaction products are visualized via thin layer chromatography (3).

Specific Activity: ~ 10,000 units/mg

Molecular Weight: 100,000 daltons

Quality Assurance: No contaminating exoglycosidase or proteolytic activity could be detected.

Quality Controls

Glycosidase Assays:

50 units of β -N-Acetyl-hexosaminidase, were incubated with 0.1 mM of fluorescently-labeled oligosaccharides and glycopeptides, in a 10 μ l reaction for 20 hours at 37°C. The reaction products were analyzed by TLC for digestion of substrate.

1 nmol of GalNAc β 1-4Gal β 1-4Glc-7-amino-4-methyl-coumarin (AMC), in 1 hour at 37°C in a total reaction volume of 10 μ l.

Unit Definition Assay: Two fold dilutions of β -N-Acetyl-hexosaminidase, are incubated with 1 nmol AMC-labeled substrate in 1X G2 Reaction Buffer in a 10 μ l reaction. The reaction mix is incubated for 1 hour at 37°C. Separation of reaction products are visualized via thin layer chromatography (3).

Specific Activity: ~ 10,000 units/mg

Molecular Weight: 100,000 daltons

Quality Assurance: No contaminating exoglycosidase or proteolytic activity could be detected.

Quality Controls

Glycosidase Assays:

50 units of β -N-Acetyl-hexosaminidase, were incubated with 0.1 mM of fluorescently-labeled oligosaccharides and glycopeptides, in a 10 μ l reaction for 20 hours at 37°C. The reaction products were analyzed by TLC for digestion of substrate.

No other glycosidase activities were detected (ND) with the following substrates:

α -Fucosidase:

Fuc α 1-2Gal β 1-4Glc-AMC Gal β 1-4 (Fuc α 1-3)GlcNAc β 1-3Gal β 1-4Glc-AMC ND

β -Galactosidase:

Galβ1-3GlcNAcβ1-4Galβ1-4Glc-AMC ND

α -Galactosidase:

 $Gal\alpha 1-3Gal\beta 1-4Gal\alpha 1-3Gal-AMC$ ND

α -Neuraminidase:

Neu5Ac α 2-3Gal β 1-3GlcNAc β 1-3Gal β 1-4Glc-AMC ND

α -Mannosidase:

Manα1-3Manβ1-4GlcNAc-AMCManα1-6Manα1-6(Manα1-3)Man-AMC

β-Glucosidase:

Glcβ1-4Glcβ1-4Glc-AMC ND

(See other side)

CERTIFICATE OF ANALYSIS

ND

α -Fucosidase:

Fuc α 1-2Gal β 1-4Glc-AMC Gal β 1-4 (Fuc α 1-3)GlcNAc β 1-3Gal β 1-4Glc-AMC ND

β-Galactosidase:

Galβ1-3GlcNAcβ1-4Galβ1-4Glc-AMC ND

α -Galactosidase:

 $\mathsf{Gal}\alpha\mathsf{1}\text{-}\mathsf{3}\mathsf{Gal}\beta\mathsf{1}\text{-}\mathsf{4}\mathsf{Gal}\alpha\mathsf{1}\text{-}\mathsf{3}\mathsf{Gal}\text{-}\mathsf{AMC} \qquad \qquad \mathsf{ND}$

α -Neuraminidase:

Neu5Ac α 2-3Gal β 1-3GlcNAc β 1-3Gal β 1-4Glc-AMC ND

α -Mannosidase:

 $\text{Man}\alpha\text{1-3Man}\beta\text{1-4GlcNAc-AMC}$

 $Man\alpha 1-6Man\alpha 1-6(Man\alpha 1-3)Man-AMC$ ND

B-Glucosidase:

Glcβ1-4Glcβ1-4Glc-AMC ND

(See other side)

CERTIFICATE OF ANALYSIS

 β -Xylosidase:

XyIβ1-4XyIβ1-4XyIβ1-4XyI-AMC ND

 β -Mannosidase:

Manβ1-4Manβ1-4Man-AMC ND

Endo F₁, F₂, H:

Dansylated invertase high mannose. ND

Endo F₂, F₃:

Dansylated fibrinogen biantennary. ND

PNGase F:

Fluoresceinated fetuin triantennary. ND

Protease Assay: After incubation of 50 units of $\beta\textsc{-N-Acetyl-hexosaminidase}_t$ with 0.2 nmol of a standard mixture of proteins in a 20 μl reaction, for 20 hours at 37°C, no proteolytic activity could be detected by SDS-PAGE.

*Note: Non-branched oligosaccharides only.

Page 2 (P0721)

β-Xvlosidase:

XyIβ1-4XyIβ1-4XyI-AMC ND

 β -Mannosidase:

Manβ1-4Manβ1-4Man-AMC ND

Endo F₁, F₂, H:

Dansylated invertase high mannose. ND

Endo F_a, F_a:

Dansylated fibrinogen biantennary. ND

PNGase F:

Fluoresceinated fetuin triantennary. ND

Protease Assay: After incubation of 50 units of β -N-Acetyl-hexosaminidase, with 0.2 nmol of a standard mixture of proteins in a 20 μ l reaction, for 20 hours at 37°C, no proteolytic activity could be detected by SDS-PAGE.

*Note: Non-branched oligosaccharides only.

References:

- 1. Robbins, P. et al. (1992) *Gene* 111, 69–76.
- 2. Guan, C. and Wong, S. New England Biolabs Inc., unpublished results.
- 3. Wong-Madden, S.T. and Landry, D. (1995) *Glycobiology* 5, 19–28.

References:

- 1. Robbins, P. et al. (1992) *Gene* 111, 69–76.
- 2. Guan, C. and Wong, S. New England Biolabs Inc., unpublished results.
- 3. Wong-Madden, S.T. and Landry, D. (1995) *Glycobiology* 5, 19–28.