

# BbvI



1-800-632-7799  
info@neb.com  
www.neb.com



R0173S 036121014101

## R0173S



300 units 2,000 U/ml Lot: 0361210

RECOMBINANT Store at -20°C Exp: 10/14

### Recognition Site:

5'...GCAGC(N)<sub>8</sub>▼...3'  
3'...CGTCG(N)<sub>12</sub>▲...5'

**Source:** An *E. coli* strain that carries the cloned BbvI gene from *Bacillus brevis* (ATCC 9999)

Supplied in: 200 mM NaCl, 10 mM Tris-HCl (pH 7.4), 0.1 mM EDTA, 1 mM dithiothreitol, 200 µg/ml BSA and 50% glycerol.

[More Units](#)

# BbvI



1-800-632-7799  
info@neb.com  
www.neb.com



R0173S 036121014101

## R0173S



300 units 2,000 U/ml Lot: 0361210

RECOMBINANT Store at -20°C Exp: 10/14

### Recognition Site:

5'...GCAGC(N)<sub>8</sub>▼...3'  
3'...CGTCG(N)<sub>12</sub>▲...5'

**Source:** An *E. coli* strain that carries the cloned BbvI gene from *Bacillus brevis* (ATCC 9999)

Supplied in: 200 mM NaCl, 10 mM Tris-HCl (pH 7.4), 0.1 mM EDTA, 1 mM dithiothreitol, 200 µg/ml BSA and 50% glycerol.

[More Units](#)

### Reagents Supplied with Enzyme:

10X NEBuffer 2

### Reaction Conditions:

1X NEBuffer 2.  
Incubate at 37°C.

### 1X NEBuffer 2:

50 mM NaCl  
10 mM Tris-HCl  
10 mM MgCl<sub>2</sub>  
1 mM dithiothreitol  
pH 7.9 @ 25°C

**Unit Definition:** One unit is defined as the amount of enzyme required to digest 1 µg of pBR322 DNA in 1 hour at 37°C in a total reaction volume of 50 µl.

**Diluent Compatibility:** Diluent Buffer B  
300 mM NaCl, 10 mM Tris-HCl, 0.1 mM EDTA,  
1 mM DTT, 500 µg/ml BSA and 50% glycerol  
(pH 7.4 @ 25°C).

### Reagents Supplied with Enzyme:

10X NEBuffer 2

### Reaction Conditions:

1X NEBuffer 2.  
Incubate at 37°C.

### 1X NEBuffer 2:

50 mM NaCl  
10 mM Tris-HCl  
10 mM MgCl<sub>2</sub>  
1 mM dithiothreitol  
pH 7.9 @ 25°C

**Unit Definition:** One unit is defined as the amount of enzyme required to digest 1 µg of pBR322 DNA in 1 hour at 37°C in a total reaction volume of 50 µl.

**Diluent Compatibility:** Diluent Buffer B  
300 mM NaCl, 10 mM Tris-HCl, 0.1 mM EDTA,  
1 mM DTT, 500 µg/ml BSA and 50% glycerol  
(pH 7.4 @ 25°C).

### Quality Control Assays

**Ligation:** After 5-fold overdigestion with BbvI, > 95% of the DNA fragments can be ligated with T4 DNA Ligase (at a 5' termini concentration of 1–2 µM) at 16°C. Of these ligated fragments, > 95% can be recut.

**16-Hour Incubation:** A 50 µl reaction containing 1 µg of pBR322 DNA and 3 units of enzyme incubated for 16 hours resulted in the same pattern of DNA bands as a reaction incubated for 1 hour with 1 unit of enzyme.

**Exonuclease Activity:** Incubation of 10 units of enzyme with 1 µg sonicated <sup>3</sup>H DNA (10<sup>5</sup> cpm/µg) for 4 hours at 37°C in 50 µl reaction buffer released < 0.1% radioactivity.

### Quality Control Assays

**Ligation:** After 5-fold overdigestion with BbvI, > 95% of the DNA fragments can be ligated with T4 DNA Ligase (at a 5' termini concentration of 1–2 µM) at 16°C. Of these ligated fragments, > 95% can be recut.

**16-Hour Incubation:** A 50 µl reaction containing 1 µg of pBR322 DNA and 3 units of enzyme incubated for 16 hours resulted in the same pattern of DNA bands as a reaction incubated for 1 hour with 1 unit of enzyme.

**Exonuclease Activity:** Incubation of 10 units of enzyme with 1 µg sonicated <sup>3</sup>H DNA (10<sup>5</sup> cpm/µg) for 4 hours at 37°C in 50 µl reaction buffer released > 0.1% radioactivity.

### Enzyme Properties

#### Activity in NEBuffers:

NEBuffer 1	100%
NEBuffer 2	<b>100%</b>
NEBuffer 3	25%
NEBuffer 4	75%

When using a buffer other than the optimal (supplied) NEBuffer, it may be necessary to add more enzyme to achieve complete digestion.

**Survival in a Reaction:** Not recommended for digest over 1 hour.

**Heat Inactivation:** 10 units of enzyme were inactivated by incubation at 65°C for 20 minutes.

**Note:** Cleaves to leave a 4-base 5' extension.

Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

= Time-Saver™ Qualified (See www.neb.com for details).

CERTIFICATE OF ANALYSIS

### Enzyme Properties

#### Activity in NEBuffers:

NEBuffer 1	100%
NEBuffer 2	<b>100%</b>
NEBuffer 3	25%
NEBuffer 4	75%

When using a buffer other than the optimal (supplied) NEBuffer, it may be necessary to add more enzyme to achieve complete digestion.

**Survival in a Reaction:** Not recommended for digest over 1 hour.

**Heat Inactivation:** 10 units of enzyme were inactivated by incubation at 65°C for 20 minutes.

**Note:** Cleaves to leave a 4-base 5' extension.

Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

= Time-Saver™ Qualified (See www.neb.com for details).

CERTIFICATE OF ANALYSIS