

# NcoI



1-800-632-7799  
info@neb.com  
www.neb.com



R0193S 030121014101

## R0193S



1,000 units    10,000 U/ml    Lot: 0301210  
RECOMBINANT    Store at -20°C    Exp: 10/14

### Recognition Site:

5' . . . C C A T G G . . . 3'  
3' . . . G G T A C C . . . 5'

**Source:** An *E. coli* strain that carries the cloned NcoI gene from *Nocardia corallina* (ATCC 19070)

Also Available In  
High Fidelity (HF™) Format

Supplied in: 50 mM KCl, 10 mM Tris-HCl (pH 7.4), 0.1 mM EDTA, 1 mM dithiothreitol, 200 µg/ml BSA and 50% glycerol.

**Reagents Supplied with Enzyme:**  
10X NEBuffer 3.

**Reaction Conditions:** 1X NEBuffer 3.  
Incubate at 37°C.

**1X NEBuffer 3:**  
100 mM NaCl  
50 mM Tris-HCl  
10 mM MgCl<sub>2</sub>  
1 mM DTT  
pH 7.9 @ 25°C

**Unit Definition:** One unit is defined as the amount of enzyme required to digest 1 µg of λ DNA in 1 hour at 37°C in a total reaction volume of 50 µl.

**Diluent Compatibility:** Diluent Buffer A  
50 mM KCl, 10 mM Tris-HCl, 0.1 mM EDTA, 1 mM dithiothreitol, 200 µg/ml BSA and 50% glycerol (pH 7.4 @ 25°C)

### Quality Control Assays

**Ligation:** After 5-fold overdigestion with NcoI, > 95% of the DNA fragments can be ligated with T4 DNA Ligase (at a 5' termini concentration of 1–2 µM) at 16°C. Of these ligated fragments, > 95% can be recut.

**16-Hour Incubation:** A 50 µl reaction containing 1 µg of DNA and 15 units of enzyme incubated for 16 hours resulted in the same pattern of DNA bands as a reaction incubated for 1 hour with 1 unit of enzyme.

**Exonuclease Activity:** Incubation of 100 units of enzyme with 1 µg sonicated <sup>3</sup>H DNA (10<sup>5</sup> cpm/µg) for 4 hours at 37°C in 50 µl reaction buffer released < 0.24% radioactivity.

**Endonuclease Activity:** Incubation of 15 units of enzyme with 1 µg φX174 RF I DNA for 4 hours at 37°C in 50 µl reaction buffer resulted in < 30% conversion to RF II.

**Blue/White Screening Assay:** This enzyme has been tested to determine the integrity of the DNA ends produced after digestion with an excess of enzyme. An appropriate vector is digested at a unique site within *lacZ<sup>α</sup>* gene with a 10-fold excess of enzyme, ligated, transformed and plated on XGal/IPTG/Amp plates. Successful expression of β-galactosidase is a function of how intact its gene remains after cloning, an intact gene gives rise to a blue colony, an interrupted gene (i.e. degraded DNA end) gives rise to a white colony. Enzymes must produce fewer than 1% white colonies to be Blue/White Certified.

### Enzyme Properties

**Activity in NEBuffers:**  
NEBuffer 1    100%  
NEBuffer 2    100%  
NEBuffer 3    100%  
NEBuffer 4    100%

**Survival in a Reaction:** A minimum of 0.25 unit is required to digest 1 µg of substrate DNA in 16 hours.

(see other side)

CERTIFICATE OF ANALYSIS

# NcoI



1-800-632-7799  
info@neb.com  
www.neb.com



R0193S 030121014101

## R0193S



1,000 units    10,000 U/ml    Lot: 0301210  
RECOMBINANT    Store at -20°C    Exp: 10/14

### Recognition Site:

5' . . . C C A T G G . . . 3'  
3' . . . G G T A C C . . . 5'

**Source:** An *E. coli* strain that carries the cloned NcoI gene from *Nocardia corallina* (ATCC 19070)

Also Available In  
High Fidelity (HF™) Format

Supplied in: 50 mM KCl, 10 mM Tris-HCl (pH 7.4), 0.1 mM EDTA, 1 mM dithiothreitol, 200 µg/ml BSA and 50% glycerol.

**Reagents Supplied with Enzyme:**  
10X NEBuffer 3.

**Reaction Conditions:** 1X NEBuffer 3.  
Incubate at 37°C.

**1X NEBuffer 3:**  
100 mM NaCl  
50 mM Tris-HCl  
10 mM MgCl<sub>2</sub>  
1 mM DTT  
pH 7.9 @ 25°C

**Unit Definition:** One unit is defined as the amount of enzyme required to digest 1 µg of λ DNA in 1 hour at 37°C in a total reaction volume of 50 µl.

**Diluent Compatibility:** Diluent Buffer A  
50 mM KCl, 10 mM Tris-HCl, 0.1 mM EDTA, 1 mM dithiothreitol, 200 µg/ml BSA and 50% glycerol (pH 7.4 @ 25°C)

### Quality Control Assays

**Ligation:** After 5-fold overdigestion with NcoI, > 95% of the DNA fragments can be ligated with T4 DNA Ligase (at a 5' termini concentration of 1–2 µM) at 16°C. Of these ligated fragments, > 95% can be recut.

**16-Hour Incubation:** A 50 µl reaction containing 1 µg of DNA and 15 units of enzyme incubated for 16 hours resulted in the same pattern of DNA bands as a reaction incubated for 1 hour with 1 unit of enzyme.

**Exonuclease Activity:** Incubation of 100 units of enzyme with 1 µg sonicated <sup>3</sup>H DNA (10<sup>5</sup> cpm/µg) for 4 hours at 37°C in 50 µl reaction buffer released < 0.24% radioactivity.

**Endonuclease Activity:** Incubation of 15 units of enzyme with 1 µg φX174 RF I DNA for 4 hours at 37°C in 50 µl reaction buffer resulted in < 30% conversion to RF II.

**Blue/White Screening Assay:** This enzyme has been tested to determine the integrity of the DNA ends produced after digestion with an excess of enzyme. An appropriate vector is digested at a unique site within *lacZ<sup>α</sup>* gene with a 10-fold excess of enzyme, ligated, transformed and plated on XGal/IPTG/Amp plates. Successful expression of β-galactosidase is a function of how intact its gene remains after cloning, an intact gene gives rise to a blue colony, an interrupted gene (i.e. degraded DNA end) gives rise to a white colony. Enzymes must produce fewer than 1% white colonies to be Blue/White Certified.

### Enzyme Properties

**Activity in NEBuffers:**  
NEBuffer 1    100%  
NEBuffer 2    100%  
NEBuffer 3    100%  
NEBuffer 4    100%

**Survival in a Reaction:** A minimum of 0.25 unit is required to digest 1 µg of substrate DNA in 16 hours.

(see other side)

CERTIFICATE OF ANALYSIS

**Heat Inactivation:** 35 units of enzyme were inactivated by incubation at 65°C for 20 minutes.


**Notes:** Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Conditions of low ionic strength, high enzyme concentration, glycerol concentration > 5%, or pH > 8.0 may result in star activity.

**Companion Products Sold Separately:**

NcoI-HF™	
#R3193S	1,000 units
#R3193L	5,000 units
#R3193T	1,000 units
#R3193M	5,000 units

NcoI-HF™ RE-Mix™	
#R5193S	50 reactions

 = Time-Saver™ Qualified (See [www.neb.com](http://www.neb.com) for details).

**Heat Inactivation:** 35 units of enzyme were inactivated by incubation at 65°C for 20 minutes.


**Notes:** Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Conditions of low ionic strength, high enzyme concentration, glycerol concentration > 5%, or pH > 8.0 may result in star activity.

**Companion Products Sold Separately:**

NcoI-HF™	
#R3193S	1,000 units
#R3193L	5,000 units
#R3193T	1,000 units
#R3193M	5,000 units

NcoI-HF™ RE-Mix™	
#R5193S	50 reactions

 = Time-Saver™ Qualified (See [www.neb.com](http://www.neb.com) for details).