

# PhoI



1-800-632-7799  
info@neb.com  
www.neb.com



R0705S 002120814081

**R0705S**

**200 units**      **2,000 U/ml**      **Lot: 0021208**  
**RECOMBINANT**    **Store at -20°C**    **Exp: 8/14**

#### Recognition Site:

5'...GG<sup>▼</sup>CC...3'  
3'...CC<sup>▲</sup>GG...5'

**Source:** An *E. coli* strain that carries the cloned PhoI gene from *Pyrococcus horikoshii* OT3 (Y. Kawarabayasi)

Supplied in: 50 mM KCl, 10 mM Tris-HCl (pH 7.4), 0.1 mM EDTA, 1 mM dithiothreitol, 200 µg/ml BSA and 50% glycerol.

**Reagents Supplied with Enzyme:**  
10X NEBuffer 3

**Reaction Conditions:** 1X NEBuffer 3.  
**Incubate at 75°C.**

**1X NEBuffer 3:**  
100 mM NaCl  
50 mM Tris-HCl  
10 mM MgCl<sub>2</sub>  
1 mM dithiothreitol  
pH 7.9 @ 25°C

**Unit Definition:** One unit is defined as the amount of enzyme required to digest 1 µg of λ DNA in 1 hour at 75°C in a total reaction volume of 50 µl.

**Diluent Compatibility:** Diluent Buffer A  
50 mM KCl, 10 mM Tris-HCl, 0.1 mM EDTA, 1 mM dithiothreitol, 200 µg/ml BSA and 50% glycerol (pH 7.4 @ 25°C).

#### Quality Control Assays

**Ligation:** After 2-fold overdigestion with PhoI, approximately 50% of the DNA fragments can be ligated with T4 DNA Ligase (at a 5' termini concentration of 1–2 µM) at 16°C. Of these ligated fragments, approximately 50% can be recut.

**16-Hour Incubation:** A 50 µl reaction containing 1 µg of λ DNA and 20 units of enzyme incubated for 16 hours resulted in the same pattern of DNA bands as a reaction incubated for 1 hour with 1 unit of enzyme.

**Exonuclease Activity:** Incubation of 20 units of enzyme with 1 µg sonicated [<sup>3</sup>H] DNA (10<sup>5</sup> cpm/µg) for 4 hours at 37°C\* in 50 µl reaction buffer released < 0.1% radioactivity.

\*This quality control was performed at 37°C to detect any *E. coli* contaminants which are not reactive at 75°C.

#### Enzyme Properties

##### Activity in NEBuffers:

NEBuffer 1    50%  
NEBuffer 2    50%  
NEBuffer 3    **100%**  
NEBuffer 4    75%

When using a buffer other than the optimal (supplied) NEBuffer, it may be necessary to add more enzyme to achieve complete digestion.

**Survival in a Reaction:** A minimum of 0.13 unit is required to digest 1 µg of substrate DNA in 16 hours.

**Heat Inactivation:** No

**Note:** PhoI is a highly thermostable restriction enzyme that can survive temperatures as high as 95°C.

Impaired by some combinations of overlapping *dcm* methylation.

Cleavage of mammalian genomic DNA is impaired by some combinations of overlapping CpG methylation.

Incubation at 37°C results in no activity.

#### Companion Products:

dam-/dcm- Competent *E. coli*  
#C2925H      20 transformation reactions  
#C2925I      24 transformation reactions

CERTIFICATE OF ANALYSIS

# PhoI



1-800-632-7799  
info@neb.com  
www.neb.com



R0705S 002120814081

**R0705S**

**200 units**      **2,000 U/ml**      **Lot: 0021208**  
**RECOMBINANT**    **Store at -20°C**    **Exp: 8/14**

#### Recognition Site:

5'...GG<sup>▼</sup>CC...3'  
3'...CC<sup>▲</sup>GG...5'

**Source:** An *E. coli* strain that carries the cloned PhoI gene from *Pyrococcus horikoshii* OT3 (Y. Kawarabayasi)

Supplied in: 50 mM KCl, 10 mM Tris-HCl (pH 7.4), 0.1 mM EDTA, 1 mM dithiothreitol, 200 µg/ml BSA and 50% glycerol.

**Reagents Supplied with Enzyme:**  
10X NEBuffer 3

**Reaction Conditions:** 1X NEBuffer 3.  
**Incubate at 75°C.**

**1X NEBuffer 3:**  
100 mM NaCl  
50 mM Tris-HCl  
10 mM MgCl<sub>2</sub>  
1 mM dithiothreitol  
pH 7.9 @ 25°C

**Unit Definition:** One unit is defined as the amount of enzyme required to digest 1 µg of λ DNA in 1 hour at 75°C in a total reaction volume of 50 µl.

**Diluent Compatibility:** Diluent Buffer A  
50 mM KCl, 10 mM Tris-HCl, 0.1 mM EDTA, 1 mM dithiothreitol, 200 µg/ml BSA and 50% glycerol (pH 7.4 @ 25°C).

#### Quality Control Assays

**Ligation:** After 2-fold overdigestion with PhoI, approximately 50% of the DNA fragments can be ligated with T4 DNA Ligase (at a 5' termini concentration of 1–2 µM) at 16°C. Of these ligated fragments, approximately 50% can be recut.

**16-Hour Incubation:** A 50 µl reaction containing 1 µg of λ DNA and 20 units of enzyme incubated for 16 hours resulted in the same pattern of DNA bands as a reaction incubated for 1 hour with 1 unit of enzyme.

**Exonuclease Activity:** Incubation of 20 units of enzyme with 1 µg sonicated [<sup>3</sup>H] DNA (10<sup>5</sup> cpm/µg) for 4 hours at 37°C\* in 50 µl reaction buffer released < 0.1% radioactivity.

\*This quality control was performed at 37°C to detect any *E. coli* contaminants which are not reactive at 75°C.

#### Enzyme Properties

##### Activity in NEBuffers:

NEBuffer 1    50%  
NEBuffer 2    50%  
NEBuffer 3    **100%**  
NEBuffer 4    75%

When using a buffer other than the optimal (supplied) NEBuffer, it may be necessary to add more enzyme to achieve complete digestion.

**Survival in a Reaction:** A minimum of 0.13 unit is required to digest 1 µg of substrate DNA in 16 hours.

**Heat Inactivation:** No

**Note:** PhoI is a highly thermostable restriction enzyme that can survive temperatures as high as 95°C.

Impaired by some combinations of overlapping *dcm* methylation.

Cleavage of mammalian genomic DNA is impaired by some combinations of overlapping CpG methylation.

Incubation at 37°C results in no activity.

#### Companion Products:

dam-/dcm- Competent *E. coli*  
#C2925H      20 transformation reactions  
#C2925I      24 transformation reactions

CERTIFICATE OF ANALYSIS