

**Total Protein Cell Lysis Buffer**

<u>Code</u>	<u>Description</u>	<u>Size</u>
M329-10ML	Total Cell Lysis Buffer <i>For Mammalian Tissue Culture Cells</i> Buffer is sufficient for 20 protein extractions of $1 \times 10^6$ cells	10 mls

**General Information:**


AMRESCO's Total Protein Cell Lysis Buffer is a non-ionic detergent lysis buffer optimized for the isolation of total proteins from adherent or non adherent cells in culture. The Total Protein Cell Lysis Buffer enables efficient solubilization of the plasma and intracellular membranes, breaks weak intermolecular bonds, and solubilizes most of the commonly studied protein antigens. Total Protein Cell Lysis Buffer reduces protein denaturation, protein complex disruption and loss of enzymatic activity since it contains less detergent than a typical RIPA lysis buffer.

**Storage/Stability:**

Store product cold, 2-8°C

**Application Disclaimer**


*For Research Use Only.  
Not for Therapeutic or Diagnostic Use.*



**Procedure:**
**Materials not supplied:**

Protease Inhibitors

Phosphate Buffered Saline (PBS) (ice cold)

 **Note:** All procedures should be performed on ice in a cold room with ice cold reagents to reduce proteolysis, dephosphorylation and denaturation.

**Add the protease inhibitor cocktail to AMRESCO's Total Protein Cell Lysis Buffer so that the final concentration of protease inhibitors is 1X.**

**Protocol below is designed for approximately 1x10<sup>6</sup> cells**

**Cell Collection:**

1. Transfer cells from tissue culture flask to an appropriately sized centrifuge tube.
2. Centrifuge at 2,000rpm for 5 minutes at 4°C.
3. Remove the media, and resuspend the cell pellet in 10ml cold 1X PBS.
4. Centrifuge at 2,000rpm for 5 minutes at 4°C.
5. Remove the supernatant and resuspend the pellet in 1ml cold 1X PBS. Transfer this resuspended pellet in PBS to a 1.5ml microcentrifuge tube.
6. Spin 1 minute at 2,000 rpm at 4°C.

**Cell Lysis:**

1. Discard the supernatant, and resuspend the pellet in 20X the packed cell volume with Total Protein Cell Lysis Buffer containing protease inhibitors. Vortex vigorously for 10 seconds.
2. Incubate the cell suspension on ice with shaking for 30 minutes with periodic vortexing.
3. Vortex vigorously for 30 seconds.
4. Centrifuge at 14,000xg for 20 minutes at 4°C.
5. Transfer the supernatant containing total cell protein into a separate microcentrifuge tube.
6. The total cell protein should be stored frozen at 20°C until needed.

**Related Products**
**Code**

M221-1ml

M222-1ml

M250-1ml

E504-100ml

E504-500ml

N655-50ml

N655-6x5ml

N673-50ml

N673-6x5ml

0260-25g

0260-50g

N182-5x10ml

K952-100ml

**Product**

Protease Inhibitor Cocktail 100X, General Use

Protease Inhibitor Cocktail 100X, General Use with EDTA

Protease Inhibitor Cocktail Mammalian Phosphate Buffered Saline (PBS), 1X Sterile Solution

SeraFree™ Cryopreservation Media (RPMI)

SeraFree™ DMEM Cryopreservation Media

Trypsin 1:300

DMSO, Ultra Pure Grade

Penicillin/Streptomycin, 100X

*Tissue Culture Grade*

ZY0484 04//09 ©Copyright 2009 by AMRESCO Inc. All Rights Reserved. AMRESCO® is a registered trademark of AMRESCO Inc.



**Corporate Headquarters**

AMRESCO Inc.  
30175 Solon Industrial Pkwy.  
Solon, OH 44139 USA  
Tel: (440) 349-1199  
Toll free: 1-800-829-2805 (US & Canada)  
Fax: (440) 349-2441  
e-mail: [info@amresco-inc.com](mailto:info@amresco-inc.com)

**To Order (USA & Canada)**

Toll free: 1-800-829-2805  
Toll free fax: 1-800-326-3733  
Fax: (440) 349-3255  
e-mail: [dom-sales@amresco-inc.com](mailto:dom-sales@amresco-inc.com)

**International Orders:**

Tel: (440) 349-1199  
Fax: (440) 349-2441  
e-mail: [int-sales@amresco-inc.com](mailto:int-sales@amresco-inc.com)

**Technical Support**

Toll free: 1-800-610-2789 (USA & Canada)  
Fax: (440) 349-0235  
e-mail: [techinquiry@amresco-inc.com](mailto:techinquiry@amresco-inc.com)

**[www.amresco-inc.com](http://www.amresco-inc.com)**

