



EZ-Vision® In-Gel

<u>Code</u>	<u>Description</u>	<u>Size</u>
N391-0.5ML	EZ-Vision® In-Gel Solution, 10,000X <i>Sufficient reagent for approximately 100 mini-gels</i>	0.5 ml
N391-5MLDRP	EZ-Vision® In-Gel, Dropper Bottle <i>Sufficient reagent for approximately 100 mini-gels</i>	5 ml
N391-15MLDRP	EZ-Vision® In-Gel, Dropper Bottle <i>Sufficient reagent for approximately 300 mini-gels</i>	15 ml
N391-SAMPLE	EZ-Vision® In-Gel, Dropper Bottle <i>Sufficient reagent for approximately 10 mini-gels</i>	0.05 ml

General Information:

EZ-Vision® In-Gel Solution, 10,000X is a non-mutagenic and non-toxic fluorescent DNA dye that is used as an in-gel (precast) stain to visualize DNA in agarose gels. EZ-Vision® In-Gel Solution is an excellent replacement for the often used but hazardous ethidium bromide for DNA applications, reducing harmful chemical exposure and eliminating hazardous waste disposal costs. As an in-gel preparation for electrophoresis, EZ-Vision® In-Gel Solution is simply added to the molten agarose prior to gel casting, allowing instant visualization of DNA bands after the run when the gel is exposed to the uv illumination of a common transilluminator. Adding EZ-Vision® In-Gel Solution to hot agarose does not reduce performance. Alternatively, EZ-Vision® In-Gel Solution can be used to stain the gel in buffer after electrophoresis.

When EZ-Vision® In-Gel is used as a pre-cast solution; the dye will run in the opposite direction of the DNA migration during electrophoresis. Therefore, the very bottom of the gel may experience a lower concentration of dye. EZ-Vision® In-Gel Solution, 10,000X is compatible with all standard uv transilluminators with a SYBR® Green (500 - 600 nm) filter recommended for documentation. EZ-Vision® In-Gel Solution is also available in two convenient dropper bottle sizes to minimize pipetting.

Storage/Stability:

Store at 2-8 °C for 1 year.

Application Disclaimer

*For research use only.
Not for therapeutic or diagnostic use.*



Procedure

For In-Gel Use

1. Dilute the EZ-Vision® In-Gel Solution 10,000X into the pre-melted agarose solution at a concentration of 1:10,000. For example, add 5 µL of the EZ-Vision® In-Gel Solution 10,000X to 50 mL of agarose gel solution. The EZ-Vision® In-Gel Solution can be added while the solution is still hot. If using the EZ-Vision® In-Gel, Dropper Bottle, add 1 - 2 drops of solution per 50 ml gel.
2. Cast the gel and allow it to solidify.
3. Load sample and run according to standard procedure.
4. After the run, remove gel and place on a UV transilluminator to immediately visualize bands. DNA bands will emit a whitish-blue fluorescence against a dark background using a standard transilluminator (254 or 302 nm). The optimal visualization filter is an SYBR® Green filter; however, an ethidium bromide filter can be used.

For Post-Stain Use

1. Run gels according to standard procedure.
2. Dilute the EZ-Vision® In-Gel Solution 10,000X into 100 mM NaCl to make a 2.5X staining solution. For example, add 25 µL to 100 mL of 100 mM NaCl (100 mM NaCl solution can be prepared by adding 10 mL of a 1 M NaCl solution to 90 mL of H₂O). If using the EZ-Vision® In-Gel, Dropper Bottle, add 1 - 2 drops of solution per 100 ml of 100 mM NaCl.
3. Place the gel in a suitable staining container and add enough 2.5X staining solution to completely submerge the gel.
4. Stain while agitating the gel at room temperature for 20-30 minutes. The optimal staining time may vary based on the gel thickness and concentration. The staining solution can be reused at least twice.
5. Destain the gel with two 10 minute exchanges of water while agitating the gel to remove additional background signal.
6. After destaining, remove gel and place on a UV transilluminator to immediately visualize bands. DNA bands will emit a whitish-blue fluorescence against a dark background using a standard transilluminator (254 or 302 nm). The optimal visualization filter is an SYBR® Green filter; however, an ethidium bromide filter can be used.

Frequently Asked Questions

Questions	Answers
Which filter is recommended for visualizing DNA stained with EZ-Vision® In-Gel Solution?	A SYBR® Green filter (500 – 600 nm) is optimal, although an ethidium bromide filter (550 – 640 nm) may also be used.
Which downstream applications are compatible with usage of EZ-Vision® In-Gel Solution?	Restriction digests, PCR, ligation and *sequencing. <i>(*In some cases, run lengths may be shorter with EZ-Vision® stained DNA compared to ethidium bromide stained DNA.)</i>
How sensitive is EZ-Vision® In-Gel staining?	EZ-Vision® In-Gel Solution can detect 6 ng DNA above 500 bp and 12 ng DNA at 50 bp.
Is migration of DNA affected when stained with EZ-Vision® In-Gel Solution?	No, EZ-Vision® In-Gel stained DNA migrates at the same rate as unstained DNA.
What is the duration of fluorescence emission upon UV exposure?	EZ-Vision® In-Gel stained DNA retains 50% of the original fluorescent intensity after 45 minutes of continuous UV exposure.
What is the shelf-life of EZ-Vision® In-Gel Solution?	EZ-Vision® In-Gel Solution, 10,000X and EZ-Vision® In-Gel, Dropper Bottle are stable for at least a year at 2-8°C when stored protected from light.
Does loading buffer need to be added to DNA samples for gels stained with EZ-Vision® In-Gel Solution?	Yes, loading buffer is necessary. There is no need, however, to add ethidium bromide or any other DNA dye.
Why can't I see my DNA?	<ol style="list-style-type: none"> 1. Wrong filter was used for detection. A SYBR® Green filter (500-600 nm) is recommended. If using an ethidium bromide filter (550-640 nm), longer exposures may be necessary. 2. A Dark Reader was used to visualize DNA. EZ-Vision® In-Gel Solution is incompatible with visualization by a Dark Reader. Use standard UV transillumination. 3. Not enough DNA was loaded on the gel. Load at

least 100 ng DNA per lane. You may need to optimize loading amounts for each sample. EZ-Vision is slightly less sensitive than ethidium bromide for smaller DNA fragments.

4. Gel running conditions were not optimized. The DNA dye in EZ-Vision® In-Gel Solution may dissociate from DNA samples with long run times. Gel running at 8 V/cm for 20 minutes is recommended.

Corporate Headquarters

AMRESCO, LLC.
6681 Cochran Road
Solon, OH 44139 USA
Tel: (440) 349-1199
Toll free: 1-800-448-4442 (US & Canada)
Fax: (440) 349-2441
e-mail: info@amresco-inc.com

To Order (USA & Canada)

www.amresco-inc.com
Toll free: 1-800-448-4442
Toll free fax: 1-800-326-3733
Fax: (440) 349-3255
e-mail: dom-sales@amresco-inc.com

International Orders:

Tel: (440) 349-1199
Fax: (440) 349-2441
e-mail: int-sales@amresco-inc.com

Technical Support

Toll free: 1-800-610-2789 (USA & Canada)
Fax: (440) 349-0235
e-mail: techinquiry@amresco-inc.com

Related Products

Code	Product
0710-100G 0710-500G	Agarose I™
J234-100G J234-500G	Agarose SFR™
E776-100G E776-500G	Agarose 3:1 HRB™
0658-1L 0658-4L 0658-20L	TBE, 10X Liquid Concentrate
K915-1.6L K915-4L K915-20L	TAE (Tris-Acetate-EDTA) Buffer, 50X
N550-300UL	Ready Ladder™ 100 bp DNA Marker
N551-600UL	Ready Ladder™ 1 kb DNA Marker
K811-50RXN	PCR DNA Marker™ With Loading Dye

ZY0521 Rev. 2 01/2012 ©Copyright 2008, 2009, 2010 by AMRESCO, LLC. All Rights Reserved. AMRESCO® is a registered trademark of AMRESCO, LLC.