



PRODUCT INFORMATION

**Thermo Scientific  
Replicator RNAi Kit**

#F-610

Lot \_                      Expiry Date \_

Store at -20°C

[www.thermoscientific.com/onebio](http://www.thermoscientific.com/onebio)

Rev.1

## CONTENTS

COMPONENTS OF THE KIT .....	4
STORAGE .....	4
1. INTRODUCTION .....	4
2. APPLICATIONS .....	4
3. METHOD OVERVIEW .....	5
4. TEMPLATE PREPARATION BY PCR .....	5
4.1 Target preparation.....	5
4.2 Primer desing .....	5
4.3 PCR.....	5
4.3.1 Cycling conditions.....	6
4.4 Purification .....	7
5. dsRNA PRODUCTION SETUP .....	8
6. PURIFICATION OF dsRNA.....	10
7. QUANTIFICATION AND STORAGE OF THE dsRNA .....	11
8. TROUBLESHOOTING .....	12
9. REFERENCES .....	13

## COMPONENTS OF THE KIT

Component	#F-610 40 rxns
Phusion High-Fidelity DNA Polymerase (2 U/ $\mu$ L)	25 $\mu$ L
10 mM dNTP mix	20 $\mu$ L
5X Phusion HF buffer	1.5 mL
Lambda control DNA (0.5 ng/ $\mu$ L)	20 $\mu$ L
Control Primer mix (25 $\mu$ M each)	20 $\mu$ L
Phi6 RNA Replicase	60 $\mu$ L
T7 RNA Polymerase	60 $\mu$ L
Pyrophosphatase	40 $\mu$ L
10X dsRNA Synthesis Buffer	250 $\mu$ L
5X NTP mix	400 $\mu$ L
50 mM MnCl <sub>2</sub>	500 $\mu$ L
8 M LiCl	2 mL

### STORAGE

The Thermo Scientific Replicator RNAi Kit is shipped in gel ice. Upon arrival, store the kit components at -20°C. The 10X dsRNA Synthesis Buffer contains DTT, which may lose its reducing capability when stored or handled inappropriately. LiCl solution can also be stored at room temperature.

### 1. INTRODUCTION

Replicator™ RNAi Kit is a highly efficient and specific tool for generating double-stranded RNA for gene silencing experiments. The kit utilizes a new technology where dsRNA is produced by combined *in vitro* transcription and replication reactions. Both reactions take place in a single incubation step. The DNA template for the dsRNA synthesis is produced by PCR using Thermo Scientific Phusion High-Fidelity DNA Polymerase. Phusion™ DNA Polymerase's high speed, fidelity and robustness ensure correct starting material for the dsRNA production and enable amplification of long templates up to 10 kb. The PCR product is transcribed into single-stranded RNA by T7 RNA Polymerase. The ssRNA is then replicated into double-stranded RNA by Phi6 RNA Replicase<sup>1-3</sup> in the same incubation step. Due to the unique combination of complementary polymerase activities, the produced dsRNA is perfectly duplexed. The reaction yield is high. The kit contains sufficient reagents for producing up to 2.5 mg of dsRNA in a total volume of 2 mL. The yield depends on the target sequence.

### 2. APPLICATIONS

dsRNA produced with this kit is a suitable substrate for enzymes such as RNase III and Dicer, which cut the dsRNA into small interfering RNAs and thereby produce a pool of siRNA targeted for a chosen sequence. Such siRNA pools have successfully been used in gene silencing experiments<sup>4-6</sup>.

### 3. METHOD OVERVIEW

The dsRNA is synthesized by combined *in vitro* transcription and replication from a DNA template.

The DNA template for the dsRNA synthesis is produced by PCR using Phusion High-Fidelity DNA Polymerase. The PCR primers are designed so that the resulting PCR product contains the target sequence flanked by T7 promoter sequence at one end and Phi6 RNA Replicase promoter sequence at the other end.

The PCR product is purified and transcribed into single-stranded RNA by T7 RNA Polymerase. This ssRNA is replicated into double-stranded RNA by Phi6 RNA Replicase in the same incubation step.

### 4. TEMPLATE PREPARATION BY PCR

#### 4.1 Target preparation

With the Replicator RNAi Kit, targets ranging between 100 bp and several kilobases in length can be amplified. The template DNA for the PCR can be either genomic DNA, plasmid DNA or cDNA produced using reverse transcriptase enzyme.

#### 4.2 Primer desing

The primers have to be designed separately for each target sequence. RNA polymerase promoter sequences needed for the dsRNA synthesis are added to both sides of the target sequence using PCR. The PCR primers are designed so that they contain RNA polymerase promoter sequences at their 5' ends. Thus in the PCR product, RNA polymerase promoter sequences flank the target sequence. In addition to promoter sequences, each primer should contain 17-22 nucleotides of target gene-specific sequence at the 3' end. The PCR amplification primers are illustrated bellow. The promoter sequences are underlined. The first nucleotides to be incorporated into the final dsRNA product at each end are marked in red. Target genespecific sequences at the 3' end of the primers are denoted by N<sub>(17-22)</sub>.

#### Amplification primers

Forward primer 5' TAATCGACTCACTATAGGGN<sub>(17-22)</sub> 3'

Reverse primer 5' GGAAAAAAN<sub>(17-22)</sub> 3'

#### 4.3 PCR

PCR reactions should be set up on ice. It is critical that the PhusionDNA Polymerase is the last component added to the PCR mixture, since the enzyme exhibits 3'→5' exonuclease activity that can degrade primers in the absence of dNTPs. Phusion DNA Polymerase produces blunt end DNA products.

The optimal amount of enzyme depends on the length of the PCR product and the amount of template used. Normally, 1 unit of Phusion DNA Polymerase per 50 µL reaction volume gives good results. Optimal enzyme amount can range from 0.5–2 units per 50 µL reaction depending on the amplicon length and difficulty. Do not exceed 2 units per 50 µL reaction. For optimal results, use 200 µM of high quality dNTPs included in the kit. The amount of template depends on the complexity of DNA. General guidelines for low complexity DNA (e.g. plasmid, lambda or BAC DNA) are: 1 pg–10 ng per 50 µL reaction volume. For high complexity genomic DNA, the amount of DNA template should be 50–200 ng per 50 µL reaction volume. If cDNA synthesis reaction mixture is used as a source of template, the volume of the template should not exceed 10% of the final PCR reaction volume. PCR and dsRNA production efficiency can be controlled using the accompanied I control DNA and control primers (see Tables 1 and 2).

Table 1. Pipetting instructions for the PCR reaction.

Component (in order)	50 µL rxn	Final conc.
H <sub>2</sub> O	add to 50 µL	
5X Phusion HF buffer	10 µL	1X
10 mM dNTPs	1 µL	200 µM each
Forward primer	X µL	0.5 µM
Reverse primer	X µL	0.5 µM
Template DNA	X µl	
Phusion DNA Polymerase	0.5 µL:	0.02 U/µL

Table 2. Pipetting instructions for the control PCR reaction.

Component	50 µL rxn
H <sub>2</sub> O	33.5 µL
5X Phusion HF buffer	10 µL
10 mM dNTPs	1 µL
Control primer mix (25 µM each)	1 µL
λ control DNA (0.5 ng/µL)	4 µL
Phusion DNA Polymerase	0.5 µL

#### 4.3.1 Cycling conditions

Due to the nature of Phusion DNA Polymerase, optimal reaction conditions may differ from standard enzyme protocols. Phusion DNA Polymerase tends to work better at elevated denaturation and annealing temperatures due to higher salt concentrations in its buffer.

Table 3. Cycling instructions.

Cycle step	Temp.	Time	Cycles
Initial denaturation	98°C	30 s	1
Denaturation	98°C	5-10 s	25-35
Annealing	45-72°C	10-30 s	
Extension	72°C	15-30 s/1 kb	
Final extension	72°C 4°C	5-10 min hold	1

Table 4. Cycling instructions for the control reaction.

Cycle step	Temp.	Time	Cycles
Initial denaturation	98°C	30 s	1
Denaturation	98°C	10 s	34
Annealing, extension	72°C	35 s	
Final extension	72°C 4°C	5 min hold	1

### Initial denaturation

30 s initial denaturation at 98°C is recommended for most templates. Some templates may require longer initial denaturation time and the length of the initial denaturation time can be extended up to 3 minutes.

### Denaturation

Keep the denaturation as short as possible. Usually 5–10 seconds at 98°C is enough for most templates.

### Primer annealing

The optimal annealing temperature for Phusion DNA Polymerase may differ significantly from that of Taq-based polymerases. For best results, calculate the  $T_m$  for the target gene-specific sequence of the primers. Always use the  $T_m$  calculator and instructions on our website ([www.thermoscientific.com/pcrwebtools](http://www.thermoscientific.com/pcrwebtools)) to determine the  $T_m$  values of primers and optimal annealing temperature.

Phusion DNA Polymerase has the ability to stabilize primer-template hybridization. As a basic rule, for primers >20 nt, anneal for 10–30 seconds at a  $T_m + 3^\circ\text{C}$  of the lower  $T_m$  primer. For primers  $\leq 20$  nt, use an annealing temperature equal to the  $T_m$  of the lower  $T_m$  primer.

If necessary, use a temperature gradient to find the optimal annealing temperature for each template-primer pair combination. The annealing gradient should extend up to the extension temperature (two-step PCR).

A 2-step protocol is recommended when primer  $T_m$  values are at least  $69^\circ\text{C}$  (> 20 nt) or  $72^\circ\text{C}$  ( $\leq 20$  nt) when calculated with our  $T_m$  calculator.

In the 2-step protocol the combined annealing/extension step should be performed at  $72^\circ\text{C}$  even when the primer  $T_m$  is  $>72^\circ\text{C}$ .

### Extension

The extension should be performed at  $72^\circ\text{C}$ . Extension time depends on amplicon length and complexity. For low complexity DNA (e.g. plasmid, lambda or BAC DNA) use extension time 15 s per 1 kb. For high complexity genomic DNA 30 s per 1 kb is recommended.

### 4.4 Purification


Following PCR, the size and purity of the PCR product should be examined on an agarose gel. The PCR product can be purified using either phenol/chloroform extraction followed by standard ethanol precipitation or commercially available DNA purification columns.

Alternatively, the PCR product can be gel-purified if necessary. Resuspend the DNA in RNase-free water or 10 mM Tris-HCl, pH 8.0–8.5. The concentration of the PCR product can be determined by spectrophotometric measurement at 260 nm or estimated from an agarose gel. The expected length of the PCR product from the control reaction is 881 bp.

## 5. dsRNA PRODUCTION SETUP

dsRNA is produced in a single reaction combining DNA transcription and replication of the ssRNA into dsRNA. The dsRNA synthesis reaction is illustrated in Figure 1.

Figure 1. dsRNA synthesis reaction.




The reaction volume can be scaled up from the 50  $\mu\text{L}$  described in Table 5. It is recommended to perform a single 20–50  $\mu\text{L}$  pilot reaction before scale up to confirm all reaction conditions are correct. Use of RNase inhibitor is not normally necessary because ssRNA produced by T7 RNA Polymerase is immediately replicated into more stable dsRNA by Phi6 RNA Replicase.

Table 5. Pipetting instructions for the dsRNA reaction.

Component (in order)	Volume/50 $\mu$ L rxn	Final conc.
RNase-free water	add to 50 $\mu$ L	
10X dsRNA Synthesis Buffer	5 $\mu$ L	1X
5X NTP mix	10 $\mu$ L	1X
DNA template	X $\mu$ L	10-40 ng/ $\mu$ L
MnCl <sub>2</sub> (50 mM)	1.5 $\mu$ L	1.5 mM
Pyrophosphatase	1 $\mu$ L	
T7 RNA Polymerase	1.5 $\mu$ L	
Phi6 RNA Replicase	1.5	

1. Thaw the reagents (10X dsRNA Synthesis Buffer, 5X NTP mix and MnCl<sub>2</sub>) at room temperature immediately before use. Keep the 10X dsRNA Synthesis Buffer and MnCl<sub>2</sub> at room temperature until use. Store the ribonucleotides and enzyme solutions on ice until use.
2. Prepare the dsRNA synthesis reaction by combining the reaction components in the following order. Reaction assembly should be done at room temperature because the dsRNA synthesis buffer contains spermidine that can precipitate nucleic acids on ice. Use 2–5  $\mu$ g of DNA template per 50  $\mu$ L reaction. Set up a separate reaction for each template of interest. It is advisable to perform a concomitant control reaction using the control PCR product to test the reaction conditions.
3. Mix the components thoroughly. Spin briefly to collect the reaction mixture at the bottom of the tube.
4. Incubate at 35°C for 2–4 hours using e.g. a thermal cycler. The reaction can also be allowed to proceed overnight. Overnight incubation doubles the yield of some dsRNAs.
5. It is recommended to examine the dsRNA product by standard agarose gel electrophoresis. Take a 5  $\mu$ L sample from the dsRNA reaction and dilute it 20- to 100-fold before the gel electrophoresis. The gel and electrophoresis buffer should be prepared with RNase-free water. Common DNA markers can be used to estimate the length of the dsRNA product with sufficient accuracy. The control reaction should yield a dsRNA of 864 bp (Figure 2). Occasionally, some ssRNA is visible after the reaction. ssRNA migrates on the gel at approximately half the size of the corresponding dsRNA.

Figure 2. Various gene-specific dsRNA fragments generated using the Replicator RNAi Kit.



PCR-amplified fragments of various sizes were used as template DNA. 0.1  $\mu$ L from each reaction mixture was loaded on an agarose gel.

M – HaeIII-digested Phi174-DNA

- 1 - dsRNA specific to 150 bp sequence of lambda DNA
- 2 - dsRNA specific to 210 bp sequence of Furin gene
- 3 – dsRNA specific to 281 bp sequence of Erk 1/2 kinase gene
- 4 - dsRNA specific to 430 bp sequence of Caspase-3 gene
- 5 - dsRNA specific to 864 bp sequence of lambda DNA

The yield of dsRNA from a 50  $\mu$ L reaction is typically 20–60  $\mu$ g, depending on the target sequence and length. The dsRNA synthesis reaction can be stored at -20°C or -80°C before purification. Storage of unpurified dsRNA may cause the formation of a white precipitate. The precipitate does not contain RNA and can be removed by centrifugation before purification.

## 6. PURIFICATION OF dsRNA

The need for purification depends on the downstream application. Double-stranded RNA can be purified from contaminating reaction components using either selective LiCl precipitation or nuclease digestion. LiCl solution is included in the kit. Both purification methods may require an additional NTP removal step (e.g. gel filtration), depending on the use of the dsRNA product.

### LiCl precipitation

LiCl precipitation selectively removes template DNA and ssRNA from the reaction. The dsRNA synthesis reaction is first precipitated with 2 M LiCl which precipitates all or most of the ssRNA present. The supernatant containing the dsRNA is then exposed to a second precipitation step using 4 M LiCl which efficiently precipitates dsRNA leaving template DNA and most NTPs in the supernatant. Keep the reactions on ice while pipetting.

1. Add 1/3 volume of 8 M LiCl to the dsRNA synthesis reaction in a 1.5 mL microcentrifuge tube and mix. The final LiCl concentration is 2 M.
2. Incubate at -20°C for 30 minutes.
3. Centrifuge at 14 000 rpm in a microcentrifuge for 20 minutes at 4°C.
4. Carefully separate the supernatant from the pellet. The dsRNA remains in the supernatant, while ssRNA is precipitated.
5. Add 1/2 volume of 8 M LiCl to the supernatant. The final LiCl concentration is 4 M. This step precipitates dsRNA.
6. Incubate at -20°C for 30 minutes.
7. Centrifuge at 14 000 rpm for 20 minutes at 4°C.

8. Remove the supernatant. The dsRNA pellet is not necessarily visible. Wash the pellet with 500  $\mu$ L of 70% (v/v) ethanol.
9. Centrifuge at 14 000 rpm for 5 minutes at 4°C.
10. Carefully remove the ethanol without disturbing the pellet.
11. Air-dry the pellet for 5–10 minutes at room temperature and resuspend the dsRNA in 50  $\mu$ L of RNase-free water or TE buffer (10 mM Tris-HCl, pH 8.0, 1 mM EDTA). Be careful not to over-dry the pellet, as it may become difficult to resuspend the dsRNA.

#### Nuclease treatment

Alternatively, the template DNA and ssRNA can be removed from the synthesis reaction by DNase (RNase-free) and RNase digestion. Perform the reactions according to the manufacturer's instructions. Make sure to choose nucleases that do not degrade double-stranded RNA.

#### Gel filtration

After LiCl precipitation or nuclease treatment, the dsRNA may be further purified using commercial gel filtration columns to ensure that all NTPs are removed from the reaction. Removing remaining NTPs will allow accurate spectrophotometric quantitation of dsRNA.

#### Analysis by gel electrophoresis

It is recommended to analyze the purified dsRNA using agarose gel electrophoresis. The size of dsRNA can be estimated comparing to common DNA markers.

### 7. QUANTIFICATION AND STORAGE OF THE dsRNA

The concentration of dsRNA can be determined spectrophotometrically at 260 nm. In order to obtain accurate quantification, the dsRNA has to be free from contaminating NTPs. Measure the absorbance at 260 nm from a 1:4 or 1:20 dilution. For accurate results, make the dilutions using 10 mM Tris-HCl, pH 7.5. The concentration of RNA can be calculated as follows: 1 A260 = 40  $\mu$ g dsRNA per ml. Store the purified RNA at -20°C or -80°C.

#### Purity of RNA

The ratio of absorbances at 260 and 280 nm (A260/A280) provides an estimate of the purity of RNA with respect to contaminants that absorb in the UV region, such as proteins. Pure RNA should exhibit an A260/A280 ratio of 1.9–2.1 in 10 mM Tris-HCl, pH 7.5. The A260/A280 ratio is influenced by pH and since water is not buffered, it is recommended to measure the absorbance of RNA sample in 10 mM Tris-HCl, pH 7.5.

### 8. TROUBLESHOOTING

PCR	
Problem	Cause and Solution
No product at all or low yield	<p>Template concentration may be too low. Use more template.</p> <p>Increase extension time.</p> <p>Increase the number of cycles.</p> <p>Decrease annealing temperature.</p> <p>Optimize enzyme concentration.</p> <p>Denaturation temperature may be too low. Optimal denaturation temperature for most templates is 98°C or higher.</p> <p>Denaturation time may be too long or too short. Optimize denaturation time.</p> <p>Check the purity and concentration of the primers.</p> <p>Check primer design.</p>
Non-specific products - High molecular weight smear	<p>Reduce enzyme concentration.</p> <p>Decrease extension time.</p> <p>Reduce total number of cycles.</p> <p>Increase annealing temperature or try a two-step protocol (see 4.3.1).</p> <p>Perform PCR setup on ice. Transfer the tubes directly to pre-heated PCR block.</p>
Non-specific products - Low molecular weight discrete bands	<p>Increase annealing temperature.</p> <p>Decrease extension time.</p> <p>Decrease enzyme concentration.</p> <p>Titrate template amount.</p> <p>Decrease primer concentration.</p> <p>Design new primers.</p> <p>Perform PCR setup on ice. Transfer the tubes directly to pre-heated PCR block.</p>
<b>dsRNA synthesis reaction</b>	
Low yield	<p>Spermidine precipitates DNA template.</p> <p>Assemble the reaction at room temperature. Spermidine may precipitate if the reaction is assembled on ice.</p> <p><b>Not enough active DTT in the buffer.</b></p> <p>The reducing capacity of DTT decreases gradually in the course of time. Add fresh DTT (10 mM).</p> <p><b>RNA is degraded.</b></p> <p>RNase contamination present in the reaction or DNA template. Be careful to avoid RNase contamination. Optionally, RNase inhibitor can also be included in the reaction.</p> <p><b>Incubation time too short.</b></p> <p>Increase the incubation time.</p> <p><b>Difficult template.</b></p> <p>Some sequences are inherently difficult templates. Design new primers for the target template.</p> <p><b>Insufficient DNA template.</b></p> <p>Check the concentration of the DNA template. Inefficient purification may leave contaminating dNTPs which cause inaccurate spectrophotometric quantitation.</p>

Problem	Cause and Solution
Multiple reaction products	<p>ssRNA. ssRNA migrates at half the size of the corresponding dsRNA on an agarose gel. ssRNA production is sequence specific and may be reduced by designing new primers for the target template.</p> <p><b>Non-specific dsRNA from non-specific PCR products.</b> dsRNA is also synthesized from non-specific PCR products. Check that the PCR reaction does not contain multiple products. Refer to PCR troubleshooting section for PCR optimization. If necessary, purify the PCR product from an agarose gel.</p>
Smear	<p>The reaction products are degraded on the gel. Make sure that the gel is prepared using RNase-free water to avoid RNase contamination.</p> <p><b>dsRNA degraded by nuclease treatment.</b> If dsRNA is purified using nucleases, make sure to choose nucleases that do not degrade dsRNA.</p>
<b>LiCl precipitation</b>	
Low yield	<p>Inefficient precipitation of dsRNA at 4 M LiCl concentration. Incubate the precipitation reaction at -20°C prior each centrifugation (excluding washing step). Skipping the incubation step may decrease the yield slightly, especially with short templates.</p> <p><b>dsRNA pellet lost.</b> Make sure not to disturb the pellet during pipetting to ensure the pellet is not lost.</p>
dsRNA not purified	<p>Inefficient precipitation of ssRNA at 2 M LiCl concentration. Incubate the precipitation reaction at -20°C prior each centrifugation. Skipping the incubation step may cause insufficient precipitation of the ssRNA.</p> <p>At some instances, all ssRNA is not removed by the 2 M LiCl precipitation step. Use RNases to remove the ssRNA.</p>

## 9. REFERENCES

1. Makeyev E.V. and Bamford D.H. EMBO J 19: 124–133, 2000.
2. Makeyev E.V. and Bamford D.H. EMBO J 19: 6275–6284, 2000.
3. Laurila M.R.L. et al. J Biol Chem 277: 17117–17124, 2002.
4. Calegari F. et al. Proc Natl Acad Sci U S A 99: 14236–14240, 2002.
5. Myers J.W. et al. Nat Biotechnol 21: 324–328, 2003.
6. Yang D. et al. Proc Natl Acad Sci U S A 99: 9942–9947, 2002.

## TECHNICAL SUPPORT

EMEA: [ts.molbio.eu@thermofisher.com](mailto:ts.molbio.eu@thermofisher.com)

Americas & APAC: [ts.molbio@thermofisher.com](mailto:ts.molbio@thermofisher.com)

## NOTICE TO PURCHASE:

The purchase price of this product includes a limited, non-transferable license under U.S. and foreign patents owned by BIO-RAD Laboratories, Inc., to use this product. No other license under these patents is conveyed expressly or by implication to the purchaser by the purchase of this product.

## PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively for research purposes and in vitro use only. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to [www.thermoscientific.com/onebio](http://www.thermoscientific.com/onebio) for Material Safety Data Sheet of the product.

© 2012 Thermo Fisher Scientific, Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.