



Instruction Manual

S.N.A.P.™ MidiPrep Kit

A Simple Nucleic Acid Prep

Catalog no. K1910-01

Version F

041702

25-0192

Overview

Introduction

The S.N.A.P.[™] MidiPrep Kit allows isolation of highly pure plasmid DNA that is suitable for transfection, manual or automated sequencing, PCR, restriction mapping, ligation, and transformation. The S.N.A.P.[™] resin is supplied in a column and preferentially binds supercoiled plasmid DNA.

Storage

Store the kit at room temperature, except for the following components which are stored at +4°C after opening the kit:

- Resuspension Buffer
 - Precipitation Salt
 - RNase A
-

Contents

The following items are included in the S.N.A.P.[™] MidiPrep Kit. Sufficient reagents are included for 20 reactions.

Note: Solutions are not interchangeable with the S.N.A.P.[™] Miniprep Kit.

Item	Composition	Quantity
S.N.A.P. [™] MidiPrep Column A (Filtering)	--	20
Column B (Binding)	--	20
Resuspension Buffer	50 mM Tris-HCl, pH 8.0, 10 mM EDTA	80 ml
RNaseA, lyophilized	--	6 mg
Lysis Buffer	0.2 M NaOH, 1% SDS	80 ml
Precipitation Salt	3 M Potassium acetate, pH 5.2	80 ml
Binding Buffer	7.5 M Guanidine-HCl	2 x 120 ml
Wash Buffer	5 M Guanidine-HCl, 50 mM MOPS, pH 7.0	100 ml
4X Final Wash	400 mM NaCl	80 ml



The Binding Buffer and Wash Buffer contain a chaotropic salt. Use gloves and protective eye wear when handling these solutions.

Using the S.N.A.P.™ MidiPrep Kit

Experimental Outline

The procedure requires the following steps:

1. Lyse cells using the Lysis Buffer
 2. Precipitate protein and genomic DNA
 3. Filter the lysate/precipitate through Column A
 4. Add Binding Buffer to the flowthrough and apply to Column B
 5. Wash the bound plasmid and dry the resin by centrifugation
 6. Elute the plasmid DNA with sterile water or TE buffer.
-

User Supplied Reagents

- 240 ml 95% ethanol for dilution of 4X Final Wash
 - 500 ml sterilized bottle for dilution of 4X Final Wash
 - Sterile 50 ml conical tubes (Falcon, Sarstedt, or Corning)
 - Centrifuge with rotor and rotor adapters for 50 ml conical tubes
 - Sterile water or TE
-


Before Starting

1. Resuspend the entire contents of the RNase A tube (6 mg) in 1 ml of Resuspension Buffer. Add the resulting solution to the remaining Resuspension Buffer. Resuspension Buffer containing RNase A should now be stored at +4°C.
 2. Add the 80 ml of 4X Final Wash to 240 ml of 95% ethanol to make 320 ml of 1X Final Wash. Store in a 500 ml sterilized bottle.
 3. Check the Lysis Buffer for a white precipitate. If present, place the buffer in a 37°C water bath for 5 minutes or until the solution clears.
 4. Use bacterial cell cultures that are 1 to 2×10^9 cells/ml ($OD_{600} < 3.0$). For cultures $>3 \times 10^9$ cells/ml, use 50 ml or less of culture per column.
-

Continued on next page

Using the S.N.A.P.™ MidiPrep Kit, Continued

Lysis and Removal of Precipitate



You will need four 50 ml conical tubes for each plasmid preparation. All centrifugations are at room temperature except as noted.

1. Centrifuge 10-100 ml of an overnight bacterial culture at $4000 \times g$ for 5-10 minutes at $+4^{\circ}\text{C}$ to pellet the cells. Pour off the medium (do not aspirate medium).
 2. Resuspend the cell pellet in 4 ml of Resuspension Buffer by vortexing or gently pipetting up and down.
 3. Add 4 ml of Lysis Solution and mix by inverting gently 5-6 times. Incubate for 3 minutes at room temperature.
 4. Add 4 ml of Precipitation Salt and invert gently 6 to 8 times. Incubate 5 minutes in ice. Invert tubes twice during incubation to ensure even formation of precipitate within the solution. Do not centrifuge.
 5. During the incubation, label two 50 ml conical tubes "A" and "B". Place one S.N.A.P.™ MidiPrep Column A (Filtering) inside the 50 ml conical tube labeled "A", and one S.N.A.P.™ MidiPrep Column B (Binding) inside the 50 ml conical tube labeled "B" (see figure).
 6. After incubation, transfer the solution from Step 4 onto Column A and centrifuge for 5 minutes at $3,000 \times g$. Discard Column A. Save the filtrate containing the plasmid DNA.
-

Plasmid Binding

1. To the filtrate add 12 ml of Binding Buffer and mix by gently inverting twice. Transfer the solution to Column B inserted in a 50 ml conical tube.
 2. Centrifuge for 2 minutes at $1,000 \times g$. The plasmid DNA is now bound to Column B. Discard the flowthrough.
 3. Add 5 ml of Wash Buffer and centrifuge for 1 minute at $2,000 \times g$. Discard flowthrough.
 4. Add 5 ml of 1X Final Wash Buffer and centrifuge for 2 minutes at $2,000 \times g$.
 5. Add 10 ml of 1X Final Wash Buffer and centrifuge for 2 minutes at $2,000 \times g$. Discard flowthrough.
 6. Centrifuge Column B at $>4,000 \times g$ for 5 minutes to dry the resin.
-

Continued on next page

Using the S.N.A.P.™ MidiPrep Kit, Continued

Plasmid Elution

1. To elute the plasmid DNA, transfer Column B to a new, sterile 50 ml conical tube and add 750 µl of sterile water or TE buffer directly to column resin.
 2. Incubate for 3 minutes at room temperature. **Note:** Do not elute in less than 750 µl. Plasmid recovery will decrease.
 3. Centrifuge for 5 minutes at >4,000 x g. The plasmid DNA is now eluted from Column B.
 4. Remove and discard the column.
-

Determining DNA Yield

The plasmid DNA is now pure and ready for all molecular biology applications. To determine the plasmid DNA yield and concentration, use a spectrophotometer and read the absorbance at 260 nm. When using a spectrophotometer, zero the machine first with water or TE at 260 nm. Use 5-10 µl of the plasmid DNA sample to determine concentration. Use the formula give below:


$$[\text{DNA}] = (A_{260}) (50^\circ\text{g/ml}) \times D$$

where D is the dilution factor. Determine the yield by multiplying the concentration by the volume of DNA.

Typical Yields

Maximum yields range from 10 to 200 µg DNA total, depending on the copy number of the plasmid and the volume of culture. Concentration of DNA will range from 10 to 275 µg/ml. It may be necessary to ethanol precipitate the DNA to concentrate it if you are using small culture volumes. An example of plasmid DNA isolated using the S.N.A.P.™ Kit is shown below.

Kb MW 1 2 3 4 MW



Legend for the gel:

Plasmid DNA was isolated from four 50 ml TOP10F' overnight cultures and electrophoresed on a 0.8% agarose gel.

Lanes 1-4: 500 ng of plasmid purified with the S.N.A.P.™ MidiPrep Kit

Troubleshooting

Troubleshooting

Review the information provided in the table below to troubleshoot your experiments.

Problem	Cause	Solution
Low yield	Incomplete lysis or too much cell lysate has clogged the column	Decrease cell culture volume used. Use up to 100 ml of an overnight culture, density 2×10^9 cells/ml ($OD_{600} < 3.0$ per ml). For more dense cultures, decrease volume of culture.
		Decrease cell culture density. Use LB medium instead of TB (Terrific Broth).
	Genomic DNA present	Perform Steps 3-5 (page 5) with gentle mixing to avoid damaging DNA. Genomic DNA must be intact to be efficiently removed.
Plasmid resistant to restriction enzyme digestion	Denatured, supercoiled plasmid present (migrates below native, supercoiled plasmid)	Do not incubate longer than 3 minutes in Lysis Buffer (Step 3, page 5).

Technical Service

Contact Us

For more information or technical assistance, please call, write, fax, or email. Additional international offices are listed on our web page (www.invitrogen.com).

United States Headquarters:

Invitrogen Corporation
1600 Faraday Avenue
Carlsbad, CA 92008 USA
Tel: 1 760 603 7200

Tel (Toll Free): 1 800 955 6288

Fax: 1 760 602 6500

E-mail:

tech_service@invitrogen.com

European Headquarters:

Invitrogen Ltd
3 Fountain Drive
Inchinnan Business Park
Paisley PA4 9RF, UK

Tel (Free Phone Orders): 0800 269 210

Tel (General Enquiries): 0800 5345 5345

Fax: +44 (0) 141 814 6287

E-mail: eurotech@invitrogen.com

MSDS Requests

To request an MSDS, go to (www.invitrogen.com).

1. On the home page, go to the left-hand column under 'Technical Resources' and select 'MSDS Requests'.
 2. Follow instructions on the page and fill out all the required fields.
 3. To request additional MSDSs, click the 'Add Another' button.
 4. All requests will be faxed unless another method is selected.
 5. When you are finished entering information, click the 'Submit' button. Your MSDS will be sent within 24 hours.
-

Continued on next page

Technical Service, Continued

Emergency Information

In the event of an emergency, customers of Invitrogen can call the 3E Company, 24 hours a day, 7 days a week for disposal or spill information. The 3E Company connects the customer with poison control or with the University of California at San Diego Medical Center doctors. 3E Company, Voice: 1-760-602-8700

Limited Warranty

Invitrogen is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, please contact our Technical Service Representatives.

Invitrogen warrants that all of its products will perform according to the specifications stated on the certificate of analysis. The company will replace, free of charge, any product that does not meet those specifications. This warranty limits Invitrogen Corporation's liability only to the cost of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. Invitrogen reserves the right to select the method(s) used to analyze a product unless Invitrogen agrees to a specified method in writing prior to acceptance of the order.

Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore Invitrogen makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Service Representatives.

Invitrogen assumes no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.

Accessory Products

Additional Products

The following products are available separately from Invitrogen.

Product	Quantity	Catalog no.
Lipofectamine™ 2000	0.75 ml	11668-027
	1.5 ml	11668-019
Lipofectamine™ Reagent	1 ml	18324-012
Lipofectin™ Reagent	1 ml	18292-011
Calcium Phosphate Transfection Kit	75 reactions	K2780-01
Plus™ Reagent	0.85 ml	11514-015
SequenceR _x Enhancer System	25 reactions	12237-012

Product Qualification

Each component of the S.N.A.P.™ MidiPrep is lot qualified for maximum performance in accordance with the S.N.A.P.™ MidiPrep protocol. A sample of 4 columns from each lot is qualified in the following manner:

Binding Capacity

300 µg of pre-purified control plasmid is applied to the column and eluted according to the protocol. The quantity of plasmid DNA eluted from the column must be >200 µg.

Kit Performance

50 ml of control plasmid is grown in TOP10F' cells and plasmid DNA is isolated using the S.N.A.P.™ MidiPrep Kit. Plasmid DNA is qualified in the following manner:

- **Visual inspection**

Four samples are run on 0.8% agarose gel for a visual inspection. Each isolated plasmid must show supercoiled plasmid DNA with no RNA contamination.

- **Restriction Digest**

Four 500 ng plasmid DNA samples are digested with *Apa*I restriction enzyme for 90 minutes. Restriction digest is run on a 0.8% agarose gel and must be >99% complete.

- **Endonuclease Activity**

Four 500 ng plasmid DNA samples are incubated with 10 mM Mg²⁺ for 4 hours and then run on a 0.8% agarose gel. Each sample must show no endonuclease contamination of supercoiled DNA present.

©1997-2002 Invitrogen Corporation. All rights reserved.



Corporate Headquarters:

*Invitrogen Corporation
1600 Faraday Avenue
Carlsbad, California 92008
Tel: 1 760 603 7200
Tel (Toll Free): 1 800 955 6288
Fax: 1 760 603 7229
Email: tech_service@invitrogen.com*

European Headquarters:

*Invitrogen Ltd
3 Fountain Drive
Inchinnan Business Park
Paisley PA4 9RF, UK
Tel (Free Phone Orders): 0800 269 210
Tel (General Enquiries): 0800 5345 5345
Fax: +44 (0) 141 814 6287
Email: eurotech@invitrogen.com*

International Offices:

*Argentina 5411 4556 0844
Australia 1 800 331 627
Austria 0800 20 1087
Belgium 0800 14894
Brazil 0800 11 0575
Canada 800 263 6236
China 10 6849 2578
Denmark 80 30 17 40*

*France 0800 23 20 79
Germany 0800 083 0902
Hong Kong 2407 8450
India 11 577 3282
Italy 02 98 22 201
Japan 03 3663 7974
The Netherlands 0800 099 3310
New Zealand 0800 600 200
Norway 00800 5456 5456*

*Spain & Portugal 900 181 461
Sweden 020 26 34 52
Switzerland 0800 848 800
Taiwan 2 2651 6156
UK 0800 838 380
For other countries see our Web site*

www.invitrogen.com

