

RNase A (Affinity Purified)

Store at -20°C .
Do not store in a frost-free freezer.

Catalog #:	AM2270	AM2271
Amount:	200 μg	1 mg
Product Description:	Affinity-purified RNase A intended for applications where the absence of DNase and other nonspecific nuclease activity is critical.	
Source:	Manufactured from Bovine pancreas by a USDA-approved source.	
Concentration:	1 mg/mL	
Unit Definition:	One unit is the amount of RNase A required to give an increase of 0.0146 A_{286} units per min in a 1 mL volume at room temperature and is equivalent to 0.1177 Kunitz Units. Assay Conditions: 100 mM Tris-acetate (pH 6.5), 1 mM EDTA, and 1 mM cyclic 2', 3'-CMP.	
Materials Not Provided:	Because of the variety of product applications, no reaction buffer is supplied.	
Storage Conditions:	Store at -20°C . Do not store in a frost-free freezer.	
Storage Buffer:	(Not included) 10 mM HEPES (pH 7.2), 20 mM NaCl, 0.1% Triton X-100, and 50% glycerol (v/v).	

USER INFORMATION

General Information:	RNase A is a ribonuclease that specifically cleaves 3' of U and C residues [1]. Ambion RNase A (Affinity Purified) is intended for critical applications when the absence of DNase and other nonspecific nuclease activities in the RNase preparation is essential.
Applications:	<p>Among other applications, RNase A (Affinity Purified) can be used to remove RNA from DNA preparations. It is active under a wide range of reaction conditions.</p> <p>For minipreps of plasmid DNA, add RNase A to 20 $\mu\text{g}/\text{mL}$ during restriction enzyme digestion reactions. For large-scale plasmid DNA preparations, add RNase A to a final concentration of 20 $\mu\text{g}/\text{mL}$ to a crude lysate that has been resuspended in TE. Incubate 20 min at 37°C, then extract with buffered phenol and chloroform/isoamyl alcohol (24:1) and precipitate using 2 M ammonium acetate and 2 volumes 100% ethanol before further purification [2].</p>
References:	<ol style="list-style-type: none"> 1. Crook EM, Mathias AP, and Rabin BR. (1960) Spectrophotometric assay of bovine pancreatic ribonuclease by the use of cytidine 2':3'-phosphate. <i>Biochem J.</i> 74:234-238. 2. Ausubel FM, Brent R, Kingston RE, Moore DD, Seidman JG, Smith JA, Struhl K, editors. (2006) <i>E. coli</i>, Plasmids, and Bacteriophages (Chapter 1). In <i>Current Protocols in Molecular Biology</i>. John Wiley & Sons, Inc. p 1.6.2 and 1.7.9.

QUALITY CONTROL

RNase A (Affinity Purified) is rigorously tested for contaminating nonspecific endonuclease, exonuclease, and protease activity.

OTHER INFORMATION

Material Safety Data Sheets:	Material Safety Data Sheets (MSDSs) can be printed or downloaded from product-specific links on our website at the following address: www.ambion.com/techlib/msds . Alternatively, e-mail your request to MSDS_Inquiry_CCRM@appliedbiosystems.com . Specify the catalog or part number(s) of the product(s), and we will e-mail the associated MSDSs unless you specify a preference for fax delivery. For customers without access to the internet or fax, our technical service department can fulfill MSDS requests placed by telephone or postal mail. (Requests for postal delivery require 1–2 weeks for processing.)
Warranty and Liability:	<p>For research use only. Not for use in diagnostic procedures.</p> <p>Applied Biosystems is committed to delivering superior product quality and performance, supported by industry-leading global service and technical support teams. Warranty information for the accompanying consumable product is available at www.ambion.com/info/warranty in "Limited Warranty for Consumables," which is subject to the exclusions, conditions, exceptions, and limitations set forth under the caption "EXCLUSIONS, CONDITIONS, EXCEPTIONS, AND LIMITATIONS" in the full warranty statement. Please contact Applied Biosystems if you have any questions about our warranties or would like information about post-warranty support.</p>
Trademarks, Patents, and Licensing:	Applied Biosystems, AB (Design), and Ambion are registered trademarks of Applied Biosystems Corporation or its subsidiaries in the US and/or certain other countries. All other trademarks are the sole property of their respective owners.
	© 2007 Ambion, Inc. All rights reserved. 4393902A