

AccuPrime™ SuperMix II

Cat. no: 12341-012
12341-020

Kit Size: 200 reactions
1,000 reactions

Store at 4°C or –20°C (non-frost-free)

Description

AccuPrime™ SuperMix II is designed for amplification of genomic DNA templates 200 bp–4 kb in size. It is **not recommended** for amplification of small genomic templates (≤ 200 bp), plasmid DNA, or cDNA templates.

AccuPrime™ SuperMix II is a specially formulated master mix for the amplification of nucleic acid templates by PCR. The mixture contains recombinant *Taq* DNA polymerase, anti-*Taq* antibodies, thermostable AccuPrime™ protein, Mg^{++} , and dNTPs at optimized concentrations. The SuperMix is supplied at 2X concentration to allow for the addition of primers and template.

Anti-*Taq* DNA polymerase antibodies inhibit polymerase activity at ambient temperatures, allowing room-temperature reaction setup and providing an automatic “hot start” in PCR (Chou *et al.*, 1992; Sharkey *et al.*, 1994). Thermostable AccuPrime™ protein enhances specific primer-template hybridization during every cycle of PCR. Antibody- and AccuPrime™ protein-mediated amplification dramatically improves PCR specificity and fidelity, providing the most robust PCR for multiplex PCR and sub-optimal primer sets.

Reagents are provided for 200 or 1,000 amplification reactions of 25 μ l each.

Storage

Store at 4°C or –20°C. No reduction of PCR performance or enzyme activity has been observed after storage for 12 months at 4°C. Note that repeated freeze-thaw cycles may reduce enzyme performance.

	<u>200-rxn kit</u>	<u>1,000-rxn kit</u>
AccuPrime™ SuperMix II	2 × 1.25 ml	12.5 ml

Part no. 12341.pps

MAN0001077

Rev. date: 11 Jun 2010

Components

40 mM Tris-HCl (pH 8.4); 100 mM KCl; 3 mM MgCl₂, 400 μM each dGTP, dATP, dTTP, dCTP; AccuPrime™ *Taq* DNA Polymerase; thermostable AccuPrime™ protein; stabilizers

Terminal Transferase Activity

Like regular *Taq* DNA polymerase, AccuPrime™ *Taq* DNA Polymerase has a nontemplate-dependent terminal transferase activity that adds a single deoxyadenosine (A) to the 3' ends of PCR products.

Guidelines for PCR

General PCR parameters and troubleshooting information are documented in Innis, et al (Innis et al., 1990). PCR reactions should be assembled in a DNA-free environment using clean, dedicated automatic pipettors and aerosol resistant barrier tips. Always keep the control DNA and other templates to be amplified isolated from the other components.

Recommendations for PCR

The protocol on the following page is designed as a starting point for PCR amplification. Optimal reaction conditions—incubation times and temperatures, primers, and template DNA—may vary.

- **Starting material:** 10 pg to 200 ng of genomic DNA (200 bp–4 kb in size).
- **Annealing temperature:** In general, the optimal annealing temperature should be 5–10° lower than the T_m of the primers used.

References

- Chou, Q., Russell, M., Birch, D., Raymond, J., and Bloch, W. (1992) Prevention of pre-PCR mis-priming and primer dimerization improves low-copy-number amplifications. *Nucl. Acids Res.*, 20, 1717-1723
- Innis, M. A., Gelfand, D. H., Sninsky, J. J., and White, T. S. (eds) (1990) *PCR Protocols: A Guide to Methods and Applications*, Academic Press, San Diego, CA
- Sharkey, D. J., Scalice, E. R., Christy, K. G., Atwood, S. M., and Daiss, J. L. (1994) Antibodies as thermolabile switches: high temperature triggering for the polymerase chain reaction. *Biotechnology*, 12, 506-509

PCR Protocol

The following general protocol is suggested as a starting point when using AccuPrime™ SuperMix II in PCR.

Due to the “hot-start” capability of AccuPrime™ SuperMix II, the reaction can be set up at room temperature.

1. Program the thermal cycler as follows:

Initial denaturation: 94°C for 2 minutes

25–35 cycles of:

Denaturation: 94°C for 15–30 seconds

Annealing: 55–60°C for 15–30 seconds

Extension: 68°C for 1 minute per kb of PCR product

2. Add the following components to a DNase/RNase-free microcentrifuge tube. The 25- μ l reaction size may be scaled up or down as needed.

<u>Component</u>	<u>Volume</u>	<u>Final Concentration</u>
AccuPrime™ SuperMix II	12.5 μ l	1X
Primer mix (10 μ M each)	0.5 μ l	0.2 μ M each
Template DNA (10 pg–200 ng)	\geq 1 μ l	as required
Autoclaved distilled water	to 25 μ l	n/a

3. Cap the tube, tap gently to mix, and centrifuge briefly to collect the contents.
4. Place the tube in the thermal cycler and run the program from Step 1.
5. After cycling, maintain the reaction at 4°C. Samples can be stored at –20°C until use.

Analyze the amplification products by agarose gel electrophoresis. We recommend using E-Gel® 1.2% gels and TrackIt™ 100 bp or 1kb Plus DNA ladders (see **Additional Products** on page 4).

Quality Control

The Certificate of Analysis (CofA) provides detailed quality control information for each product. The CofA is available on our website at

www.invitrogen.com/cofa, and is searchable by product lot number, which is printed on each box.

Additional Products

<u>Product</u>	<u>Amount</u>	<u>Catalog no.</u>
E-Ge ^l ® 1.2% Starter Pak	6 gels plus PowerBase™	G6000-01
TrackIt™ 100 bp DNA Ladder	100 applications	10488-058
TrackIt™ 1kb Plus DNA Ladder	100 applications	10488-085

Limited Use Label License No. 1: Thermostable Polymerases

Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the U.S.: 5,789,224, 5,618,711, 6,127,155. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim, no right to perform any patented method, and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

Limited Use Label License No. 5: Invitrogen Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For products that are subject to multiple limited use label licenses, the terms of the most restrictive limited use label license shall control. Life Technologies Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Life Technologies Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Life Technologies is willing to accept return of the product with a full refund. For information about purchasing a license to use this product or the technology embedded in it for any use other than for research use please contact Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008; Phone (760) 603-7200 or e-mail: outlicensing@lifetech.com.

Limited Use Label License No. 14: Direct Inhibition by Anti-Polymerase Antibodies

Licensed to Life Technologies Corporation, under U.S. Patent Nos. 5,338,671; 5,587,287, and foreign equivalents for use in research only.

Limited Use Label License No. 33: AccuPrime™ Reagent

This product is sold under licensing arrangements with Stratagene. The purchase price of this product includes limited, nontransferable rights under U.S. Patents Nos. 5,449,603; 5,605,824; 5,646,019; and 5,773,257 and foreign equivalents owned by Stratagene to use only this amount of the product to practice the claims in said patents solely for activities of the purchaser within the field of research. Further information on purchasing licenses under the above patents may be obtained by contacting the Director of Business Development, Stratagene, 11011 North Torrey Pines Road, La Jolla, California 92037.

©2010 Life Technologies Corporation. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.

The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.