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GoTaq[®] Probe qPCR Master Mix

Cat. # A6101, A6102

www.promega.com

GoTaq[®] Probe qPCR Master Mix

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- Applied Biosystems 7500 and 7500 FAST Real-Time PCR System
- Bio-Rad CFX96 Real-Time PCR Detection System
- Bio-Rad DNA Engine Opticon® and Opticon® 2 Real Time PCR Detection Systems
- Bio-Rad/MJ Research Chromo4™ Real-Time Detector
- Cepheid SmartCycler® system
- Corbett Rotor-Gene™ 3000 and 6000 Real-Time Rotary Analyzer
- Eppendorf Mastercycler® ep realplex Real-Time PCR System
- Roche LightCycler® 480 Real-Time PCR System
- Stratagene Mx3000P® and Mx3005P® Real-Time PCR Systems
- Stratagene Mx4000® Multiplex Quantitative PCR System

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- Applied Biosystems ABI PRISM® 7000 and 7700 Sequence Detection System
- Applied Biosystems 7300 and 7900HT Real-Time PCR System
- Applied Biosystems GeneAmp® 5700 Thermal Cycler
- Applied Biosystems StepOne™ and StepOnePlus™ Real-Time PCR Systems

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1. Bustin, S.A. et al. (2009) The MIQE guidelines: Minimum information for publication of quantitative real-time PCR experiments. *Clin. Chem.* 55, 611–22.
2. Dorak, M.T. (2009) Glossary of real-time PCR terms. This can be viewed online at: www.dorak.info/genetics/glosrt.html
3. Fleige, S. and Pfaffl, M.W. (2006) RNA integrity and the effect on the real-time qRT-PCR performance. *Mol. Aspects Med.* 27, 126–39.
4. Lefever, S. et al. (2009) RDML: Structured language and reporting guidelines for real-time quantitative PCR data. *Nucleic Acids Res.* 37, 2065–9.
5. Livak, K.J. and Schmittgen, T.D. (2001) Analysis of relative gene expression data using real-time quantitative PCR and the 2(-delta delta C(T)) method. *Methods* 25, 402–8.