



Promega

Technical Manual

pCAT[®]3 Reporter Vectors

INSTRUCTIONS FOR USE OF PRODUCTS E1851, E1861, E1871
AND E1881.



www.promega.com

PRINTED IN USA.
Revised 3/08



AF9TM036 0308TM036

Part# TM036

pCAT[®]3 Reporter Vectors

All technical literature is available on the Internet at: www.promega.com/tbs/
 Please visit the web site to verify that you are using the most current version of this
 Technical Manual. Please contact Promega Technical Services if you have questions on use
 of this system. E-mail: techserv@promega.com

1. Description.....	2
2. Product Components and Storage Conditions	2
3. General Considerations	3
A. Common Structural Elements of the pCAT [®] 3 Reporter Vectors	3
B. Advantages of the pCAT [®] 3 Vectors	3
C. Distinguishing Features of the pCAT [®] 3 Reporter Vectors	5
D. CAT Assay	10
E. Mapping Genetic Elements Located Within DNA Fragments.....	10
F. Site-Specific Mutagenesis of the DNA Fragments	10
4. Cloning Methods	12
A. Cloning Strategies.....	12
B. Preparation of pCAT [®] 3 Vectors and Insert DNA for Cloning.....	12
C. Transformation Protocols for pCAT [®] 3 Vectors	12
D. Isolation of Plasmid DNA.....	13
5. Transfection of Mammalian Cells.....	13
6. Assay of CAT Activity	13
7. Generation of Nested Deletions	14
8. Generation of Single-Stranded DNA and Site-Specific Mutations	14
A. Production of Single-Stranded DNA.....	14
B. Generation of Site-Specific Mutations.....	15
9. Sequencing of pCAT [®] 3 Reporter Vectors.....	15
10. References	16
11. Appendix	18
A. pCAT [®] 3-Basic Vector Restriction Sites.....	18
B. pCAT [®] 3-Enhancer Vector Restriction Sites	20
C. pCAT [®] 3-Promoter Vector Restriction Sites	22
D. pCAT [®] 3-Control Vector Restriction Sites	24
E. Related Products.....	26

1. Description

The pCAT[®]3 Reporter Vectors provide a basis for the quantitative analysis of factors that potentially regulate mammalian gene expression. These factors may be *cis*-acting, such as promoters and enhancers, or *trans*-acting, such as various DNA-binding factors. The backbone of the pCAT[®]3 Reporter Vectors (1) is similar to the pGL3 Luciferase Vectors, with the exception of an intron located 5' of the chloramphenicol acetyltransferase (CAT) gene. The vector backbones were designed to increase expression of the reporter gene for easier monitoring of transcriptional activity in transfected eukaryotic cells. In addition, the pCAT[®]3 Reporter Vectors contain numerous features that aid in the structural characterization of putative regulatory sequences.

For peer-reviewed articles that cite use of the pCAT[®]3 Vectors, visit:
www.promega.com/citations/

2. Product Components and Storage Conditions

<u>Product</u>	<u>Size</u>	<u>Cat.#</u>
pCAT [®] 3-Control Vector	20µg	E1851
pCAT [®] 3-Promoter Vector	20µg	E1861
pCAT [®] 3-Basic Vector	20µg	E1871
pCAT [®] 3-Enhancer Vector	20µg	E1881

Information on related products, including the CAT Assay System, is provided in Section 11.E.

Storage Conditions: Store the Vector DNA at -20°C.

3. General Considerations

3.A. Common Structural Elements of the pCAT[®]3 Reporter Vectors

Except for the inclusion of promoters and enhancers, the four pCAT[®]3 Reporter Vectors are structurally identical. The distinguishing features of each vector are summarized in Section 3.C. The pCAT[®]3 Vectors each contain a high-copy-number replicon for maintenance in *E. coli*, an ampicillin-resistance gene for selection, and a filamentous phage origin of replication (f1 ori) for single-stranded DNA (ssDNA) production. Restriction sites for insertion of DNA fragments are located upstream and downstream of the CAT gene.

3.B. Advantages of the pCAT[®]3 Vectors

The pCAT[®]3 Reporter Vector family provides significant advances over the first generation of pCAT[®] Reporter Vectors (Table 1). The pCAT[®]3 Reporter Vectors contain a redesigned vector backbone, which increases CAT expression, improves in vivo vector stability, provides greater flexibility in performing genetic manipulations and facilitates transfer of cloned fragments between the pCAT[®]3 Vectors and pGL3 Luciferase Reporter Vectors. The modifications in the reporter vectors result in CAT expression levels higher than those obtained with the first generation of pCAT[®] Reporter Vectors, while maintaining relatively low background CAT expression.

The increase in the CAT expression observed with the pCAT[®]3 Vectors provides greater sensitivity. It may now be possible to obtain measurable CAT expression in cell types that are difficult to transfect or when studying weak promoter elements. However, relative expression profiles may vary between cell types (1). Therefore, it is important to include the appropriate control vectors in all experiments.

Table 1. Characteristics of the pCAT®3 Vectors.

Feature	Purpose of Feature
Poly(A) signal derived from the late SV40 poly(A) signal.	Late SV40 poly(A) signal is more efficient than the early SV40 poly(A) (2).
Synthetic poly(A) transcriptional pause site located 5' of the multiple cloning site.	Reduces background CAT expression while avoiding possible recombination between two SV40 poly(A) sequences in the same plasmid (3,4).
Chimeric intron located 5' of the CAT gene.*	The SV40 small-t antigen intron in the first generation of pCAT® Vectors reduced expression when placed 3' of the CAT gene due to cryptic splicing (5-7).
Kozak consensus sequence located at the 5' end of the CAT gene.	Provides for optimal translation efficiency (8).
Multiple cloning region.	Increases convenience and provides compatibility with the pGL3 Vectors.
f1 ori site in the vector backbone.	Provides the ability to generate ssDNA for sequencing or mutagenesis.
Unique NcoI site in the Kozak sequence.	Facilitates subcloning.
Unique XbaI site downstream of the CAT gene.	Facilitates subcloning.
An EcoRI site has been removed from the CAT gene, resulting in the amino acid change Phe73→Leu73.	Facilitates subcloning.
<p>*Note: The chimeric intron is composed of the donor site from the first intron of the human β-globin gene and the branch and acceptor site from the intron of an immunoglobulin gene (9). The sequences of the donor and acceptor sites, along with the branchpoint site, have been optimized to match the consensus sequence (10).</p> <p>While the chimeric intron placed 5' of the CAT gene increases CAT expression, its location 5' of the gene may also result in spurious transcription in some cell lines. We have removed some of the known potential regulatory sequences from the intron; however, any unidentified regulatory sequences still present within the vector can lead to increased background CAT expression. The user should recognize this possibility and use the proper experimental controls. If the intron is not required for increased expression, it is conveniently flanked by HindIII sites for easy removal from the pCAT®3 Vectors.</p>	

3.C. Distinguishing Features of the pCAT[®]3 Reporter Vectors

Maps of the pCAT[®]3-Basic, Enhancer, Promoter and Control Vectors are shown in Figures 1–4. Sequence accession numbers and information on restriction enzyme sites present in these vectors are provided in Section 11.

pCAT[®]3-Basic Vector

The pCAT[®]3-Basic Vector lacks eukaryotic promoter and enhancer sequences, allowing maximum flexibility in cloning putative regulatory sequences. Expression of CAT activity in cells transfected with this plasmid depends on insertion and proper orientation of a functional promoter upstream from the intron and the CAT gene. Potential enhancer elements can also be inserted upstream of the promoter or in the BamHI or Sall sites downstream of the CAT transcription unit.

pCAT[®]3-Enhancer Vector

The pCAT[®]3-Enhancer Vector contains an SV40 enhancer located downstream of the CAT gene and the poly(A) signal. This aids in the verification of functional promoter elements because the presence of an enhancer will often result in transcription of the CAT gene at higher levels.


pCAT[®]3-Promoter Vector

The pCAT[®]3-Promoter Vector contains an SV40 promoter upstream of the intron and the CAT gene. DNA fragments containing putative enhancer elements can be inserted either upstream or downstream of the promoter-CAT transcriptional unit.

pCAT[®]3-Control Vector

The pCAT[®]3-Control Vector contains SV40 promoter and enhancer sequences, resulting in strong expression of CAT in many mammalian cell types. The pCAT[®]3-Control Vector is useful for monitoring transfection efficiency, and provides a convenient standard for comparing promoter and enhancer activities expressed by pCAT[®]3 recombinants.

Note: The specific transcriptional characteristics of the pCAT[®]3 Vectors will vary for different cell types. This may be particularly true for COS cells, which contain the SV40 large-T antigen. The SV40 large-T antigen promotes replication from the SV40 origin, which is found in the promoter of the pCAT[®]3-Promoter and pCAT[®]3-Control Vectors. The combination of large-T antigen and SV40 origin will result in a higher copy number of these vectors in COS cells, which in turn may result in increased expression of the reporter gene compared to other cell and vector combinations.




1131VA08_5A

Figure 1. The pCAT®3-Basic Vector circle map. -[^]-, position of intron; CAT, cDNA encoding the chloramphenicol acetyltransferase gene; Amp^r, gene conferring ampicillin resistance in *E. coli*; f1 ori, origin of replication derived from filamentous phage; ori, origin of plasmid replication in *E. coli*. Arrows within CAT and the Amp^r genes indicate the direction of transcription; the arrow in f1 ori indicates the direction of ssDNA strand synthesis. Restriction sites shown in parentheses are not unique sites.

Sequence reference points:

SV40 Promoter	(none)
SV40 Promoter-directed transcriptional start sites	(none)
SV40 Enhancer	(none)
SV40 late poly(A) region	981-1202
CAT gene	291-947
chimeric intron	103-235
upstream poly(A) region	3867-4020
multiple cloning region	1-58
RVprimer3 binding site	3969-3988
RVprimer4 binding site	1289-1270
β -lactamase gene (Amp ^r)	3149-2292
f1 origin	3282-3736
ColE1-derived plasmid replication origin	1527




1129VA08_5A

Figure 2. The pCAT[®]3-Enhancer Vector circle map. -[^]-, position of intron; CAT, cDNA encoding the chloramphenicol acetyltransferase gene; Amp^r, gene conferring ampicillin resistance in *E. coli*; f1 ori, origin of replication derived from filamentous phage; ori, origin of plasmid replication in *E. coli*. Arrows within CAT and the Amp^r genes indicate the direction of transcription; the arrow in f1 ori indicates the direction of ssDNA strand synthesis. Restriction sites shown in parentheses are not unique sites.

Sequence reference points:

SV40 Promoter	(none)
SV40 Promoter-directed transcriptional start sites	(none)
SV40 Enhancer	1222-1458
SV40 late poly(A) region	981-1202
CAT gene	291-947
chimeric intron	102-234
upstream poly(A) region	4113-4266
multiple cloning region	1-58
RVprimer3 binding site	4215-4234
RVprimer4 binding site	1535-1516
β-lactamase gene (Amp ^r)	3395-2538
f1 origin	3528-3982
ColE1-derived plasmid replication origin	1773




1132VA08_5A

Figure 3. The pCAT[®]3-Promoter Vector circle map. -[^]-, position of intron; CAT, cDNA encoding the chloramphenicol acetyltransferase gene; Amp^r, gene conferring ampicillin resistance in *E. coli*; f1 ori, origin of replication derived from filamentous phage; ori, origin of plasmid replication in *E. coli*. Arrows within CAT and the Amp^r genes indicate the direction of transcription; the arrow in f1 ori indicates the direction of ssDNA strand synthesis. Restriction sites shown in parentheses are not unique sites.

Sequence reference points:

SV40 Promoter	48-250
SV40 Promoter-directed transcriptional start sites	185, 191, 196
SV40 Enhancer	(none)
SV40 late poly(A) region	1173-1394
CAT gene	483-1139
chimeric intron	295-427
upstream poly(A) region	4059-4212
multiple cloning region	1-41
RVprimer3 binding site	4161-4180
RVprimer4 binding site	1481-1462
β-lactamase gene (Amp ^r)	3341-2484
f1 origin	3474-3928
ColE1-derived plasmid replication origin	1719



1130VA08_5A

Figure 4. The pCAT®3-Control Vector circle map. -[^]-, position of intron; CAT, cDNA encoding the chloramphenicol acetyltransferase gene; Amp^r, gene conferring ampicillin resistance in *E. coli*; f1 ori, origin of replication derived from filamentous phage; ori, origin of plasmid replication in *E. coli*. Arrows within CAT and the Amp^r genes indicate the direction of transcription; the arrow in f1 ori indicates the direction of ssDNA strand synthesis. Restriction sites shown in parentheses are not unique sites.

Sequence reference points:

SV40 Promoter	48-250
SV40 Promoter-directed transcriptional start sites	185, 191, 196
SV40 Enhancer	1414-1650
SV40 late poly(A) region	1173-1394
CAT gene	483-1139
chimeric intron	295-427
upstream poly(A) region	4305-4458
multiple cloning region	1-41
RVprimer3 binding site	4407-4426
RVprimer4 binding site	1727-1708
β -lactamase gene (Amp ^r)	3587-2730
f1 origin	3720-4174
ColE1-derived plasmid replication origin	1965

3.D. CAT Assay

Chloramphenicol acetyltransferase (CAT), encoded by a bacterial drug-resistance gene, inactivates chloramphenicol by acetylating the drug at one or both of its two hydroxyl groups (11). This gene is not found in eukaryotes, and therefore eukaryotic cells contain no background CAT activity. This characteristic, along with assay sensitivity, has made the CAT gene a good reporter for studies of mammalian gene expression (12,13).

3.E. Mapping Genetic Elements Located Within DNA Fragments

The locations of functional elements within a DNA fragment are often determined by making a set of unidirectional nested deletions following the method of Henikoff (14) and then assaying for changes in biological activity. This method takes advantage of the unique properties of Exonuclease III (ExoIII), which will digest 5' overhangs but not 3' overhangs or α -phosphorothioate nucleotide filled-in overhangs. Nested deletions of an insert DNA can be made directly in the pCAT[®]3 family of reporter vectors using this method, eliminating the need for subcloning steps. The multiple cloning site of the pCAT[®]3 Vectors contains upstream KpnI and SacI restriction sites, which can be used to generate the 3' overhangs resistant to ExoIII (Figures 1-5). After treatment with ExoIII, S1 Nuclease is added to remove the resulting ssDNA overhangs, Klenow Fragment is added to flush the ends, and the ends are ligated to circularize the vectors. Deletion clones can be screened by gel electrophoresis of miniprep DNA, and the precise deletion endpoints within the promoter region can be determined by DNA sequencing (see Section 9).

3.F. Site-Specific Mutagenesis of the DNA Fragments

Once a DNA fragment's biological activity has been identified, site-specific mutagenesis may be used to further define the sequences associated with the activity. To provide the ssDNA template necessary for some mutagenesis reactions, the pCAT[®]3 Reporter Vectors contain an origin of replication derived from filamentous phage. This allows single-stranded plasmid DNA to be produced and secreted in phage-like particles from *E. coli* infected with the appropriate helper phage.

Mutagenesis is performed using oligonucleotides that are complementary to the ssDNA but contain the desired changes to the nucleic acid sequence (15). The oligonucleotides are hybridized to the template DNA, and double-stranded DNA is synthesized using a DNA polymerase. After amplification of the DNA in *E. coli*, the mutations may be verified by DNA sequencing.




Figure 5. pCAT[®]3 Vector multiple cloning regions. Upstream and downstream cloning sites and the locations of sequencing primer (RVprimer3 and RVprimer4) binding sites are shown. The large primer arrows indicate the direction of sequencing. The positions of the promoter (in the pCAT[®]3-Promoter and pCAT[®]3-Control Vectors) and the enhancer (in the pCAT[®]3-Enhancer and pCAT[®]3-Control Vectors) are shown as insertions into the sequence of the pCAT[®]3-Basic Vector. (Note that the promoter replaces four bases [AAGT] of the pCAT[®]3-Basic Vector.) The sequence shown is of the DNA strand generated from the fl ori.

4. Cloning Methods

4.A. Cloning Strategies

The restriction sites for XhoI and Sall have compatible ends, as do BglII and BamHI. Therefore, cloning into the XhoI or BglII sites upstream of CAT or the downstream Sall or BamHI sites allows easy positioning of DNA inserts either upstream or downstream of the CAT reporter gene. Thus, positional effects of a putative genetic element may be readily tested. Cloning fragments into a single site generally will yield both possible orientations relative to the reporter gene, making these effects also readily testable. It should be noted that the area encompassing the XhoI restriction site in the multiple cloning region exhibits considerable secondary structure when the plasmid is in a supercoiled configuration and, as a result, is resistant to digestion. We therefore recommend gel purifying the vector after digestion with XhoI to avoid an excessive number of background colonies.

The other upstream restriction sites may be used for cloning. However, note that some of the sites are required for generation of nested deletions (see Section 7). Specifically, the KpnI or SacI site is needed to generate a 3' overhang upstream of the insert.

4.B. Preparation of pCAT[®]3 Vectors and Insert DNA for Cloning

The fragment and vector DNA should be digested with restriction enzymes that generate compatible ends for cloning. In most cases, the ends of the DNA fragment may require modification, either by using synthetic linkers, by using PCR primers containing recognition sites for the desired restriction enzymes, or by filling in restriction site overhangs. It may be advantageous to treat the vector DNA with Calf Intestinal Alkaline Phosphatase (CIAP) to remove 5'-phosphate groups, thus preventing recircularization of vector without insert. Sufficient DNA should be prepared to perform control reactions for digestion, ligation and transformation steps.

To ensure capture of the correct insert DNA, the desired restriction fragment can be purified by electrophoresis on an agarose gel and then recovered from the gel using the Wizard[®] SV Gel and PCR Clean-Up System (Cat.# A9281), or an equivalent method. Alternatively, unfractionated restriction fragments can be cloned into the target plasmid, and the desired recombinant identified by gel electrophoresis of plasmid DNA.

4.C. Transformation Protocols for pCAT[®]3 Vectors

Because the CAT Reporter Vectors are supplied as modified DNA, *E. coli* hosts may be either restriction⁺ or restriction⁻. Use of a *recA* host such as JM109 is preferred to prevent undesirable recombination between the insert and the host chromosomal DNA. A strain that has an F' episome is required for ssDNA production.

Grow JM109 on minimal plates (M-9) supplemented with thiamine HCl prior to preparation of competent cells and transformation. This selects for the F' episome.

4D. Isolation of Plasmid DNA

The Wizard® Plus SV Minipreps DNA Purification System (Cat.# A1340, A1470) may be used for small-scale preparation of plasmid DNA for screening clones. DNA suitable for transfection may be purified using the PureYield™ Plasmid Midipreps System (Cat.# A2492, A2495).

5. Transfection of Mammalian Cells

Transfection of DNA into eukaryotic cells may be mediated by cationic liposomes (16), calcium phosphate (13,17), DEAE-dextran (13,18), or electroporation (17). Transfection systems based on cationic lipids (TransFast™ Transfection Reagent, Transfectam® Reagent and Tfx™ Reagents) and calcium phosphate (Profection® Mammalian Transfection Systems) are available from Promega. For more information on these transfection reagents, please request the *TransFast™ Transfection Reagent Technical Bulletin* (#TB260), the *Transfectam® Reagent Technical Bulletin* (#TB116), the *Tfx™-Reagents Technical Bulletin* (#TB216) or the *Profection® Mammalian Transfection Systems Technical Manual* (#TM012). All of these documents are available on our web site at: www.promega.com/tbs/

6. Assay of CAT Activity

For CAT transient expression assays, cell extracts are typically prepared 48–72 hours post-transfection. CAT activity may then be monitored by two alternative methods using the CAT Enzyme Assay System with Reporter Lysis Buffer (Cat.# E1000). The most rapid, sensitive, and convenient of these is based on liquid scintillation counting (LSC) of CAT reaction products. Cell extracts are incubated in a reaction mix containing ¹⁴C- or ³H-labeled chloramphenicol and n-Butyryl Coenzyme A. CAT transfers the n-butyryl moiety of the cofactor to chloramphenicol. For the LSC assay, the reaction products are extracted with a small volume of xylene. The n-butyryl chloramphenicol partitions mainly into the xylene phase, while unmodified chloramphenicol remains predominantly in the aqueous phase (19). The xylene phase is mixed with scintillant and counted in a scintillation counter. This assay can be completed in as little as 2–3 hours, is linear for nearly three orders of magnitude and can detect as little as 3×10^{-4} units of CAT. CAT activity can also be analyzed using thin layer chromatography (TLC). This method is more time-consuming than LSC but allows visual confirmation of the data. For more information, see the *CAT Enzyme Assay System with Reporter Lysis Buffer Technical Bulletin* #TB084 (available at www.promega.com/tbs/).

7. Generation of Nested Deletions

Unidirectional deletions of any inserted DNA can be made using a procedure developed by Henikoff (14) in which Exonuclease III (ExoIII) is used to specifically digest insert DNA from a 5'-protruding or blunt-end restriction site. In the pCAT[®]3 Reporter Vectors, these 5' overhangs or blunt ends are supplied by digesting the plasmid with BglII, MluI, NheI, XhoI or XmaI. When the plasmids are cut with KpnI or SacI, which yield 3' overhangs, the ExoIII will be unable to digest in the other direction.

The uniform rate of enzyme digestion allows deletions of various lengths to be made simply by removing timed aliquots from the reaction. Given that small deletions (less than 500 bases) are probably desired, we recommend performing the reaction at a lower temperature (between 4–16°C). Samples from the ExoIII reaction are removed at timed intervals to tubes containing S1 nuclease, which removes the remaining single-stranded tails. The low pH and the presence of zinc cations in the S1 buffer effectively inhibit further digestion by ExoIII. After neutralization and heat inactivation of the S1 nuclease, Klenow is added to flush the ends, and the ends are ligated to circularize the vectors. The ligation mixtures are used directly to transform competent cells. Each successive timepoint yields a collection of subclones containing clustered deletions extending further into the original insert. For a more detailed protocol, please see the *Erase-a-Base[®] System Technical Manual* (#TM006; www.promega.com/tbs/).

8. Generation of Single-Stranded DNA and Site-Specific Mutations

8.A. Production of Single-Stranded DNA

To generate single-stranded DNA (ssDNA) from the pCAT[®]3 Vectors, bacterial cells containing pCAT[®]3 Vectors are infected with an appropriate helper phage. The plasmid then enters the f1 replication mode, and the resulting ssDNA is exported from the cell as an encapsulated phage-like particle. The single-stranded plasmid DNA is purified from the supernatant by simple precipitation and extraction procedures (20–22).

8.B. Generation of Site-Specific Mutations

Site-specific mutagenesis, as developed by Hutchinson *et al.* (15), is accomplished by hybridizing a synthetic oligonucleotide that is complementary to the single-stranded template except for a region of mismatch near the center. This region contains the desired nucleotide change or changes. Following hybridization to the single-stranded target DNA, the oligonucleotide is extended with DNA polymerase to create a double-stranded structure. The nick is then sealed with DNA ligase, and the duplex structure is transformed into an *E. coli* host. Theoretically, the yield of mutants using the Hutchinson procedure should be 50% (due to semi-conservative replication). In practice, however, the mutant yield may be much lower, often only a few percent or less. This is presumably due to factors such as incomplete *in vitro* polymerization, primer displacement by the DNA polymerase used in the fill-in reaction, and *in vivo* host-directed mismatch repair mechanisms, which favor repair of the unmethylated newly synthesized DNA strand. Because of the low mutant yield, methods have been developed to increase the mutation frequency (23).

9. Sequencing of pCAT[®]3 Reporter Vectors

It may be desirable to sequence the DNA inserted into the CAT[®] 3 Reporter Vectors. Two examples of such applications are to determine the exact position of generated deletions (Section 7) and to confirm production of a site-specific mutation (Section 8). Two primers are available for sequencing the pCAT[®]3 Vectors: RVprimer3 (Reporter Vector Primer 3) for sequencing clockwise across the upstream cloning sites and RVprimer4 for sequencing counterclockwise across the BamHI and Sall cloning sites downstream of CAT.

RVprimer3 5'-CTAGCAAATAGGCTGTCCC-3'

RVprimer4 5'-GACGATAGTCATGCCCCGCG-3'

RVprimer3 is especially useful for identifying positions of nested deletions.

Note: Both primers can be used for dsDNA sequencing, but only the RVprimer4 may be used for ssDNA sequencing.

10. References

1. Groszkreutz, D.J. *et al.* (1996) CAT reporter systems: New pCAT[®]3 Reporter Vectors and antibodies provide increased expression and detection capabilities. *Promega Notes*, **55**, 2-9.
2. Carswell, S. and Alwine, J.C. (1989) Efficiency of utilization of the simian virus 40 late polyadenylation site: Effects of upstream sequences. *Mol. Cell. Biol.* **9**, 4248-58.
3. Levitt, N. *et al.* (1989) Definition of an efficient synthetic poly(A) site. *Genes Dev.* **3**, 1019-25.
4. Enriquez-Harris, P. *et al.* (1991) A pause site for RNA polymerase II is associated with termination of transcription. *EMBO J.* **10**, 1833-42.
5. Brondyk, B. (1994) pCI and pSI Mammalian Expression Vectors. *Promega Notes* **49**, 7-11.
6. Evans, M.J. and Scarpulla, R.C. (1989) Introns in the 3'-untranslated region can inhibit chimeric CAT and beta-galactosidase gene expression. *Gene* **84**, 135-42.
7. Huang, M.T.F. and Gorman, C.M. (1990) The simian virus 40 small-t intron, present in many common expression vectors, leads to aberrant splicing. *Mol. Cell. Biol.* **10**, 1805-10.
8. Kozak, M. (1989) The scanning model for translation: An update. *J. Cell Biol.* **108**, 229-41.
9. Bothwell, A.L.M. *et al.* (1981) Heavy chain variable region contribution to the NP_b family of antibodies: somatic mutation evident in a gamma 2a variable region. *Cell* **24**, 625-37.
10. Senapathy, P., Shapiro, M.B. and Harris, N.L. (1990) Splice junctions, branch point sites, and exons: Sequence statistics, identification, and applications to genome project. *Methods Enzymol.* **183**, 252-78.
11. Shaw, W.V. (1975) Chloramphenicol acetyltransferase from chloramphenicol-resistant bacteria. *Methods Enzymol.* **43**, 737-55.
12. Gorman, C.M., Moffat, L.F. and Howard, B.H. (1982) Recombinant genomes which express chloramphenicol acetyltransferase in mammalian cells. *Mol. Cell. Biol.* **2**, 1044-51.
13. Cullen, B.R. (1987) Use of eukaryotic expression technology in the functional analysis of cloned genes. *Methods Enzymol.* **152**, 684-704.
14. Henikoff, S. (1987) Unidirectional digestion with exonuclease III in DNA sequence analysis. *Methods Enzymol.* **155**, 156-65.
15. Hutchison, C.A. *et al.* (1978) Mutagenesis at a specific position in a DNA sequence. *J. Biol. Chem.* **253**, 6551-60.
16. Schenborn, E., Goiffon, V. and Oler, J. (1998) An efficient new transfection reagent for eukaryotic cells: TransFast™ Transfection Reagent. *Promega Notes* **65**, 2-6.
17. Ausubel, F.M. *et al.* (1988) *Current Protocols in Molecular Biology*, John Wiley and Sons, NY.
18. Rosenthal, N. (1987) Identification of regulatory elements of cloned genes with functional assays. *Methods Enzymol.* **152**, 704-20.

19. Seed, B. and Sheen, J.Y. (1988) A simple phase-extraction assay for chloramphenicol acyltransferase activity. *Gene* **67**, 271-7.
20. Dotto, G.P., Enea, V. and Zinder, N.D. (1981) Functional analysis of bacteriophage f1 intergenic region. *Virology* **114**, 463-73.
21. Dotto, G.P. and Zinder, N.D. (1983) The morphogenetic signal of bacteriophage f1. *Virology* **130**, 252-6.
22. Dotto, G.P., Huriuchi, K. and Zinder, N.D. (1984) The functional origin of bacteriophage f1 DNA replication. Its signals and domains. *J. Mol. Biol.* **172**, 507-21.
23. Kunkel, T.A. (1985) Rapid and efficient site-specific mutagenesis without phenotypic selection. *Proc. Natl. Acad. Sci. USA* **82**, 488-92.

11. Appendix

11.A. pCAT[®]3-Basic Vector Restriction Sites

The following restriction enzyme tables were constructed using DNASTAR[®] sequence analysis software. Please note that we have not verified this information by restriction digestion with each enzyme listed. The location given specifies the 3'-end of the cut DNA (the base to the left of the cut site). For more information on the cut sites of these enzymes or if you identify a discrepancy, please contact your local Promega Branch or Distributor. In the U.S., contact Promega Technical Services at 800-356-9526. The pCAT[®]3-Basic Vector sequence is available in the GenBank[®] database (GenBank[®]/EMBL Accession Number **U57024**) and at www.promega.com/vectors/

Table 2. Restriction Enzymes That Cut the pCAT[®]3-Basic Vector 1-5 Times.

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
AccB7I	1	736	DraIII	1	3514
AccI	1	1220	DrdI	2	1577, 3558
AccIII	1	500	DsaI	1	289
Acc65I	1	1	EagI	5	60, 269, 964, 968, 3860
AcyI	1	2899	EarI	3	1353, 3157, 3795
AflII	2	66, 263	EclHfKI	1	2362
AflIII	2	15, 1469	Eco47III	1	1345
Alw26I	5	128, 153, 728, 2423, 3199	Eco52I	5	60, 269, 964, 968, 3860
Alw44I	2	1783, 3029	EcoICRI	1	9
AlwNI	1	1885	FokI	5	196, 484, 2328, 2509, 2796
AspHI	4	11, 1787, 2948, 3033	FseI	1	970
AvaI	2	26, 32	FspI	2	2584, 3757
AvaII	2	2500, 2722	HaeII	4	1347, 1717, 3356, 3364
BalI	1	771	HgaI	4	1580, 2158, 2888, 3289
BamHI	1	1213	HincII	2	1111, 1221
BanII	3	11, 33, 3440	HindII	2	1111, 1221
BbsI	2	174, 1298	HindIII	2	53, 276
BglI	2	2482, 3750	HpaI	1	1111
BglII	1	36	Hsp92I	1	2899
BsaI	2	128, 2423	KpnI	1	5
BsaAI	1	3511	MluI	1	15
BsaBI	1	1212	MspAII	5	404, 848, 1811, 2056, 2997
BsaHI	1	2899	NaeI	3	968, 1339, 3408
BsaMI	4	497, 904, 1032, 1125	NciI	5	27, 28, 1849, 2545, 2896
BsmI	4	497, 904, 1032, 1125	NcoI	1	289
BspHI	2	2189, 3197	NgoMIV	3	966, 1337, 3406
BspMI	2	90, 3990	NheI	1	21
BssSI	2	1642, 3026			
Bst98I	2	66, 263			
BstZI	5	60, 269, 964, 968, 3860			
Cfr10I	4	966, 1337, 2442, 3406			
ClaI	3	1206, 3918, 4022			

Note: The enzymes listed in boldface type are available from Promega.

Table 2. Restriction Enzymes That Cut the pCAT[®]3-Basic Vector 1-5 Times (continued).

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
NotI	3	60, 269, 3860	SmaI	1	28
NspI	1	1473	SrfI	1	28
PflMI	1	736	SspI	4	816, 3166, 3719, 3834
PshAI	1	1284	StyI	1	289
PspAI	1	26	TfiI	3	854, 1444, 3915
PstI	1	76	VspI	1	2534
PvuI	2	2732, 3778	XbaI	1	951
PvuII	1	404	XhoI	1	32
SacI	1	11	XmaI	1	26
Sall	1	1219	XmnI	1	2961
ScaI	3	921, 2842, 3925			
SinI	2	2500, 2722			

Table 3. Restriction Enzymes That Do Not Cut the pCAT[®]3-Basic Vector.

AatII	Bpu1102I	DraII	NdeI	RsrII	StuI
AgeI	Bsp120I	Eco72I	NruI	SacII	Swal
ApaI	BsrGI	Eco81I	NsiI	SfiI	Tth111I
AscI	BssHIII	EcoNI	Pacl	Sgfi	Xcml
AvrII	Bst1107I	EcoRI	PinAI	SgrAI	
BbeI	BstEII	EcoRV	PmeI	SnaBI	
BbrPI	BstXI	EheI	PmlI	SpeI	
BbuI	Bsu36I	I-PpoI	Ppu10I	SphI	
BclI	CspI	KasI	PpuMI	SpII	
BlpI	Csp45I	NarI	Psp5II	Sse8387I	

Table 4. Restriction Enzymes That Cut the pCAT[®]3-Basic Vector 6 or More Times.

AcI	BsrSI	DraI	Hsp92II	MspI	ScrFI
AluI	Bst71I	EaeI	MaeI	NdeII	SfaNI
BanI	BstOI	Fnu4HI	MaeII	NlaIII	TaqI
BbvI	BstUI	HaeIII	MaeIII	NlaIV	Tru9I
BsaOI	CfoI	HhaI	MboI	PleI	XhoII
BsaJI	DdeI	HinfI	MboII	RsaI	
Bsp1286I	DpnI	HpaII	MnlI	Sau3AI	
BsrI	DpnII	HphI	MseI	Sau96I	

Note: The enzymes listed in boldface type are available from Promega.

11.B. pCAT[®]3-Enhancer Vector Restriction Sites

The following restriction enzyme tables were constructed using DNASTAR[®] sequence analysis software. Please note that we have not verified this information by restriction digestion with each enzyme listed. The location given specifies the 3'-end of the cut DNA (the base to the left of the cut site). For more information on the cut sites of these enzymes or if you identify a discrepancy, please contact your local Promega Branch or Distributor. In the U.S., contact Promega Technical Services at 800-356-9526. The pCAT[®]3-Enhancer Vector sequence is available in the GenBank[®] database (GenBank[®]/EMBL Accession Number U57026) and at www.promega.com/vectors/

Table 5. Restriction Enzymes That Cut the pCAT[®]3-Enhancer Vector 1-5 Times.

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
AccB7I	1	736	Cfr10I	4	966, 1583, 2688, 3652
AccI	1	1466	ClaI	3	1206, 4164, 4268
AccIII	1	500	DraIII	1	3760
Acc65I	1	1	DrdI	2	1823, 3804
AcyI	1	3145	DsaI	1	289
AflII	1	263	EagI	5	60, 269, 964, 968, 4106
AflIII	2	15, 1715	EarI	3	1599, 3403, 4041
Alw26I	5	127, 152, 728, 2669, 3445	EclHKI	1	2608
Alw44I	2	2029, 3275	Eco47III	1	1591
AlwNI	1	2131	Eco52I	5	60, 269, 964, 968, 4106
AspHI	4	11, 2033, 3194, 3279	EcoICRI	1	9
AvaI	2	26, 32	FseI	1	970
AvaII	2	2746, 2968	FspI	2	2830, 4003
BalI	1	771	HaeII	4	1593, 1963, 3602, 3610
BamHI	1	1459	HgaI	4	1826, 2404, 3134, 3535
BanII	3	11, 33, 3686	HincII	2	1111, 1467
BbsI	2	173, 1544	HindII	2	1111, 1467
BbuI	2	1317, 1389	HindIII	2	53, 276
BglI	2	2728, 3996	HpaI	1	1111
BglII	1	36	Hsp92I	1	3145
BsaI	2	127, 2669	KpnI	1	5
BsaAI	1	3757	MluI	1	15
BsaBI	1	1212	NaeI	3	968, 1585, 3654
BsaHI	1	3145	NciI	5	27, 28, 2095, 2791, 3142
BsaMI	4	497, 904, 1032, 1125	NcoI	1	289
BsmI	4	497, 904, 1032, 1125	NgoMIV	3	966, 1583, 3652
BspHI	2	2435, 3443	NheI	1	21
BspMI	2	89, 4236	NotI	3	60, 269, 4106
BssSI	2	1888, 3272	NsiI	2	1315, 1387
Bst98I	1	263	NspI	3	1317, 1389, 1719
BstZI	5	60, 269, 964, 968, 4106			

Note: The enzymes listed in boldface type are available from Promega.

Table 5. Restriction Enzymes That Cut the pCAT[®]3-Enhancer Vector 1-5 Times (continued).

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
PfI	1	736	SphI	2	1317, 1389
Ppu10I	2	1311, 1383	SrfI	1	28
PshAI	1	1530	SspI	4	816, 3412, 3965, 4080
PspAI	1	26	StyI	1	289
PstI	1	75	TfiI	3	854, 1690, 4161
PvuI	2	2978, 4024	VspI	1	2780
PvuII	1	404	XbaI	1	951
SacI	1	11	XhoI	1	32
SalI	1	1465	XmaI	1	26
ScaI	3	921, 3088, 4171	XmnI	1	3207
SinI	2	2746, 2968			
SmaI	1	28			

Table 6. Restriction Enzymes That Do Not Cut the pCAT[®]3-Enhancer Vector.

AatII	Bpu1102I	Csp45I	KasI	Psp5II	Sse8387I
AgeI	Bsp120I	DraII	NarI	RsrII	StuI
ApaI	BsrGI	Eco72I	NdeI	SacII	SwaI
AscI	BssHIII	Eco81I	NruI	SfiI	Tth111I
AvrII	Bst1107I	EcoNI	PacI	SgfI	XcmI
BbeI	BstEII	EcoRI	PinAI	SgrAI	
BbrPI	BstXI	EcoRV	PmeI	SnaBI	
BclI	Bsu36I	EheI	PmlI	SpeI	
BlpI	CspI	I-PpoI	PpuMI	SplI	

Table 7. Restriction Enzymes That Cut the pCAT[®]3-Enhancer Vector 6 or More Times.

AcI	Bst71I	Fnu4HI	MaeII	NlaIII	Tru9I
AluI	BstOI	FokI	MaeIII	NlaIV	XhoII
BanI	BstUI	HaeIII	MboI	PleI	
BbvI	CfoI	HhaI	MboII	RsaI	
BsaOI	DdeI	Hinfi	MnlI	Sau3AI	
BsaJI	DpnI	HpaII	MseI	Sau96I	
Bsp1286I	DpnII	HphI	MspI	ScrFI	
BsrI	DraI	Hsp92II	MspA1I	SfaNI	
BsrSI	EaeI	MaeI	NdeII	TaqI	

Note: The enzymes listed in boldface type are available from Promega.

11.C. pCAT[®]3-Promoter Vector Restriction Sites

The following restriction enzyme tables were constructed using DNASTAR[®] sequence analysis software. Please note that we have not verified this information by restriction digestion with each enzyme listed. The location given specifies the 3'-end of the cut DNA (the base to the left of the cut site). For more information on the cut sites of these enzymes or if you identify a discrepancy, please contact your local Promega Branch or Distributor. In the U.S., contact Promega Technical Services at 800-356-9526. The pCAT[®]3-Promoter Vector sequence is available in the GenBank[®] database (GenBank[®]/EMBL Accession Number U57027) and at www.promega.com/vectors/

Table 8. Restriction Enzymes That Cut the pCAT[®]3-Promoter Vector 1-5 Times.

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
AccB7I	1	928	Cfr10I	4	1158, 1529, 2634, 3598
AccI	1	1412	ClaI	3	1398, 4110, 4214
AccIII	1	692	DraIII	1	3706
Acc65I	1	1	DrdI	2	1769, 3750
AcyI	1	3091	DsaI	1	481
AflII	2	258, 455	EagI	5	252, 461, 1156, 1160, 4052
AflIII	2	15, 1661	EarI	3	1545, 3349, 3987
Alw26I	5	320, 345, 920, 2615, 3391	EclHfKI	1	2554
Alw44I	2	1975, 3221	Eco47III	1	1537
AlwNI	1	2077	Eco52I	5	252, 461, 1156, 1160, 4052
AspHI	4	11, 1979, 3140, 3225	EcoICRI	1	9
AvaI	2	26, 32	FseI	1	1162
AvaII	2	2692, 2914	FspI	2	2776, 3949
AvrII	1	229	HaeII	4	1539, 1909, 3548, 3556
BalI	1	963	HgaI	4	1772, 2350, 3080, 3481
BamHI	1	1405	HincII	2	1303, 1413
BanII	3	11, 33, 3632	HindII	2	1303, 1413
BbsI	2	366, 1490	HindIII	2	245, 468
BglI	3	182, 2674, 3942	HpaI	1	1303
BglII	1	36	Hsp92I	1	3091
BsaI	2	320, 2615	KpnI	1	5
BsaAI	1	3703	MluI	1	15
BsaBI	2	48, 1404	MspAII	5	596, 1040, 2003, 2248, 3189
BsaHI	1	3091	NaeI	3	1160, 1531, 3600
BsaMI	4	689, 1096, 1224, 1317	NciI	5	27, 28, 2041, 2737, 3088
BsmI	4	689, 1096, 1224, 1317	NcoI	1	481
BspHI	2	2381, 3389	NgoMIV	3	1158, 1529, 3598
BspMI	2	282, 4182	NheI	1	21
BssSI	2	1834, 3218			
Bst98I	2	258, 455			
BstZI	5	252, 461, 1156, 1160, 4052			

Note: The enzymes listed in boldface type are available from Promega.

Table 8. Restriction Enzymes That Cut the pCAT[®]3-Promoter Vector 1-5 Times (continued).

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
NotI	3	252, 461, 4052	SmaI	1	28
NspI	1	1665	SrfI	1	28
PflMI	1	928	SspI	4	1008, 3358, 3911, 4026
PshAI	1	1476	StuI	1	228
PspAI	1	26	StyI	2	229, 481
PstI	1	268	TfiI	3	1046, 1636, 4107
PvuI	2	2924, 3970	VspI	1	2726
PvuII	1	596	XbaI	1	1143
SacI	1	11	XhoI	1	32
Sall	1	1411	XmaI	1	26
ScaI	3	1113, 3034, 4117	XmnI	1	3153
SfiI	1	182			
SinI	2	2692, 2914			

Table 9. Restriction Enzymes That Do Not Cut the pCAT[®]3-Promoter Vector.

AatII	Bpu1102I	Csp45I	KasI	PpuMI	SpII
AgeI	Bsp120I	DraII	NarI	Psp5II	Sse8387I
ApaI	BsrGI	Eco72I	NdeI	RsrII	Swal
AscI	BssHII	Eco81I	NruI	SacII	Tth111I
BbeI	Bst1107I	EcoNI	NsiI	SgfI	XcmI
BbrPI	BstEII	EcoRI	PacI	SgrAI	
BbuI	BstXI	EcoRV	PinAI	SnaBI	
BclI	Bsu36I	EheI	PmeI	SpeI	
BlpI	CspI	I-PpoI	PmlI	SphI	

Table 10. Restriction Enzymes That Cut the pCAT[®]3-Promoter Vector 6 or More Times.

AcI	BsrSI	DraI	HphI	MseI	Sau96I
AluI	Bst7II	EaeI	Hsp92II	MspI	ScrFI
BanI	BstOI	Fnu4HI	MaeI	NdeII	SfaNI
BbvI	BstUI	FokI	MaeII	NlaIII	TaqI
BsaOI	CfoI	HaeIII	MaeIII	NlaIV	Tru9I
BsaJI	DdeI	HhaI	MboI	PleI	XhoII
Bsp1286I	DpnI	HinfI	MboII	RsaI	
BsrI	DpnII	HpaII	MnlI	Sau3AI	

Note: The enzymes listed in boldface type are available from Promega.

11.D. pCAT[®]3-Control Vector Restriction Sites

The following restriction enzyme tables were constructed using DNASTAR[®] sequence analysis software. Please note that we have not verified this information by restriction digestion with each enzyme listed. The location given specifies the 3'-end of the cut DNA (the base to the left of the cut site). For more information on the cut sites of these enzymes or if you identify a discrepancy, please contact your local Promega Branch or Distributor. In the U.S., contact Promega Technical Services at 800-356-9526. The pCAT[®]3-Control Vector sequence is available in the GenBank[®] database (GenBank[®]/EMBL Accession Number U57025) and at www.promega.com/vectors/

Table 11. Restriction Enzymes That Cut the pCAT[®]3-Control Vector 1-5 Times.

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
AccB7I	1	928	Cfr10I	4	1158, 1775, 2880, 3844
AccI	1	1658	ClaI	3	1398, 4356, 4460
AccIII	1	692	DraIII	1	3952
Acc65I	1	1	DrdI	2	2015, 3996
AcyI	1	3337	DsaI	1	481
AflII	2	258, 455	EagI	5	252, 461, 1156, 1160, 4298
AflIII	2	15, 1907	EarI	3	1791, 3595, 4233
Alw26I	5	320, 345, 920, 2861, 3637	EclHKI	1	2800
Alw44I	2	2221, 3467	Eco47III	1	1783
AlwNI	1	2323	Eco52I	5	252, 461, 1156, 1160, 4298
AspHI	4	11, 2225, 3386, 3471	EcoICRI	1	9
AvaI	2	26, 32	FseI	1	1162
AvaII	2	2938, 3160	FspI	2	3022, 4195
AvrII	1	229	HaeII	4	1785, 2155, 3794, 3802
BalI	1	963	HgaI	4	2018, 2596, 3326, 3727
BamHI	1	1651	HincII	2	1303, 1659
BanII	3	11, 33, 3878	HindII	2	1303, 1659
BbsI	2	366, 1736	HindIII	2	245, 468
BbuI	2	1509, 1581	HpaI	1	1303
BglI	3	182, 2920, 4188	Hsp92I	1	3337
BglII	1	36	KpnI	1	5
BsaI	2	320, 2861	MluI	1	15
BsaAI	1	3949	NaeI	3	1160, 1777, 3846
BsaBI	2	48, 1404	NciI	5	27, 28, 2287, 2983, 3334
BsaHI	1	3337	NcoI	1	481
BsaMI	4	689, 1096, 1224, 1317	NgoMIV	3	1158, 1775, 3844
BsmI	4	689, 1096, 1224, 1317	NheI	1	21
BspHI	2	2627, 3635	NotI	3	252, 461, 4298
BspMI	2	282, 4428	NsiI	2	1507, 1579
BssSI	2	2080, 3464	NspI	3	1509, 1581, 1911
Bst98I	2	258, 455	PfMI	1	928
BstZI	5	252, 461, 1156, 1160, 4298			

Note: The enzymes listed in boldface type are available from Promega.

Table 11. Restriction Enzymes That Cut the pCAT[®]3-Control Vector 1-5 Times (continued).

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
Ppu10I	2	1503, 1575	SphI	2	1509, 1581
PshAI	1	1722	SrfI	1	28
PspAI	1	26	SspI	4	1008, 3604, 4157, 4272
PstI	1	268	StuI	1	228
PvuI	2	3170, 4216	StyI	2	229, 481
PvuII	1	596	TfiI	3	1046, 1882, 4353
SacI	1	11	VspI	1	2972
SalI	1	1657	XbaI	1	1143
ScaI	3	1113, 3280, 4363	XhoI	1	32
SfiI	1	182	XmaI	1	26
SinI	2	2938, 3160	XmnI	1	3399
SmaI	1	28			

Table 12. Restriction Enzymes That Do Not Cut the pCAT[®]3-Control Vector.

AatII	Bpu1102I	CspI	EheI	PmeI	SnaBI
AgeI	Bsp120I	Csp45I	I-PpoI	PmlI	SpeI
ApaI	BsrGI	DraII	KasI	PpuMI	SplI
AscI	BssHIII	Eco72I	NarI	Psp5II	Sse8387I
BbeI	Bst1107I	Eco81I	NdeI	RsrII	SwaI
BbrPI	BstEII	EcoNI	NruI	SacII	Tth111I
BclI	BstXI	EcoRI	PacI	SgfI	XcmI
BlpI	Bsu36I	EcoRV	PinAI	SgrAI	

Table 13. Restriction Enzymes That Cut the pCAT[®]3-Control Vector 6 or More Times.

AciI	BsrSI	DraI	HphI	MseI	Sau3AI
AluI	Bst71I	EaeI	Hsp92II	MspI	Sau96I
BanI	BstOI	Fnu4HI	MaeI	MspAII	ScrFI
BbvI	BstUI	FokI	MaeII	NdeII	SfaNI
BsaOI	CfoI	HaeIII	MaeIII	NlaIII	TaqI
BsaJI	DdeI	HhaI	MboI	NlaIV	Tru9I
Bsp1286I	DpnI	HinfI	MboII	PleI	XhoII
BsrI	DpnII	HpaII	MnlI	RsaI	

Note: The enzymes listed in boldface type are available from Promega.

11.E. Related Products

CAT Assay Reagents

Product	Size	Cat.#
CAT Enzyme Assay System With Reporter Lysis Buffer	50 reactions	E1000
Reporter Lysis Buffer, 5X	30ml	E3971
Chloramphenicol Acetyltransferase (CAT)	100u	E1051

Reporter Vector Sequencing Primers

Product	Size	Cat.#
RVprimer3 (clockwise)	2µg	E4481
RVprimer4 (counterclockwise)	2µg	E4491

Mammalian Transfection Systems

Product	Size	Cat.#
TransFast™ Transfection Reagent	1.2mg	E2431
Tfx™-50 Reagent	2.1mg	E1811
Tfx™-20 Reagent	4.8mg	E2391
ProFection® Mammalian Transfection System Calcium Phosphate	40 reactions	E1200

Competent Cells

Product	Size	Cat.#
JM109 Competent Cells, >10 ⁸ cfu/µg*	1ml	L2001
JM109 Competent Cells, >10 ⁷ cfu/µg	1ml	L1001

*For Laboratory Use. Each 1ml order is provided as 5 × 200µl aliquots. Also included with each order is pGEM®-3Z Competent Cell Control DNA.

pGL3 Luciferase Reporter Vectors

Please visit www.promega.com to see a complete listing of our reporter vectors.

Product	Size	Cat.#
pGL3-Control Vector	20µg	E1741
pGL3-Enhancer Vector	20µg	E1771
pGL3-Promoter Vector	20µg	E1761
pGL3-Basic Vector	20µg	E1751

pGL4 Luciferase Reporter Vectors

Vector	Multiple Cloning Region	Reporter Gene	Protein Sequence	Reporter Gene		Mam-alian Selectable Marker	Cat.#
				Promoter/Element	Degradation Response		
pGL4.10[luc2]	Yes	<i>luc2^A</i>	No	No	No	No	E6651
pGL4.11[luc2P]	Yes	"	hPEST	No	No	No	E6661
pGL4.12[luc2CP]	Yes	"	hCL1-hPEST	No	No	No	E6671
pGL4.13[luc2/SV40]	No	"	No	SV40	No	No	E6681
pGL4.14[luc2/Hygro]	Yes	"	No	No	Hygro	No	E6691
pGL4.15[luc2P/Hygro]	Yes	"	hPEST	No	Hygro	No	E6701
pGL4.16[luc2CP/Hygro]	Yes	"	hCL1-hPEST	No	Hygro	No	E6711
pGL4.17[luc2/Neo]	Yes	"	No	No	Neo	No	E6721
pGL4.18[luc2P/Neo]	Yes	"	hPEST	No	Neo	No	E6731
pGL4.19[luc2CP/Neo]	Yes	"	hCL1-hPEST	No	Neo	No	E6741
pGL4.20[luc2/Puro]	Yes	"	No	No	Puro	No	E6751
pGL4.21[luc2P/Puro]	Yes	"	hPEST	No	Puro	No	E6761
pGL4.22[luc2CP/Puro]	Yes	"	hCL1-hPEST	No	Puro	No	E6771
pGL4.23[luc2/minP]	Yes	"	No	minP	No	No	E6691
pGL4.24[luc2P/minP]	Yes	"	hPEST	"	No	No	E6701
pGL4.25[luc2CP/minP]	Yes	"	hCL1-hPEST	"	No	No	E6711
pGL4.26[luc2/minP/Hygro]	Yes	"	No	"	Hygro	No	E6721
pGL4.27[luc2P/minP/Hygro]	Yes	"	hPEST	"	Hygro	No	E6731
pGL4.28[luc2CP/minP/Hygro]	Yes	"	hCL1-hPEST	"	Hygro	No	E6741
pGL4.29[luc2P/CRE/Hygro]	No	"	hPEST	CRE	Hygro	No	E6751
pGL4.30[luc2P/NFAT-RE/Hygro]	No	"	hPEST	NFAT-RE	Hygro	No	E6761
pGL4.31[luc2P/GAL4UAS/Hygro]	No	"	hPEST	GAL4UAS	Hygro	No	E6771
pGL4.70[hRluc]	Yes	<i>hRluc^B</i>	No	No	No	No	E6881
pGL4.71[hRlucP]	Yes	"	hPEST	No	No	No	E6891
pGL4.72[hRlucCP]	Yes	"	hCL1-hPEST	No	No	No	E6901
pGL4.73[hRluc/SV40]	No	"	No	SV40	No	No	E6911
pGL4.74[hRluc/TK]	No	"	No	HSV-TK	No	No	E6921
pGL4.75[hRluc/CMV]	No	"	No	CMV	No	No	E6931
pGL4.76[hRluc/Hygro]	Yes	"	No	No	Hygro	No	E6941
pGL4.77[hRlucP/Hygro]	Yes	"	hPEST	No	Hygro	No	E6951
pGL4.78[hRlucCP/Hygro]	Yes	"	hCL1-hPEST	No	Hygro	No	E6961
pGL4.79[hRluc/Neo]	Yes	"	No	No	Neo	No	E6971
pGL4.80[hRlucP/Neo]	Yes	"	hPEST	No	Neo	No	E6981
pGL4.81[hRlucCP/Neo]	Yes	"	hCL1-hPEST	No	Neo	No	E6991
pGL4.82[hRluc/Puro]	Yes	"	No	No	Puro	No	E7501
pGL4.83[hRlucP/Puro]	Yes	"	hPEST	No	Puro	No	E7511
pGL4.84[hRlucCP/Puro]	Yes	"	hCL1-hPEST	No	Puro	No	E7521

^A*luc2* = synthetic firefly luciferase gene. ^B*hRluc* = synthetic *Renilla* luciferase gene.

DNA Purification Products

Product	Size	Cat.#
Wizard® SV Gel and PCR Clean-Up System*	50 preps	A9281
	250 preps	A9282
Wizard® <i>Plus</i> SV Minipreps DNA Purification System*	50 preps	A1330
	250 preps	A1460
PureYield™ Plasmid Midipreps System	50 preps	A2492
	100 preps	A2495

*For Laboratory Use

© 1995-2008 Promega Corporation. All Rights Reserved.

Erase-a-Base, pCAT, ProFection, Transfectam and Wizard are registered trademarks of Promega Corporation. PureYield, TransFast and Tfx are trademarks of Promega Corporation.

Bacto is a registered trademark of Difco Laboratories, Detroit, Michigan. DNASTAR is a registered trademark of DNASTAR, Inc. GenBank is a registered trademark of the US Dept of Health and Human Services.

Products may be covered by pending or issued patents or may have certain limitations. Please visit our Web site for more information.

All prices and specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.