



Promega

Technical Bulletin

pGEM[®]-5Zf(+) Vector

INSTRUCTIONS FOR USE OF PRODUCT P2241.



www.promega.com

PRINTED IN USA.
Revised 11/06



AF9TB047 1106TB047

Part# TB047

pGEM[®]-5Zf(+) Vector

All technical literature is available on the Internet at: www.promega.com/tbs/
 Please visit the web site to verify that you are using the most current version of this
 Technical Bulletin. Please contact Promega Technical Services if you have questions on use
 of this system. E-mail: techserv@promega.com

| | |
|---|---|
| I. Description..... | 1 |
| II. Product Components and Storage Conditions | 2 |
| III. pGEM [®] -5Zf(+) Vector Multiple Cloning Region and Map | 2 |
| IV. pGEM [®] -5Zf(+) Vector Restriction Enzyme Sites..... | 4 |
| V. Related Products | 6 |
| VI. Reference | 7 |

I. Description

The pGEM[®]-5Zf(+) Vector is a derivative of the pGEM[®]-3Zf(+) Vector and contains the origin of replication of the filamentous phage f1. The plasmid serves as a standard cloning vector, as a template for in vitro transcription, and can be used for the production of circular ssDNA. pGEM[®]-5Zf(+) contains T7 and SP6 RNA polymerase promoters flanking a multiple cloning region within the α -peptide coding region of β -galactosidase (1). Insertional inactivation of the α -peptide allows recombinant clones to be directly identified by color screening on indicator plates. The multiple cloning region contains unique restriction sites for ApaI, AatII, SphI, NcoI, SacII, EcoRV, SpeI, NotI, PstI, Sall, NdeI, SacI, BstXI and NsiI. This arrangement is designed specifically for generating unidirectional deletions with the Erase-a-Base[®] System. The polylinker contains restriction enzyme sites that produce 5' overhangs or blunt ends (sensitive to Exonuclease III), flanked on both sides by blocks of restriction sites that generate 3' overhangs (resistant to Exonuclease III).

For induction of ssDNA, bacterial cells containing pGEM[®]-5Zf(+) recombinants are infected with an appropriate helper phage. The plasmid then enters the f1 replication mode, and the resulting ssDNA is exported from the cell as an encapsidated virus-like particle. The sequence of the ssDNA rescued upon infection with helper phage is complementary to the sequence shown in Figure 1. The exported ssDNA can be used for mutagenesis in vitro or can be sequenced using the T7 Promoter Primer (Cat.# Q5021) and pUC/M13 Forward Primer (Cat.# Q5391, Q5601).

Promega vectors sequences are available online at: www.promega.com/vectors/
 and from the GenBank[®] database.

II. Product Components and Storage Conditions

| Product | Size | Cat. # |
|----------------------------------|------|--------|
| pGEM [®] -5Zf(+) Vector | 20µg | P2241 |

The pGEM[®]-5Zf(+) Vector is supplied with a glycerol stock of bacterial strain JM109. The JM109 cells do not contain the vector and are not competent cells.

Storage Conditions: Store the pGEM[®]-5Zf(+) Vector at -20°C and the glycerol stock of JM109 cells at -70°C.

III. pGEM[®]-5Zf(+) Vector Multiple Cloning Region and Map




Figure 1. pGEM[®]-5Zf(+) Vector promoter and multiple cloning region sequence. The sequence shown corresponds to RNA synthesized by T7 RNA polymerase and is complementary to RNA synthesized by SP6 RNA polymerase. The strand shown is complementary to the ssDNA produced by this vector.




Figure 2. pGEM®-5Zf(+) Vector map.

! The pGEM®-5Zf(+) and pGEM®-5Zf(-) Vectors are identical except for the orientation of the f1 origin.

pGEM®-5Zf(+) Vector sequence reference points:

| | |
|---|--------------------|
| T7 RNA Polymerase transcription initiation site | 1 |
| Multiple cloning region | 10-113 |
| SP6 RNA Polymerase promoter (-17 to +3) | 124-143 |
| SP6 RNA Polymerase transcription initiation site | 126 |
| Binding site of pUC/M13 Reverse Sequencing Primer | 161-177 |
| <i>lacZ</i> start codon | 165 |
| <i>lac</i> operon sequences | 151-380; 2821-2981 |
| <i>lac</i> operator | 185-201 |
| β -lactamase coding region | 1322-2182 |
| phage f1 region | 2365-2820 |
| binding site of pUC/M13 Forward Sequencing Primer | 2941-2957 |
| T7 RNA Polymerase promoter (-17 to +3) | 2984-3 |

! Use the T7 or pUC/M13 Forward Primer to sequence ssDNA produced by the pGEM®-5Zf(+) Vector.

Specialized applications of the pGEM®-5Zf(+) Vector:

- used with the Erase-a-Base® System (For information, please request Technical Manual #TM006.)
- ssDNA production
- blue/white screening for recombinants
- transcription in vitro from dual-opposed promoters (For protocol information, please request the Riboprobe® in vitro Transcription Systems Technical Manual #TM016.)

IV. pGEM[®]-5Zf(+) Vector Restriction Enzyme Sites

The following restriction enzyme tables were constructed using DNASTAR[®] sequence analysis software. Please note that we have not verified this information by restriction digestion with each enzyme listed. The location given specifies the 3' end of the cut DNA (the base to the left of the cut site). For more information on the cut sites of these enzymes, or if you identify a discrepancy, please contact your local Promega Branch or Distributor. In the U.S., contact Promega Technical Services at 800-356-9526. Vector sequences are also available in the GenBank[®] database (GenBank[®]/EMBL Accession Number X65308) and on the Internet at: www.promega.com/vectors/

Table 1. Restriction Enzymes That Cut the pGEM[®]-5Zf(+) Vector Between 1 and 5 Times.

| Enzyme | # of Sites | Location | Enzyme | # of Sites | Location |
|----------------|------------|-----------------------------|----------------|------------|--------------------------------|
| AatII | 1 | 20 | DrdI | 2 | 610, 2544 |
| AccI | 1 | 76 | DsaI | 2 | 37, 43 |
| AcyI | 2 | 17, 1932 | EagI | 2 | 31, 62 |
| AflIII | 2 | 99, 502 | EarI | 3 | 386, 2190, 2878 |
| Alw26I | 2 | 1456, 2232 | EclHKI | 1 | 1395 |
| Alw44I | 2 | 816, 2062 | Eco52I | 2 | 31, 62 |
| AlwNI | 1 | 918 | EcoICRI | 1 | 92 |
| Apal | 1 | 14 | EcoRV | 1 | 51 |
| AspHI | 4 | 94, 820, 1981, 2066 | FokI | 5 | 119, 1361, 1542, 1829, 2916 |
| AvaII | 2 | 1533, 1755 | FspI | 2 | 1617, 2840 |
| BanI | 3 | 246, 1343, 2626 | HaeII | 4 | 380, 750, 2740, 2748 |
| BanII | 3 | 14, 94, 2664 | HgaI | 4 | 613, 1191, 1921, 2806 |
| BbuI | 1 | 26 | HincII | 1 | 77 |
| BglI | 3 | 39, 1515, 2833 | HindIII | 1 | 77 |
| BsaI | 1 | 1456 | Hsp92I | 2 | 17, 1932 |
| BsaAI | 1 | 2589 | MaeI | 5 | 56, 997, 1250, 1585, 2740 |
| BsaHI | 2 | 17, 1932 | MluI | 1 | 99 |
| BsaJI | 5 | 37, 43, 241, 662, 2936 | NaeI | 1 | 2692 |
| Bsp120I | 1 | 10 | NciI | 4 | 30, 882, 1578, 1929 |
| BspHI | 2 | 1222, 2230 | NcoI | 1 | 37 |
| BspMI | 1 | 62 | NdeI | 1 | 82 |
| BssSI | 2 | 675, 2059 | NgoMIV | 1 | 2690 |
| BstOI | 5 | 242, 530, 651, 664, 2937 | NotI | 1 | 62 |
| BstXI | 1 | 103 | NsiI | 1 | 112 |
| BstZI | 2 | 31, 62 | NspI | 2 | 26, 506 |
| Cfr10I | 2 | 1475, 2690 | Ppu10I | 1 | 108 |
| DdeI | 4 | 777, 1186, 1352, 1892 | PstI | 1 | 73 |
| DraI | 3 | 1261, 1280, 1972 | PvuI | 2 | 1765, 2861 |
| DraIII | 1 | 2589 | PvuII | 2 | 326, 2890 |

Table 1. Restriction Enzymes That Cut the pGEM®-5Zf(+) Vector Between 1 and 5 Times (continued).

| Enzyme | # of Sites | Location | Enzyme | # of Sites | Location |
|--------|------------|------------|----------|------------|------------------------|
| RsaI | 1 | 1875 | Sse8387I | 1 | 73 |
| SacI | 1 | 94 | SspI | 2 | 2199, 2381 |
| SacII | 1 | 46 | StyI | 1 | 37 |
| Sall | 1 | 75 | TaqI | 4 | 76, 602, 2046, 2622 |
| ScaI | 1 | 1875 | TfiI | 2 | 337, 477 |
| SfiI | 1 | 39 | VspI | 3 | 273, 332, 1567 |
| SinI | 2 | 1533, 1755 | XmnI | 1 | 1994 |
| SpeI | 1 | 55 | | | |
| SphI | 1 | 26 | | | |

Table 2. Restriction Enzymes That Do Not Cut the pGEM®-5Zf(+) Vector.

| | | | | | |
|---------------|---------------|-----------------|----------------|-------------|----------------|
| AccB7I | BbsI | Bst98I | EheI | PfiMI | SnaBI |
| AccIII | BclI | BstEII | FseI | PinAI | SplI |
| Acc65I | BglIII | Bsu36I | HindIII | PmeI | SrfI |
| AflIII | BlpI | ClaI | HpaI | PmlI | StuI |
| AgeI | Bpu1102I | CspI | I-PpoI | PpuMI | Swal |
| AscI | BsaBI | Csp45I | KasI | PshAI | Tth111I |
| AvaI | BsaMI | Drall | KpnI | Psp5II | XbaI |
| AvrII | BsmI | Eco47III | NarI | PspAI | XcmI |
| BalI | BsrBRI | Eco72I | NheI | RsrII | XhoI |
| BamHI | BsrGI | Eco8II | NruI | SgfI | XmaI |
| BbeI | BssHII | EcoNI | PacI | SgrAI | |
| BbrPI | Bst1107I | EcoRI | PaeR7I | SmaI | |

Table 3. Restriction Enzymes That Cut the pGEM®-5Zf(+) Vector 6 or More Times.

| | | | | | |
|-----------------|-------------|----------------|---------------|---------------|--------------|
| AcI | Bst7II | HaeIII | MaellI | NdeII | SfaNI |
| AluI | BstUI | HhaI | MboI | NlaIII | Tru9I |
| BbvI | CfoI | HinfI | MboII | NlaIV | XhoII |
| BsaOI | DpnI | HpaII | MnlI | PleI | |
| Bsp1286I | DpnII | HphI | MseI | Sau3AI | |
| BsrI | EaeI | Hsp92II | MspI | Sau96I | |
| BsrSI | Fnu4HI | MaeII | MspAII | ScrFI | |

Note: The enzymes listed in boldface type are available from Promega.

V. Related Products

Vectors

| Product | Size | Cat. # |
|-----------------------------------|------|--------|
| pGEM [®] -3Z Vector | 20μg | P2151 |
| pGEM [®] -4Z Vector | 20μg | P2161 |
| pGEM [®] -3Zf(+) Vector | 20μg | P2271 |
| pGEM [®] -3Zf(-) Vector | 20μg | P2261 |
| pGEM [®] -5Zf(-) Vector | 20μg | P2351 |
| pGEM [®] -7Zf(+) Vector | 20μg | P2251 |
| pGEM [®] -7Zf(-) Vector | 20μg | P2371 |
| pGEM [®] -9Zf(-) Vector | 20μg | P2391 |
| pGEM [®] -11Zf(+) Vector | 20μg | P2411 |
| pGEM [®] -11Zf(-) Vector | 20μg | P2421 |
| pGEM [®] -13Zf(+) Vector | 20μg | P2541 |

All pGEM[®] Vectors are provided with a glycerol stock of bacterial strain JM109. The JM109 cells do not contain vector and are not competent cells.

| Product | Size | Cat. # |
|----------------------|------|--------|
| pSP64 Poly(A) Vector | 20μg | P1241 |
| pSP72 Vector | 20μg | P2191 |
| pSP73 Vector | 20μg | P2221 |

Competent Cells

| Product | Size | Cat. # |
|---|-----------------|--------|
| Single Step (KRX) Competent Cells | 5 × 200μl | L3001 |
| | 20 × 50μl | L3002 |
| JM109 Competent Cells, >10 ⁸ cfu/μg* | 1ml (5 × 200μl) | L2001 |
| JM109 Competent Cells, >10 ⁷ cfu/μg | 1ml (5 × 200μl) | L1001 |
| HB101 Competent Cells, >10 ⁸ cfu/μg | 1ml (5 × 200μl) | L2011 |
| HB101 Competent Cells, >10 ⁷ cfu/μg | 1ml (5 × 200μl) | L1011 |

*For Laboratory Use.

Riboprobe[®] in vitro Transcription Systems

| Product | Size | Cat. # |
|-------------------------------------|----------|--------|
| Riboprobe [®] System – SP6 | 1 system | P1420 |
| Riboprobe [®] System – T3 | 1 system | P1430 |
| Riboprobe [®] System – T7 | 1 system | P1440 |

For Laboratory Use.

Sequencing Primers

| Product | Size | Cat. # |
|---|----------|--------|
| SP6 Promoter Primer | 2µg | Q5011 |
| T7 Promoter Primer | 2µg | Q5021 |
| pUC/M13 Primer, Reverse (17mer) | 2µg | Q5401 |
| pUC/M13 Primer, Forward (17mer) | 2µg | Q5391 |
| pUC/M13 Primer, Forward (24mer) | 2µg | Q5601 |
| pUC/M13 Primer, Reverse (22mer) | 2µg | Q5421 |
| Erase-a-Base® System (minus vectors and bacterial strain) | 1 system | E5750 |

RiboMAX™ Large Scale RNA Production Systems

| Product | Size | Cat.# |
|--|----------|-------|
| RiboMAX™ Large Scale RNA Production System – SP6 | 1 system | P1280 |
| RiboMAX™ Large Scale RNA Production System – T7 | 1 system | P1300 |

For Laboratory Use.

VI. Reference

1. Yanisch-Perron, C., Vieira, J., Messing, J. (1985) Improved M13 phage cloning vectors and host strains: Nucleotide sequences of the M13mp18 and pUC19 vectors. *Gene* **33**, 103-19.

© 1988-2006 Promega Corporation. All Rights Reserved.

Erase-a-Base, pGEM and Riboprobe are registered trademarks of Promega Corporation. RiboMAX is a trademark of Promega Corporation.

DNASTAR is a registered trademark of DNASTAR, Inc. GenBank is a registered trademark of the U.S. Department of Health and Human Services.

Products may be covered by pending or issued patents or may have certain limitations. Please visit our Web site for more information.

All prices and specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.