



Promega

Technical Bulletin

pGEM[®]-9Zf(-) Vector

INSTRUCTIONS FOR USE OF PRODUCT P2391.



www.promega.com

PRINTED IN USA.
Revised 5/07



AF9TB070 0507TB070

Part# TB070

pGEM[®]-9Zf(-) Vector

All technical literature is available on the Internet at: www.promega.com/tbs/
Please visit the web site to verify that you are using the most current version of this
Technical Bulletin. Please contact Promega Technical Services if you have questions on use
of this system. E-mail: techserv@promega.com

I. Description	1
II. Product Components and Storage Conditions	2
III. pGEM[®]-9Zf(-) Vector Multiple Cloning Region and Map	2
IV. pGEM[®]-9Zf(-) Vector Restriction Sites	4
V. Related Products	6
VI. Reference	7

I. Description

The pGEM[®]-9Zf(-) Vector is a recombinant plasmid designed to provide a versatile range of cloning strategies, efficient synthesis of RNA *in vitro* and the production of single-stranded DNA. The plasmid contains T7 and SP6 RNA polymerase promoters flanking a multiple cloning region within the α -peptide coding region of β -galactosidase (1). Insertional inactivation of the α -peptide allows recombinant clones to be directly identified by color screening on indicator plates. The multiple cloning region is unique and includes restriction sites for NsiI, SpeI, HindIII, XbaI, EcoRI, SalI and SacI.

For induction of ssDNA, bacterial cells containing pGEM[®]-9Zf(-) recombinants are infected with an appropriate helper phage. The plasmid then enters the f1 replication mode, and the resulting ssDNA is exported from the cell as an encapsidated virus-like particle. The sequence of the ssDNA produced is identical to the sequence shown in Figure 1. The exported ssDNA can be used for mutagenesis *in vitro* or can be sequenced using the Promega T7 Promoter Primer or pUC/M13 Reverse Primer.

Promega vector sequences are available online at: www.promega.com/vectors/ and are also available from the GenBank[®] database.

II. Product Components and Storage Conditions

Product	Size	Cat.#
pGEM [®] -9Zf(-) Vector	20μg	P2391

The pGEM[®]-9Zf(-) Vector is provided with a glycerol stock of bacterial strain JM109. The JM109 cells do not contain the vector and are not competent.

Storage Conditions: Store the pGEM[®]-9Zf(-) Vector at -20°C and the glycerol stock of JM109 cells at -70°C.

III. pGEM[®]-9Zf(-) Vector Multiple Cloning Region and Map





Figure 1. pGEM[®]-9Zf(-) Vector promoter and multiple cloning region sequence. The sequence shown corresponds to RNA synthesized by SP6 RNA polymerase and is complementary to RNA synthesized by T7 RNA polymerase. The strand shown is the same as the ssDNA strand produced by this vector.



0288VA05_4A

Figure 2. pGEM®-9Zf(-) Vector map and sequence reference points.

pGEM®-9Zf(-) Vector sequence reference points:

SP6 RNA polymerase transcription initiation site	1
T7 RNA polymerase promoter (-17 to +3)	53-72
T7 RNA polymerase transcription initiation site	55
<i>lac</i> operon sequences	93-321, 2712-2868
binding site of pUC/M13 Reverse Sequencing Primer	102-118
<i>lacZ</i> start codon	106
<i>lac</i> operator	126-142
β-lactamase (Amp ^r) coding region	1263-2123
phage f1 region	2255-2710
binding site of pUC/M13 Forward Sequencing Primer	2831-2847
multiple cloning region	2876-87
SP6 RNA polymerase promoter (-17 to +3)	2896-3

Note: Use the T7 or pUC/M13 Reverse Primer to sequence ssDNA produced by the pGEM®-9Zf(-) Vector.

Specialized applications of the pGEM®-9Zf(-) Vector:

- Allows excision of insert containing SP6 and T7 promoters
- ssDNA production
- Blue/white screening for recombinants
- Transcription *in vitro* from dual-opposed promoters (For protocol information please request the *Riboprobe® in vitro Transcription Systems Technical Manual*, #TM016.)

IV. pGEM[®]-9Zf(-) Vector Restriction Sites

The following restriction enzyme tables were constructed using DNASTAR[®] sequence analysis software. Please note that we have not verified this information by restriction digestion with each enzyme listed. The location given specifies the 3' end of the cut DNA (the base to the left of the cut site). For more information on the cut sites of these enzymes, or if you identify a discrepancy, please contact your local Promega Branch or Distributor. In the U.S., contact Promega Technical Services at 800-356-9526. Vector sequences are also available through GenBank[®] (GenBank[®]/EMBL Accession Number X65312) and on the Internet at: www.promega.com/vectors/

Table 1. Restriction Enzymes That Cut the pGEM[®]-9Zf(-) Vector Between 1 and 5 Times.

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
AccI	1	43	HaeII	4	321, 691, 2330, 2338
AflIII	1	443	HgaI	4	554, 1132, 1862, 2263
Alw26I	2	1397, 2173	HincII	1	44
Alw44I	2	757, 2003	HindIII	1	18
AlwNI	1	859	Hsp92I	1	1873
AvaII	2	1474, 1696	NaeI	1	2382
BalI	1	85	NciI	3	823, 1519, 1870
BanI	3	187, 1284, 2444	NgoMIV	1	2380
BanII	2	52, 2414	NotI	1	2887
BglI	3	82, 1456, 2723	NsiI	1	9
BsaI	1	1397	NspI	1	447
BsaAI	1	2485	PvuI	2	1706, 2751
BsaII	3	182, 603, 2826	PvuII	2	267, 2780
Bsp1286I	5	52, 761, 1922, 2007, 2414	RsaI	1	1816
BspHI	2	1163, 2171	SacI	1	52
BssSI	2	616, 2000	SalI	1	42
BstOI	5	183, 471, 592, 605, 2827	ScaI	1	1816
BstZI	1	2887	SfaNI	5	16, 540, 1592, 1783, 2032
DdeI	4	718, 1127, 1293, 1833	SfiI	1	82
DraI	3	1202, 1221, 1913	SinI	2	1474, 1696
DraIII	1	2488	SpeI	1	12
DrdI	2	551, 2532	SspI	2	2140, 2693
EaeI	5	83, 282, 1724, 2860, 2887	TaqI	4	43, 543, 1987, 2450
EarI	3	327, 2131, 2768	TfiI	2	278, 418
EclHKI	1	1336	Tth111I	1	2879
EcoICRI	1	50	VspI	3	214, 273, 1508
EcoRI	1	36	XbaI	1	27
FokI	4	1302, 1483, 1770, 2806	XmnI	1	1935
FspI	2	1558, 2730			

Table 2. Restriction Enzymes That Do Not Cut the pGEM[®]-9Zf(-) Vector.

AatII	BclI	Csp45I	NdeI	SnaBI
Acc65I	BglII	CspI	NheI	SphI
AccB7I	BlpI	Eco47III	NruI	StuI
AccIII	BsaBI	EcoNI	Pacl	StyI
AgeI	BsaMI	EcoRV	PmeI	Swal
ApaI	BspMI	FseI	PmlI	XcmI
AscI	BsrGI	HpaI	PpuMI	XhoI
AvaI	BssHII	I-PpoI	PshAI	XmaI
AvrII	Bst98I	KasI	PstI	
BamHI	BstEII	KpnI	SacII	
BbeI	BstXI	MluI	SgfI	
BbsI	Bsu36I	NarI	SgrAI	
BbuI	Clal	NcoI	SmaI	

Table 3. Restriction Enzymes That Cut the pGEM[®]-9Zf(-) Vector 6 or More Times.

AcII	DpnI	HphI	MspAII	Sau96I
AluI	Fnu4HI	Hsp92II	MspI	ScrFI
BbvI	HaeIII	MaeIII	NdeII	Tru9I
BsrSI	HhaI	MboI	NlaIV	XhoII
BstUI	Hinfl	MboII	PleI	
CfoI	HpaII	MnlI	Sau3AI	

V. Related Products

Product	Size	Cat. #
pGEM [®] -3Z Vector	20µg	P2151
pGEM [®] -4Z Vector	20µg	P2161
pGEM [®] -3Zf(+) Vector	20µg	P2271
pGEM [®] -3Zf(-) Vector	20µg	P2261
pGEM [®] -5Zf(+) Vector	20µg	P2241
pGEM [®] -5Zf(-) Vector	20µg	P2351
pGEM [®] -7Zf(+) Vector	20µg	P2251
pGEM [®] -7Zf(-) Vector	20µg	P2371
pGEM [®] -11Zf(+) Vector	20µg	P2411
pGEM [®] -11Zf(-) Vector	20µg	P2421
pGEM [®] -13Zf(+) Vector	20µg	P2541

All pGEM[®] Vectors are provided with a glycerol stock of bacterial strain JM109. The JM109 cells do not contain the vector and are not competent.

Product	Size	Cat.#
pSP64 Poly(A) Vector	20µg	P1241
pSP72 Vector	20µg	P2191
pSP73 Vector	20µg	P2221

Product	Size	Cat.#
SP6 Promoter Primer	2µg	Q5011
T7 Promoter Primer	2µg	Q5021
pUC/M13 Primer, Reverse (17mer)	2µg	Q5401
pUC/M13 Primer, Forward (17mer)	2µg	Q5391
pUC/M13 Primer, Forward (24mer)	2µg	Q5601
pUC/M13 Primer, Reverse (22mer)	2µg	Q5421

VI. Reference

1. Yanish-Perron, C. *et al.* (1985) Improved M13 phage cloning vectors and host strains: Nucleotide sequences of the M13mp18 and pUC19 vectors. *Gene* **33**, 103-19.

©1988-2007 Promega Corporation. All Rights Reserved.

pGEM and Riboprobe are registered trademarks of Promega Corporation.

DNASTAR is a registered trademark of DNASTAR, Inc. GenBank is a registered trademark of the U.S. Dept. of Health and Human Services.

Products may be covered by pending or issued patents or may have certain limitations. Please visit our Web site for more information.

All prices and specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.