



Promega

Technical Bulletin

pSP-*luc*+NF Fusion Vector

INSTRUCTIONS FOR USE OF PRODUCT E4471.



www.promega.com

pSP-*luc*+NF Fusion Vector

All technical literature is available on the Internet at: www.promega.com/tbs/
Please visit the web site to verify that you are using the most current version of this
Technical Bulletin. Please contact Promega Technical Services if you have questions on use
of this system. E-mail: techserv@promega.com

1. Description.....	1
2. Product Components and Storage Conditions	5
3. Related Products	6
4. References	6

1. Description

The pSP-*luc*+NF Fusion Vector^(a,b,c) is a luciferase cassette vector containing the engineered firefly luciferase gene, *luc*+NF. The *luc*+NF gene is related to the *luc*+ gene found in the pGL3 family of eukaryotic reporter vectors (1) but has been further modified for maximum flexibility in constructing N-terminal fusions (NF) with luciferase. The pSP-*luc*+NF Fusion Vector is not itself intended for the expression of luciferase in eukaryotic cells because it does not contain eukaryotic promoters, enhancers or polyadenylation signals.

The *luc*+NF gene is positioned downstream of an SP6 promoter and ribosome binding site. An opposing T7 promoter is located immediately downstream of *luc*+NF. Thus, the pSP-*luc*+NF Fusion Vector provides a convenient template for the *in vitro* synthesis of both sense and antisense luciferase transcripts for studies involving *in situ* hybridization, RNA processing, RNA transfection or coupled *in vitro* transcription/translation and protein folding. Multiple cloning regions containing recognition sequences for commonly used restriction enzymes are positioned at the 5' and 3' ends of *luc*+NF to provide maximum flexibility in cloning. Subcloning *luc*+NF into expression vectors provides a useful genetic reporter with exceptional sensitivity. Luciferase enzymatic activity can be assayed most efficiently using the Luciferase Assay System.

Luciferase is a 61kDa monomeric protein that does not require post-translational modifications for enzymatic activity. Thus, it can function as a genetic reporter immediately upon translation (2-4). Luciferase synthesized by *in vitro* translation can be labeled with ³⁵S, as the protein contains 4 cysteine and 14 methionine residues. To ensure full enzymatic activity of luciferase, no more than 5 codons can be deleted from either the 5'- or 3'-end of the coding region.

Properties of the *luc*+NF Gene

The numerous modifications that distinguish *luc*+NF from the native luciferase gene are detailed in Table 1. Briefly, the *luc*+NF gene (and the related *luc*+ gene) contains four general classes of modifications: 1) The C-terminal tripeptide has been mutated to eliminate peroxisome targeting of the expressed protein; 2) codon usage was improved for expression in plant and animal cells; 3) two potential sites of N-glycosylation were removed; and 4) several DNA sequence changes were made to disrupt extended palindromes, remove internal restriction sites and eliminate consensus sequences recognized by genetic regulatory binding proteins (5), helping to ensure that the reporter gene itself is unaffected by spurious host transcriptional signals.

In addition, the *luc*+NF gene contains a unique BstEII restriction site inserted immediately downstream of the luciferase ATG translation start codon. The BstEII site allows a cloned insert to be positioned immediately downstream of the *luc*+NF initiation codon. If the cloned insert already contains an initiation codon, the *luc*+NF initiation sequences can be removed using BstEII and another restriction enzyme that cuts in the upstream multiple cloning region. This eliminates the possibility of expressing low levels of functional luciferase through internal initiation in the luciferase portion of the gene fusion.

As a result of the insertion of the BstEII site, the glutamic acid at position 2 of the native protein has been replaced with a valine-threonine at positions 2 and 3 of the *luc*+NF protein. Because of the altered N-terminal amino acid sequence, luciferase encoded by *luc*+NF may exhibit slightly lower enzymatic activity than that encoded by *luc*+.

Note: The numbering scheme for amino acids is relative to the first amino acid, Met #1, of luciferase (*luc*). The numbering of nucleotides is relative to the first base, "A", of the luciferase open reading frame. All numberings for Table 1 refer to the amino acid or nucleic acid positions for the native luciferase gene. Because of the addition of the BstEII site at position 4 of the native nucleic acid sequence to create *luc*+NF, all of the nucleic acid numbering will be shifted by 3 nucleotides and all amino acid numbering will be shifted by 1 residue for the *luc*+NF gene.

Table 1. Summary of Modifications in the *luc*+NF Gene.

Purpose of Modification	Sequence Modification in <i>luc</i>+NF					
Introduce NcoI and BstEII sites for the construction of N-terminal fusions with <i>luc</i> +NF.	<i>luc</i> :	AAA	Met ₁ ATG	Glu ₂ GAA	Asp ₃ GAC	
	<i>luc</i> +NF:	TCC	Met ₁ ATG	Val ₂ GTC	Thr ₃ ACC	Asp ₄ GAC
			NcoI		BstEII	
		#50	XbaI			
Remove internal XbaI site; disrupt extended palindrome.	<i>luc</i> :		CTCTAGAGG			
	<i>luc</i> +NF:		CGCTGGAAG			
Eliminate potential glycosylation and ATF sites.	<i>luc</i> :	Asn ₅₁	Ile ₅₂ ATC	Thr ₅₃ ACG	TACGCGGAA	
			Asp ₅₁	Ile ₅₂ Thr ₅₃		
	<i>luc</i> +NF:	GAC	ATC	ACT	TACGCTGAG	
Eliminate potential glycosylation sites.	<i>luc</i> :	Asn ₁₂₀	Ile ₁₂₁ ATT	Ser ₁₂₂ TCG		
			Gly ₁₂₀	Ile ₁₂₁ Ser ₁₂₂		
	<i>luc</i> +NF:	GGC	ATT	TCG		
		#376				
Remove potential TGT-3 site; improve codon usage.	<i>luc</i> :		GTAGTGTTTGTT			
	<i>luc</i> +NF:		GTTGGTTCGTT			
		#429				
Improve codon usage.	<i>luc</i> :		ATTACCAATAATCCAG			
	<i>luc</i> +NF:		GCTCCAATCATCCAA			
		#549				
Improve codon usage.	<i>luc</i> :		ACCAGAGTCCCTTTGATCGTGACAAA			
	<i>luc</i> +NF:		GCCAGAGTCCCTCGATAGGGACAAG			
		#586	EcoRI			
Remove internal EcoRI site, improve codon usage.	<i>luc</i> :		ATAATGAATTCC			
	<i>luc</i> +NF:		ATCATGAACTCC			
		#611	BstEII			
Remove internal BstEII site and potential AP-2 and LF-A1 sites; improve codon usage.	<i>luc</i> :		GGTTACCTAAGGGTGTGGCCCTTCCG			
	<i>luc</i> +NF:		GTCGTGCCTAAAGGTGTCGCTCTGCT			
		#649				
Remove potential AP-1 site.	<i>luc</i> :		TGCGTCAG			
	<i>luc</i> +NF:		TGCGTGAG			

Table 1. Summary of Modifications in the *luc*+NF Gene (continued).

Purpose of Modification	Sequence Modification in <i>luc</i>+NF
	#823
Improve codon usage.	<i>luc</i> : TTACGATCCCTTCAGGATTACAAA <i>luc</i> +NF: CTGAGGAGCCTTCAGGATTACAAG
	#859
Improve codon usage.	<i>luc</i> : TTGCTAGTACCAACCCTATTTTCA <i>luc</i> +NF: CTGCTGGTGCCAACCCTATTCICC
	#948
Eliminate internal palindrome; improve codon usage.	<i>luc</i> : GGGCGCACCTCTTTTCGAAA <i>luc</i> +NF: TGGCGCTCCCCTCTCTAAG
	#987
Improve codon usage.	<i>luc</i> : AAAACGCTTCCATCTTCCAGGGATACGA <i>luc</i> +NF: CAAGAGGTTCCATCTGCCAGGTATCAGG
	#1161
Eliminate potential AP-1 site; improve codon usage.	<i>luc</i> : GAGAGGCGAATTATGTGTGTCAGAGGA <i>luc</i> +NF: AAGAGGCGAACTGTGTGTGAGAGGT
	#1305
Eliminate palindrome structure; improve codon usage.	<i>luc</i> : AGTTGACCGCTTGAAGTCTTTAATTAATAATAC <i>luc</i> +NF: CGTTGACCGCCTGAAGTCTCTGATTAAGTAC
	#1336 EcoRV
Remove internal EcoRV site; improve codon usage.	<i>luc</i> : AAAGGATATCAGGTGGCC <i>luc</i> +NF: AAAGGCTATCAGGTGGCT
	#1368 ClaI
Remove internal ClaI site; improve codon usage.	<i>luc</i> : ATCGATATTGTTA <i>luc</i> +NF: ATCCATCTTGCTC
	#1403
Remove potential SP-1 site and AP-2 sites.	<i>luc</i> : CGGGCGTGGC <i>luc</i> +NF: CAGGTGTCCG
	Gly ₅₄₇ Lys ₅₄₈ Ser ₅₄₉ Lys ₅₅₀ Leu ₅₅₁ stop
Remove peroxisome targeting sequence.	<i>luc</i> : GGA AAG TCC AAA TTG TAA Gly ₅₄₇ Lys ₅₄₈ Ile ₅₄₉ Ala ₅₅₀ Val ₅₅₁ stop <i>luc</i> +NF: GGA AAG ATC GCC GTG TAA

2. Product Components and Storage Conditions

Product	Size	Cat.#
pSP- <i>luc</i> +NF Fusion Vector	20µg	E4471

Storage Conditions: Store the pSP-*luc*+NF Fusion Vector at -20°C.




Figure 1. pSP-*luc*+NF Fusion Vector circle map and sequence reference points.
(rbs = ribosome binding site)

Sequence reference points:

SP6 RNA polymerase initiation site	1
upstream multiple cloning region	11-53
luciferase coding region	45-1700
downstream multiple cloning region	1702-1731
T7 RNA polymerase promoter (-17 to +3)	1740-1759
T7 RNA polymerase initiation site	1742
ColE1-derived replication origin	2014
β-lactamase (Amp ^r) coding region	3636-2776
SP6 RNA polymerase promoter (-17 to +3)	4087-3




Figure 2. pSP-*luc*+NF Fusion Vector promoter and multiple cloning region sequence. The sequence shown corresponds to RNA synthesized by SP6 RNA polymerase and corresponds to the sequence of *luc*+NF mRNA. Note that cleavage at the ClaI site is blocked when the pSP-*luc*+NF Fusion Vector is propagated in host strains containing *dam* methylase activity.

3. Related Products

Product	Size	Cat.#
Luciferase Assay System (with Cell Culture Lysis Reagent)	100 assays	E1500
Luciferase Assay System with Reporter Lysis Buffer	100 assays	E4030
Luciferase Assay Reagent 10-Pack	1,000 assays	E1501
Luciferase Assay Reagent (bulk liquid reagent)	100ml	E1483
pGL3-Control Vector	20µg	E1741
pGL3-Enhancer Vector	20µg	E1771
pGL3-Promoter Vector	20µg	E1761
pGL3-Basic Vector	20µg	E1751

4. References

1. Sherf, B. and Wood, K.V. (1994) Increased expression and convenience with the new pGL3 Luciferase Reporter Vectors. *Promega Notes* **49**, 14-21.
2. Ow, D. *et al.* (1986) Transient and stable expression of the firefly luciferase gene in plant cells and transgenic plants. *Science* **234**, 856-9.
3. deWet, J.R. *et al.* (1987) Firefly luciferase gene: Structure and expression in mammalian cells. *Mol. Cell Biol.* **7**, 725-37.
4. Wood, K.V. (1990) Firefly luciferase: A new tool for molecular biologists. *Promega Notes* **28**, 1-3.
5. Faisst, S. and Meyer, S. (1992) Compilation of vertebrate-encoded transcription factors. *Nucleic Acids Res.* **20**, 3-26.

ⓂREAD THIS FIRST BEFORE OPENING PRODUCT

For research use only. The terms of the limited license conveyed with the purchase of this product are as follows: Researchers may use this product in their own research and they may transfer derivatives to others for such research use provided that at the time of transfer a copy of this label license is given to the recipients and the recipients agree to be bound by the conditions of this label license. Researchers shall have no right to modify or otherwise create variations of the nucleotide sequence of the luciferase gene except that Researchers may: (1) clone heterologous DNA sequences at either or both ends of said luciferase gene so as to create fused gene sequences provided that the coding sequence of the resulting luciferase gene has no more than four deoxynucleotides missing at the affected terminus when compared to the intact luciferase gene sequence, and (2) insert and remove nucleic acid sequences in furtherance of splicing research predicated on the inactivation or reconstitution of the luminescent activity of the encoded luciferase. In addition, Researchers must do one of the following: (1) use luminescent assay reagents purchased from Promega Corporation for all determinations of luminescence activity resulting from the research use of this product and its derivatives; or (2) contact Promega to obtain a license for the use of the product and its derivatives. No other use or transfer of this product or its derivatives is authorized without the express written consent of Promega including, without limitation, Commercial Use. Commercial Use means any and all uses of this product and derivatives by a party for monetary or other consideration and may include but is not limited to use in: (1) product manufacture; and (2) to provide a service, information or data; and/or resale of the product or its derivatives, whether or not such product or derivatives are resold for use in research. With respect to such Commercial Use, or any diagnostic, therapeutic or prophylactic uses, please contact Promega for supply and licensing information. If the purchaser is not willing to accept the conditions of this limited use statement, Promega is willing to accept the return of the unopened product and provide the purchaser with a full refund. However, in the event the product is opened, then the purchaser agrees to be bound by the conditions of this limited use statement. The above license relates to Promega patents and/or patent applications on improvements to the luciferase gene.

ⓂU.S. Pat. No. 5,670,356.

ⓂThe method of recombinant expression of *Coleoptera* luciferase is covered by U.S. Pat. Nos. 5,583,024, 5,674,713 and 5,700,673. A license (from Promega for research reagent products and from The Regents of the University of California for all other fields) is needed for any commercial sale of nucleic acid contained within or derived from this product.

© 1994-2009 Promega Corporation. All Rights Reserved.

Products may be covered by pending or issued patents or may have certain limitations. Please visit our Web site for more information.

All prices and specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.