

FreeStyle™ 293 Expression Medium

Protocol Outline

- A. Thaw cells.
- B. Passage cells every 2–3 days.

FreeStyle™ 293-F Cell Culturing Protocol

-  See page 3 to view a typical procedure for subculturing.

Scaling Up FreeStyle™ 293-F Cell Culture

You can scale up FreeStyle™ 293-F cultures in spinner flasks or bioreactors. Determine the optimal spinner or impeller speed and seeding density for your culture system.









If the split ratio of cells to fresh media is less than 1:2, you may need to spin down the cell suspension and resuspend in fresh, pre-warmed FreeStyle™ 293 Expression Medium prior to inoculating the spinner or bioreactor culture.

At high stirring speeds (i.e. >130 rpm) and/or depending on the impeller design, you may need to supplement the FreeStyle™ 293 Expression Medium with additional Pluronic® F-68 (2.5–5 mL/L of 10% Pluronic® F-68) to avoid shear stress in the culture.

Adapting Other 293 Cells to FreeStyle™ 293 Expression Medium

Cryopreserving FreeStyle™ 293-F Cells






Limited Product Warranty and Disclaimer Details

	Catalog Number	Size
 Package Contents	12338-018	1000 mL
	12338-026	6 × 1000 mL
	12338-001	10 L
	12338-002	20 L
 Storage Conditions	<ul style="list-style-type: none"> ▪ Store at 4°C for a 12-month shelf life. ▪ Protect from light. 	
	<ul style="list-style-type: none"> ▪ FreeStyle™ 293-F Cells ▪ 125-mL polycarbonate, disposable, sterile, vent-cap Erlenmeyer shaker flask or other appropriate vessel for culturing suspension cells ▪ Orbital shaker in temperature and CO₂ controlled incubator ▪ Reagents and equipment to determine cell viability (e.g., hemocytometer with trypan blue or cell counter) 	
 Required Materials	<ul style="list-style-type: none"> ▪ FreeStyle™ 293-F Cells ▪ 125-mL polycarbonate, disposable, sterile, vent-cap Erlenmeyer shaker flask or other appropriate vessel for culturing suspension cells ▪ Orbital shaker in temperature and CO₂ controlled incubator ▪ Reagents and equipment to determine cell viability (e.g., hemocytometer with trypan blue or cell counter) 	
	<ul style="list-style-type: none"> ▪ FreeStyle™ 293-F Cells ▪ 125-mL polycarbonate, disposable, sterile, vent-cap Erlenmeyer shaker flask or other appropriate vessel for culturing suspension cells ▪ Orbital shaker in temperature and CO₂ controlled incubator ▪ Reagents and equipment to determine cell viability (e.g., hemocytometer with trypan blue or cell counter) 	
 Timing	Thawing and Recovery: 2–3 days	
	Subculturing: Every 2–3 days	
 Selection Guide	Protein Expression Systems	
	Go online to view related products.	
 Product Description	<ul style="list-style-type: none"> ▪ FreeStyle™ 293 Expression Medium is a chemically defined and serum-free medium, specifically developed to support the growth and transfection of FreeStyle™ 293-F Cells under suspension culture conditions. ▪ This medium does not contain any proteins, hydrolysates, or components of animal origin. 	
	<ul style="list-style-type: none"> ▪ FreeStyle™ 293 Expression Medium contains GlutaMAX™-I supplement and does not require supplementation with L-glutamine or GlutaMAX™-I supplement. ▪ Subculture FreeStyle™ 293-F Cells when they reach a density of approximately 1–3 × 10⁶ viable cells/mL, typically every 2–3 days. Split the FreeStyle™ 293-F culture to 0.2–0.5 × 10⁶ cells/mL. ▪ Keep cell densities between 1–3 × 10⁶ cells/mL of culture for best performance. 	
 Important Guidelines	<ul style="list-style-type: none"> ▪ FreeStyle™ 293 Expression Medium contains GlutaMAX™-I supplement and does not require supplementation with L-glutamine or GlutaMAX™-I supplement. ▪ Subculture FreeStyle™ 293-F Cells when they reach a density of approximately 1–3 × 10⁶ viable cells/mL, typically every 2–3 days. Split the FreeStyle™ 293-F culture to 0.2–0.5 × 10⁶ cells/mL. ▪ Keep cell densities between 1–3 × 10⁶ cells/mL of culture for best performance. 	
	<ul style="list-style-type: none"> ▪ FreeStyle™ 293 Expression Medium contains GlutaMAX™-I supplement and does not require supplementation with L-glutamine or GlutaMAX™-I supplement. ▪ Subculture FreeStyle™ 293-F Cells when they reach a density of approximately 1–3 × 10⁶ viable cells/mL, typically every 2–3 days. Split the FreeStyle™ 293-F culture to 0.2–0.5 × 10⁶ cells/mL. ▪ Keep cell densities between 1–3 × 10⁶ cells/mL of culture for best performance. 	
 Online Resources	Visit our product page for additional information and protocols. For support, visit www.lifetechnologies.com/support .	
	Visit our product page for additional information and protocols. For support, visit www.lifetechnologies.com/support .	

For Research Use Only. Not for use in diagnostic procedures.

Thawing and Passaging FreeStyle™ 293-F Cells in FreeStyle™ 293 Medium

Follow the procedure below to recover and subculture FreeStyle™ 293-F Cells in FreeStyle™ 293 Expression Medium.

	Timeline	Steps	Procedure Details		
Day 1	1 	Thaw cells	Rapidly thaw the cells in a water bath, decontaminate the vial using 70% ethanol, and open the cryovial in a class II biological cabinet.		
	2 	Add cells to medium	Add cells to 29 mL of pre-warmed medium in 125-mL shake flask.		
	3 	Count cells and determine viability	Within 1–2 hours post-thaw, count cells and determine viability. Use hemocytometer and trypan blue exclusion method or automated cell counter. Cell density should be approximately 0.3×10^6 cells/mL and cell viability >90%.		
	4 	Incubate	Temperature 37°C	Humidified Atmosphere 8% CO ₂ in air	Orbital Shaker Platform 125 rpm
Days 3–4	5 	Subculture cells	<p>First passage: When cell density reaches $>1 \times 10^6$ cells/mL at $\geq 90\%$ viability (typically 2–3 days post-thaw), split cells to 0.3×10^6 cells/mL in FreeStyle™ 293 Expression Medium.</p> <p>Subsequent passages: Every 2–3 days, cells should reach $1\text{--}3 \times 10^6$. Split to $0.2\text{--}0.5 \times 10^6$ cells/mL. Do not grow above 3×10^6 cells/mL.</p> <p>We recommend using a 125- or 250-mL flask containing 30 or 60 mL of medium, respectively.</p>		