
# Technical Data Sheet **Purified Mouse Anti-Human IFN-y**

Product Information	
Material Number:	554549
Size:	0.1 mg
Concentration:	0.5 mg/ml
Clone:	4S.B3
Immunogen:	Partially purified human IFN- $\gamma$ from supernatants of human PBMC stimulated with Staphylococcus aureus
Isotype:	Mouse IgG1, κ
Reactivity:	QC Testing: Human
Storage Buffer:	Aqueous buffered solution containing ≤0.09% sodium azide.

## Description

The 4S.B3 antibody reacts with human interferon- $\gamma$  (IFN- $\gamma$ ). The immunogen used to generate this hybridoma was partially purified human IFN- $\gamma$  obtained from supernatants of human PBMC stimulated with *Staphylococcus aureus*. This is a neutralizing antibody.



Expression of IFN-y by stimulated human peripheral blood mononuclear cells (PBMC). Human PBMC were stimulated for 6 hr with PMA (50 ng/ml; Sigma) and calcium ionophore A23187 (500 ng/ml; Sigma) in the presence of 2 µM BD GolgiStop<sup>™</sup> (Cat. No. 554724). The PBMC were stained with PE-Cy5-anti-CD3 (PE-Cy5 UCH71, Cat. No. 555334), fixed, permeabilized, and subsequently stained with 0.125 µg of PE-mouse anti-human IFN-γ antibody (PE-4S.B3, Cat. No. 554552) by using Pharmingen's staining protocol (left panel). The binding of PE-4S.B3 was blocked by preincubation of cells with unlabeled 4S.B3 antibody (5 µg, right panel). The quadrant markers for the bivariate dot plot were set based on the autofluorescence controls and verified using the unlabeled antibody and ligand blocking controls.

# **Preparation and Storage**

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography. Store undiluted at 4° C.

## **Application Notes**

Intracellular block/flow cytometry	Routinely Tested
Western blot	Reported
ELISA Detection	Reported

# **Recommended Assay Procedure:**

**Blocking Control for Intracellular Staining:** The purified 4S.B3 antibody (Cat. No. 554549) can be used as a blocking control to demonstrate specificity of IFN-γ staining by directly conjugated clone 4S.B3. To perform this control, the fixed/permeabilized cells (~1 million) can be

## **BD Biosciences**

Jnited States	Canada	Europe	Japan	Asia Pacific	Latin America/Caribbean
377.232.8995	888.259.0187	32.53.720.550	0120.8555.90	65.6861.0633	55.11.5185.9995
or country-spe	ecific contact inf	formation, visit <b>v</b>	www.bdbioscien	ces.com/how_to	_order/
ise of our produc product or as a co vritten authoriza for Research Use	ts. Purchase does n mponent of anoth tion of Becton Dick Only. Not for use in	ot include or carry er product. Any us kinson and Compan diagnostic or there	any right to resell or	r transfer this produce than the permitte d. Not for resale.	lations that may occur with the ct either as a stand-alone ed use without the express 07 BD

incubated with 1-10  $\mu$ g of purified 4S.B3 antibody (Cat. No. 554549) for 20 minutes at 4°C, prior to staining with the directly conjugated antibody (e.g., 0.1 -0.5  $\mu$ g mAb/1 million cells) (right panel). The intracellular cytokine staining technique and use of blocking controls are described in detail by C. Prussin and D. Metcalfe. For specific methodology, please visit our web site, www.bdbiosciences.com, and go to the protocols section or the chapter on intracellular staining in the Immune Function Handbook.

**ELISA Detection:** The biotinylated 4S.B3 antibody (Cat. No. 554550) is useful as a detection antibody in a sandwich ELISA for measuring human IFN- $\gamma$  protein levels. Biotinylated 4S.B3 antibody can be paired with the purified NIB42 antibody (Cat. No. 551221) as the capture antibody, with recombinant human IFN- $\gamma$  (Cat. No. 554616) as the standard. For testing human IFN- $\gamma$  complex in biological fluids like serum or plasma, our human IFN- $\gamma$  specific OptEIA<sup>TM</sup> sandwich ELISA set (Cat.No. 555142I) and OptEIA<sup>TM</sup> sandwich ELISA kit (Cat. No. 550612) are recommended.

Western Blot: The 4S.B3 antibody has been found useful for Western blotting.

# **Suggested Companion Products**

Size	Clone (none)	
0.7 ml		
100 tests	UCHT1	
100 tests	UCHT1	
0.5 mg	4S.B3	
1.0 mg	NIB42	
25 µg	(none)	
2 plates	(none)	
20 tests	(none)	
	25 μg 2 plates	

## **Product Notices**

1. Since applications vary, each investigator should titrate the reagent to obtain optimal results.

- 2. Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols.
- 3. Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.

#### References

Meager A. Characterization of interferons and immunoassays. In: Clemens MJ, Morris AG, Gearing AJH, ed. Lymphockies and Interferons. A Practical Approach. Oxford: IRL Press Ltd; 1987:105-127.(Biology)

Meager A, Parti S, Barwick S, Spragg J, O'Hagan K. Detection of hybridomas secreting monoclonal antibodies to human gamma interferon using a rapid screening technique and specificity of certain monoclonal antibodies to gamma interferon. J Interferon Res. 1984; 4(4):619-625 (Biology)

Prussin C, Metcalfe DD. Detection of intracytoplasmic cytokine using flow cytometry and directly conjugated anti-cytokine antibodies. J Immunol Methods. 1995; 188(1):117-128.(Methodology)