

Technical Data Sheet

Anti-Mouse IgG1 Magnetic Particles - DM

Product Information

Material Number:	557983
Size:	5.0 ml
Clone:	A85-1
Reactivity:	QC Testing: Cells (such as rat splenocytes) labelled with mouse IgG1 antibodies
Storage Buffer:	Aqueous buffered solution containing BSA and ≤0.09% sodium azide.

Description

BD IMag™ Anti-Mouse IgG1 Magnetic Particles -DM are magnetic nanoparticles that have monoclonal antibody conjugated to their surfaces. The A85-1 clone has been reported to react specifically with mouse IgG1 of Igh-Ca and Igh-Cb haplotypes, but reportedly does not react with other Ig isotypes. These particles are optimized for the positive selection or depletion of leukocyte subpopulations labelled with mouse IgG1 antibodies using the BD IMagnet™.

Leukocytes are labelled with BD IMag™ Anti-Mouse IgG1 Magnetic Particles - DM according to the Magnetic Labeling and Separation Protocol. In brief, cells are labeled with a mouse IgG1 antibody that recognizes the subpopulation of interest. After washing away excess antibody, BD IMag™ Anti-Mouse IgG1 Magnetic Particles -DM are added to the cell suspension and bind the mouse IgG1 antibody on the cells. This labeled cell suspension is then placed within the magnetic field of the BD IMagnet™ (Cat.No. 552311). Positive selection or depletion is then performed. Labeled cells migrate toward the magnet (positive fraction), leaving the unlabeled cells in suspension so they can be drawn off (depleted or negative fraction). The tube is then removed from the magnetic field for resuspension of the positive fraction. The magnetic separation steps are diagrammed in the Depletion and Positive Selection Flow Charts. After the positive fraction has been washed, the small size of the magnetic particles allows the positive fraction to be further evaluated in downstream applications such as flow cytometry.




Figure 1. Depletion of human T and B lymphocytes and NK cells. Human PBMC were stained with APC mouse anti-human CD4 mAb RPA-T4 (Cat. No. 555349), CD8 mAb RPA-T8 (Cat. No. 555369), CD19 mAb H1B19 (Cat. No. 555415), and CD56 mAb B159 (Cat. No. 555518), then labeled with BD IMag™ Anti-Mouse IgG1 Magnetic Particles - DM. After labeling, the cells were separated using the BD IMagnet™, and the depleted and positive fractions were collected as described in the Magnetic Labeling and Separation Protocol. Please refer to the Depletion Flow Chart to identify the separated cell populations represented in this figure. Unseparated PBMC (left panel), the Final Depleted Fraction (middle panel), and the Positive Fraction (right panel) were analyzed by flow cytometry. Nonviable cells were eliminated from analysis by staining with propidium iodide, and mononuclear cells were identified by scatter profile. The percentage of non-T-B-NK leukocytes in each sample is given. Flow cytometry was performed on a BD FACSCalibur™ flow cytometry system.

Preparation and Storage

Store undiluted at 4° C.

BD Biosciences

www.bdbiosciences.com

United States Canada Europe Japan Asia Pacific Latin America/Caribbean
877.232.8995 888.259.0187 32.53.720.550 0120.8555.90 65.6861.0633 55.11.5185.9995

For country-specific contact information, visit www.bdbiosciences.com/how_to_order/

Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton Dickinson and Company is strictly prohibited.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.

BD, BD Logo and all other trademarks are the property of Becton, Dickinson and Company. ©2007 BD






Figure 2. Positive selection of rat $\alpha\beta$ TCR-expressing T lymphocytes. Lewis rat splenocytes were stained with FITC mouse anti-rat $\alpha\beta$ TCR mAb R73 (Cat. No. 554913), and then labeled with BD IMag™ Anti-Mouse IgG1 Magnetic Particles - DM. After labeling, the cells were separated using the BD IMagnet™, and the negative and positive fractions were collected as described in the Magnetic Labeling and Separation Protocol. Please refer to the Positive Selection Flow Chart to identify the separated cell populations represented in this figure. Unseparated Spleen (left panel), the Negative Fraction (middle panel), and the Positive Fraction (right panel) were analyzed by flow cytometry. Nonviable cells were eliminated from analysis by staining with propidium iodide, and all viable leukocytes are displayed. The percentage of $\alpha\beta$ TCR+ T lymphocytes in each sample is given. The Negative Fraction includes B lymphocytes and T lymphocytes expressing the $\gamma\delta$ TCR. Flow cytometry was performed on a BD FACSCalibur™ flow cytometry system.

Application Notes

Application

Cell separation	Routinely Tested
-----------------	------------------

Recommended Assay Procedure:

Magnetic Labeling and Separation Protocol

- Prepare the following buffers and place on ice.
 - Cell-staining buffer: Phosphate Buffered Saline, 3% heat inactivated fetal calf serum, 0.1% sodium azide.
 - 1X BD IMag™ buffer: Dilute BD IMag™ Buffer (10X) (Cat. no. 552362) 1:10 with sterile distilled water or alternatively, prepare Phosphate Buffered Saline, supplemented with 0.5% BSA, 2 mM EDTA, and 0.09% sodium azide.*
- Prepare a single-cell suspension from the lymphoid tissue of interest or prepare peripheral blood mononuclear cells (PBMC) from anti-coagulated blood, preferably by density gradient centrifugation using the appropriate density Ficoll-Paque™ solution. Remove clumps of cells and/or debris by passing the suspended cells through a 70- μ m nylon cell strainer.
- Count the cells, and resuspend them in cell-staining buffer at a concentration of 2×10^7 cells/ml.
- Add the mouse IgG1 antibody (or cocktail of mouse IgG1 antibodies) at the appropriate concentration, and incubate on ice for 15 minutes.†
- Wash the labeled cells with an excess volume of 1X BD IMag™ buffer, and carefully aspirate ALL the supernatant. For depletions, proceed with Step 6. For positive selections, proceed with Step 17.

Depletions:

- Vortex the BD IMag™ Anti-Mouse IgG1 Magnetic Particles - DM thoroughly, and add 50 μ l of particles for every 1×10^7 total cells.
- MIX THOROUGHLY. Refrigerate rat or mouse leukocytes for 30 minutes at 6°C -12°C. Incubate human PBMC at room temperature for 30 minutes.†
- Bring the labeling volume up to $2-8 \times 10^7$ cells/ml with 1X BD IMag™ buffer or culture medium.*
- Transfer the labeled cells to a 12 x 75 mm round-bottom test tube (eg, BD Falcon™, Cat. No. 352058), maximum volume added not to exceed 1.0 ml. Place this positive-fraction tube on the BD IMagnet™ (horizontal position) for 6-8 minutes.
 - For greater volume, transfer the cells to a 17 x 100 mm round-bottom test tube (eg, BD Falcon™, Cat. no. 352057), maximum volume added not to exceed 3.0 ml. Place this positive-fraction tube on the BD IMagnet™ (vertical position) for 8 minutes.
- With the tube on the BD IMagnet™ and using a glass Pasteur pipette, carefully aspirate the supernatant (depleted fraction) and place in a new tube.
- Remove the positive-fraction tube from the BD IMagnet™, and add 1X BD IMag™ buffer (or medium) to the same volume as in Step 9. Resuspend the positive fraction well by pipetting up and down 10-15 times and place back on the BD IMagnet™ for 6-8 minutes.
 - 17 x 100 mm tube: Place on the BD IMagnet™ for 8 minutes.
- Using a new Pasteur pipette, carefully aspirate the supernatant and combine with the depleted fraction from Step 11 above.
- Repeat Steps 11 and 12. The Combined Depleted Fraction contains cells with no bound antibodies or magnetic particles. These cells are ready for downstream applications, or they can be further enriched by proceeding to Step 15.
- The positive-fraction cells remaining in the original tube can be resuspended in an appropriate buffer or culture medium for downstream applications, including flow cytometry.
- To increase the purity of the Combined Depleted Fraction, place the tube on the BD IMagnet™ for another 6-8 minutes.
 - 17 x 100 mm tube: Place on the BD IMagnet™ for 8 minutes.
- Carefully aspirate the supernatant and place in a new tube. This is the Final Depleted Fraction. The cells are ready to be processed for downstream applications.


Positive Selections:

17. Vortex the BD IMag™ Anti-Mouse IgG1 Magnetic Particles - DM thoroughly, and add 10-50 µl of particles for every 1 x 10⁷ total cells.
18. MIX THOROUGHLY. Refrigerate rat or mouse leukocytes for 30 minutes at 6°C -12°C. Incubate human PBMC at room temperature for 30 minutes.†
19. Bring the labeling volume up to 2-8 x 10⁷ cells/ml with 1X BD IMag™ buffer.
20. Immediately place the tube onto the BD IMagnet™ and incubate for 6-8 minutes.
21. With the tube on the BD IMagnet™, carefully aspirate the supernatant. This supernatant is considered the Negative Fraction.
22. Remove the tube from the BD IMagnet™, and add 1X BD IMag™ buffer to the same volume as in Step 19. Gently resuspend the cells by pipetting up and down, and return the tube to the BD IMagnet™ for another 2-4 minutes.
23. With the tube on the BD IMagnet™, carefully remove the supernatant.
24. Repeat Steps 22 and 23.
25. After the final wash step, remove the tube from the BD IMagnet™. Resuspend the Positive Fraction in an appropriate buffer or culture medium, and proceed with desired downstream application(s), including flow cytometry.

NOTES:

* For depletion of rat leukocytes, tissue culture medium usually results in a slight increase in viability and recovery, when compared to IMag buffer, without reducing cell purity. Because applications can vary, researchers are encouraged to run a trial comparison of culture media and IMag buffer to demonstrate that there are no adverse effects.

† Avoid non-specific labeling by working quickly and adhering to recommended incubation times.




Suggested Companion Products

Catalog Number	Name	Size	Clone
552362	Buffer (10X)	100 ml	(none)
552311	Cell Separation Magnet	each	(none)

POSITIVE SELECTION FLOW CHART

(The circled numbers correspond to the steps of the Protocol on the following page.)



Product Notices

1. Since applications vary, each investigator should titrate the reagent to obtain optimal results.
2. BD IMag™ particles are prepared from carboxy-functionalized magnetic particles which are manufactured by Skold Technology.
3. Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.
4. Source of all serum proteins is from USDA inspected abattoirs located in the United States.
5. Ficoll-Paque is a trademark of Amersham Biosciences Limited.
6. Please refer to www.bdbiosciences.com/pharming/en/protocols for technical protocols.