

Technical Data Sheet

Mouse NK Cell Separation Set - DM




Product Information

Material Number:	558004
Component:	51-9002979
Description:	PE Rat Anti-Mouse CD49b
Size:	10 ml (1 ea)
Clone Name:	DX5
Isotype:	Rat (LEW) IgM, κ
Storage Buffer:	Aqueous buffered solution containing BSA and $\leq 0.09\%$ sodium azide.
Component:	51-9002978
Description:	Anti-PE Particles 2 - DM
Size:	10 ml (1 ea)
Clone Name:	E31-1459
Isotype:	Mouse IgG1, κ
Storage Buffer:	Aqueous buffered solution containing BSA and $\leq 0.09\%$ sodium azide.

Description

The BD IMag™ Mouse NK Cell Separation Set - DM is optimized for the positive selection or depletion of NK and NK-T cells using the BD IMagnet™ (Cat. No. 552311). The PE rat anti-mouse CD49b antibody (clone DX5) recognizes CD49b, which has been reported to be expressed at high density on most NK cells and at low density on subsets of NK-T cells in many mouse strains (e.g., A/J, AKR, BALB/c, C3H/HeJ, C57BL/6, C57BL/10, C57BR, C57L, C58, CBA/Ca, CBA/J, DBA/1, DBA/2, NOD, SJL, SWR, 129/J). The BD IMag™ Anti-PE Particles 2 - DM are magnetic nanoparticles that have anti-PE monoclonal antibody covalently conjugated to their surfaces. This Separation Set has been optimized for use with the BD IMagnet™ and it contains sufficient reagents to label 2×10^9 leukocytes.

In summary, the PE rat anti-mouse CD49b antibody stains NK and NK-T cells. After washing away excess antibody, BD IMag™ Anti-PE Particles 2 - DM are added to the cell suspension and bind the cells bearing the PE-conjugated antibody. The tube containing this labeled cell suspension is then placed within the magnetic field of the BD IMagnet™. Positive selection or depletion is then performed. Labeled cells migrate toward the magnet (positive fraction), leaving the unlabeled cells in suspension so they can be drawn off (depleted or negative fraction). The tube is then removed from the magnetic field for resuspension of the positive fraction. The selections are repeated twice to increase the purity of the positive fraction. For clarification of the procedure, the magnetic separation steps are diagrammed in the Positive Selection Flow Chart. The small size of the BD IMag™ particles allows the positive fraction to be further evaluated in downstream applications such as flow cytometry.



Positive selection of Mouse NK cells. BALB/c splenocytes were stained with PE rat anti-mouse CD49b (clone DX5) and FITC hamster anti-mouse CD3e (clone 145-2C11) (Cat. No. 553061) monoclonal antibodies, and then labeled with BD IMag™ Anti-PE Particles 2 - DM. After labeling, the cells were separated using the BD IMagnet™, and the negative (CD49b-) and positive (CD49b+) fractions were collected as described in the Magnetic Labeling and Separation Protocol. Please refer to the Positive Selection Flow Chart to identify the separated cell populations represented in this figure. Unseparated splenocytes (left panel), the negative fraction (middle panel), and the positive fraction (right panel) were analyzed by flow cytometry. Nonviable cells were eliminated from analysis by staining with propidium iodide, and all viable leukocytes are displayed. The combined percentage of NK and NK-T cells (CD3e-CD49b[bright] and CD3e+ CD49b[dim], respectively) in each sample is given. Flow cytometry was performed on a BD FACSCalibur™ instrument.

BD Biosciences

bdbiosciences.com

United States	Canada	Europe	Japan	Asia Pacific	Latin America/Caribbean
877.232.8995	888.259.0187	32.53.720.550	0120.8555.90	65.6861.0633	55.11.5185.9995

For country-specific contact information, visit bdbiosciences.com/how_to_order/

Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton Dickinson and Company is strictly prohibited.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.

BD, BD Logo and all other trademarks are the property of Becton, Dickinson and Company. ©2006 BD



BD


BD Biosciences

Preparation and Storage

Store undiluted at 4° C.

POSITIVE SELECTION FLOW CHART

(The circled numbers correspond to the steps of the following Protocol.)



Application Notes

Application

Cell separation	Routinely Tested
-----------------	------------------

Recommended Assay Procedure:

MAGNETIC LABELING AND SEPARATION PROTOCOL

1. Dilute BD IMag™ Buffer (10X) (Cat. No. 552362) 1:10 with sterile distilled water or prepare 1X BD IMag™ buffer by supplementing Phosphate Buffered Saline with 0.5% BSA, 2 mM EDTA, and 0.09% sodium azide. Place on ice.
2. Aseptically prepare a single-cell suspension from the peripheral lymphoid tissue of interest. Remove clumps of cells and/or debris by passing the suspended cells through a 70-µm nylon cell strainer. NK cells are fragile and therefore should be prepared in tissue culture medium [e.g., Dulbecco's Minimum Essential Medium (DMEM) supplemented with 10% fetal bovine serum and L-glutamine].
3. Count the cells. If the concentration is between 1 x 10⁷ and 2 x 10⁷ cells/ml, then proceed to Step 4. If cells are more dilute than 1 x 10⁷ cells/ml, then spin down the cells and resuspend them in tissue culture medium at a concentration of 2 x 10⁷ cells/ml.
4. **Optional:** Add BD Mouse Fc Block™ purified anti-mouse CD16/CD32 mAb 2.4G2 (Cat. No. 553141) at 0.25 µg per 1 x 10⁶ cells, and incubate on ice for 15 minutes.*
5. Add the BD IMag™ PE Rat Anti-Mouse CD49b mAb (clone DX5) at 5 µl per 1 x 10⁶ cells, and refrigerate for 15 minutes at 6-12°C.†
6. Wash the labeled cells with an excess volume of 1X BD IMag™ buffer, and carefully aspirate ALL the supernatant.
7. Vortex the BD IMag™ Anti-PE Particles 2 - DM thoroughly, and add 5 µl of particles for every 1 x 10⁶ total cells.
8. **Mix thoroughly.** Refrigerate for 30 minutes at 6-12°C.‡
9. Bring the labeling volume up to 2 to 8 x 10⁷ cells/ml with 1X BD IMag™ buffer.
10. Immediately place the tube onto the BD IMagnet™ and incubate for 6 to 8 minutes.
11. With the tube on the BD IMagnet™, carefully aspirate the supernatant. This supernatant is considered the Negative (or NK cell-depleted) Fraction.
12. Remove the tube from the BD IMagnet™, and add 1X BD IMag™ buffer to the same volume as in Step 9. Gently resuspend the cells by pipetting up and down, and return the tube to the BD IMagnet™ for another 2 to 4 minutes. Longer incubation time will increase the percentage of NK-T cells in the positive fraction.
13. With the tube on the BD IMagnet™, carefully remove the supernatant (wash fraction) and discard.

14. Repeat Steps 12 and 13.

15. After the final wash step, remove the tube from the BD IMagnet™. Resuspend the Positive Fraction in an appropriate buffer or culture medium, and proceed with desired downstream application(s), including flow cytometry.

Notes:

* The use of BD Mouse Fc Block™ purified anti-mouse CD16/CD32 mAb 2.4G2 in step 4 can increase the purity and recovery of the NK cells by up to 5%. Please note that this results in enriched NK cells that may have purified anti-mouse CD16/CD32 mAb bound to their surface, which may affect the function of those NK cells.

† Avoid nonspecific labeling by working quickly and adhering to recommended incubation times.

Suggested Companion Products

<u>Catalog Number</u>	<u>Name</u>	<u>Size</u>	<u>Clone</u>
552362	Buffer (10X)	100 ml	(none)
552311	Cell Separation Magnet	each	(none)
553141	Purified Rat Anti-Mouse CD16/CD32 (Mouse BD Fc Block™)	0.1 mg	2.4G2
553061	FITC Hamster Anti-Mouse CD3e	0.1 mg	145-2C11

Product Notices

1. Since applications vary, each investigator should titrate the reagent to obtain optimal results.
2. BD IMag™ particles are prepared from carboxy-functionalized magnetic particles which are manufactured by Skold Technology.
3. Please refer to www.bdbiosciences.com/pharming/en/protocols for technical protocols.
4. Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.
5. Source of all serum proteins is from USDA inspected abattoirs located in the United States.

References

Arase H, Saito T, Phillips JH, Lanier LL. Cutting edge: the mouse NK cell-associated antigen recognized by DX5 monoclonal antibody is CD49b (alpha 2 integrin, very late antigen-2). *J Immunol.* 2001; 167(3):1141-1144.(Biology)

Gadue P, Stein PL. NK T cell precursors exhibit differential cytokine regulation and require Itk for efficient maturation. *J Immunol.* 2002; 169(5):2397-2406. (Biology)

Karnbach C, Daws MR, Niemi EC, Nakamura MC. Immune rejection of a large sarcoma following cyclophosphamide and IL-12 treatment requires both NK and NK T cells and is associated with the induction of a novel NK T cell population. *J Immunol.* 2001; 167(5):2569-2576.(Biology)

Lee IF, Qin H, Trudeau J, Dutz J, Tan R. Regulation of autoimmune diabetes by complete Freund's adjuvant is mediated by NK cells. *J Immunol.* 2004; 172(2):937-942.(Biology)

Ortaldo JR, Winkler-Pickett R, Mason AT, Mason LH. The Ly-49 family: regulation of cytotoxicity and cytokine production in murine CD3+ cells. *J Immunol.* 1998; 160(1):1158-1165.(Biology)