

pHrodo™ Red and Green BioParticles® Conjugates for Phagocytosis

Table 1 Contents and storage

Material	Cat. no.	Amount	Ex/Em Maxima	Storage*
pHrodo™ Red <i>E. coli</i> BioParticles® Conjugate	P35361	5 vials each containing 2 mg lyophilized product	560/585 nm	<ul style="list-style-type: none"> • ≤ -20°C • Desiccate • Protect from light
pHrodo™ Green <i>E. coli</i> BioParticles® Conjugate	P35366		509/533 nm	
pHrodo™ Red Zymosan A BioParticles® Conjugate	P35364	5 vials each containing 1 mg lyophilized product	560/585 nm	
pHrodo™ Green Zymosan A BioParticles® Conjugate	P35365		509/533 nm	
pHrodo™ Red <i>S. aureus</i> BioParticles® Conjugate	A10010	5 vials each containing 2 mg lyophilized product	560/585 nm	
pHrodo™ Green <i>S. aureus</i> BioParticles® Conjugate	P35367		509/533 nm	
* When stored as directed, the product is stable for at least 6 months.				
Number of assays: Sufficient for 100 assays when using described protocol.				

Introduction

pHrodo™ BioParticles® Conjugates are novel, no-wash fluorogenic reagents developed for quantitative measurements of phagocytosis and its regulation by drugs, genetic, or environmental factors. With an optional no-cell background subtraction method, a large and specific signal is obtained from cells that ingest the particles, providing a specific index of phagocytosis with a variety of pretreatments or conditions. The unique pHrodo™-based system measures phagocytic activity based on acidification of the particles as they are ingested, eliminating the wash and quenching steps that are necessary with nonfluorogenic indicators of bacterial uptake.¹⁻² To achieve this, the particles are conjugated to pHrodo™ dye, a novel, fluorogenic dye that dramatically increases in fluorescence as the pH of its surroundings becomes more acidic (Figure 1, page 2). We have included sufficient pHrodo™ BioParticles® Conjugate for ~100 wells in a 96-well format, with step-by-step instructions for performing this assay in a fluorescence microplate reader. With proper changes to the protocol, other plate formats can easily be adapted to platforms such as traditional imaging, high-content screening (HCS), high-throughput screening (HTS), fluorescent microplate reader, and flow cytometry. With proper settings, these reagents can also be adapted for bectop instruments, such as FLoid® Cell Imaging Station, Tali® Image-based Cytometer, and Attune® Acoustic Focusing Cytometer.

The methodology for this reagent's use was developed using adherent RAW and MMM (J774A.1) murine macrophage cells,² but can be adapted for use with other adherent cells, primary cells, or even cells in suspension.³ Cells assayed for phagocytic activity with pHrodo™ BioParticles® conjugates may also be fixed with standard 2–4% paraformaldehyde solutions for later analysis, preserving differences in signal between control and experimental samples with high fidelity for up to 48 hours. pHrodo™ BioParticles® conjugate preparations are also amenable to opsonization (Cat. nos. E2870, S2860), which can greatly enhance their uptake and signal strength in the assay.

Figure 1 The fluorescence emission spectra of pHrodo™ Red and pHrodo™ Green conjugates




Figure 2 Workflow for pHrodo™ BioParticles® conjugates





Figure 3 Effect of different concentrations of Cytochalasin D in MMM cells detected with **(A)** pHrodo™ Red E. coli BioParticles® or **(B)** pHrodo™ Green E. coli BioParticles® conjugates in RAW cells while using a fluorescent microplate reader (shown here without the background subtraction, with N=3, ± SEM).



Before You Begin

Materials Recommended but Not Provided

- Macrophage RAW and MMM cells cultured in cell growth medium. Other cell lines may also be used, if preferred.
- Uptake Buffer; we recommend using Live Cell Imaging Solution (Cat. no. A14291DJ) for best results. Alternatively, you may use any other appropriate buffer at pH 7.4.
- 96-well microplate or any other plate, with proper installment capable of detecting fluorescence according to emission/excitation maxima given in Table 1 (page 1).
- *Optional:* Opti-MEM® culture medium (Cat. no. 31985-062).
- *Optional:* Stock solutions of experimental effector compounds for testing their influence on phagocytosis (for example, cytochalasin D inhibits phagocytosis by inhibiting actin cytoskeletal rearrangements).

Experimental Protocol

Phagocytosis Assay Protocol

The following protocol describes an experimental test of phagocytic function with appropriate controls in a 96-well plate; however, other plate formats can be adopted with proper adjustments to the protocol. Cellular auto fluorescence background is determined with cells plated free of pHrodo™ BioParticles® Conjugates (but otherwise under control and experimental conditions), and reagent background fluorescence is determined using wells that contain the pHrodo™ BioParticles™ Conjugates but no cells. This protocol describes the use of one vial of particles, prepared at 1 mg/mL in the buffered saline solution of your choice. To minimize background fluorescence from non-ingested pHrodo™ BioParticles® Conjugate, we strongly recommend controlling the extracellular pH by adding Live Cell Imaging Solution (Cat. no. A14291DJ) for best results. Alternatively, you may use any other appropriate buffer at pH 7.4.

Assay Controls

To minimize experimental errors, we recommend making measurements from a minimum of three replicates of positive control, experimental, and no-cell control samples, though the numbers of experimental and control wells can be adjusted as required to meet the needs of the particular study.

Amount of BioParticles® A single vial of pHrodo™ BioParticles® conjugate dilutes to 2 mL of volume that will be used in the assay, which is distributed across 20 wells. 100 µL of this suspension is used per sample well, including no-cell background controls. The average fluorescence value of these no-cell background control wells is subtracted from all cell-containing wells at the end of the assay to yield a cell-specific, net phagocytosis signal. Note that a large specific signal can readily be obtained with or without background subtraction, as shown in Figure 3, page 3.

Preparing the Cells

- 1.1 Subculture the RAW or MMM macrophage cells (or preferred cell type) in complete medium for 3–4 days in advance of performing the assay.
- 1.2 On the day of the assay, harvest the cells from the culture plates and centrifuge the suspension. Resuspend the pellet in Opti-MEM® medium or preferred culture medium at 10⁶ cells/mL. Scale your culture to aim for ~2 × 10⁶ cells per vial of pHrodo™ BioParticles® Conjugate. Alternatively, cells can be plated into the 96-well plate a day or more in advance, with the aim of having 100,000 viable cells per well on the day of the assay.

If you are using cells other than RAW or MMM, you may need to determine optimal cell culture conditions and densities for your specific cell type. In general, better signals in the plate reader are obtained with maximal cell densities.

- 1.3 Plate the cells into a 96-well plate at 100,000 cells/well using 100 µL per well. We recommend plating your positive control and experimental wells in triplicate or greater. Be sure to leave one well empty of cells for every positive control well, so that a no-cell control background subtraction may be performed.

For example, plate four columns of four wells, leaving the fifth column of four empty for the no-cell control. Note that higher background fluorescence levels may be seen with acidic poly-D-lysine coated microplates.

- 1.4 Add 100 µL of Opti-MEM® medium or complete culture medium to the wells left aside for the no-cell background determination.
- 1.5 Cover the loaded microplate and allow the cells to settle and adhere to the microplate for at least one hour in a humidified incubator with 5% CO₂ at 37°C.
- 1.6 Prepare the experimental wells by adding the experimental phagocytosis effector (e.g., Cytochalasin D) at the desired concentrations, taking care to add vehicle controls to untreated wells. Note that the time and concentration of experimental effector pretreatment may vary greatly with the agent or treatment under study.

Preparing the BioParticles® Conjugate

- 2.1 Thaw one vial each of the pHrodo™ BioParticles® fluorescent particles for every 20 wells to be tested. This number includes the no-cell control wells that will receive fluorescent particles, but no cells. Pipette 2 mL Uptake Buffer such as the Live Cell Imaging Solution (Cat. no. A14291DJ) into the vial containing lyophilized product and briefly vortex the solution to completely resuspend the particles so that you have 1 mg/mL for pHrodo™ *E. coli* or *S. aureus* BioParticles®, or 0.5 mg/mL for pHrodo™ Zymosan A BioParticles®.
- 2.2 Transfer the suspension into a clean glass tube and sonicate for 5 minutes, until all the fluorescent particles are homogeneously dispersed.

Adding the Fluorescent Particles

- 3.1 After the cells have adhered and the phagocytosis effectors have been added, remove the culture medium from each of the microplate wells by vacuum aspiration.
- 3.2 Quickly replace the culture medium with 100 μ L of the prepared pHrodo™ BioParticles® suspension from step 2.2, adding it to the positive control, experimental, and no-cell background subtraction wells. Experimental effector solutions may be prepared ahead of time with separate vials of pHrodo™ BioParticles® suspension to keep them present throughout the assay.
- 3.3 Cover and transfer the microplate to an incubator warmed to 37°C for 1–2 hours to allow phagocytosis and acidification to reach its maximum. Although aseptic techniques have been used to produce pHrodo™ BioParticles® conjugates, these products are not sterile and incubation of more than 4 hours is not recommended. Do not use an elevated CO₂ cell culture incubator unless the Uptake Buffer in use has a bicarbonate buffering system, because elevated CO₂ levels will artificially acidify the buffer and elevate the background fluorescence.

Fluorescence Measurements and Results

- 4.1 Scan or image all experimental, control, and no-cell control wells of the microplate in the fluorescence plate reader or any other appropriate fluorescent instrument using the appropriate settings according to the excitation and emission maxima given in Table 1, page 1.
- 4.2 If using a fluorescence plate reader, calculate the net phagocytosis and the response to the phagocytosis effector agent. Net phagocytosis is calculated by subtracting the average fluorescence intensity of the no-cell negative-control wells from all positive-control and experimental wells. The mean and standard deviation of the net positive control and net experimental wells should then be calculated. The phagocytosis response to the experimental effector (% Effect) can then be calculated as a fraction of the net positive control phagocytosis as follows:

$$\% \text{ Effect} = \frac{\text{Net experimental phagocytosis} \times 100\%}{\text{Net positive control phagocytosis}}$$

References

1. J Immunol Methods 60, 115 (1983); 2. J Immunol Methods 162, 1 (1993); 3. J Biol Chem 273, 14813 (1998); 4. Methods Enzymol 58, 141 (1979).

Product List Current prices may be obtained from our website or from our Customer Service Department.

Cat #	Product Name	Unit Size
A10010	pHrodo™ Red <i>S. aureus</i> BioParticles® conjugate for phagocytosis	5 × 2 mg
P35361	pHrodo™ Red <i>E. coli</i> BioParticles® conjugate for phagocytosis	5 × 2 mg
P35364	pHrodo™ Red Zymosan A BioParticles® conjugate *for phagocytosis*	5 × 1 mg
P35365	pHrodo™ Green Zymosan A BioParticles® conjugate *for phagocytosis*	5 × 1 mg
P35366	pHrodo™ Green <i>E. coli</i> BioParticles® conjugate *for phagocytosis*	5 × 2 mg
P35367	pHrodo™ Green <i>S. aureus</i> BioParticles® conjugate *for phagocytosis*	5 × 2 mg
Related Products		
P10361	pHrodo™ Red dextran, 10,000 MW *for endocytosis*	0.5 mg
P35368	pHrodo™ Green dextran, 10,000 MW *for endocytosis*	0.5 mg
P36600	pHrodo™ Red, succinimidyl ester (pHrodo™ Red, SE)	1 mg
P35370	pHrodo™ Green C ₂ -maleimide	1 mg
P35371	pHrodo™ Red C ₂ -maleimide	1 mg
A10025	pHrodo™ Red <i>E. coli</i> BioParticles® Phagocytosis Kit *for flow cytometry* *100 tests*	1 kit
A10026	pHrodo™ Red Phagocytosis Particle Labeling Kit *for flow cytometry* *100 tests*	1 kit
P35362	pHrodo™ Red Avidin *Fluorogenic pH sensor*	1 mg
P35363	pHrodo™ Red Microscale Labeling Kit *Fluorogenic pH sensor* *3 labelings*	1 kit
A14291DJ	Live Cell Imaging Solution	500 mL
R37602	Image-iT® Fixation/Permeabilization Kit	1 kit
R37603	BackDrop™ Background Suppressor *for live cells*	1 kit
R37605	NucBlue™ Live Cell Stain *Hoechst 33342 special formulation*	1 kit
R37606	NucBlue™ Fixed Cell Stain *DAPI special formulation*	1 kit
E2870	<i>Escherichia coli</i> BioParticles® opsonizing reagent	1 unit
S2860	<i>Staphylococcus aureus</i> BioParticles® opsonizing reagent	1 unit

Related Platforms



Attune® Acoustic Focusing Cytometer
(Cat. no. 4469120)



Tali® Image-based Cytometer
(Cat. no. T10796)



FLoid® Cell Imaging Station
(Cat. no. 4471136)

Purchaser Notification

Corporate Headquarters

5791 Van Allen Way
Carlsbad, CA 92008
USA
Phone: +1 760 603 7200
Fax: +1 760 602 6500
Email: techsupport@lifetech.com

European Headquarters

Inchinnan Business Park
3 Fountain Drive
Paisley PA4 9RF
UK
Phone: +44 141 814 6100
Toll-Free Phone: 0800 269 210
Toll-Free Tech: 0800 838 380
Fax: +44 141 814 6260
Tech Fax: +44 141 814 6117
Email: euroinfo@invitrogen.com
Email Tech: eurotech@invitrogen.com

Japanese Headquarters

LOOP-X Bldg. 6F
3-9-15, Kaigan
Minato-ku, Tokyo 108-0022
Japan
Phone: +81 3 5730 6509
Fax: +81 3 5730 6519
Email: jpinfo@invitrogen.com

Additional international offices are listed at
www.lifetechnologies.com

These high-quality reagents and materials must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Read the Safety Data Sheet provided for each product; other regulatory considerations may apply.

Obtaining Support

For the latest services and support information for all locations, go to www.lifetechnologies.com.

At the website, you can:

- Access worldwide telephone and fax numbers to contact Technical Support and Sales facilities
- Search through frequently asked questions (FAQs)
- Submit a question directly to Technical Support (techsupport@lifetech.com)
- Search for user documents, SDSs, vector maps and sequences, application notes, formulations, handbooks, certificates of analysis, citations, and other product support documents
- Obtain information about customer training
- Download software updates and patches

SDS

Safety Data Sheets (SDSs) are available at www.lifetechnologies.com/sds.

Certificate of Analysis

The Certificate of Analysis provides detailed quality control and product qualification information for each product. Certificates of Analysis are available on our website. Go to www.lifetechnologies.com/support and search for the Certificate of Analysis by product lot number, which is printed on the product packaging (tube, pouch, or box).

Limited Product Warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

Disclaimer

LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

Limited Use Label License: Research Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial services of any kind, including, without limitation, reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, California 92008.

The trademarks mentioned herein are the property of Life Technologies Corporation and/or its affiliate(s) or their respective owners.

©2013 Life Technologies Corporation. All rights reserved.