

Ammonium Chloride Solution



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Catalog # 07800	100 mL
07850	500 mL

FOR RESEARCH USE ONLY. NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES.

Product Description

Ammonium Chloride Solution is recommended for lysis of red blood cells (RBCs) in preparations of human umbilical cord blood, and human or mouse peripheral blood, spleen, or bone marrow cells.

This RBC lysis solution is buffered and optimized for gentle lysis of erythrocytes, with minimal effect on leukocytes. This Ammonium Chloride Solution does not contain a fixative agent, therefore leukocytes are viable following RBC lysis.

Properties

- Storage:** Store at -20°C.
- Shelf Life:** Stable until expiry date (EXP) on label.
- Contains:**
- 0.8% NH₄Cl
 - 0.1 mM EDTA in water
- Buffered with KHCO₃ to achieve a final pH of 7.2 - 7.6.

Handling / Directions For Use

1. Thaw Ammonium Chloride Solution overnight at 2 - 8°C.
NOTE: If not used immediately store at 2 - 8°C for up to 2 weeks.
2. Add appropriate volume of Ammonium Chloride Solution to the sample being processed (see below).
3. Thoroughly mix the cell suspension by inverting the tube and place on ice for 10 minutes.
4. Wash the cells twice in the appropriate medium prior to use.

HUMAN BONE MARROW

Add Ammonium Chloride Solution to the sample at a volume:volume ratio of 4:1.
For example, add 4 mL of Ammonium Chloride Solution to 1 mL of sample.

HUMAN PERIPHERAL BLOOD, CORD BLOOD

Add Ammonium Chloride Solution to the sample at a volume:volume ratio of 9:1.
For example, add 9 mL of Ammonium Chloride Solution to 1 mL of sample.

MOUSE BONE MARROW, SPLEEN, PERIPHERAL BLOOD

Add Ammonium Chloride Solution to the sample at a volume:volume ratio of 9:1.
For example, add 9 mL of Ammonium Chloride Solution to 1 mL of sample.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485 MEDICAL DEVICE STANDARDS.

Copyright © 2015 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Inc. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.