



pEF1/V5-His A, B, and C

Catalog no. V920-20

Rev. date: 7 June 2010

Manual part no. 28-0184

MAN0000664

User Manual

Table of Contents

Kit Contents and Storage.....	iv
Introduction	1
Product Overview	1
Methods	2
Cloning into pEF1/V5-His A, B, and C.....	2
Transformation and Transfection.....	6
Appendix	9
Human EF-1 α Promoter	9
pEF1/V5-His Vector	10
pEF1/V5-His/ <i>lacZ</i>	12
Accessory Products	13
Technical Support.....	15
Purchaser Notification	16
References.....	17

Kit Contents and Storage

Shipping and Storage

pEF1/V5-His vectors are shipped on wet ice. Upon receipt, store vectors at -20°C.

Kit Contents

All vectors are supplied as detailed below. **Store the vectors at -20°C.**

Vector	Composition	Amount
pEF1/V5-His A, B, and C	40 µL of 0.5 µg/µL vector in 10 mM Tris-HCl, 1 mM EDTA, pH 8.0	20 µg
pEF1/V5-His / <i>lacZ</i>	40 µL of 0.5 µg/µL vector in 10 mM Tris-HCl, 1 mM EDTA, pH 8.0	20 µg

Intended Use

For research use only. Not intended for human or animal diagnostic or therapeutic uses.

Introduction

Product Overview

Description of the System

pEF1/V5-His A, B, and C are 6.2 kb vectors derived from pcDNA™3.1/V5-His and designed for overproduction of recombinant proteins in mammalian cell lines. Features of the vectors allow purification and detection of expressed proteins (see pages 10–11 for more information). High-level stable and transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

- Human elongation factor 1 α -subunit promoter (hEF-1 α) for high-level expression across a broad range of species and cell types (Goldman *et al.*, 1996; Mizushima and Nagata, 1990) (see page 9 for more information).
- Three reading frames to facilitate in-frame cloning with a C-terminal tag encoding the V5 epitope and a polyhistidine metal-binding peptide.
- Episomal replication in cell lines that are latently infected with SV40 or that express the SV40 large T antigen (*e.g.* COS7).

The control plasmid, pEF1/V5-His/*lacZ*, is the pEF1/V5-His A vector with a 3.1 kb fragment containing the β -galactosidase gene cloned in frame with the C-terminal peptide (see page 12). It is included for use as a positive control for transfection, expression, purification, and detection in the cell line of choice.

Experimental Outline

Use the following outline to clone and express your gene of interest in pEF1/V5-His.

1. Consult the multiple cloning sites described on pages 3–5 to determine which vector (A, B, or C) should be used to clone your gene in frame with the C-terminal V5 epitope and polyhistidine tag.
 2. Ligate your insert into the appropriate vector and transform into *E. coli*. Select transformants on 50 to 100 $\mu\text{g}/\text{mL}$ ampicillin.
 3. Analyze your transformants for the presence of insert by restriction digestion.
 4. Select a transformant with the correct restriction pattern and confirm that your gene is in frame with the C-terminal peptide by sequencing.
 5. Transfect your construct into the cell line of choice.
 6. Test for expression of your recombinant gene by western blot analysis or functional assay. If you do not have an antibody to your protein, Invitrogen offers the Anti-V5 Antibody or the Anti-His(C-term) Antibody to detect your recombinant protein. See page 14 for more information.
 7. To purify your recombinant protein, you may use metal-chelating resin such as ProBond™. ProBond™ resin is available separately.
-

Methods

Cloning into pEF1/V5-His A, B, and C

General Molecular Biology Techniques

For help with DNA ligations, *E. coli* transformations, restriction enzyme analysis, purification of single-stranded DNA, DNA sequencing, and DNA biochemistry, refer to *Molecular Cloning: A Laboratory Manual* (Sambrook *et al.*, 1989) or *Current Protocols in Molecular Biology* (Ausubel *et al.*, 1994).

E. coli Strain

Many *E. coli* strains are suitable for the growth of this vector including TOP10F', DH5 α F', JM109, and INV α F'. We recommend that you propagate vectors containing inserts in *E. coli* strains that are recombination deficient (*recA*) and endonuclease A deficient (*endA*).

Maintaining pEF1/V5-His

To propagate and maintain the pEF1/V5-His vectors, use a small amount of the supplied 0.5 μ g/ μ L stock solution in TE, pH 8.0 to transform a *recA*, *endA* *E. coli* strain like TOP10F', DH5 α , JM109, or equivalent. Select transformants on LB plates containing 50–100 μ g/mL ampicillin. Be sure to prepare a glycerol stock of each plasmid for long term storage.

Cloning Considerations

Your insert should contain a Kozak consensus sequence with an ATG initiation codon for proper initiation of translation (Kozak, 1987; Kozak 1990). An example of a Kozak consensus sequence is provided below. Other sequences are possible, but the G or A at position -3 and the G at position +4 (shown in bold) illustrates the most commonly occurring sequence with strong consensus. Replacing one of the two bases at these positions provides moderate consensus, while having neither results in weak consensus. The ATG initiation codon is shown underlined.

(G/A)NNATGG

If you wish to express your protein WITHOUT the C-terminal peptide, be sure to include a stop codon.

Continued on next page

Cloning into pEF1/V5-His A, B, and C, Continued

Multiple Cloning Site of Version A

Below is the multiple cloning site for pEF1/V5-His A. Restriction sites are labeled to indicate the cleavage site. The boxed nucleotides indicate the variable region. **Note that there is a stop codon between the *Spe* I site and the *Bst*X I site.** The multiple cloning site has been confirmed by sequencing and functional testing. For more information on the hEF-1 α promoter, see page 9.

```

3' end of hEF-1 $\alpha$  Intron 1
1579  GTTTGGATCT TGGTTCATTC TCAAGCCTCA GACAGTGGTT CAAAGTTTTT TTCTTCCATT TCAGGTGTCG TGAGGAATTA
5' end of hEF-1 $\alpha$  Exon 2

1659  GCTTGGTACT AATACGACTC ACTATAGGGA GACCCAAGCT GGCTAGGTAA GCT TGG TAC CGA GCT CGG ATC CAC
T7 promoter/priming site Kpn I BamH I Spe I
Trp Tyr Arg Ala Arg Ile His

1733  TAG TCC AGT GTG GTG GAA TTC TGC AGA TAT CCA GCA CAG TGG CGG CCG CTC GAG TCT AGA GGG CCC
BstX I EcoR I EcoR V BstX I Not I Xba I
*** Ser Ser Val Val Glu Phe Cys Arg Tyr Pro Ala Gln Trp Arg Pro Leu Glu Ser Arg Gly Pro

1799  TTC GAA GGT AAG CCT ATC CCT AAC CCT CTC CTC GGT CTC GAT TCT ACG CGT ACC GGT CAT CAT CAC
BstB I V5 epitope Polyhistidine
Phe Glu Gly Lys Pro Ile Pro Asn Pro Leu Leu Gly Leu Asp Ser Thr Arg Thr Gly His His His

1865  CAT CAC CAT TGA GT TTAAACCCGC TGATCAGCCT CGACTGTGCC TTCTAGTTGC CAGCCATCTG TTGTTTGCCC
region Pme I BGH Reverse priming site
His His His ***

1939  CTCCCCGTG CCTTCCTGA CCCTGGAAGG TGCCACTCCC ACTGTCCTTT CTAATAAAA TGAGGAAATT GCATGCGATT
BGH polyadenylation signal

2019  GTCTGAGTAG GTGTCATTCT ATTCTGGGGG GTGGGGTGGG GCAGGACAGC AAGGGGGAGG ATTGGGAAGA CAATAGCAGG

```

Continued on next page

Cloning into pEF1/V5-His A, B, and C, Continued

Multiple Cloning Site of Version B

Below is the multiple cloning site for pEF1/V5-His B. Restriction sites are labeled to indicate the cleavage site. The boxed nucleotides indicate the variable region. The multiple cloning site has been confirmed by sequencing and functional testing. For more information on the hEF-1 α promoter, see page 9.

3' end of hEF-1 α Intron 1
|
5' end of hEF-1 α Exon 2

1579 GTTTGGATCT TGGTTCATTC TCAAGCCTCA GACAGTGGTT CAAAGTTTTT TTCTTCCATT TCAGGTGTCG TGAGGAATTA

Kpn I BamH I Spe I

1659 GCTTGGTACT AATACGACTC ACTATAGGGA GACCCAAGCT GGCTAGGTAA G CTT GGT ACC GAG CTC GGA TCC ACT
Leu Gly Thr Glu Leu Gly Ser Thr

T7 promoter/priming site

BstX I EcoR I EcoR V BstX I Not I Xba I

1734 AGT CCA GTG TGG TGG AAT TCT GCA GAT ATC CAG CAC AGT GGC GGC CGC TCG AGT CTA GAG GGC CCG
Ser Pro Val Trp Trp Asn Ser Ala Asp Ile Gln His Ser Gly Gly Arg Ser Ser Leu Glu Gly Pro

BstB I V5 epitope Polyhistidine

1800 CGG TTC GAA GGT AAG CCT ATC CCT AAC CCT CTC CTC GGT CTC GAT TCT ACG CGT ACC GGT CAT CAT
Arg Phe Glu Gly Lys Pro Ile Pro Asn Pro Leu Leu Gly Leu Asp Ser Thr Arg Thr Gly His His

region Pme I BGH Reverse priming site

1866 CAC CAT CAC CAT TGA GTTTAAAC CCGCTGATCA GCCTCGACTG TGCCTTCTAG TTGCCAGCCA TCTGTTGTTT
His His His His ***

BGH polyadenylation signal

1939 GCCCCTCCC CGTGCCTTCC TTGACCCTGG AAGGTGCCAC TCCCCTGTC CTTTCCTAAT AAAATGAGGA AATTGCATCG

2019 CATTGTCTGA GTAGGTGTCA TTCTATTCTG GGGGGTGGGG TGGGGCAGGA CAGCAAGGGG GAGGATTGGG AAGACAATAG

Continued on next page

Cloning into pEF1/V5-His A, B, and C, Continued

Multiple Cloning Site of Version C

Below is the multiple cloning site for pEF1/V5-His C. Restriction sites are labeled to indicate the cleavage site. The boxed nucleotides indicate the variable region. The multiple cloning site has been confirmed by sequencing and functional testing. For more information on the hEF-1 α promoter, see page 9.

```

                                                    3' end of hEF-1 $\alpha$  Intron 1
                                                    |
1579  GTTTGGATCT TGGTTCATTC TCAAGCCTCA GACAGTGGTT CAAAGTTTTT TTCTTCCATT TCAGGTGTCG TGAGGAATTA
                                                    5' end of hEF-1 $\alpha$  Exon 2

1659  GCTTGGTACT AATACGACTC ACTATAGGGA GACCCAAGCT GGCTAGGTAA GC TTG GTA CCG AGC TCG GAT CCA CTA
      T7 promoter/priming site
      Kpn I      BamH I      Spe I
      Leu Val Pro Ser Ser Asp Pro Leu

1735  GTC CAG TGT GGT GGA ATT CTG CAG ATA TCC AGC ACA GTG GCG GCC GCT CGA GGT CAC CCA TTC GAA
      BstX I  EcoR I      EcoR V      BstX I  Not I      BstE II      BstB I
      Val Gln Cys Gly Gly Ile Leu Gln Ile Ser Ser Thr Val Ala Ala Ala Arg Gly His Pro Phe Glu

1801  GGT AAG CCT ATC CCT AAC CCT CTC CTC GGT CTC GAT TCT ACG CGT ACC GGT CAT CAT CAC CAT CAC
      V5 epitope
      Polyhistidine region
      Gly Lys Pro Ile Pro Asn Pro Leu Leu Gly Leu Asp Ser Thr Arg Thr Gly His His His His His

1867  CAT TGA GTTTAA ACCCGCTGAT CAGCCTCGAC TGTGCCTTCT AGTTGCCAGC CATCTGTTGT TTGCCCTCC
      Pme I      BGH Reverse priming site
      His ***

1939  CCCGTGCCTT CCTTGACCCT GGAAGGTGCC ACTCCCACTG TCCTTTCCTA ATAAAATGAG GAAATTGCAT CGCATTGTCT
      BGH polyadenylation signal

2019  GAGTAGGTGT CATTCTATTC TGGGGGGTGG GGTGGGGCAG GACAGCAAGG GGGAGGATTG GGAAGACAAT AGCAGGCATG
  
```

Transformation and Transfection

E. coli **Transformation**

Transform your ligation mixtures into a competent *recA*, *endA* *E. coli* strain (e.g., TOP10F', DH5 α) and select on LB plates containing 50–100 $\mu\text{g}/\text{mL}$ ampicillin. Select 10–20 clones and analyze for the presence and orientation of your insert.



We recommend that you sequence your construct with the T7 Forward and BGH Reverse primers sequences to confirm that your gene is fused in frame with the V5 epitope and the C-terminal polyhistidine tag. Refer to the diagrams on pages 3–5 for location and sequence of primer binding sites.

Plasmid **Preparation**

Plasmid DNA for transfection into eukaryotic cells must be very clean and free from phenol and sodium chloride. Contaminants will kill the cells and salt will interfere with lipid complexing decreasing transfection efficiency. We recommend isolating plasmid DNA using the PureLink™ HiPure Miniprep Kit or the PureLink™ HiPure Midiprep Kit (see page 13 for ordering information).

Methods of **Transfection**

For established cell lines (e.g. HeLa), consult original references or the supplier of your cell line for the optimal method of transfection. It is recommended that you follow exactly the protocol for your cell line. Pay particular attention to medium requirements, when to pass the cells, and at what dilution to split the cells. Further information is provided in *Current Protocols in Molecular Biology*.

Methods for transfection include calcium phosphate (Chen and Okayama, 1987; Wigler *et al.*, 1977), lipid-mediated (Felgner *et al.*, 1989; Felgner and Ringold, 1989), and electroporation (Chu *et al.*, 1987; Shigekawa and Dower, 1988). Invitrogen offers the Lipofectamine™ 2000 Reagent for mammalian transfection.

Positive Control

pEF1/V5-His/*lacZ* is provided as a positive control vector for mammalian transfection and expression (see page 12). It may be used to optimize transfection conditions for your cell line. The gene encoding β -galactosidase is expressed in mammalian cells under the hEF-1 α promoter. A successful transfection will result in β -galactosidase expression that can be easily assayed (see next page).

Continued on next page

Transformation and Transfection, Continued

Assay for β -galactosidase Activity

You may assay for β -galactosidase expression by activity assay using cell-free lysates (Miller, 1972) or by staining the cells for activity. Invitrogen offers the β -Gal Assay Kit and the β -Gal Staining Kit for fast and easy detection of β -galactosidase expression.

Detecting Fusion Proteins

A number of antibodies are available from Invitrogen that can be used to detect expression of your fusion protein from pEF1/V5-His (see page 13).

Neomycin (Geneticin[®]) Activity

Geneticin[®] blocks protein synthesis in mammalian cells by interfering with ribosomal function. It is an aminoglycoside, similar in structure to neomycin, gentamycin, and kanamycin. Expression of the bacterial aminoglycoside phosphotransferase gene (APH), derived from Tn5, in mammalian cells results in detoxification of Geneticin[®] (Southern and Berg, 1982).

Geneticin[®] Selection Guidelines

Geneticin[®] is available from Invitrogen (see page 13). Use as follows:

- Prepare Geneticin[®] in a buffered solution (*e.g.*, 100 mM HEPES, pH 7.3).
- Use 100 to 1,000 $\mu\text{g}/\text{mL}$ of Geneticin[®] in complete medium.
- Calculate concentration based on the amount of active drug (check the lot label).
- Test varying concentrations of Geneticin[®] on your cell line to determine the concentration that kills your cells (kill curve). Cells differ in their susceptibility to Geneticin[®].

Cells will divide once or twice in the presence of lethal doses of Geneticin[®], so the effects of the drug take several days to become apparent. Complete selection can take from 3 to 6 weeks of growth in selective medium.

Continued on next page

Transformation and Transfection, Continued

Preparing Cells for Lysis

Use the procedure below to prepare cells for lysis prior to purification of your protein on ProBond™. You will need 5×10^6 to 1×10^7 cells for purification of your protein on a 2 mL ProBond™ column (see the ProBond™ Purification manual).

1. Seed cells in five T-75 flasks or 2 to 3 T-175 flasks.
 2. Grow the cells in selective medium until they are 80–90% confluent.
 3. Harvest the cells by treating with trypsin-EDTA for 2 to 5 minutes or by scraping the cells in PBS.
 4. Inactivate the trypsin by diluting with fresh medium (if necessary) and transfer the cells to a sterile microcentrifuge tube.
 5. Centrifuge the cells at $240 \times g$ for 5 minutes. You may lyse the cells immediately or freeze in liquid nitrogen and store at -80°C until needed.
-

Lysis of Cells

If you are using ProBond™ resin, refer to the ProBond™ Purification manual for details about sample preparation for chromatography. If you are using another metal-chelating resin, refer to the manufacturer's instruction for recommendations on sample preparation.

Appendix

Human EF-1 α Promoter

Description

The diagram below shows all the features of the EF-1 α promoter used in pEF1/V5-His vectors (Mizushima and Nagata, 1990). Features are marked as per Uetsuki, *et al.*, 1989.

```


      ┌── 5' end of human EF-1 $\alpha$  promoter
461  GGAGTGCCTC GTGAGGCTCC GGTGCCCGTC AGTGGGCAGA GCGCACATCG CCCACAGTCC
521  CCGAGAAGTT GGGGGGAGGG GTCGGCAATT GAACCGGTGC CTAGAGAAGG TGGCGCGGGG
581  TAAACTGGGA AAGTGATGTC GTGTACTGGC TCCGCCTTTT TCCCAGGGT GGGGAGAAC
      TATA box                               Start of Transcription
641  CGTATATAAG TGCAGTAGTC GCCGTGAACG TTCTTTTTCG CAACGGGTTT GCCGCCAGAA
      ┌── 5' end of Intron 1
701  CACAGGTAAG TGCCGTGTGT GTTCCCGCG GGCCTGGCCT CTTTACGGGT TATGGCCCTT
761  GCGTGCCTTG AATTACTTCC ACCTGGCTGC AGTACGTGAT TCTTGATCCC GAGCTTCGGG
821  TTGGAAGTGG GTGGGAGAGT TCGAGGCCCT GCGCTTAAGG AGCCCTTCG CCTCGTGCTT
881  GAGTTGAGGC CTGGCCTGGG CGCTGGGGCC GCCGCGTGGC AATCTGGTGG CACCTTCGCG
941  CCTGTCTCGC TGCTTTCGAT AAGTCTCTAG CCATTTAAAA TTTTGTATGA CCTGCTGCGA
1001 CGCTTTTTTT CTGGCAAGAT AGTCTTGTA ATGCGGGCCA AGATCTGCAC ACTGGTATTT
1061 CGGTTTTTGG GCGCGCGGGC GCGCACGGGG CCGTGCGTCC CCAGGCACA TGTTCGGCGA
      Sp 1
1121 GCGGGGCCT GCGAGCGCG CCACCGAGAA TCGGACGGGG GTAGTCTCAA GCTGGCCGGC
      Sp 1                               Sp 1
1181 CTGCTCTGGT GCCTGGCCTC GCGCCCGCGT GTATCGCCCC GCCCTGGGCG GCAAGGCTGG
1241 CCCGGTCGGC ACCAGTTGCG TGAGCGGAAA GATGGCCGCT TCCCAGCCCT GCTGCAGGGA
1301 GCTCAAAATG GAGGACGCGG CGCTCGGGAG AGCGGGCGGG TGAGTCACCC ACACAAAGGA
      Sp 1
1361 AAAGGCCTT TCCGTCCTCA GCCGTCGCTT CATGTGACTC CACGGAGTAC CGGGCGCCGT
1421 CCAGGCACCT CGATTAGTTC TCGAGCTTTT GGAGTACGTC GTCTTTAGGT TGGGGGAGG
1481 GGTTTTATGC GATGGAGTTT CCCACACTG AGTGGGTGGA GACTGAAGTT AGGCCAGCTT
1541 GGCACCTGAT GTAATTCTCC TTGGAATTG CCCTTTTGA GTTTGGATCT TGGTTCATTC
1601 TCAAGCCTCA GACAGTGGTT CAAAGTTTTT TTCTTCCATT TCAGGTGTCG TGA...
      3' end of Intron 1 ─┘
      5' end of Exon 2

```

pEF1/V5-His Vector

Map of pEF1/V5-His

The figure below summarizes the features of the pEF1/V5-His vectors. The sequences for pEF1/V5-His A, B, and C are available for downloading from www.invitrogen.com or by contacting **Technical Support** (page 15).



Comments for pEF1/V5-His A 6174 nucleotides

- EF-1 α promoter: bases 474-1651
- T7 promoter/priming site: bases 1668-1687
- Multiple cloning site: bases 1713-1881
- V5 epitope: bases 1805-1846
- Polyhistidine tag: bases 1856-1873
- BGH reverse priming site: bases 1896-1913
- BGH polyadenylation sequence: bases 1895-2122
- f1 origin of replication: bases 2172-2600
- SV40 promoter and origin: bases 2628-2935
- Neomycin resistance gene (ORF): bases 3010-3804
- SV40 polyadenylation sequence: bases 3980-4110
- pUC origin: bases 4493-5155 (opposite strand)
- Ampicillin resistance gene (ORF): bases 5300-6160 (opposite strand)

* Frame-dependent variations.
In Version C, there is a unique BstE II site, but no Xba I site.

Continued on next page

pEF1/V5-His Vector, Continued

Features of pEF1/V5-His

pEF1/V5-His A (6174 bp), pEF1/V5-His B (6178 bp), and pEF1/V5-His C (6170 bp) contain the following elements. All features have been functionally tested.

Feature	Benefit
Human elongation factor 1 α (hEF-1 α) promoter	Allows overexpression of your recombinant protein in a broad range of mammalian cell types (Goldman <i>et al.</i> , 1996; Mizushima and Nagata, 1990).
T7 promoter/priming site	Allows for <i>in vitro</i> transcription in the sense orientation and sequencing through the insert.
Multiple cloning site in three reading frames	Allows insertion of your gene and facilitates cloning in frame with the V5 epitope and polyhistidine C-terminal tag.
V5 epitope (Gly-Lys-Pro-Ile-Pro-Asn-Pro-Leu-Leu-Gly-Leu-Asp-Ser-Thr)	Allows detection of your recombinant protein with the Anti-V5 Antibody or Anti-V5-HRP Antibody (Southern <i>et al.</i> , 1991) (see page 14 for ordering information).
C-terminal polyhistidine tag	Allows purification of your recombinant protein on metal-chelating resin such as ProBond™. In addition, the C-terminal polyhistidine tag is the epitope for the Anti-His(C-term) Antibody and the Anti-His.(C-term)-HRP Antibody (see page 14 for ordering).
BGH reverse priming site	Allows sequencing through the insert.
Bovine growth hormone (BGH) polyadenylation signal	Efficient transcription termination and polyadenylation of mRNA (Goodwin and Rottman, 1992).
f1 origin	Allows rescue of single-stranded DNA.
SV40 early promoter and origin	Allows efficient, high-level expression of the neomycin resistance gene and episomal replication in cells expressing the SV40 large T antigen.
Neomycin (Geneticin®) resistance gene	Selection of stable transfectants in mammalian cells (Southern and Berg, 1982).
SV40 polyadenylation signal	Efficient transcription termination and polyadenylation of mRNA.
pUC origin	High-copy number replication and growth in <i>E. coli</i> .
Ampicillin resistance gene (β -lactamase)	Selection of vector in <i>E. coli</i> .


pEF1/V5-His/lacZ

Description

pEF1/V5-His/lacZ is a 9186 bp control vector containing the gene for β -galactosidase. pEF1/V5-His A was digested with *Kpn* I and *Pme* I. A 3.1 kb *Kpn* I-*Pme* I fragment containing the β -galactosidase gene was then ligated into the digested pEF1/V5-His A vector in frame with the C-terminal peptide.

Map of Control Vector

The figure below summarizes the features of the pEF1/V5-His/lacZ vector. The nucleotide sequence for pEF1/V5-His/lacZ is available for downloading from www.invitrogen.com or by contacting **Technical Support** (page 15).



Comments for pEF1/V5-His/lacZ 9186 nucleotides

- EF-1 α promoter: bases 470-1653
- LacZ portion of the fusion: 1721-4777
- V5 epitope: bases 4805-4846
- Polyhistidine tag: bases 4856-4873
- BGH reverse priming site: bases 4896-4913
- BGH polyadenylation sequence: bases 4899-5126
- f1 origin of replication: bases 5172-5600
- SV40 promoter and origin: bases 5654-5936
- Neomycin resistance gene (ORF): bases 6011-6805
- SV40 polyadenylation sequence: bases 6981-7111
- pUC origin: bases 7494-8167 (opposite strand)
- Ampicillin resistance gene (ORF): bases 8312-9172 (opposite strand)

Accessory Products

Introduction

The following products may be used with the pEF1/V5-His vectors. For details, visit www.invitrogen.com or contact **Technical Support** (page 15).

Item	Amount	Catalog no.
ProBond™ Purification System	6 × 2 mL precharged, prepacked ProBond™ resin columns and buffers for native and denaturing purification	K850-01
ProBond™ Resin	50 mL	R801-01
	150 mL	R801-15
Electrocomp™ TOP10F'	5 × 80 µL	C665-55
One Shot® TOP10F' Chemically Competent <i>E. coli</i>	20 × 50 µL	C3030-03
PureLink™ HiPure Plasmid Miniprep Kit	100 preps	K2100-03
PureLink™ HiPure Plasmid Midiprep Kit	25 preps	K2100-04
β-Gal Assay Kit	80 mL	K1455-01
β-Gal Staining Kit	1 kit	K1465-01
Geneticin®	1 gram	11811-023
	5 grams	11811-031
	25 grams	11811-098
Lipofectamine™ 2000 Reagent	0.75 mL	11668-027

Primers

For your convenience, Invitrogen offers a custom primer synthesis service. Visit www.invitrogen.com for more details.

Continued on next page

Accessory Products, Continued

Antibodies

If you do not have an antibody specific to your protein, Invitrogen offers the Anti-V5, or Anti-His(C-term) antibodies to detect your recombinant fusion protein. Horseradish peroxidase (HRP)- and alkaline phosphatase (AP)-conjugated antibodies are available for convenient one-step detection.

Antibody	Epitope	Catalog no.
Anti-V5	Detects a 14 amino acid epitope derived from the P and V proteins of the paramyxovirus, SV5 (Southern et al., 1991): GKPIP NPLLGLDST	R960-25
Anti-V5-HRP		R961-25
Anti-V5-AP		R962-25
Anti-His(C-term)	Detects the C-terminal polyhistidine tag (requires the free carboxyl group for detection) (Lindner et al., 1997): HHHHHH-COOH	R930-25
Anti-His(C-term)-HRP		R931-25
Anti-His(C-term)-AP		R932-25

Technical Support

Web Resources



Visit the Invitrogen website at www.invitrogen.com for:

- Technical resources, including manuals, vector maps and sequences, application notes, SDSs, FAQs, formulations, citations, handbooks, etc.
 - Complete technical support contact information
 - Access to the Invitrogen Online Catalog
 - Additional product information and special offers
-

Contact Us

For more information or technical assistance, call, write, fax, or email. Additional international offices are listed on our website (www.invitrogen.com).

Corporate Headquarters:

5791 Van Allen Way
Carlsbad, CA 92008 USA
Tel: 1 760 603 7200
Tel (Toll Free): 1 800 955 6288
Fax: 1 760 602 6500
E-mail: tech_support@invitrogen.com

Japanese Headquarters:

LOOP-X Bldg. 6F
3-9-15, Kaigan
Minato-ku, Tokyo 108-0022
Tel: 81 3 5730 6509
Fax: 81 3 5730 6519
E-mail: jpinfo@invitrogen.com

European Headquarters:

Inchinnan Business Park
3 Fountain Drive
Paisley PA4 9RF, UK
Tel: +44 (0) 141 814 6100
Tech Fax: +44 (0) 141 814 6117
E-mail: eurotech@invitrogen.com

SDS

Safety Data Sheets (SDSs) are available on our website at www.invitrogen.com/sds.

Certificate of Analysis

The Certificate of Analysis provides detailed quality control and product qualification information for each product. Certificates of Analysis are available on our website. Go to www.invitrogen.com/support and search for the Certificate of Analysis by product lot number, which is printed on the box.

Limited Warranty

Invitrogen (a part of Life Technologies Corporation) is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, contact our Technical Support Representatives. All Invitrogen products are warranted to perform according to specifications stated on the certificate of analysis. The Company will replace, free of charge, any product that does not meet those specifications. This warranty limits the Company's liability to only the price of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. The Company reserves the right to select the method(s) used to analyze a product unless the Company agrees to a specified method in writing prior to acceptance of the order.

Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore the Company makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Support Representatives.

Life Technologies Corporation shall have no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.

Purchaser Notification

**Limited Use Label
License No. 22:
Vectors and
Clones Encoding
Histidine Hexamer**

This product is licensed under U.S. Patent Nos. 5,284,933 and 5,310,663 and foreign equivalents from Hoffmann-LaRoche, Inc., Nutley, NJ and/or Hoffmann-LaRoche Ltd., Basel, Switzerland and is provided only for use in research. Information about licenses for commercial use is available from QIAGEN GmbH, Max-Volmer-Str. 4, D-40724 Hilden, Germany.

**Limited Use Label
License No. 60:
EF-1 α Promoter**

EF-1 α promoter products are sold under license for research purposes only. The use of this product for any commercial purpose, including but not limited to, use in any study for the purpose of a filing of a new drug application, requires a license from: Mochida Pharmaceutical Co., Ltd., 7, Yotsuya 1-Chome, Shinjuku-Ku, Tokyo 160, Japan. Tel: 81-3-3225-5451; Fax: 81-3-3225-6091.

References

- Ausubel, F. M., Brent, R., Kingston, R. E., Moore, D. D., Seidman, J. G., Smith, J. A., and Struhl, K. (1994). *Current Protocols in Molecular Biology* (New York: Greene Publishing Associates and Wiley-Interscience).
- Chen, C., and Okayama, H. (1987). High-Efficiency Transformation of Mammalian Cells by Plasmid DNA. *Molec. Cell. Biol.* *7*, 2745-2752.
- Chu, G., Hayakawa, H., and Berg, P. (1987). Electroporation for the Efficient Transfection of Mammalian Cells with DNA. *Nucleic Acids Res.* *15*, 1311-1326.
- Felgner, P. L., Holm, M., and Chan, H. (1989). Cationic Liposome Mediated Transfection. *Proc. West. Pharmacol. Soc.* *32*, 115-121.
- Felgner, P. L., and Ringold, G. M. (1989). Cationic Liposome-Mediated Transfection. *Nature* *337*, 387-388.
- Goldman, L. A., Cutrone, E. C., Kotenko, S. V., Krause, C. D., and Langer, J. A. (1996). Modifications of Vectors pEF-BOS, pcDNA1, and pcDNA3 Result in Improved Convenience and Expression. *BioTechniques* *21*, 1013-1015.
- Goodwin, E. C., and Rottman, F. M. (1992). The 3'-Flanking Sequence of the Bovine Growth Hormone Gene Contains Novel Elements Required for Efficient and Accurate Polyadenylation. *J. Biol. Chem.* *267*, 16330-16334.
- Kozak, M. (1987). An Analysis of 5'-Noncoding Sequences from 699 Vertebrate Messenger RNAs. *Nucleic Acids Res.* *15*, 8125-8148.
- Kozak, M. (1990). Downstream Secondary Structure Facilitates Recognition of Initiator Codons by Eukaryotic Ribosomes. *Proc. Natl. Acad. Sci. USA* *87*, 8301-8305.
- Miller, J. H. (1972). *Experiments in Molecular Genetics* (Cold Spring Harbor, New York: Cold Spring Harbor Laboratory).
- Mizushima, S., and Nagata, S. (1990). pEF-BOS, a Powerful Mammalian Expression Vector. *Nucleic Acids Res.* *18*, 5322.
- Sambrook, J., Fritsch, E. F., and Maniatis, T. (1989). *Molecular Cloning: A Laboratory Manual*, Second Edition (Plainview, New York: Cold Spring Harbor Laboratory Press).
- Shigekawa, K., and Dower, W. J. (1988). Electroporation of Eukaryotes and Prokaryotes: A General Approach to the Introduction of Macromolecules into Cells. *BioTechniques* *6*, 742-751.
- Southern, J. A., Young, D. F., Heaney, F., Baumgartner, W., and Randall, R. E. (1991). Identification of an Epitope on the P and V Proteins of Simian Virus 5 That Distinguishes Between Two Isolates with Different Biological Characteristics. *J. Gen. Virol.* *72*, 1551-1557.
- Southern, P. J., and Berg, P. (1982). Transformation of Mammalian Cells to Antibiotic Resistance with a Bacterial Gene Under Control of the SV40 Early Region Promoter. *J. Molec. Appl. Gen.* *1*, 327-339.
- Uetsuki, T., Naito, A., Nagata, S., and Kaziro, Y. (1989). Isolation and Characterization of the Human Chromosomal Gene for Polypeptide Chain Elongation Factor-1 α . *J. Biol. Chem.* *264*, 5791-5798.
- Wigler, M., Silverstein, S., Lee, L.-S., Pellicer, A., Cheng, Y.-C., and Axel, R. (1977). Transfer of Purified Herpes Virus Thymidine Kinase Gene to Cultured Mouse Cells. *Cell* *11*, 223-232.

Notes



Corporate Headquarters

5791 Van Allen Way

Carlsbad, CA 92008

T: 1 760 603 7200

F: 1 760 602 6500

E: tech_support@invitrogen.com

For country-specific contact information, visit our web site at www.invitrogen.com