

USER GUIDE

invitrogen™
by *life* technologies™

pVAX1™

Catalog no. V260-20

Revision date 2 March 2012

Publication Part number 25-0256

MAN0000094

**For Research Use Only. Not intended for any animal or human
therapeutic or diagnostic use.**

life
technologies™

Contents

Kit Contents and Storage.....	iv
Introduction.....	1
Product Overview	1
Methods.....	2
Cloning into pVAX1™	2
Transient Transfection	5
Appendix	6
Map of pVAX1™ Vector	6
Features of pVAX1™ Vector	7
Map of pVAX1™ /lacZ	8
Accessory Products	9
Technical Support.....	10
Purchaser Notification	11
References.....	12

Kit Contents and Storage

Shipping and Storage

pVAX1™ and pVAX1™/*lacZ* are shipped at room temperature. Upon receipt, store vectors at -20°C.

Kit Contents

20 µg each of pVAX1™ and pVAX1™/*lacZ*, are supplied at 0.5 µg/µL in 10 mM Tris-HCl, 1 mM EDTA, pH 8.0 in a total volume of 40 µL.

For research use only. Not intended for any human or animal therapeutic or diagnostic use.

Introduction

Product Overview

Description of the System

pVAX1™ is a 3.0 kb plasmid vector designed for use in the development of DNA vaccines. The vector was constructed to be consistent with the Food and Drug Administration (FDA) document, “Points to Consider on Plasmid DNA Vaccines for Preventive Infectious Disease Indications”, published December 22, 1996 (see FDA “Points to Consider” below). Features of the vector allow high-copy number replication in *E. coli* and high-level transient expression of the protein of interest in most mammalian cells (see page 6–7 for more information). The vector contains the following elements:

- Human cytomegalovirus immediate-early (CMV) promoter* for high-level expression in a wide range of mammalian cells
- Bovine growth hormone (BGH) polyadenylation signal for efficient transcription termination and polyadenylation of mRNA
- Kanamycin resistance gene for selection in *E. coli*

The control plasmid, pVAX1™/lacZ, is included for use as a positive control for transfection and expression in the cell line of choice.

Experimental Outline

Use the following outline to clone and express your gene of interest in pVAX1™.

1. Consult the multiple cloning site described on page 3 to design a strategy to clone your gene into pVAX1™.
 2. Ligate your insert into the appropriate vector and transform into *E. coli*. Select transformants on LB plates containing 50 µg/mL kanamycin.
 3. Analyze your transformants for the presence of insert by restriction digestion.
 4. Select a transformant with the correct restriction pattern and use sequencing to confirm that your gene is cloned in the proper orientation.
 5. Transfect your construct into the mammalian cell line of choice and test for transient expression of your protein of interest.
-

FDA “Points to Consider”

pVAX1™ was constructed by modifying the vector, pcDNA™3.1, to accommodate the following considerations put forth by the FDA Center for Biologics Evaluation and Research (CBER) in its document, “Points to Consider on Plasmid DNA Vaccines for Preventive Infectious Diseases Indications” (Docket no. 96N-0400).

- Sequences not necessary for replication in *E. coli* or for expression of recombinant protein in mammalian cells were removed to limit DNA sequences with possible homology to the human genome and to minimize the possibility of chromosomal integration.
 - The kanamycin resistance gene was substituted for the ampicillin resistance gene because aminoglycoside antibiotics are less likely to elicit allergic responses in humans.
-

Methods

Cloning into pVAX1™

Introduction

A diagram is provided on the next page to help you design a cloning strategy for ligating your gene of interest into pVAX1™. General considerations for transformation are listed below.

General Molecular Biology Techniques

For help with DNA ligation, *E. coli* transformation, restriction enzyme analysis, DNA sequencing, and DNA biochemistry, refer to *Molecular Cloning: A Laboratory Manual* (Sambrook *et al.*, 1989) or *Current Protocols in Molecular Biology* (Ausubel *et al.*, 1994) (See **References**, page 12).

E. coli strain for Transformation

Many *E. coli* strains are suitable for the propagation of this vector, including TOP10, DH5α™, and DH10. We recommend that you propagate vectors containing inserts in *E. coli* strains that are recombination deficient (*recA*¹) and endonuclease A deficient (*endA*).

For your convenience, TOP10 is available as chemically competent or electrocompetent cells for purchase (see page 9 for ordering information).

Transformation Method

You may use any method of your choice for transformation. Chemical transformation is the most convenient for most researchers. Electroporation is the most efficient and the method of choice for large plasmids.

Propagating pVAX1™

To propagate and maintain the pVAX1™ vector, use the supplied 0.5 µg/µL stock solution in TE, pH 8.0 to transform a *recA*¹, *endA* *E. coli* strain like TOP10, DH5α™, or equivalent. Select transformants on LB plates containing 50 µg/mL kanamycin. Be sure to prepare a glycerol stock of your plasmid-containing *E. coli* strain for long-term storage (see page 4).

Continued on next page

Cloning into pVAX1™, Continued

Cloning Considerations


pVAX1™ is a nonfusion vector. Your insert must contain a Kozak translation initiation sequence and an initiation codon (ATG) for proper initiation of translation (Kozak, 1987; Kozak, 1991; Kozak, 1990). An example of a Kozak consensus sequence is provided below. Note that other sequences are possible (see references above), but the A at position -3 and the G at position +4 are the most critical (shown in bold). The ATG initiation codon is shown underlined.

ANNATGG

Your insert must also contain a stop codon for proper termination of your gene. Note that the *Xba* I site contains an internal stop codon (TCTAGA).

Multiple Cloning Site of pVAX1™

Below is the multiple cloning site for pVAX1™. Restriction sites are labeled to indicate the cleavage site. The multiple cloning site has been confirmed by sequencing and functional testing. The vector sequence may be downloaded from www.lifetechnologies.com/support or from Technical Support (page 10).



*Note that there are two *Pme* I sites and two *BstX* I sites in the polylinker.

Transformation of Ligation Mixture

Transform your ligation mixture into a competent *recA*¹, *endA* *E. coli* strain (e.g. TOP10, DH5α™) and select on LB plates containing 50 µg/mL kanamycin. Select 10–20 clones and analyze for the presence and orientation of your insert.



We recommend that you sequence your construct with the T7 Forward and BGH Reverse primers to confirm that your gene is cloned in the proper orientation for expression and that it contains an ATG and a stop codon. See the diagram above for the sequences and location of the priming sites. For your convenience, we offer the T7 Promoter Primer and BGH Reverse Primer for purchase (see page 9 for ordering) as well as a custom primer synthesis service.

Continued on next page

Cloning into pVAX1™, Continued

Preparing a Glycerol Stock

Once you have identified the correct clone, purify the colony and make a glycerol stock for long-term storage. It is also a good idea to keep a DNA stock of your plasmid at -20°C .

1. Streak the original colony out for single colonies on an LB plate containing $50\ \mu\text{g}/\text{mL}$ kanamycin. Incubate the plate at 37°C overnight.
 2. Isolate a single colony and inoculate into 1–2 mL of LB containing $50\ \mu\text{g}/\text{mL}$ kanamycin.
 3. Grow the culture to mid-log phase ($\text{OD}_{600} = 0.5\text{--}0.7$).
 4. Mix 0.85 mL of culture with 0.15 mL of sterile glycerol and transfer to a cryovial.
 5. Store at -80°C .
-

Transient Transfection

Introduction

Once you have verified that your gene is cloned in the correct orientation and contains an initiation ATG and a stop codon, then you are ready to transiently transfect your mammalian cell line of choice to check for protein expression. We recommend that you include the positive control vector and a mock transfection (negative control) to evaluate your results.

Plasmid Preparation

Plasmid DNA for transfection into eukaryotic cells must be very clean and free from phenol and sodium chloride. Contaminants will kill the cells and salt will interfere with lipids, decreasing transfection efficiency. We recommend isolating plasmid DNA using the PureLink[®] HiPure Miniprep Kit or the PureLink[®] HiPure Midiprep Kit (page 9 for ordering information) or CsCl gradient centrifugation.

Methods of Transfection

For established cell lines (e.g. HeLa), consult original references or the supplier of your cell line for the optimal method of transfection. It is recommended that you follow exactly the protocol for your cell line. Pay particular attention to medium requirements, when to pass the cells, and at what dilution to split the cells. Further information is provided in *Current Protocols in Molecular Biology* (Ausubel et al., 1994).

Methods for transfection include calcium phosphate (Chen and Okayama, 1987; Wigler et al., 1977), lipid-mediated (Felgner et al., 1987; Felgner et al., 1989), and electroporation (Chu et al., 1987; Shigekawa and Dower, 1988). For high efficiency transfection in a broad range of mammalian cells, use Lipofectamine[®] 2000 Reagent available for purchase (page 9). For more information on Lipofectamine[®] 2000 and other transfection reagents, visit www.lifetechnologies.com/support or contact Technical Support (see page 10).

Positive Control

pVAX1[™]/*lacZ* is provided as a positive control vector for mammalian transfection and expression (see page 8). It may be used to optimize transfection conditions for your cell line. The gene encoding β -galactosidase is expressed in mammalian cells under the control of the CMV promoter. A successful transfection will result in β -galactosidase expression that can be easily assayed (see below).

Assay for β -galactosidase Activity


You may assay for β -galactosidase expression by activity assay using cell-free lysates (Miller, 1972) or by staining the cells for activity. We offer the β -Gal Assay Kit and the β -Gal Staining Kit for fast and easy detection of β -galactosidase expression (see page 9 for ordering).

Appendix

Map of pVAX1™ Vector

Map of pVAX1™

The figure below summarizes the features of the pVAX1™ vector. The sequence for pVAX1™ is available for downloading from www.lifetechnologies.com or from Technical Support (see page 10).



Comments for pVAX1™:
2999 bp

CMV promoter: bases 33-620
T7 promoter/priming site: bases 664-683
Multiple cloning site: bases 696-811
BGH reverse priming site: bases 823-840
BGH polyadenylation signal: bases 829-1053
Kanamycin resistance gene: bases 1226-2020
pUC origin: bases 2320-2993

Continued on next page

Features of pVAX1™ Vector

Features of pVAX1™

pVAX1™ (2999 bp) contains the following elements. All features have been functionally tested.

Feature	Benefit
Human cytomegalovirus (CMV) immediate-early promoter/enhancer	Permits efficient, high-level expression of your recombinant protein (Andersson <i>et al.</i> , 1989; Boshart <i>et al.</i> , 1985; Nelson <i>et al.</i> , 1987)
T7 promoter/priming site	Allows for <i>in vitro</i> transcription in the sense orientation and sequencing through the insert
Multiple cloning site	Allows insertion of your gene and facilitates cloning
BGH reverse priming site	Permits sequencing through the insert
Bovine growth hormone (BGH) polyadenylation signal	Efficient transcription termination and polyadenylation of mRNA (Goodwin and Rottman, 1992)
Kanamycin resistance gene	Selection of vector in <i>E. coli</i> (Davies and Smith, 1978)
pUC origin (pUC-derived)	High-copy number replication and growth in <i>E. coli</i> (Bolivar <i>et al.</i> , 1977; Bolivar <i>et al.</i> , 1977)


Map of pVAX1™//lacZ

Description

pVAX1™//lacZ is a 6050 bp control vector containing the gene for β-galactosidase. The vector was constructed by cloning a 3.1 kb *Pst* I-*Xba* I fragment containing the *lacZ* gene into the *Pst* I-*Xba* I site of pVAX1™.

Map of Control Vector

The figure below summarizes the features of the pVAX1™//lacZ vector. The nucleotide sequence for pVAX1™//lacZ is available for downloading from www.lifetechnologies.com or by contacting Technical Support. See page 10 for more information.



Comments for pVAX1™//lacZ:
6050 bp

CMV promoter: bases 33-620

T7 promoter/priming site: bases 664-683

LacZ ORF: bases 773-3829

BGH reverse priming site: bases 3874-3891

BGH polyadenylation signal: bases 3880-4104

Kanamycin resistance gene: bases 4277-5071

pUC origin: bases 5371-6044

Accessory Products

Additional Products

The following additional products may be used with the pVAX1™ vectors. For more information, visit www.lifetechnologies.com or contact Technical Support (see page 10).

Item	Quantity	Cat. no.
One Shot® TOP10 Chemically Competent <i>E. coli</i>	10 reactions	C4040-10
One Shot® TOP10 Electrocomp™ <i>E. coli</i>	20 reactions	C4040-50
TOP10 Electrocomp™ <i>E. coli</i>	20 reactions	C664-55
MAX Efficiency® DH10B™ Competent Cells	5 × 0.2 mL	18297-010
T7 promoter primer	2 µg	N560-02
BGH reverse primer	2 µg	N575-02
PureLink® HiPure Plasmid Miniprep Kit	100 preps	K2100-03
PureLink® HiPure Plasmid Midiprep Kit	25 preps	K2100-04
Lipofectamine® 2000 Reagent	1.5 mL	11668-019
β-Gal Assay Kit	1 kit	K1455-01
β-Gal Staining Kit	1 kit	K1465-01

Technical Support

Obtaining support For the latest services and support information for all locations, go to www.lifetechnologies.com/support.

At the website, you can:

- Access worldwide telephone and fax numbers to contact Technical Support and Sales facilities
- Search through frequently asked questions (FAQs)
- Submit a question directly to Technical Support (techsupport@lifetech.com)
- Search for user documents, SDSs, vector maps and sequences, application notes, formulations, handbooks, certificates of analysis, citations, and other product support documents
- Obtain information about customer training
- Download software updates and patches

Safety Data Sheets (SDS)

Safety Data Sheets (SDSs) are available at www.lifetechnologies.com/support.

Certificate of Analysis

The Certificate of Analysis provides detailed quality control and product qualification information for each product. Certificates of Analysis are available on our website. Go to www.lifetechnologies.com/support and search for the Certificate of Analysis by product lot number, which is printed on the box.

Limited warranty

Life Technologies Corporation is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about a Life Technologies product or service, contact our Technical Support Representatives. All Life Technologies products are warranted to perform according to specifications stated on the certificate of analysis. The Company will replace, free of charge, any product that does not meet those specifications. This warranty limits the Company's liability to only the price of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. The Company reserves the right to select the method(s) used to analyze a product unless the Company agrees to a specified method in writing prior to acceptance of the order.

Life Technologies makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore the Company makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, report it to our Technical Support Representatives.

Life Technologies Corporation shall have no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.

Purchaser Notification

**Limited Use Label
License No:
Research Use
Only**

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

References

- Andersson, S., Davis, D. L., Dahlbäck, H., Jörnvall, H., and Russell, D. W. (1989). Cloning, Structure, and Expression of the Mitochondrial Cytochrome P-450 Sterol 26-Hydroxylase, a Bile Acid Biosynthetic Enzyme. *J. Biol. Chem.* *264*, 8222-8229.
- Ausubel, F. M., Brent, R., Kingston, R. E., Moore, D. D., Seidman, J. G., Smith, J. A., and Struhl, K. (1994). *Current Protocols in Molecular Biology* (New York: Greene Publishing Associates and Wiley-Interscience).
- Bolivar, F., Rodriguez, R. L., Betlach, M. C., and Boyer, H. W. (1977). Construction and Characterization of New Cloning Vehicles. I. Ampicillin-Resistant Derivatives of the Plasmid pMB9. *Gene* *2*, 75-93.
- Bolivar, F., Rodriguez, R. L., Greene, P. J., Betlach, M. C., Heyneker, H. L., and Boyer, H. W. (1977). Construction and Characterization of New Cloning Vehicles. II. A Multipurpose Cloning System. *Gene* *2*, 95-113.
- Boshart, M., Weber, F., Jahn, G., Dorsch-Häsler, K., Fleckenstein, B., and Schaffner, W. (1985). A Very Strong Enhancer is Located Upstream of an Immediate Early Gene of Human Cytomegalovirus. *Cell* *41*, 521-530.
- Chen, C., and Okayama, H. (1987). High-Efficiency Transformation of Mammalian Cells by Plasmid DNA. *Mol. Cell. Biol.* *7*, 2745-2752.
- Chu, G., Hayakawa, H., and Berg, P. (1987). Electroporation for the Efficient Transfection of Mammalian Cells with DNA. *Nuc. Acids Res.* *15*, 1311-1326.
- Davies, J., and Smith, D. I. (1978). Plasmid-Determined Resistance to Antimicrobial Agents. *Ann. Rev. Microbiol.* *32*, 469-518.
- Felgner, P. L., Gadek, T. R., Holm, M., Roman, R., Chan, H. W., Wenz, M., Northrop, J. P., Ringold, G. M., and Danielsen, M. (1987). Lipofectin: A Highly Efficient, Lipid-mediated DNA-transfection Procedure. *Proc. Natl. Acad. Sci. USA* *84*, 7413-7417.
- Felgner, P. L., Holm, M., and Chan, H. (1989). Cationic Liposome Mediated Transfection. *Proc. West. Pharmacol. Soc.* *32*, 115-121.
- Goodwin, E. C., and Rottman, F. M. (1992). The 3'-Flanking Sequence of the Bovine Growth Hormone Gene Contains Novel Elements Required for Efficient and Accurate Polyadenylation. *J. Biol. Chem.* *267*, 16330-16334.
- Kozak, M. (1987). An Analysis of 5'-Noncoding Sequences from 699 Vertebrate Messenger RNAs. *Nuc. Acids Res.* *15*, 8125-8148.
- Kozak, M. (1991). An Analysis of Vertebrate mRNA Sequences: Intimations of Translational Control. *J. Cell Biol.* *115*, 887-903.
- Kozak, M. (1990). Downstream Secondary Structure Facilitates Recognition of Initiator Codons by Eukaryotic Ribosomes. *Proc. Natl. Acad. Sci. USA* *87*, 8301-8305.
- Miller, J. H. (1972). *Experiments in Molecular Genetics* (Cold Spring Harbor, New York: Cold Spring Harbor Laboratory).

Continued on next page

References, Continued

- Nelson, J. A., Reynolds-Kohler, C., and Smith, B. A. (1987). Negative and Positive Regulation by a Short Segment in the 5'-Flanking Region of the Human Cytomegalovirus Major Immediate-Early Gene. *Mol. Cell. Biol.* 7, 4125-4129.
- Sambrook, J., Fritsch, E. F., and Maniatis, T. (1989). *Molecular Cloning: A Laboratory Manual*, Second Edition (Plainview, New York: Cold Spring Harbor Laboratory Press).
- Shigekawa, K., and Dower, W. J. (1988). Electroporation of Eukaryotes and Prokaryotes: A General Approach to the Introduction of Macromolecules into Cells. *BioTechniques* 6, 742-751.
- Wigler, M., Silverstein, S., Lee, L.-S., Pellicer, A., Cheng, Y.-C., and Axel, R. (1977). Transfer of Purified Herpes Virus Thymidine Kinase Gene to Cultured Mouse Cells. *Cell* 11, 223-232.
-

©2012 Life Technologies Corporation. All rights reserved.

The trademarks mentioned herein are the properties of Life Technologies Corporation or their respective owners.

Notes

