

BacMam ERK2 [pThr185/pTyr187] Cellular Assay

Catalog Number: A12888

Literature Part Number: A12888.PIS

Literature Lot Number: V1

Revision date: 10 January 2011

FAST FACTS

For first-time BacMam users, we recommend using cells like U-2 OS (ATCC®, HTB-96™), which can be transduced exceptionally well.


BacMam Enhancer Solution (PV5835) is **not** required for this assay using U-2 OS cells. However, it is **highly recommended when using CHO cells** (see quick protocol below).

Optimal Virus %: We recommend that you perform a titration of the BacMam ERK2 Reagent to determine the optimal percentage of virus for your transduction in your cell background of interest or when you receive a new lot of virus. Select the lowest percentage of BacMam Reagent that yields the largest assay window (response ratio). See example below.

Kit Components	SKU#	Amount	Storage	Handling
BacMam ERK2 Reagent	A12889	2 × 25 mL	4°C	<ul style="list-style-type: none"> Do not freeze Minimize exposure to ambient light Use sterile technique Aliquot to minimize handling,
LanthaScreen® Tb-anti-ERK [pThr185/pTyr187] Antibody	PV5833*	10 µg	-20°C	Aliquot if necessary to avoid multiple freeze/thaw cycles
BacMam Enhancer Solution (1000X)	PV5835	150 µL	-20°C	Aliquot if necessary to avoid multiple freeze/thaw cycles
6X LanthaScreen® Cellular Assay Lysis Buffer	A12891	6 mL	4°C	Supplement with inhibitors and antibody

*To order additional antibody, use catalog no. PV5269 (25 µg amount).

Titration of BacMam ERK2 Reagent in GeneBLazer® 5HT1A-Gα15-NFAT-bla CHO-K1 cells (Invitrogen catalog no. K1712)
(Detection of ERK2 phosphorylation induced by 5-HT)



Additional Materials Required, but not provided	Source	Part #
Positive Control GPCR agonist We recommend: 5-HT (if using CHO cells expressing 5-HT1A)	Sigma	H9523
Positive Control Growth Factor We recommend: EGF (if using U-2 OS)	Invitrogen	PHG0314
Cell Line of Interest	Various	Various
DMSO	Fluka	41647
Protease Inhibitor	Sigma	P8340
Phosphatase Inhibitor	Sigma	P0044
White tissue culture-treated, 384-well assay plates	Corning	3570
Fluorescence plate reader with top-read and TR-FRET capability	www.invitrogen.com/instrumentsetup for details	
Optional: Clear-bottom, tissue-culture treated, 384-well plates	Corning	3712

Detailed Protocols and Additional Assay Performance Data Available

Please visit www.invitrogen.com and search for A12888 to download the full protocol and validation packet for this assay.

Protocols and validation packets are located under the "How to Use" tab on the product page.

Validation Packets include assay performance under variable experimental conditions.

Technical Support

For additional assistance in running this BacMam-enabled Cellular Assay, please contact our technical support team at drugdiscoverytech@invitrogen.com or 760-603-7200 (enter 3 for “know your party’s extension”, then enter 40266).

Quick Assay Protocol for GPCR-mediated ERK2 phosphorylation in CHO cells

The following protocol is developed for a CHO cell line over-expressing the serotonin receptor HTR1A. Conditions such as plating density, stimulant, or stimulation time may need to be optimized for CHO cell lines expressing different GPCRs. Phosphorylation of ERK at [Thr185/Tyr187] can be induced in CHO cells with either 10% serum or 100 nM PMA.

		Cell-Free Control Wells	Unstimulated Control Wells	Stimulated Control Wells	Test Compound Wells
BacMam Transduction	Step 1 Grow and Transduce Cells	<ul style="list-style-type: none"> Grow cells in Growth Medium* to 90–100% confluence (~1–2 × 10⁵ cells/cm²). Perform five 3-fold serial dilutions of BacMam reagent in PBS with Ca²⁺ and Mg²⁺. Remove media from cells and add the serial diluted BacMam virus. Incubate at room temperature protected from light for 3–4 hrs. 			
	Step 2 Add Enhancer and Incubate Cells	Remove virus and add Growth Medium* plus 1X BacMam Enhancer Solution. Incubate the plate at 37°C/5% CO ₂ for 20–24 hours.			
LanthaScreen® Assay	Step 3 Harvest Cells	Harvest cells, wash once with Assay Medium** and resuspend in Assay Medium at 0.75 × 10 ⁶ cells/mL (you need to determine the optimal cell density for your cell line).			
	Step 4 Plate Cells	20 µL/well Assay Media only	20 µL transduced cells/well (about 15,000 cells/well), quick spin of the plate		
	Step 5 Serum-starve and Incubate Cells	Incubate the plate at 37°C/5% CO ₂ for 4 hours			
	Step 6 Prepare Complete 6X Lysis Buffer	To 1 mL 6X Lysis Buffer, add 30 µL of 100x protease inhibitor, 30 µL 100x phosphatase inhibitor, and Tb-anti-ERK2 [pThr185/pTyr187] Antibody to 12 nM.			
	Step 7 Add Agonist	10 µL/well of 0.3% DMSO in Assay Media	10 µL/well of 3X 5-HT (900 nM) in Assay Media	10 µL/well of 3X Test Compound in Assay Media	
	Step 8 Stimulate Cells	Incubate the plate at 37°C/5% CO ₂ for 6–8 minutes			
	Step 8 Add Lysis Buffer (including Tb-Ab)	Add 6 µL/ well of Complete 6X Lysis Buffer to each well			
	Step 9 Cell Lysis/ Assay Equilibration	Incubate plate for ~3 hours at room temperature in the dark			
Step 10 Read Plate and Analyze Data	See LanthaScreen® Detection in the online protocol—Excitation filter: 337 nm (30 nm bandwidth); Emission filters: 490 nm (10 nm bandwidth) and 520 nm (25 nm bandwidth)				

* **Growth Media:** D-MEM Media (Invitrogen 10569) with 10% dFBS, 10 mM HEPES, 0.1 mM NEAA, and 100 U/mL Penicillin/100 µg/mL Streptomycin

** **Assay Media:** Opti-MEM® I (Invitrogen 11058) with 0.1% cdFBS (or dialyzed FBS), 0.1 mM NEAA, 1 mM Sodium Pyruvate, and 100 U/mL Penicillin/100 µg/mL Streptomycin

Quick Assay Protocol for Growth Factor-mediated ERK2 phosphorylation in U-2 OS cells

The cell harvesting and plating densities, growth medium, and assay medium must be optimized for your particular cell line(s). The following protocol is developed for U-2 OS cells. It may be applied to many other cell types such as HEK293, HeLa and A549. Conditions may need to be optimized for different cell types.

		Cell-Free Control Wells	Unstimulated Control Wells	Stimulated Control Wells	Test Compound Wells
BacMam Transduction	Step 1 Grow, Harvest Cells and Transduce	<ul style="list-style-type: none"> Grow cells in Growth Medium* to 90-100% confluence ($\sim 0.8\text{--}1.2 \times 10^5$ cells/cm²). Harvest and wash cells and resuspend in Assay Medium** at 8.4×10^5 cells/mL. Perform five 3-fold serial dilutions of BacMam reagent in Assay Medium. Add 0.4 mL undiluted or serial diluted BacMam reagent to 1 mL cells to generate a virus titration range of 0.1% to 30% (v/v) final concentration. 			
	Step 2 Plate Cells/virus mixture	Add 20 μ L/well Assay Media only	Add 20 μ L cells and BacMam mixture per well ($\sim 12,000$ cells/well), and quick spin the plate		
	Step 3 Incubate Cells	Incubate the plate at 37°C/5% CO ₂ for 20–24 hours			
LanthaScreen® Assay	Step 4 Prepare Complete 6X Lysis Buffer	To 1 mL 6X Lysis Buffer, add 30 μ L of 100x protease inhibitor, 30 μ L 100x phosphatase inhibitor, and Tb-anti-ERK2 [pThr185/pTyr187] Antibody to 12 nM			
	Step 5 Add media or ligand	Add 10 μ L/well of Assay Media with 0.3% DMSO	Add 10 μ L/well of 3X Agonist in Assay Media with 0.3% DMSO (for EGF, use ~ 600 ng/mL)	Add 10 μ L/well of 3X Test Compound in Assay Media (0.3% DMSO)	
	Step 6 Stimulate Cells	Incubate the plate at 37°C/5% CO ₂ for 6–8 minutes			
	Step 7 Add Lysis Buffer (including Tb-Ab)	Add 6 μ L/well of Complete 6X Lysis Buffer to each well			
	Step 8 Cell Lysis/Assay Equilibration	Incubate plate for ~ 3 hours at room temperature in the dark			
	Step 9 Read Plate and Analyze Data	See LanthaScreen® Detection in the online protocol—Excitation filter: 337 nm (30 nm bandwidth); Emission filters: 490 nm (10 nm bandwidth) and 520 nm (25 nm bandwidth)			

* **Growth Media for U-2 OS Cells:** McCoy's 5A Media with 10% dFBS, 10 mM HEPES, 0.1 mM NEAA, 1 mM Sodium Pyruvate, and 100 U/mL Penicillin/ 100 μ g/mL Streptomycin

** **Assay Media:** Opti-MEM® I (Invitrogen 11058) with 0.1% cdFBS (or dialyzed FBS), 0.1 mM NEAA, 1 mM Sodium Pyruvate, and 100 U/mL Penicillin/ 100 μ g/mL Streptomycin

We would like your feedback

Please let us know how you like this Product Information Sheet.

Visit www.invitrogen.com/feedbacksurvey to take a quick 5-minute survey about your experience and get a free gift upon completion.*

All customers who complete this survey will receive a free gift.

*Available while supplies last. Offer void where prohibited by federal, state, or local laws or regulation or agency/institutional policy.

The performance of this product is guaranteed for six months from the date of purchase if stored and handled properly.

The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners. ATCC and HTB-96 are trademarks of ATCC.

©2011 Life Technologies. All rights reserved.

Purchaser Notification

Limited Use Label License No. 21: Bac-to-Bac® and Bac-to-Bac® HT

This product is sold under patent license from Monsanto for research purposes only and no license for commercial use is included. Requests for licenses for commercial manufacture or use should be directed to Director, Monsanto Corporate Research, 800 N. Lindbergh, St. Louis, Missouri 63167.

Limited Use Label License No. 127: GFP with Heterologous Promoter

This product is sold under license from Columbia University. Rights to use this product are limited to research use only, and expressly exclude the right to manufacture, use, sell or lease this product for use for measuring the level of toxicity for chemical agents and environmental samples in cells and transgenic animals. No other rights are conveyed. Not for human use or use in diagnostic or therapeutic procedures. Inquiry into the availability of a license to broader rights or the use of this product for commercial purposes should be directed to Columbia Innovation Enterprise, Columbia University, Engineering Terrace-Suite 363, New York, New York 10027.

Limited Use Label License No. 168: Lanthanide Chelates

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use, or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, California 92008. Phone (760) 603-7200. Fax (760) 602-6500.

Limited Use Label License No. 198: Fluorescent Proteins and Stable Cell Lines Expressing Such Proteins

The purchase of this product conveys to the buyer the nontransferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for profit entity). No rights are conveyed to modify or clone the gene encoding GFP contained in this product. The buyer cannot sell or otherwise transfer (a) this product, (b) its components, or (c) materials made by the employment of this product or its components to a third party or otherwise use this product or its components or materials made by the employment of this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the employment of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Life Technologies Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that none of this product, or any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Life Technologies Corporation is willing to accept return of the product with a full refund. For information on purchasing a license to use this product for purposes other than those permitted above, contact Licensing Department, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, California 92008. Phone (760) 603-7200 or outlicensing@lifetech.com.

Limited Use Label License No. 272: Humanized GFP

This product is sold for research use only. Not for therapeutic or diagnostic use in humans.

Limited Use Label License No. 306: Baculovirus Vectors

This product is for research use only by those researchers in laboratories of academic, government, industrial and/or clinical institutions engaged in the investigation of biological or biochemical processes, or research and development of biological products. This product is not to be used in the manufacture, use or sale of human or animal diagnostic, therapeutic or prophylactic products.

Limited Use Label License No. 308: WPRE Element

This product contains the Woodchuck Post-transcriptional Regulatory Element ("WPRE") which is the subject of intellectual property owned by The Salk Institute for Biological Studies, and licensed to Life Technologies Corporation. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; and/or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. In addition, any use of WPRE outside of this product or the product's authorized use requires a separate license from the Salk Institute. Life Technologies will not assert a claim against the buyer of infringement of patents owned by Life Technologies and claiming this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product or for a Commercial Purpose. If the purchaser is not willing to accept the limitations of this limited use statement, Life Technologies is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, California 92008, Phone (760) 603-7200. Fax (760) 602-6500, or The Salk Institute for Biological Studies, 10010 North Torrey Pines Road, La Jolla, CA 92037, Attn.: Office of Technology Management, Phone: (858) 453-4100 extension 1275, Fax: (858) 546-8093.

Limited Use Label License No. 332: BacMam Virus Use

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer solely in accordance with the accompanying product literature or manual. Purchase of this product does not convey a license to expand, amplify, or otherwise propagate the provided viral particles or to otherwise modify or alter the virus by any means.

Limited Use Label License No. 358: Research Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.