

Jump-In™ GripTite™ HEK293 Retargeting Kit

Catalog Number: A14150

Literature Part Number: A14150PIS (MAN0005284)

Literature Lot Number: V01.00

Revision date: 5 October 2011

Kit components	Amount/Composition	Storage
Jump-In™ GripTite™ HEK293 Cells	2 vials (~3 × 10 ⁶ cells/vial in Freezing Medium*)	Liquid nitrogen
pJTI™ R4 Int (integrase vector)	200 µL at 0.5 µg/µL in TE buffer, pH 8.0**	-20°C
pJTI™ R4 DEST CMV pA (destination vector)	200 µL at 0.5 µg/µL in TE buffer, pH 8.0	-20°C

* Recovery™ Cell Culture Freezing Medium; **TE buffer, pH 8.0: 10 mM Tris-HCl, 1 mM EDTA, pH 8.0

FAST FACTS

- For additional materials required but not provided, visit www.lifetechnologies.com and search for A14150 to download the full detailed protocol.
- For first-time users, we recommend following the detailed protocol available online.
- After the initial thaw and passage, Jump-In™ GripTite™ HEK293 cells usually double in about 24 hours.
- **We highly recommend that** you include all the suggested negative and positive controls in your transfection (retargeting) experiments. The controls provide a good indicator for the principle success of the retargeting reaction and are crucial for any troubleshooting.
- **Positive control:** pJTI™ R4 EXP CMV EmGFP pA (positive control vector, Cat. no. A14146) + pJTI™ R4 Int (integrase vector, provided in this kit).

Experiment Outline

The table below describes the major steps required for retargeting the Jump-In™ GripTite™ HEK293 cell line.

Table 1 Retargeting experiment workflow

Step	Action
1	Thaw and expand the Jump-In™ GripTite™ HEK293 cells
2	Create an entry clone by cloning your gene of interest into a Gateway® entry vector
3	Generate a retargeting construct by performing an LR recombination reaction between the entry clone and pJTI™ R4 DEST CMV pA (i.e., the destination vector)
4	Cotransfect your retargeting construct and the integrase vector into the Jump-In™ GripTite™ HEK293 cells
5	Select for retargeted Jump-In™ GripTite™ HEK293 cells in Selection Medium containing Blasticidin
6	Confirm retargeting of the Jump-In™ GripTite™ HEK293 cells by PCR
7	Characterize the retargeted clones



Note that this product information sheet offers instructions and guidelines for thawing and propagating Jump-In™ GripTite™ HEK293 cells, and provides only an overview of retargeting experiments. For detailed instructions for creating a retargeting construct, transfecting (retargeting) the Jump-In™ GripTite™ HEK293 cell line, and selecting and characterizing the retargeted clones, refer to the full detailed protocol provided online.

Technical Support

For assistance, contact our technical support team at drugdiscoverytech@lifetech.com or 760-603-7200 (enter 3 for “know your party’s extension”, then enter 40266).

Jump-In™ GripTite™ HEK293 Cell Culture

Table 2 Media used in culturing Jump-In™ GripTite™ HEK293 cells

Component	Thawing Medium	Growth Medium	Retargeting Selection Medium	Catalog number
D-MEM with GlutaMAX™-I (high glucose)	90%	90%	90%	10569-010
Dialyzed FBS (<i>Do not substitute!</i>)	10%	10%	10%	26400-036
MEM Non-Essential Amino Acids Solution	0.1 mM	0.1 mM	0.1 mM	11140-050
HEPES Buffer (pH 7.3)	25 mM	25 mM	25 mM	15630-080
Penicillin(antibiotic)	100 U/mL	100 U/mL	100 U/mL	15140-122
Streptomycin (antibiotic)	100 µg/mL	100 µg/mL	100 µg/mL	15140-122
Geneticin®	—	600 µg/mL	600 µg/mL	11811-031
Hygromycin B	—	200 µg/mL	—	10687-010
Blasticidin	—	—	10 µg/mL	A11139-02

Thawing Jump-In™ GripTite™ HEK293 Cells

1. Rapidly thaw the cells with gentle agitation in a 37°C water bath.
2. Exchange media by transferring the thawed cells into 10 mL of Thawing Medium in a sterile 15-mL tube, centrifuge at 200 × g for 5 minutes, and resuspend the cells in 1 mL of fresh Thawing Medium.
3. Transfer the cells to a T-75 tissue culture flask containing 20 mL of pre-equilibrated Thawing Medium and place the flask in a humidified 37°C/5% CO₂ incubator.
4. At first passage, switch to Growth Medium.

Propagating Jump-In™ GripTite™ HEK293 Cells

1. Aspirate medium from growing cells, rinse once in PBS, and add the appropriate amount of 0.05% Trypsin/EDTA (3 mL for a 100-mm dish, 5 mL for a T-75 flask, 10 mL for a T-175 flask).
2. Add an equal volume of Growth Medium to inactivate the 0.05% Trypsin/EDTA.
3. Verify under a microscope that cells have detached and clumps have completely dispersed.
4. Determine the viable cell number using a hemocytometer or a cell counter. Cell number and viability can be quickly and conveniently determined using the Countess® Automated Cell Counter. We recommend determining cell health frequently to ensure optimal performance in retargeting experiments.
5. Centrifuge the cells at 200 × g for 5 minutes and resuspend then in Growth Medium.
6. Seed fresh culture vessel containing pre-warmed Growth Medium at the appropriate cell density. We recommend a split ratio of 1:3 to 1:10.

Important: Do not allow the cells to reach confluence.

Purchaser Notification

Limited Use Label License: Research Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

Limited Use Label License No. 19: Gateway® Cloning Products

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for profit entity). The purchase of this product does not convey a license under any method claims in the foregoing patents or patent applications, or to use this product with any recombination sites other than those purchased from Life Technologies Corporation or its authorized distributor. The right to use methods claimed in the foregoing patents or patent applications with this product for research purposes only can only be acquired by the use of Clonase™ purchased from Life Technologies Corporation or its authorized distributors. The buyer cannot modify the recombination sequence(s) contained in this product for any purpose. The buyer cannot sell or otherwise transfer (a) this product, (b) its components, or (c) materials made by the employment of this product or its components to a third party or otherwise use this product or its components or materials made by the employment of this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the employment of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Notwithstanding the preceding, any buyer who is employed in an academic or government institution may transfer materials made with this product to a third party who has a license from Life Technologies under the patents identified above to distribute such materials. Transfer of such materials and/or information to collaborators does not convey rights to practice any methods claimed in the foregoing patents or patent applications. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Life Technologies Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that none of (i) this product, (ii) any of its components, or (iii) a method claim of the foregoing patents, was used in the manufacture of such product. Life Technologies Corporation will not assert a claim against the buyer of infringement of the above patents based upon the use of this product to manufacture a protein for sale, provided that no method claim in the above patents was used in the manufacture of such protein. If the purchaser is not willing to accept the limitations of this limited use statement, Life Technologies is willing to accept return of the product with a full refund. For information on purchasing a license to use this product for purposes other than those permitted above, contact Licensing Department, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, California 92008. Phone (760) 603-7200.

LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

The performance of this product is guaranteed for six months from the date of purchase if stored and handled properly.

The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

©2011 Life Technologies Corporation. All rights reserved. Reproduction forbidden without permission.

For research use only. Not intended for human or animal therapeutic or diagnostic use.