

Table of Contents

Kit Contents and Handling.....	1
Overview.....	2
LanthaScreen® Cellular Assay (Terbium-based TR-FRET detection)	2
BacMam Technology	2
Workflow for the BacMam Histone H3K4me2 Cellular Assay.....	2
Before Starting.....	3
Materials Required but Not Provided.....	3
U-2 OS Cell Culture Reagents (optional, if using U-2 OS cells).....	3
Guidelines for Optimizing BacMam Histone H3K4me2 Cellular Assays.....	3
Assay Protocol.....	5
Quick Reference Protocol Example (U-2 OS Cells): 384-well Assay Plate Format.....	5
Detailed Protocol.....	6
Terbium TR-FRET Detection.....	8
Instruments and Filters	8
Reading the Assay Plate and Data Analysis.....	8
Testing Terbium TR-FRET Detection Using the Instrument Controls.....	8
Appendix.....	9
Alternative Transduction Protocol A (For difficult-to-transduce cells).....	9
Alternative Transduction Protocol B (For HeLa, MCF7, T47D, MDA-MB-231 and others)	11
Troubleshooting Guide	12
Purchaser Notification	13

Kit Contents and Handling

Component	Part no.	Amount	Storage	Handling
LanthaScreen® Tb-anti-Histone H3K4me2 Antibody	A14139	5 µg	-20°C	<ul style="list-style-type: none"> Protect from light Avoid multiple freeze/thaw cycles
LanthaScreen® 6X Cellular Assay Lysis Buffer	A12891	6 mL	4°C	On the day of assay, supplement with protease inhibitor cocktail* and antibody
BacMam Histone H3 Reagent	A12894	25 mL	4°C	<ul style="list-style-type: none"> DO NOT FREEZE Use sterile technique Avoid extended exposure to ambient room light
Instrument Control Terbium TR-FRET Kit				
Low Instrument Control, 1 mL	A14138	1 kit	4°C	Protect from light (do not vortex)
High Instrument Control, 1 mL				

* See **Materials Required but Not Provided**, page 3.

Overview

BacMam Cellular Assays use the BacMam gene delivery system in conjunction with LanthaScreen® Cellular Assays to measure post-translational modifications of a target substrate. The combination of the two technologies provides a fast, convenient, and robust method for interrogating specific signal transduction events in a cell background of choice.

LanthaScreen® Cellular Assay (Terbium-based TR-FRET detection)

LanthaScreen® Cellular Assays are HTS-compatible immunoassays that are used for interrogating target-specific post-translational modifications in a cell-based format. Target proteins are expressed as fusions with the green fluorescent protein (GFP) in living cells, and modification-specific antibodies labeled with Terbium (Tb) are used to detect post-translational modifications of the target protein in a time-resolved fluorescence resonance energy transfer (TR-FRET) format.

The use of GFP as a FRET acceptor circumvents the need for complex antigen-capturing reagents, thereby providing a high-throughput alternative to commonly used analytical methods such as Western blot and ELISA.

For more information, visit www.invitrogen.com/lanthascreeencellular.

BacMam Technology

While GFP-Histone H3 can be delivered to cells via multiple methods (e.g., stable cell line generation, transient transfection, electroporation, retroviral transduction), BacMam technology is a convenient approach that uses a modified baculovirus to efficiently deliver and transiently express genes (in this case, GFP-Histone H3) in mammalian cells. BacMam viruses are non-replicating in mammalian cells, rendering them safe as research reagents.

This technology has several advantages over traditional transient methods for heterologous gene expression, including:

- High transduction efficiency across a broad range of cell types, including primary and stem cells
- Little-to-no observable cytopathic effects
- Reproducible and titratable target gene expression
- Compatibility with simultaneous delivery of multiple genes.

For more information on BacMam, visit www.invitrogen.com/bacmam.

Workflow for the BacMam Histone H3K4me2 Cellular Assay

Day 1: Cells are transduced with BacMam GFP-Histone H3 reagent and plated onto a 384-well assay plate. Cells are left untreated or treated with compound for 20 to 24 hours.

Day 2: The expression of GFP-Histone H3 fusion protein in the nucleus can be visualized by fluorescence microscopy. Cells are lysed in the presence of a Tb-labeled anti-Histone H3K4me2 specific antibody, and the level of K4 di-methylation on GFP-Histone H3 is measured on a TR-FRET-compatible plate reader. Little or no TR-FRET generally indicates little or no modification whereas high TR-FRET indicates high K4me2 level.




Figure 1 Representative Assay Workflow. Cells are mixed with BacMam Reagent encoding the GFP-tagged Histone H3 protein and are plated onto a 384-well assay plate. Cells are left untreated or treated with compound for 20 to 24 hours. Cells are then lysed in the presence of a terbium-anti-Histone H3K4me2 antibody and TR-FRET is detected using a fluorescence microplate reader with standard TR-FRET settings.

For Technical Support for this or other Drug Discovery Products, dial 760-603-7200, option 3, extension 40266.

Before Starting

Materials Required but Not Provided

Materials	Recommended Source	Cat. no.
Cell Line of Interest (refer to www.invitrogen.com/bacmam for more information regarding BacMam compatible cell types)	Various	Various
Assay Medium (commonly used growth media are compatible with the terbium TR-FRET readout; where possible, avoid media containing phenol-red because its presence in the assay can interfere with the TR-FRET signal and lead to reduced assay performance)	Various	Various
Protease Inhibitor Cocktail	Sigma	P8340
Assay Plates (white opaque plates) Tissue culture-treated 384-well assay plates	Corning	3570
Fluorescence plate reader with top-read and TR-FRET capability	www.invitrogen.com/instrumentsetup	
<i>Optional:</i> Clear-bottom, tissue culture-treated 384-well assay plates setup in parallel to the TR-FRET assay plates for visualizing GFP-Histone H3 expression	Corning	3712

U-2 OS Cell Culture Reagents (optional, if using U-2 OS cells)

U-2 OS cells transduce well with BacMam viruses, and we recommend them as a control cell line for detecting H3K4me2.

Media/Reagents	Recommended source	Part no.
U-2 OS cells	ATCC	HTB-96
McCoy's 5A Medium	Invitrogen	16600-108
Fetal Bovine Serum (dialyzed)	Invitrogen	26400-036
Nonessential amino acids (NEAA)	Invitrogen	11140-050
Sodium Pyruvate	Invitrogen	11360-070
Penicillin/Streptomycin (antibiotic)	Invitrogen	15140-122
HEPES Buffer Solution (1 M)	Invitrogen	15630-080
Dulbecco's Phosphate-buffered saline (PBS) without Ca ²⁺ and Mg ²⁺	Invitrogen	14190-136
Trypsin/EDTA	Invitrogen	25300-062

Guidelines for Optimizing BacMam Histone H3K4me2 Cellular Assays

First-time Terbium TR-FRET users

Prior to setting up an assay, we strongly recommend that you check your plate reader setup for terbium-based TR-FRET detection using the Terbium TR-FRET Instrument Controls (provided in this kit) as outlined on page 7. For more information about your specific instrument and to purchase filters, visit www.invitrogen.com/instrumentsetup. For technical assistance, contact Drug Discovery Technical Support at drugdiscoverytech@lifetech.com or call 760-603-7200, option 3, extension 40266).

For Technical Support for this or other Drug Discovery Products, dial 760-603-7200, option 3, extension 40266.

Working with BacMam Histone H3 reagent

- For first time users of BacMam Reagent, we recommend including a control cell line which transduces exceptionally well, such as U-2 OS (ATCC® number: HTB-96) (see protocol on page 5).
- Since assay performance depends upon effective expression of GFP-Histone H3, we strongly recommend that you perform a titration of the BacMam GFP-Histone H3 Reagent to determine the optimal percentage of virus (v/v) for transduction in your cell background of interest. For initial testing of BacMam GFP-Histone H3 in a given cell-type, we recommend that you test ~30%, 20%, 10%, 5%, 2.5%, and 1.25% (v/v) final dilutions. Select the lowest percentage of BacMam Reagent that yields the largest assay window. See **example below**.




Figure 2 Detection of Histone H3K4me2 in U-2 OS cells transduced with various concentrations of BacMam Histone H3 Reagent. Emission Ratios are plotted against the final BacMam reagent concentration. In this case, 10% BacMam Histone H3 generated the highest emission ratio.

- Whereas many cell types can be transduced efficiently using the detailed protocol described here for U-2 OS, some challenging cell types (e.g., CHO, HeLa, MCF7 etc.) may require alternative protocols described in the **Appendix**.
- For best results, use healthy, well-maintained cells when performing BacMam transductions. Note that for many cell-types, cryopreserved cells can also be transduced immediately following thaw.
- To minimize day-to-day variability, be sure to use the same growth conditions (e.g., similar harvest density, similar passage numbers).
- BacMam Enhancer (Cat. no. PV5835, supplied at 1,000X) can improve the transduction efficiency of difficult-to-transduce cell lines (see Alternate Transduction Protocol A, page 9). However, it is not required for transducing U-2 OS cells (see protocol for U-2 OS, page 5). When using the Enhancer, we recommend testing at least two concentrations of it (e.g., 0.5X and 1X) in comparison to a control without Enhancer to identify the best concentration for your cell-type of interest. To minimize the effect of the Enhancer solution on your assay of interest, we recommend removing the Enhancer 20 to 24 hours prior to the LanthaScreen® assay.

Important assay parameters for optimization

- Confluence of cells at harvest for assay set-up may impact results, such as the assay window.
- Cell plating density (i.e., cell number per well in the assay plate)
- Expression level of GFP-fusion protein (i.e., titration of BacMam GFP-Histone H3 reagent)
- Compound concentrations and/or treatment times
- Antibody equilibration time

Application Note

Visit www.lifetechnologies.com and search for catalog number A14164 to download the application note for the BacMam Histone H3K4me2 Cellular Assay. The application note is located under the “How to Use” tab on the product page, and contains more information about H3K4me2 assay optimization and applications.

For Technical Support for this or other Drug Discovery Products, dial 760-603-7200, option 3, extension 40266.

Assay Protocol

In the following protocol, the cells are incubated with the BacMam GFP-Histone H3 Reagent in a 384-well assay plate for 20–24 hours. Two alternative transduction protocols provided in the **Appendix** on page 9 may be used with difficult-to-transduce cell lines (**Protocol A**) or other easy-to-transduce cells such as HeLa, T47D, MDA-MB-231 and MCF7 (**Protocol B**).

The cell harvesting and plating densities, growth medium, and assay medium must be optimized for your particular cell line(s). The following protocol was developed for U-2 OS cells.

Quick Reference Protocol Example (U-2 OS Cells): 384-well Assay Plate Format

This quick reference protocol is designed for experienced users using U-2 OS cells, with testing performed in the presence of various concentrations of BacMam Histone H3 reagent. Conditions may need to be optimized for different cell types. For the detailed protocol, see page 6.

		Non-transduced Wells	Transduced Wells
BacMam Transduction	Step 1 Grow, harvest and plate cells	<ul style="list-style-type: none"> Grow cells in Growth Medium* to 60–95% confluence ($\sim 0.5 \times 10^5$ to 1.0×10^5 cells/cm²). Harvest cells and resuspend in Growth Medium at 3.75×10^5 cells/mL. Plate 20 μL/well of cell suspension (about 7,500 cells/well) onto a 384-well assay plate (and optionally a separate plate with clear-bottom for GFP imaging later) Quick spin the plate at $30 \times g$ for 1 minute (if performing the experiment manually) 	
	Step 2** Add BacMam Reagent	Add 5 μ L/well of Growth Medium	Add 5 μ L/well of BacMam GFP-Histone H3 reagent (undiluted or diluted with growth medium to result in different concentrations of the BacMam)
	Step 3 Incubate Cells/BacMam	<ul style="list-style-type: none"> Quick spin the plate at $30 \times g$ for 1 minute (if performing the experiment manually) Incubate the plate at 37°C/5% CO₂ for 20–24 hours (allows for GFP-Histone H3 expression) 	
LanthaScreen® Histone H3K4me2 Assay	Step 4 (Optional) GFP Imaging	If desired, observe and image GFP-Histone H3 expression under a fluorescence microscope using standard FITC filter sets if cells/virus were plated on a separate plate with clear-bottom.	
	Step 5 Prepare Complete 6X Lysis Buffer	To 1 mL of 6X Lysis Buffer, add 30 μ L of 100X protease inhibitor cocktail, and LanthaScreen® Tb-anti-Histone H3K4me2 Antibody to 12 nM. Scale the volume needed to the number of wells $\times 5 \mu$ L/well $\times 1.2$ to ensure extra buffer.	
	Step 6 Add Lysis Buffer (including Tb-Ab)	<ul style="list-style-type: none"> Add 5 μL/well of Complete 6X Lysis Buffer (including Tb-Ab and protease inhibitor) Quick spin the plate at $30 \times g$ for 1 minute (if performing the experiment manually) Incubate plate for ~ 2 to 3 hours at room temperature in the dark 	
	Step 7 Read Plate and Analyze Data	See Terbium TR-FRET Detection on page 8	

* **Growth Medium** for U-2 OS Cells: McCoy's 5A Media supplemented with 10% dFBS, 10 mM HEPES, 0.1 mM NEAA, 1 mM Sodium Pyruvate, and 100 U/mL Penicillin/100 μ g/mL Streptomycin

** Once the optimal BacMam concentration is determined, BacMam reagent can be added to the cell suspension in **Step 1** to the optimal concentration (v/v). Cells/virus mixture can then be plated onto the 384-well assay plate at 20 μ L/well (7,500 cells/well). For **inhibitory compound treatment**, add 5 μ L/well of the 5X compound in Growth Medium and then incubate for 20 to 24 hours prior to **Step 5** and the addition of Lysis Buffer.

Detailed Protocol

For first-time Terbium TR-FRET users, we strongly recommend testing your instrument setup using the set of HIGH and LOW instrument controls provided in the kit prior to setting up the assay. See **Terbium TR-FRET Detection** on page 8 for details.

In this protocol, cells are incubated with virus for 20 to 24 hours at the time of plating onto the assay plate.

Day 1: BacMam Transduction and/or Compound Incubation in a 384-well assay plate

1. Begin with U-2 OS cells in Growth Medium (McCoy's 5A Medium supplemented with 10% dFBS, 10 mM HEPES, 0.1 mM NEAA, 1 mM Sodium Pyruvate, and 100 U/mL Penicillin/100 µg/mL Streptomycin) at 60–95% confluency ($\sim 0.5 \times 10^5$ to 1.0×10^5 cells/cm²). Confluence of cells may impact results, such as the assay window.

Note: For many cell types, such as U-2 OS, a cell seeding density of $\sim 30,000$ cells/cm² for 2–3 days with a harvest density of $\sim 0.5 \times 10^5$ to 1.0×10^5 cells/cm² is optimal. To minimize day-to-day variability, use consistent growth conditions (e.g., similar harvest density, similar passage numbers).

Note: For many cell-types, cryopreserved cells can be transduced immediately following thaw to save culturing time.

2. Harvest the cells by trypsinization and resuspend them in Growth Medium at 3.75×10^5 cells/mL.

Note: The number of cells per well can affect the assay window and should be optimized for your cell background of interest. We recommend starting with 5,000 –20,000 cells per well seeded in 20 µL of growth medium for 384-well format.

3. Plate 20 µL/well of cells onto a 384-well assay plate (about 7,500 cells/well for U-2 OS).
4. *Optional:* Plate the cells in a separate clear-bottom plate for image analysis of GFP expression on Day 2.
5. If the experiment is performed manually, we recommend quickly spinning the assay plates at $30 \times g$ for 1 minute after plating the cells.
6. Prepare BacMam Reagent dilution: for initial testing, we recommend preparing four 2-fold dilutions of the BacMam reagent. (e.g., add 100 µL of BacMam Reagent to 100 µL Growth Medium and mix, then take 100 µL to mix with 100 µL Growth Medium and so on).
7. Add 5 µL of Growth Medium to the Non-transduced Control Wells and 5 µL of undiluted or diluted BacMam Reagent from Step 6 to the Transduced Wells (so that the final virus concentration is 20%, 10%, 5%, 2.5%, and 1.25%).
8. If the experiment is performed manually, we recommend quickly spinning the assay plates at $30 \times g$ for 1 minute.
9. Place cells in a humidified 37°C/5% CO₂ incubator for 20–24 hours to allow for the transduction and expression of the GFP fusion protein.

Note: Once the optimal BacMam concentration is determined, BacMam reagent can be added to the cell suspension in **Step 2** to the optimal concentration (v/v). Cells/virus mixture can then be plated onto the 384-well assay plate at 20 µL/well (7,500 cells/well). For **inhibitory compound treatment**, add 5 µL/well of the 5X compound in growth medium and then incubate for 20 to 24 hrs.

Day 2: LanthaScreen® Cellular Assay

1. *Optional:* After 20 to 24 hours of incubation (Step 9, previous page), analyze GFP expression levels in the clear-bottom plate by fluorescence microscopy using standard FITC filter sets.
2. Prepare **Complete 6X Lysis Buffer** by adding protease inhibitor cocktail at a 1:33 dilution (e.g., 30 µL of 100X stock protease inhibitor per 1 mL of LanthaScreen® 6X Cellular Assay Lysis Buffer) and adding Tb-anti-Histone H3K4me2 antibody at a concentration of 12 nM. Mix by pipetting up and down gently. Store the Lysis Buffer on ice until use (prepare fresh on day of assay).

Note: We recommend scaling the volume of Complete 6X Lysis Buffer as follows:

Number of assay wells to be lysed × 5 µL/well × 1.2 scaling factor

(Scaling factor ensures some extra Lysis Buffer to offset any loss to pipette dispensing tips, dead volumes, etc.; the scaling factor may need to be adjusted based on your dispensing setup).

3. Add 5 µL/well of Complete 6X Lysis Buffer. Cover the plate.
4. *Optional:* Separately add 20 µL/well (in 384-well format) of the HIGH and LOW instrument controls to the blank assay plate wells. We recommend plating a minimum of 3 replicates of each control.
5. If the experiment is performed manually, we recommend quickly spinning the assay plates at 30 × g for 1 minute after adding the Lysis Buffer.
6. Incubate the covered plate at room temperature in the dark for ~2–3 hours or for the appropriate antibody equilibration time. The equilibration time can be optimized for your cell line of interest (typically 60 minutes to several hours).

Note: The assay plate may also be stored at 4°C overnight prior to reading. Allow the plate to warm to room temperature prior to reading. (Longer storage time at 4°C is also possible if evaporation is minimized).

7. Proceed to reading the plate, as described in the next section.

Terbium TR-FRET Detection

Instruments and Filters

Terbium TR-FRET-based cellular assays can be performed on a variety of plate readers, such as the PE Envision. For more information on your particular instrument, refer to www.invitrogen.com/instrumentsetup or contact Invitrogen Discovery Sciences technical support (drugdiscoverytech@lifetech.com or 760-603-7200, option 3, extension 40266).

Reading the Assay Plate and Data Analysis

All measurements should be taken at room temperature from the top of the wells, **with the plate lid or plate seal removed**.

1. Let the assay plate warm to room temperature before reading, if necessary.
2. Set the fluorescence plate reader to top/time-resolved read mode. Allow the lamp in the plate reader to warm up for at least 10 minutes before taking measurements.
3. Remove the lid or plate seal and read the plate using the appropriate filters. Note that filter bandwidths are critical and cannot be approximated. For instrument-specific setup details, refer to www.invitrogen.com/instrumentsetup.

Note: We do not recommend using monochromator-based instruments without adjustable bandwidth, because the sensitivity of these instruments is not sufficient to adequately detect the TR-FRET signal.

4. Calculate the acceptor/donor Emission Ratio (520 nm for acceptor and 490 nm or 495 nm for donor) for each well, by dividing the acceptor emission values by the donor emission values. Do **not** average the 520 nm and 490 nm or 495 nm reading and then take the ratio.
5. *Optional:* Convert the data to Assay Window by dividing each emission ratio value by the value from Non-transduced control wells (or cells not expressing any GFP-Histone H3).

Testing Terbium TR-FRET Detection Using the Instrument Controls

To test your instrument set-up for performing terbium-based TR-FRET cellular assays:

1. Separately add 20 μ L/well (for 384-well format) of the HIGH and LOW instrument controls to blank assay plate wells. We recommend plating a minimum of 3 replicates of each control.
2. Read the plate, as described in the next section.
3. Calculate the acceptor/donor Emission Ratio (GFP/Tb) for each well, by dividing the acceptor emission values (GFP) by the donor emission values (Tb).
4. Determine the HIGH/LOW ratio by dividing the average Emission Ratio from the HIGH control wells by the average Emission ratio from the LOW control wells. This value should be between 2–4, depending on the specific plate reader used.

Appendix

The following alternative transduction protocols may be used with difficult-to-transduce cell lines (**Protocol A**) or other easy-to-transduce cells such as HeLa, T47D, MDA-MB-231 and MCF7 (**Protocol B**).

In short, **Protocol A** requires that the cells be allowed to adhere to the tissue culture flasks prior to the transduction with the BacMam Reagent (this is the longest protocol, but it generally allows for a higher transduction efficiency for difficult cell types). **Protocol B** requires that cells are incubated with the virus at the time of adhering to the tissue-culture plate, usually 24 hours prior to re-plating onto the assay plate.

Alternative Transduction Protocol A (For difficult-to-transduce cells)

In this protocol, cells are allowed to adhere to the tissue-culture flasks before transduction with BacMam Reagent.

Day 1: Adhere Cells

1. Begin with cells grown to complete confluency in normal tissue-culture flasks. Confluency of the cells may impact results, such as the assay window.
2. Trypsinize and harvest the adherent cells as recommended by the cell line manufacturer.
3. Plate the desired number of cells in Growth Medium and allow them to adhere (typically 16–24 hours).

Tip: For many cell types (with a doubling time of approximately 24 hours), a seeding density of approximately 2×10^4 – 4×10^4 cells/cm² results in 50–80% confluency 24 hours after seeding. This has proven optimal for transducing cell lines such as CHO. It may be necessary to optimize the cell density for specific cell backgrounds.

Day 2: Transduce Cells

1. Determine the volume of BacMam Reagent necessary to cover the adhered cells in the tissue culture flask. We recommend using ~1 mL of BacMam solution (diluted as described in the next step) for every 10 cm² of flask surface area.
2. Prepare a dilution of BacMam Reagent (v/v) in Dulbecco's Phosphate Buffered Saline (dPBS) containing Ca²⁺ and Mg²⁺ (Cat. no. 14040-133). We recommend testing a range of v/v dilutions of BacMam Reagent and using 30%, 20%, 10%, 3%, and 1% (v/v) as a starting point (e.g., add 100 µL of BacMam Reagent to 900 µL of dPBS for a 10% v/v dilution).
3. Gently wash the cells once with dPBS containing Ca²⁺ and Mg²⁺.
4. Remove dPBS from Step 3, and gently add the solution of the diluted BacMam Reagent from Step 2 to the cells. Incubate the cells at room temperature (20–25°C) for 3–4 hours with gentle rocking, protected from light.
5. Aspirate the transduction solution from the cell culture dish.
6. Add an appropriate volume of complete cell culture growth medium or growth medium containing the BacMam Enhancer Solution.

Note: BacMam Enhancer Solution (Cat. no. PV5835) can improve the transduction efficiency with difficult-to-transduce cell lines. When using the Enhancer, we recommend testing at least two concentrations of it (e.g., 0.5X and 1X) in comparison to a control without Enhancer to identify the best concentration for your cell-type (i.e., the concentration that yields good GFP-Histone H3 expression with little or no detectable Enhancer-associated toxicity).

7. Incubate cells for 20–24 hours in optimal growth conditions (e.g., a humidified 37°C/5% CO₂ incubator).

Day 3: Harvest & Plate Cells

1. *Optional:* Analyze GFP expression levels by fluorescence microscopy using standard FITC filter sets.
2. Trypsinize and harvest the transduced cells. Be careful not to over-trypsinize the cells because this can result in poor viability and a decreased assay window.
3. Resuspend the harvested cells in Growth Medium with serum to inactivate the trypsin. Centrifuge the cells at $200 \times g$ for 5 minutes. Aspirate the Growth Medium, resuspend the cell pellet in Assay Medium, and plate the cells.

Tip: The number of cells per well can affect the assay window and should be optimized. We recommend starting with 7,500–20,000 cells per well seeded in a 384-well format.

4. Plate 20 μL /well of cells onto a 384-well assay plate
5. *Optional:* Plate cells in a parallel clear-bottom plate for image analysis of GFP expression.
6. If the experiment is performed manually, we recommend quickly spinning the assay plates at $30 \times g$ for 1 minute after plating cells.
7. If you are testing compound effect on the H3K4me2 level, add 5 μL /well of 5X compound in Growth Medium to the Compound-treated wells, and 5 μL /well of Growth Medium to the No-compound wells. If the compound is dissolved in DMSO, make sure that the amount of DMSO per well remains consistent in all the wells.
8. If the experiment is performed manually, we recommend quickly spinning the assay plates at $30 \times g$ for 1 minute after compound addition.
9. Incubate plates for 20–24 hours in a humidified incubator at $37^\circ\text{C}/5\% \text{CO}_2$.

Day 4. LanthaScreen® Cellular Assay

1. Proceed to the **LanthaScreen® Cellular Assay**, page 7.

Alternative Transduction Protocol B (For HeLa, MCF7, T47D, MDA-MB-231 and others)

Day 1: BacMam Transduction

1. Begin with cells (such as HeLa) in Growth Medium* grown to 60–90% confluency ($\sim 0.2 \times 10^5$ to 0.8×10^5 cells/cm²). Confluency of the cells may impact results, such as the assay window.
2. Trypsinize and harvest cells as recommended by the cell line manufacturer.
3. Prepare Growth Medium containing 0.75X BacMam Enhancer Solution by adding 7.5 μ L of the 1,000X Enhancer Solution to 10 mL of Growth Medium.
4. Resuspend cells in Growth Medium containing the BacMam Enhancer Solution at 5×10^5 cells/mL.
5. Plate 2 mL ($\sim 1 \times 10^6$ cells) cell suspension onto each well of a 6-well plate.
6. Immediately after seeding the cells, add the desired amount of BacMam Reagent to the cells. For initial optimization, we recommend testing ~30%, 20%, 10%, 3%, 1%, and 0.3% v/v dilutions of the BacMam Reagent (e.g., for 6-well plate format containing ~ 2 mL of cell suspension per well, add 300 μ L of the BacMam Reagent for an ~10% v/v dilution, and then add 700 μ L of Growth Medium so that the final total volume is 3 mL/well).
7. Place cells in a humidified 37°C/5% CO₂ incubator for 20–24 hours to allow for the transduction and expression of the GFP fusion protein.

Day 2: Cell Harvest and Plating onto a 384-well Assay Plate

1. *Optional:* Analyze GFP expression levels by fluorescence microscopy using standard FITC filter sets.
2. Trypsinize and harvest the transduced cells. Be careful not to over-trypsinize the cells because this can result in poor viability and a decreased assay window.
3. Resuspend the harvested cells in Growth Medium with serum to inactivate the trypsin.
4. Count the cells and then centrifuge them at $200 \times g$ for 5 minutes. Aspirate the Growth Medium and resuspend the cell pellet with Growth Medium (without the BacMam Enhancer solution) at a density of 0.5×10^6 cells/mL.

Tip: The number of cells per well can affect the assay window and should be optimized for your cell background of interest. We recommend starting with 7,500–20,000 cells per well seeded in 20 μ L of growth medium for 384-well format.

5. Plate 20 μ L/well of cells onto a 384-well assay plate.
6. *Optional:* Plate cells in a parallel clear-bottom plate for image analysis of GFP expression.
7. If the experiment is performed manually, we recommend quickly spinning the assay plates at $30 \times g$ for 1 minute after plating cells.
8. If you are testing compound effect on the H3K4me2 level, add 5 μ L/well of 5X compound in Growth Medium to the Compound-treated wells and 5 μ L/well of Growth Medium to the No-compound wells. If the compound is dissolved in DMSO, make sure that the amount of DMSO per well remains consistent in all the wells.
9. If the experiment is performed manually, we recommend quickly spinning the assay plates at $30 \times g$ for 1 minute after compound addition.
10. Incubate plates for 20–24 hours in a humidified incubator at 37°C/5% CO₂.

Day 3. LanthaScreen® Cellular Assay

1. Proceed to the **LanthaScreen Cellular Assay**, page 7.

Troubleshooting Guide

Observation	Potential Solutions
Weak or no expression of GFP-fusion in the cell line of interest in a clear-bottom assay plate.	Confirm that your fluorescence microscope is configured appropriately for detection of GFP/FITC.
	Perform a virus titration to find the optimal virus concentration for your cell background.
	Confirm that the BacMam Reagent is not contaminated. Contaminated reagent looks cloudy or contains chunks of mold.
	For first-time users, we recommend the standard transduction protocol using U-2 OS cells (page 6).
	If the standard protocol works for U-2 OS cells but not for your cells, try the Alternative Transduction Protocol A or B in the Appendix (pages 9 and 11, respectively).
>50% expression of GFP-fusion is observed, but weak or no TR-FRET signal over background is detected.	Confirm that the fluorescence plate reader is configured appropriately for Terbium TR-FRET detection. Filter bandwidth requirements are exact. For more information about your specific instrument and to purchase filters, visit www.invitrogen.com/instrumentsetup . Contact Life Technologies Discovery Sciences Technical Support at 760-603-7200, option 3, extension 40266 or drugdiscoverytech@lifetech.com for more information.
	For first-time users, we recommend following the standard transduction protocol using U-2 OS cells (page 6).
	Perform a cell density experiment to find out the optimal cell harvesting density for your cell line of interest.
	Image the cells in clear-bottom microtiter plates. Ensure that the cells have adhered to the bottom of the plate and are not expressing very high levels of GFP. Dimly green cells are desirable. Excessive expression of the GFP-fusion may be deleterious to cell health.
Day-to-day fluctuations in assay window are observed.	Be sure to use identical growth conditions for culturing your cells (e.g., same harvest density).

For Technical Support for this or other Drug Discovery Products, dial 760-603-7200, option 3, extension 40266.

Purchaser Notification

Limited Use Label License: Research Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

This product is sold under license from Columbia University. Rights to use this product are limited to research use only. No other rights are conveyed. Inquiry into the availability of a license to broader rights or the use of this product for commercial purposes should be directed to Columbia Innovation Enterprise, Columbia University, Engineering Terrace-Suite 363, New York, New York 10027.

THIS PRODUCT IS SOLD UNDER PATENT LICENSE FROM MONSANTO FOR RESEARCH PURPOSES ONLY AND NO LICENSE FOR COMMERCIAL USE IS INCLUDED. REQUESTS FOR LICENSES FOR COMMERCIAL MANUFACTURE OR USE SHOULD BE DIRECTED TO DIRECTOR, MONSANTO CORPORATE RESEARCH, 800 N. Lindbergh, St. Louis, Missouri 63167.

Limited Use Label License No. 308: WPRE Element

This product contains the Woodchuck Post-transcriptional Regulatory Element ("WPRE") which is the subject of intellectual property owned by The Salk Institute for Biological Studies, and licensed to Life Technologies Corporation. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; and/or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. In addition, any use of WPRE outside of this product or the product's authorized use requires a separate license from the Salk Institute. Life Technologies will not assert a claim against the buyer of infringement of patents owned by Life Technologies and claiming this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product or for a Commercial Purpose. If the purchaser is not willing to accept the limitations of this limited use statement, Life Technologies is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, California 92008, Phone (760) 603-7200. Fax (760) 602-6500, or The Salk Institute for Biological Studies, 10010 North Torrey Pines Road, La Jolla, CA 92037, Attn.: Office of Technology Management, Phone: (858) 453-4100 extension 1275, Fax: (858) 546-8093.

Limited Use Label License No. 332: BacMam Virus Use

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer solely in accordance with the accompanying product literature or manual. Purchase of this product does not convey a license to expand, amplify, or otherwise propagate the provided viral particles or to otherwise modify or alter the virus by any means.

LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

© 2011 Life Technologies Corporation. All rights reserved. Reproduction forbidden without permission.

The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

For research use only. Not intended for human or animal therapeutic or diagnostic use.