

**BacMam Histone H3K4me2 Cellular Assay**

Catalog Number: A14164

Literature Part Number: A14164.PIS (MAN0005183)

Literature Lot Number: V1

Revision date: 16 August 2011


**FAST FACTS**

For first-time BacMam users, we recommend using cells like U-2 OS following the detailed protocol available online.

**Optimal Virus %:** We recommend performing a titration of the BacMam Histone H3 Reagent to determine the optimal percentage of virus for transduction in your cell background of interest or when you receive a new lot of virus. Select the lowest percentage of BacMam Reagent that yields the largest assay window (response ratio). See example below.

Component	Part no.	Amount	Storage	Handling
LanthaScreen® Tb-anti-Histone H3K4me2 Antibody	A14139	5 µg	-20°C	<ul style="list-style-type: none"> <li>Protect from light</li> <li>Avoid multiple freeze/thaw cycles</li> </ul>
LanthaScreen® 6X Cellular Assay Lysis Buffer	A12891	6 mL	4°C	<ul style="list-style-type: none"> <li>On the day of assay, supplement with protease inhibitor cocktail and antibody</li> </ul>
BacMam Histone H3 Reagent	A12894	25 mL	4°C	<ul style="list-style-type: none"> <li>DO NOT FREEZE</li> <li>Use sterile technique</li> <li>Avoid extended exposure to ambient room light</li> </ul>
Instrument Control Terbium TR-FRET Kit Low Instrument Control, 1 mL High Instrument Control, 1 mL	A14138	1 kit	4°C	<ul style="list-style-type: none"> <li>Protect from light (do not vortex)</li> </ul>

**Titration of BacMam Histone H3 Reagent in U-2 OS cells (Detection of Histone H3K4me2)**



Additional materials required, but not provided	Source	Part no.
Cell Line of Interest	Various	Various
Protease Inhibitor	Sigma	P8340
White tissue culture-treated, 384-well assay plates	Corning	3570
Fluorescence plate reader with top-read and TR-FRET capability	<a href="http://www.invitrogen.com/instrumentsetup">www.invitrogen.com/instrumentsetup</a> for details	
Optional: Clear-bottom, tissue-culture treated, 384-well plates	Corning	3712

**Detailed Protocols and Additional Assay Performance Data Available**

Visit [www.lifetechnologies.com](http://www.lifetechnologies.com) and search for A14164 to download the full detailed protocol and application note for this assay. Protocols and application note are located under the "How to Use" tab on the product page. Application Notes include assay performance under variable experimental conditions.

**Technical Support**

For additional assistance in running this BacMam-enabled Cellular Assay, contact our technical support team at [drugdiscoverytech@lifetech.com](mailto:drugdiscoverytech@lifetech.com) or 760-603-7200 (enter 3 for "know your party's extension", then enter 40266).

## Quick Reference Protocol for Transduction and LanthaScreen® Cellular Assay using U-2 OS Cells

This quick reference protocol is designed for experienced users using U-2 OS cells, with testing performed in the presence of various concentrations of BacMam Histone H3 reagent. Conditions may need to be optimized for different cell types. For a detailed protocol, see the protocol on our web site.

		Non-transduced Wells	Transduced Wells
BacMam Transduction	<b>Step 1</b> Grow, harvest and plate cells	<ul style="list-style-type: none"> <li>Grow cells in Growth Medium* to 60–95% confluency (<math>\sim 0.5 \times 10^5</math> to <math>1.0 \times 10^5</math> cells/cm<sup>2</sup>).</li> <li>Harvest cells and resuspend in Growth Medium at <math>3.75 \times 10^5</math> cells/mL.</li> <li>Plate 20 <math>\mu</math>L/well cell suspension (about 7,500 cells/well) onto a 384-well assay plate (and optionally a separate plate with clear-bottom for GFP imaging later)</li> <li>Quick spin the plate at <math>30 \times g</math> for 1 minute (if performing the experiment manually)</li> </ul>	
	<b>Step 2**</b> Add BacMam Reagent	Add 5 $\mu$ L/well Growth Medium	Add 5 $\mu$ L/well BacMam GFP-Histone H3 reagent (undiluted or diluted with growth medium to result in different concentrations of BacMam)
	<b>Step 3</b> Incubate Cells/BacMam	<ul style="list-style-type: none"> <li>Quick spin the plate at <math>30 \times g</math> for 1 minute (if performing the experiment manually)</li> <li>Incubate the plate at 37°C/5% CO<sub>2</sub> for 20–24 hours (allows for GFP-Histone H3 expression)</li> </ul>	
LanthaScreen® Histone H3K4me2 Assay	<b>Step 4 (Optional)</b> GFP Imaging	If desired, observe and image GFP-Histone H3 expression under a fluorescence microscope using standard FITC filter sets if cells/virus were plated on a separate plate with clear-bottom.	
	<b>Step 5</b> Prepare Complete 6X Lysis Buffer	For 1 mL of 6X Lysis Buffer, add 30 $\mu$ L 100X protease inhibitor cocktail, and LanthaScreen® Tb-anti-Histone H3K4me2 Antibody to 12 nM. Scale volume needed to the number of wells $\times 5 \mu$ L/well $\times 1.2$ to ensure extra buffer.	
	<b>Step 6</b> Add Lysis Buffer (including Tb-Ab)	<ul style="list-style-type: none"> <li>Add 5 <math>\mu</math>L/ well of Complete 6X Lysis Buffer (including Tb-Ab and protease inhibitor)</li> <li>Quick spin the plate at <math>30 \times g</math> for 1 minute (if performing the experiment manually)</li> <li>Incubate plate for ~2 to 3 hours at room temperature in the dark</li> </ul>	
	<b>Step 7</b> Read Plate and Analyze Data	See <b>Terbium TR-FRET Detection</b> on page 7 in the detailed protocol (available online)	

\***Growth Media** for U-2 OS Cells: McCoy's 5A Media supplemented with 10% dFBS, 10 mM HEPES, 0.1 mM NEAA, 1 mM Sodium Pyruvate, and 100 U/mL Penicillin/100  $\mu$ g/mL Streptomycin

\*\* Once the optimal BacMam concentration is determined, BacMam reagent can be added to the cell suspension in **Step 1** to the optimal concentration (v/v). Cells/virus mixture can then be plated onto the 384-well assay plate at 20  $\mu$ L/well (7,500 cells/well). For **inhibitory compound treatment**, add 5  $\mu$ L/well of the 5X compound in growth medium, and then incubate for 20 to 24 hours prior to **Step 5** and the addition of lysis buffer.

This product is subject to one or more limited use label licenses.

Refer to [www.lifetechnologies.com](http://www.lifetechnologies.com) for the corresponding limited use label licenses.

The performance of this product is guaranteed for six months from the date of purchase if stored and handled properly.

The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

©2011 Life Technologies Corporation. All rights reserved.

---

## Purchaser Notification

### Limited Use Label License: Research Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com) or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

This product is sold under license from Columbia University. Rights to use this product are limited to research use only. No other rights are conveyed. Inquiry into the availability of a license to broader rights or the use of this product for commercial purposes should be directed to Columbia Innovation Enterprise, Columbia University, Engineering Terrace-Suite 363, New York, New York 10027.

THIS PRODUCT IS SOLD UNDER PATENT LICENSE FROM MONSANTO FOR RESEARCH PURPOSES ONLY AND NO LICENSE FOR COMMERCIAL USE IS INCLUDED. REQUESTS FOR LICENSES FOR COMMERCIAL MANUFACTURE OR USE SHOULD BE DIRECTED TO DIRECTOR, MONSANTO CORPORATE RESEARCH, 800 N. Lindbergh, St. Louis, Missouri 63167.

### Limited Use Label License No. 308: WPRE Element

This product contains the Woodchuck Post-transcriptional Regulatory Element ("WPRE") which is the subject of intellectual property owned by The Salk Institute for Biological Studies, and licensed to Life Technologies Corporation. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; and/or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. In addition, any use of WPRE outside of this product or the product's authorized use requires a separate license from the Salk Institute. Life Technologies will not assert a claim against the buyer of infringement of patents owned by Life Technologies and claiming this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product or for a Commercial Purpose. If the purchaser is not willing to accept the limitations of this limited use statement, Life Technologies is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, California 92008, Phone (760) 603-7200. Fax (760) 602-6500, or The Salk Institute for Biological Studies, 10010 North Torrey Pines Road, La Jolla, CA 92037, Attn.: Office of Technology Management, Phone: (858) 453-4100 extension 1275, Fax: (858) 546-8093.

### Limited Use Label License No. 332: BacMam Virus Use

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer solely in accordance with the accompanying product literature or manual. Purchase of this product does not convey a license to expand, amplify, or otherwise propagate the provided viral particles or to otherwise modify or alter the virus by any means.

LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

© 2011 Life Technologies Corporation. All rights reserved. Reproduction forbidden without permission.

The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

For research use only. Not intended for human or animal therapeutic or diagnostic use.