



by Thermo Fisher Scientific

Expi293F™ Cells

Cat. nos.	Size	Conc. 1×10^7 cells/vial
A14527	1 vial	Store in liquid nitrogen
A14528	6 × 1 vial	
Pub. Part no. 100014707	Pub. no. MAN0006283	Rev. A.0

Description

The Expi293F™ cell line is a component of the Expi293™ Expression System. They are a variant of the 293 cell line that has been adapted to high-density suspension growth in Expi293™ Expression Medium. The 293 cell line is a permanent line established from primary embryonic human kidney transformed with sheared human adenovirus type 5 DNA (1,2). The E1A adenovirus gene expressed in 293 cells participates in transactivation of some viral promoters, allowing these cells to produce very high levels of protein.

The Expi293F™ cell line exhibits the following characteristics:

- Prepared from Master Cell Bank cultures derived from parental 293-F cells that were selected for high transfection efficiency and protein expression.
- Adapted to high density, serum-free, suspension growth in Expi293™ Expression Medium.
- Demonstrates high transfection efficiencies with ExpiFectamine™ 293 Reagent.
- Suspension cultures may be transfected in Expi293™ Expression Medium without any requirement for media change.
- Designed for scalable transient transfection and protein expression at culture volumes ranging from less than 1 mL (i.e., multi-well plates) to tens of liters.



As with other human cell lines, when working with Expi293F™ Cells, handle as potentially biohazardous material under at least Biosafety Level-2 (BL-2) containment. For more information on BL-2 guidelines, refer to *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed., published by the Centers for Disease Control, which is available for downloading at: www.cdc.gov/od/ohs/biosfty/bmb15/bmb15toc.htm.

Preparing Expi293™ Expression Medium

Expi293™ Expression Medium (available separately, Cat. no. A14351) is a complete, serum-free, protein-free, and chemically defined medium that does not require any supplementation. Pre-warm the medium to 37°C prior to use.

Thawing Expi293F™ Cells

1. Remove the vial of cells from liquid nitrogen and place it in a 37°C water bath for 1 to 2 minutes to thaw the cells rapidly with gentle agitation. Do **not** submerge the vial in the water.
2. Just before the cells are completely thawed, decontaminate the vial by wiping it with 70% ethanol before opening it in a Class II biological safety cabinet.
3. Transfer the entire contents of the cryovial into a 125-mL polycarbonate, disposable, sterile Erlenmeyer shaker flask containing 29 mL of pre-warmed Expi293™ Expression Medium.
4. Incubate the cells in a 37°C incubator with humidified atmosphere of 8% CO₂ in air on an orbital shaker platform rotating at 125 rpm.
Note: If not using a vent-cap flask, loosen the cap of the flask a quarter turn from snug to allow oxygenation/aeration of the culture.
5. Once the culture has reached $>1 \times 10^6$ viable cells/mL (typically 2–4 days), aseptically take a sample of the cell suspension, and determine viable and total cell counts.
6. Subculture the Expi293F™ Cells by seeding shaker flasks at 3×10^5 viable cells/mL in pre-warmed Expi293™ Expression Medium. We generally use a 125-mL or 250-mL polycarbonate, disposable, sterile, Erlenmeyer flasks containing 30 mL or 60 mL total working volume of cell suspension, respectively.

Guidelines for passaging Expi293F™ Cells

- Subculture the Expi293F™ Cells a minimum of three times to allow them to recover from thawing before using them in transfection experiments.
- Subculture Expi293F™ Cells when they reach a density of approximately 3×10^6 – 5×10^6 viable cells/mL, typically every 3–4 days.
- Split the Expi293F™ culture to 0.3×10^6 – 0.5×10^6 cells.
- When maintaining Expi293F™ Cells, we generally use a 125-mL or 250-mL polycarbonate, disposable, sterile Erlenmeyer flask containing 25–40 mL or 50–80 mL total working volume of cell suspension, respectively.
- Glass flasks without baffles may be used, but thorough cleaning after each use is essential to avoid potential toxicity, which is more problematic in serum-free cultures.

Passaging Expi293F™ Cells

1. Aseptically take a small aliquot of the cell suspension and determine the viable and total cell counts manually using a hemocytometer or electronically using an automated cell counter.
2. Using the cell density determined in Step 1, calculate the split ratio needed to seed the new shaker flask at 0.3×10^6 – 0.5×10^6 viable cells/mL.
3. Dilute the cells in fresh, pre-warmed Expi293™ Expression Medium to give a final cell density of 0.3×10^6 – 0.5×10^6 viable cells/mL in the desired final volume.
4. Incubate flasks in a 37°C incubator containing a humidified atmosphere of 8% CO₂ in air on an orbital shaker platform rotating at 125 rpm.
5. Repeat Steps 1–4 as necessary to maintain or expand cells.

References

1. Graham, F.L., *et al.* (1977) *J. Gen. Virol.* 36, 59–74.
2. Harrison, T., *et al.* (1977) *Virology* 77, 319–329.

Technical support

For additional product and technical information, such as Safety Data Sheets (SDS), Certificates of Analysis, etc., visit our website at www.lifetechnologies.com. For further assistance, email our Technical Support team at techsupport@lifetech.com.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

Disclaimer

TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

LIMITED USE LABEL LICENSE No. 561: Internal Research and Development Use of Mammalian and Insect Cell Lines

Notice to Purchaser: The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product and progeny, to perform internal research and development for the sole benefit of the purchaser. The purchase of this product does not grant the purchaser any additional rights, including, without limitation: (a) the right to use the product or its components in commercial applications of any kind, including bioproduction (e.g. use of the product to manufacture therapeutic agents or diagnostic test components, such as proteins, antibodies, viral particles or virus-like particles), quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration; (b) the right to use the product or its components as a therapeutic agent or diagnostic test component; (c) the right to use the product or its components to produce material for use in human clinical trials; or (d) the right to transfer or resell the product or its components in any form, progeny or derivative. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, California 92008.

© 2015 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries.

For support visit www.lifetechnologies.com/support or email techsupport@lifetech.com

www.lifetechnologies.com