

# pJTI<sup>™</sup> R4 EXP CMV-T0 EmGFP pA Vector

Catalog Number: A15005

Doc. Part no. A15005PIS

Pub. no. MAN0007119

Rev. 2.0

Component	Amount	Storage
pJTI <sup>™</sup> R4 EXP CMV-TO EmGFP pA vector	100 µg in TE buffer, pH 8.0*	-20°C

\* TE Buffer, pH 8.0: 10 mM Tris-HCl, 1 mM EDTA, pH 8.0

### Description

The pJTI<sup> $^{\text{TM}}$ </sup> R4 EXP CMV-TO EmGFP pA vector is a positive control vector for assessing the success of a retargeting reaction in a Jump-In<sup> $^{\text{TM}}</sup> T-REx<sup><math>^{\text{TM}}</sup>$ </sup> parental cell line. When co-transfected with the integrase vector (pJTI<sup> $^{\text{TM}}$ </sup> R4 Int vector, included in the Jump-In<sup> $^{\text{TM}}</sup> T-REx<sup><math>^{\text{TM}}</sup>$ </sup> parental kits) into a Jump-In<sup> $^{\text{TM}}</sup> T-REx<sup><math>^{\text{TM}}</sup>$ </sup> parental cell line containing the genomic R4 site, the Emerald Green Fluorescent Protein (EmGFP) is expressed under the control of the TET repressor. After induction with doxycycline, successfully integrated cells will fluoresce green.</sup></sup></sup>

**FAST FACTS** 

•	For detailed information on the use and propagation of pJTI <sup>™</sup> R4 vectors, refer to the user guide for your Jump-In <sup>™</sup> T-REx <sup>™</sup>
parental cell line available on our website. Visit www.lifetechnologies.com to search for the part number of the Jur	
T-REx <sup>™</sup> parental kit and select manuals.	

• We highly recommend including the pJTI<sup>™</sup> R4 EXP CMV-TO EmGFP pA vector in your Jump-In<sup>™</sup> T-REx<sup>™</sup> retargeting experiment as a positive control along with negative controls (no plasmid DNA, no integrase vector) to visualize the results and optimize the retargeting conditions.

## Map of pJTI<sup>™</sup> R4 EXP CMV-TO EmGFP pA vector

The pJTI<sup> $^{\text{M}}$ </sup> R4 EXP CMV-TO EmGFP pA vector (5799 base pairs) is designed for inducible expression of EmGFP after retargeting into the genomic R4 site of a Jump-In<sup> $^{\text{M}}</sup>$  T-REx<sup> $^{\text{M}}</sup> platform cell line. The vector is intended for use as a positive control for retargeting experiments. Co-transfection with the integrase vector pJTI<sup><math>^{\text{M}}$ </sup> R4 Int into Jump-In<sup> $^{\text{M}}</sup> T-REx<sup><math>^{\text{M}}</sup>$ </sup> parental cell lines containing the genomic R4 site results in Geneticin<sup>®</sup> resistance for selected retargeted cells. The complete sequence of the pJTI<sup> $^{\text{M}</sup></sup> R4 EXP CMV-TO EmGFP pA vector is available for downloading at www.lifetechnologies.com or by contacting Technical Support.</sup></sup>$ </sup></sup>



#### **Technical Support**

For assistance, contact our technical support team at **drugdiscoverytech@lifetech.com** or 760-603-7200 (enter 3 for "know your party's extension", then enter 40266).

For Research Use Only. Not for use in diagnostic procedures.

## Features of pJTI<sup>™</sup> R4 EXP CMV-TO EmGFP pA vector

The pJTI<sup>™</sup> R4 EXP CMV-TO EmGFP pA vector contains the following elements. All features have been functionally tested.

Feature	Benefit
2X TetO <sub>2</sub>	Tetracycline operator sequence (TO) consisting of two tandem 19 nucleotide repeats that serve as a binding site for Tet repressor homodimers
CMV promoter	Human cytomegalovirus immediate-early (CMV) promoter/enhancer for high-level expression in a wide range of mammalian cells
EF1α promoter	Drives the high level expression of the Geneticin <sup>®</sup> resistance gene following the retargeting into the genomic R4 site of the parental Jump-In <sup>™</sup> T-REx <sup>™</sup> cell line
R4 attB site	Allows R4 integrase-mediated integration into the complementary genomic R4 <i>att</i> P site in the parental Jump-In <sup>™</sup> cell line
Ampicillin resistance	Allows the selection and propagation of the vector in <i>E. coli</i>
HSV TK polyA signal	The Herpes Simplex Virus thymidine kinase polyadenylation signal for proper termination and processing of the recombinant transcript
EmGFP	Emerald Green Fluorescent Protein for fluorescent detection of retargeted cells

## Propagating pJTI<sup>™</sup> R4 EXP CMV-T0 EmGFP pA vector

To propagate and maintain the pJTI<sup>T</sup> R4 EXP CMV-TO EmGFP pA vector, we recommend using 10 ng of the vector to transform a *recA*, *endA E. coli* strain such as TOP10F', DH5 $\alpha^{T}$ -T1<sup>R</sup>, TOP10, or equivalent. Select transformants on LB plates containing 50–100 µg/mL of ampicillin. Be sure to prepare a glycerol stock of a transformant containing the plasmid for long-term storage.

### Limited Product Warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

#### **Purchaser Notification**

#### Limited Use Label License: Research Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

#### Limited Use Label License No. 127: GFP with Heterologous Promoter

Notice to Purchaser: This product is sold under license from Columbia University. Rights to use this product are limited to research use only, and expressly exclude the right to manufacture, use, sell or lease this product for use for measuring the level of toxicity for chemical agents and environmental samples in cells and transgenic animals. No other rights are conveyed. Not for human use or use in diagnostic or therapeutic procedures. Inquiry into the availability of a license to broader rights or the use of this product for commercial purposes should be directed to Columbia Innovation Enterprise, Columbia University, Engineering Terrace-Suite 363, New York, New York 10027.

#### Limited Use Label License No. 328: Phi C31 Recombinase Technology

Notice to Purchaser: This product contains site specific integration technology licensed from Stanford University to be used solely for the purchaser's internal research. The purchase of this product does not grant the purchaser any additional rights including, without limitation, the right to use this product for the production of therapeutic agents or diagnostics test components. For information on obtaining additional rights, please contact **outlicensing@lifetech.com** or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

#### Disclaimer

LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

The performance of this product is guaranteed for six months from the date of purchase if stored and handled properly.

© 2013 Life Technologies Corporation. All rights reserved.

The trademarks mentioned herein are the property of Life Technologies Corporation and/or their affiliate(s) or their respective owners.