



Bacterial Alkaline Phosphatase

Cat. No. 18011-015

Size: 2,500 units

Conc.: 150 U/μl

Store at -20°C (not frost-free).

Description:

Bacterial Alkaline Phosphatase (BAP) is a phosphomonoesterase that hydrolyzes 3' and 5' phosphates from DNA and RNA. It is suitable for removing 5' phosphates prior to end-labeling and for dephosphorylating vectors prior to insert ligation. BAP is active at 65°C for at least 1 h. It is sensitive to inhibition by micromolar amounts of inorganic phosphate. The enzyme is purified from *E. coli* C90.

Components:

18011-015 BAP
Y02290 10X Dephosphorylation Buffer

Unit Definition:

One unit hydrolyzes 1 nmol of ATP in 30 min at 37°C.

Storage Buffer:

10 mM Tris-HCl (pH 8.0)
120 mM NaCl
50% (v/v) glycerol

10X Dephosphorylation Buffer:

100 mM Tris-HCl (pH 8.0)
Store buffer at 4°C or -20°C.

Quality Control:

This product has passed the following quality control assays: absence of detectable endodeoxyribonuclease, exodeoxyribonuclease and ribonuclease activities; performance in dephosphorylating blunt-ended DNA. The enclosed buffers were assayed with the enzyme and met quality control specifications.

Doc. Rev. 021202

This product is distributed for laboratory research only. CAUTION: Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established.

For technical questions about this product, call the Invitrogen Tech-LineSM U.S.A. 800 955 6288

Traditional Protocol:

This protocol dephosphorylates 1 pmol of 5'-DNA termini from purified DNA. DNA dephosphorylated by this method is suitable for cloning or for labeling by T4 polynucleotide kinase using the Forward Reaction:

1. Determine the mass of DNA required for 1 pmol of the type of DNA 5' end.
2. To a 1.5-ml microcentrifuge tube, add 4 μ l of 10X Dephosphorylation Buffer and 1 pmol of DNA ends.
3. Add autoclaved, distilled water to 39 μ l.
4. Dilute BAP in dilution buffer such that 1 μ l contains 70 units of enzyme.
5. Incubate the reaction at 65°C for 1 h.
6. Inactivate/remove the BAP according to the protocol described below.

Simplified Protocol:

This protocol allows for the dephosphorylation of DNA directly in restriction endonuclease buffer in the presence of the restriction endonuclease. This is a convenient way of preparing DNA for cloning.

1. Restriction endonuclease digest the vector DNA. (**NOTE:** Heat inactivation of the restriction endonuclease and subsequent purification of the vector DNA are not necessary.)
2. Add 1 μ l of BAP (150 units) to the restriction endonuclease digest.
3. Incubate the reaction at 65°C for 1 h.
4. Inactivate/remove the BAP according to the protocol described below.

Inactivation/Removal of Bacterial Alkaline Phosphatase (Organic Extraction):

Add an equal volume of phenol:chloroform:isoamyl alcohol (25:24:1). Vortex thoroughly and centrifuge at 14,000 \times g at room temperature for 5 min. Carefully remove the upper, aqueous phase and transfer it to a fresh microcentrifuge tube. Add 0.1 volume of 3 M sodium acetate. Vortex. Add 2.5 volumes of 100% EtOH. (**NOTE:** Do not substitute NH₄OAc for NaOAc because NH₄ ions inhibit T4 polynucleotide kinase.) Vortex the mixture thoroughly and centrifuge at 14,000 \times g at room temperature for 5 min.