

## Taq DNA Polymerase PCR Buffer

Cat. No. 18067-017

Size: 2 x 1 ml

Store at  $-20^{\circ}\text{C}$ .

### Description

For use with *Taq* DNA Polymerase (Cat. No. 18038). The PCR Buffer is supplied as a 10X concentrate and should be diluted 1:10 in the final reaction (*e.g.*, use 5  $\mu\text{l}$  in a 50- $\mu\text{l}$  PCR reaction).

Buffer Composition (10X): 200 mM Tris-HCl (pH 8.4), 500 mM KCl.

### Components

10X PCR Buffer minus Mg

50 mM Magnesium Chloride

### Quality Control

No detectable contaminating activity is observed in the endodeoxyribonuclease and ribonuclease assays.

### Guidelines and Recommendations

Critical parameters and troubleshooting information for PCR are well documented (1).

PCR reactions should be assembled in a DNA-free environment. Use of "clean" dedicated automatic pipettes and aerosol-resistant barrier tips are recommended.

Always keep the control DNA and other templates to be amplified isolated from the other components.

## Basic PCR Protocol

The following protocol serves as a general guideline and a starting point for any PCR amplification. Optimal reaction conditions (incubation times and temperatures, concentration of *Taq* DNA Polymerase, primers, MgCl<sub>2</sub>, and template DNA) may vary and need to be evaluated by the customer.

1. Add the following components to a DNase/RNase-free 0.5-ml microcentrifuge tube sitting on ice. Scale the reaction volumes as needed. Prepare a master mix for multiple reactions, to minimize reagent loss and to enable accurate pipetting.

<u>Components</u>	<u>Volume</u>	<u>Final Concentration</u>
10X PCR buffer minus Mg	10 $\mu$ l	1X
10 mM dNTP mixture	2 $\mu$ l	0.2 mM each
50 mM MgCl <sub>2</sub>	3 $\mu$ l	1.5 mM
Primer mix (10 $\mu$ M each)	5 $\mu$ l	0.5 $\mu$ M each
Template DNA	1-20 $\mu$ l	----
<i>Taq</i> DNA Polymerase (5 U/ $\mu$ l)	0.5 $\mu$ l	2.5 units
Autoclaved distilled water to	100 $\mu$ l	

2. Mix contents of tube and overlay with 50  $\mu$ l of mineral or silicone oil.
3. Cap tubes and centrifuge briefly to collect the contents to the bottom.
4. Incubate tubes in a thermal cycler at 94°C for 3 minutes to completely denature the template.
5. Perform 25-35 cycles of PCR amplification as follows:
  - Denature 94°C for 45 seconds
  - Anneal 55°C for 30 seconds
  - Extend 72°C for 1 minute, 30 seconds
6. Incubate for an additional 10 minutes at 72°C and maintain the reaction at 4°C. The samples can be stored at -20° until use.
7. Analyze the amplification products by agarose gel electrophoresis and visualize by ethidium bromide staining.

**“Hot Start” PCR Protocol**

In the "hot-start" method, the addition of *Taq* DNA Polymerase is withheld until the reaction temperature is at 80°C, to ensure high specificity of the products being synthesized.

1. Add all components as in the Basic PCR Protocol, except for the *Taq* DNA Polymerase.
2. Mix contents of tube and overlay with 50 µl of mineral or silicone oil.
3. Cap tubes and centrifuge briefly to collect the contents to the bottom.
4. Incubate tubes in a thermal cycler at 94°C for 3 minutes to completely denature the template.
5. After denaturation at 94°C, maintain the reaction at 80°C.
6. Add 0.5 µl of *Taq* DNA Polymerase (2.5 U) to each reaction. Be certain to add the enzyme beneath the layer of oil.
7. Continue with 25-35 cycles of denaturation, annealing and extension as in the Basic PCR Protocol.

**Reference**

1. Innis, M.A., Gelfand, D.H., Sninsky, J.J., and White, T.J., eds. (1990) *PCR Protocols: A Guide to Methods and Applications*, Academic Press, San Diego, CA.

**Limited Use Label License No. 4: Products for PCR that include no rights to perform PCR**

This product is compatible for use in the Polymerase Chain Reaction (PCR) process claimed in patents owned by Roche Molecular Systems, Inc. and F. Hoffmann-La Roche Ltd. No license under these patents is conveyed expressly, by implication, by estoppel or otherwise to the purchaser by the purchase of this product. A license to use these patented processes for certain research and development activities accompanies the purchase of certain reagents of Applied Biosystems and other licensed suppliers when used in conjunction with an authorized thermal cycler, or is available from Applied Biosystems. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

©2002–2007 Invitrogen Corporation. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.