

## SuperScript® III First-Strand Synthesis SuperMix

**Cat. No: 18080-400**

**Size: 50 reactions**

**Store at –20°C**

### Description

The SuperScript® III First-Strand Synthesis System SuperMix is an optimized SuperMix formulation for first-strand cDNA synthesis from purified poly(A)<sup>+</sup> or total RNA. RNA targets from 100 bp to >12 kb can be detected with this system. The amount of starting material can vary from 0.1 pg to 5 µg of total RNA.

The kit includes SuperScript® III/RNaseOUT™ Enzyme Mix, 2X First-Strand Reaction Mix, and Annealing Buffer. SuperScript® III Reverse Transcriptase. Included in the Enzyme Mix, is a version of M-MLV RT that has been engineered to reduce RNase H activity and provide increased thermal stability (1, 2). The enzyme can be used to synthesize cDNA at a temperature range of 45–55°C, providing increased specificity, higher yields of cDNA, and more full-length product than other reverse transcriptases. Because SuperScript® III RT is not significantly inhibited by ribosomal and transfer RNA, it can be used to synthesize cDNA from total RNA.

RNaseOUT™ Recombinant RNase Inhibitor is included in the Enzyme Mix to safeguard against degradation of target RNA due to ribonuclease contamination.

The 2X First-Strand Reaction Mix includes 10 mM MgCl<sub>2</sub> and 1 mM of each dNTP in a buffer formulation that has been optimized for first-strand synthesis of cDNA. The Annealing Buffer is used in the initial template-primer annealing step. Separate tubes of oligo(dT)<sub>20</sub> and random hexamers are also provided.

cDNA synthesis can be performed using either total RNA or poly(A)<sup>+</sup>-selected RNA primed with oligo(dT), random primers, or a gene-specific primer. Recommendations for PCR amplification of the first-strand cDNA are provided on page 3.

### Kit Contents

<u>Component</u>	<u>Amount</u>
SuperScript® III/RNaseOUT™ Enzyme Mix	100 µL
2X First-Strand Reaction Mix (contains 10 mM MgCl <sub>2</sub> , and 1 mM each dNTP)	500 µL
Annealing Buffer	50 µL
Oligo(dT) <sub>20</sub> (50 µM)	50 µL
Random hexamers (50 ng/µL)	50 µL

### Related Products

<u>Product</u>	<u>Amount</u>	<u>Catalog No.</u>
Platinum® <i>Taq</i> DNA Polymerase	100 units	10966-018
	250 units	10966-026
	500 units	10966-034
Platinum® <i>Taq</i> DNA Polymerase High Fidelity	100 units	11304-011
	500 units	11304-029
Platinum® <i>Pfx</i> DNA Polymerase	100 units	11708-013
	250 units	11708-021
	500 units	11708-039
PureLink™ RNA Mini Kit	50 rxns	12183-018A
TRIZOL® Reagent	100 mL	15596-026
	200 mL	15596-018
FastTrack® 2.0 mRNA Isolation Kit	6 reactions	K1593-02
DNase I, Amplification Grade	100 units	18068-015
Custom Primers	to order, visit <a href="http://www.invitrogen.com">www.invitrogen.com</a>	

### Intended Use

For research use only. Not intended for any animal or human therapeutic or diagnostic use.

### Product Qualification

The Certificate of Analysis provides detailed quality control information for each product. Certificates of Analysis are available on our website at [www.invitrogen.com/support](http://www.invitrogen.com/support).

## Recommendations and Guidelines

### Required Materials

- 0.2-mL, nuclease-free thin-walled PCR tubes
- Thermal cycler, preheated to 65°C
- Ice
- RNase/DNase-free water
- Microcentrifuge
- Vortex mixer

### RNA

- High-quality, intact RNA is essential for full-length, high-quality cDNA synthesis. This kit is designed for use with 0.1 pg to 5 µg of total RNA or 0.1 pg to 500 ng of poly(A)<sup>+</sup> RNA.
- RNaseOUT™ Recombinant RNase Inhibitor is included in the Enzyme Mix to safeguard against degradation of target RNA due to ribonuclease contamination of the RNA preparation.
- To isolate total RNA, we recommend TRIzol® Reagent (Cat. nos. 15596-026/-018) or the PureLink™ RNA Mini-Kit (Cat. no. 12183-018A). Isolation of poly(A)<sup>+</sup> RNA is typically not necessary, although it may improve the yield of specific cDNAs; we recommend the FastTrack® 2.0 mRNA Isolation Kit (Cat. no. K1593-02).
- Small amounts of genomic DNA in the RNA preparation may be amplified along with the target cDNA. If your application requires removal of all genomic DNA from your RNA preparation, we recommend using DNase I, Amplification Grade (Cat. no. 18068-015). DNase I, Amplification Grade, has been extensively purified to remove trace ribonuclease activities commonly associated with other “RNase-free” enzyme preparations, and does not require the addition of placental RNase inhibitor.

### cDNA Synthesis Reaction Conditions

- For difficult or high GC-content templates, use a 55°C cDNA synthesis temperature.
- Keep all components, reaction mixes, and samples on ice. After preparation of the annealing mixture, transfer the mixture to the preheated thermal cycler.

### Primers

First-strand cDNA synthesis can be primed using random hexamers, oligo(dT), or gene-specific primers (GSPs):

- Random hexamers are the most nonspecific priming method, and are typically used for difficult or high GC-content mRNA. Using random hexamers, all RNAs in a population are templates for first-strand cDNA synthesis, and PCR primers confer specificity during PCR.

**Note:** For most RT-PCR applications, 50 ng of random hexamers per 5 µg of total RNA is adequate. Increasing hexamers to 250 ng per 5 µg of RNA may increase yield of small PCR products (<500 bp), but may decrease the yield of longer PCR products and full-length transcripts.

- Oligo(dT), a more specific priming method, is used to hybridize to 3' poly(A) tails, which are found in the vast majority of eukaryotic mRNAs. Since poly(A)<sup>+</sup> RNA constitutes approximately 1% to 2% of total RNA, the amount and complexity of cDNA is considerably less than with random hexamers. We recommend using the oligo(dT)<sub>20</sub> provided in the kit.

**Note:** Oligo(dT) is recommended over random hexamers or GSPs for new mRNA targets. Oligo(dT) produces an RT-PCR product more consistently than random hexamers or GSPs.

- The most specific priming method uses a gene-specific primer (GSP) for the sequence of interest. First-strand synthesis can be primed with the PCR primer that hybridizes nearest to the 3' terminus of the mRNA. Note that some GSPs fail to prime cDNA synthesis even though they work in PCR on DNA templates. If gene-specific priming fails in RT-PCR, repeat first-strand synthesis using oligo(dT) as the primer.

### Guidelines for PCR

The first-strand cDNA from the synthesis reaction may be amplified directly using PCR. We recommend using 10% of the first-strand reaction (2 µL) for PCR. However, for some targets, increasing the amount of first-strand reaction to up to 10 µL may result in increased product yield.

We recommend the following DNA polymerases (for ordering information, see page 1):

- **Platinum® Taq DNA Polymerase** provides automatic hot-start conditions for increased specificity and sensitivity. It is recommended for targets up to 4 kb.
- **Platinum® Taq DNA Polymerase High Fidelity** provides increased fidelity and higher yields for targets up to 15 kb.
- **Platinum® Pfx DNA Polymerase** possesses a proofreading 3' to 5' exonuclease activity and provides maximum fidelity for PCR. It is recommended for targets up to 12 kb.

Consult the product documentation provided with each DNA polymerase for recommended protocols and optimization guidelines. Documentation is also available on our website at [www.invitrogen.com](http://www.invitrogen.com).

## First-Strand cDNA Synthesis

The following procedure is designed to convert 0.1 pg to 5 µg of total RNA or 0.1 pg to 500 ng of poly(A)<sup>+</sup> RNA into first-strand cDNA:

- Mix and briefly centrifuge each component before use. Preheat the thermal cycler to 65°C.
- Combine the following in a 0.2-mL thin-walled PCR tube on ice:
 

<u>Component</u>	<u>Amount</u>
up to 5 µg total RNA	<i>n</i> µL
Primer (50 µM oligo(dT) <sub>20</sub> , or 2 µM gene-specific primer, or 50 ng/µL random hexamers)	1 µL
Annealing Buffer	1 µL
RNase/DNase-free water	to 8 µL
- Incubate in a thermal cycler at 65°C for 5 minutes, and then immediately place on ice for at least 1 minute. Collect the contents of the tube by brief centrifugation.
- Add the following to the tube on ice:
 

2X First-Strand Reaction Mix	10 µL
SuperScript <sup>®</sup> III/RNaseOUT <sup>™</sup> Enzyme Mix	2 µL
- Vortex the sample briefly to mix, and collect by brief centrifugation. Incubate as follows:
 

Oligo(dT) <sub>20</sub> or GSP primed:	50 minutes at 50°C
Random hexamer primed:	5–10 minutes at 25°C, followed by 50 minutes at 50°C
- Terminate the reactions at 85°C for 5 minutes. Chill on ice.

Store the cDNA synthesis reaction at –20°C, or proceed directly to PCR.

## Troubleshooting Guide

Problem	Probable Solution
<b>No or faint bands after analysis of amplified products</b>	
Procedural error in first-strand cDNA synthesis	Repeat the procedure, being careful to follow each step. Be careful to include the Annealing Buffer when adding primers and template for optimal yield.
RNase contamination	Maintain aseptic conditions to prevent RNase contamination. RNaseOUT <sup>™</sup> is included in the Enzyme Mix to inhibit RNases.
Polysaccharide coprecipitation of RNA	Precipitate RNA with lithium chloride to remove polysaccharides, as described in Sambrook <i>et al.</i>
Target mRNA contains strong transcriptional pauses	Use random hexamers instead of oligo(dT) in the first-strand reaction. Increase the reaction temperature to 55°C. Prewarm the 2X First-Strand Reaction Mix to 55°C and hold the primer-template mixture at 55°C, then add the Enzyme Mix. Perform cDNA synthesis at 55°C. Use PCR primers closer to the 3' terminus of the target cDNA.
Too little first-strand product was used in PCR	Use up to 10 µL of the first-strand reaction.
GSP was used for first-strand synthesis	Try another GSP or switch to oligo(dT). Make sure the GSP is the antisense sequence.
Inhibitors of RT present	Remove inhibitors by ethanol precipitation of mRNA preparation before the first-strand reaction. Include a 70% (v/v) ethanol wash of the mRNA pellet. <b>Note:</b> Inhibitors of RT include sodium dodecyl sulfate (SDS), EDTA, guanidinium salts, formamide, sodium pyrophosphate, and spermidine.
<b>Unexpected bands after electrophoretic analysis</b>	
Contamination by genomic DNA	Pretreat RNA with DNase I, Amplification Grade (Cat. no. 18068-015), as described in the DNase I documentation. Design primers that anneal to sequence in exons on both sides of an intron or at the exon/exon boundary of the mRNA to differentiate between amplified cDNA and potential contaminating genomic DNA.
Nonspecific annealing of primers	Vary the annealing conditions. Use Platinum <sup>®</sup> Taq DNA Polymerase for automatic hot-start PCR. Optimize magnesium concentration for each template and primer combination.
Primers formed dimers	Design primers without complementary sequences at the 3' ends.

## References

- Berger, S.L. and Kimmel, A.R. (1987) *Methods Enzymol* 152, 316.
- Bracete, A.M., Mertz, L.M., Fox, D.K. (1999) *Focus*<sup>®</sup> 21, 38.
- Chomczynski, P. (1993) *Biotechniques* Vol. 15, 532.
- Chomczynski, P. and Sacchi, N. (1987) *Anal. Biochem.* 162, 156.
- Compton, T. (1990) in *PCR Protocols: A Guide to Methods and Applications* (Innis, M., Gelfand, D., Sninsky, J., and White, T., eds.), p. 39, Academic Press, Inc.
- Frohman, M.A., Dush, M.K, and Martin, G.R. (1988) *Proc. Nat. Acad. Sci USA* 85, 8998.
- Gerard, G.F. (1994) *Focus*<sup>®</sup> 16, 102.
- Sambrook J., Fritsch, E.F., and Maniatis, T. (1989) *Molecular Cloning: A Laboratory Manual*, Cold Spring Harbor Laboratory Press, Cold Spring Harbor.
- Simms, D., Cizdziel, P.E., and Chomczynski, P. (1993) *Focus*<sup>®</sup> 15, 99.
- Westfall, B., Sitaraman, K., Solus, J., Hughes, J., and Rashtchian, A. (1997) *Focus*<sup>®</sup> 19, 46.
- Westfall, B., Sitaraman, K., Lee, J., Borman, J. and Rashtchian, A. (1999) *Focus*<sup>®</sup> 21, 49.
- Takagi, M., Nishioka, M., Kakahara, H., Kitabayashi, M., Inoue, H., Kawakami, B., Oka, M., and Imanaka, T. (1997) *Appl. Environ. Microbiol.* 63, 4504.
- Sitaraman, K., Darfler, M., and Westfall, B. (1999) *Focus*<sup>®</sup> 21, 10.
- Nathan, M., Mertz, L., Fox, D. (1995) *Focus*<sup>®</sup> 17, 78.
- Schwabe, W., Lee, J.E., Nathan, M., Xu, R.H., Sitaraman, K., Smith, M., Potter, R.J., Rosenthal, K., Rashtchian, A., Gerard, G.F. (1998) *Focus*<sup>®</sup> 20, 30.

## Limited Use Label License No:5 Invitrogen Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For products that are subject to multiple limited use label licenses, the terms of the most restrictive limited use label license shall control. Life Technologies Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Life Technologies Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Life Technologies is willing to accept return of the product with a full refund. For information about purchasing a license to use this product or the technology embedded in it for any use other than for research use please contact Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008 ; Phone (760) 603-7200 or e-mail: [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com)