



Instruction Manual

MultiShot™ StripWell Mach1™-T1^R Chemically Competent *E. coli*

Catalog no. C8696-01

Version C

13 April 2004

25-0651

Table of Contents

| | |
|-----------------------------------|----------|
| Table of Contents..... | iii |
| Important Information..... | v |
| Methods..... | 1 |
| Overview | 1 |
| Transforming Competent Cells..... | 2 |
| Appendix..... | 4 |
| Technical Service..... | 4 |
| Purchaser Notification | 6 |

Important Information

Shipping and Storage

The MultiShot™ StripWell Mach1™-T1^R Chemically Competent *E. coli* kit is shipped on dry ice. Upon receipt, store the kit at -80°C.

Contents

The MultiShot™ StripWell Mach1™-T1^R Chemically Competent *E. coli* kit contains the following reagents. Fifty microliters of chemically competent cells are supplied per well with a transformation efficiency of $>1 \times 10^9$ cfu/ μ g DNA.

| Item | Composition | Amount |
|--|--|----------------------|
| S.O.C. Medium (store at room temperature or +4°C) | 2% Tryptone 0.5% Yeast Extract 10 mM NaCl 2.5 mM KCl 10 mM MgCl ₂ 10 mM MgSO ₄ 20 mM glucose | 2 x 15 ml |
| Mach1™-T1 ^R Cells | --- | 1 96-stripwell plate |
| pUC19 Control DNA | 10 pg/ μ l in 5 mM Tris-HCl, 0.5 mM EDTA, pH 8 | 50 μ l |

Genotype

F⁻ ϕ 80(*lacZ*) Δ M15 Δ *lacX74* *hsdR*(r_K⁻m_K⁺) Δ *recA1398* *endA1* *tonA*

Product Qualification

- Competent cells are tested for transformation efficiency using the control plasmid included in the kit and following the procedure on page 2. Transformation efficiency should be $>1 \times 10^9$ cfu/ μ g of pUC19 DNA.
 - Untransformed cells are tested for the appropriate antibiotic sensitivity, the absence of phage contamination, and resistance to phage T5 (a standard test that demonstrates resistance to phage T1).
-

continued on next page

Important Information, continued

**Information for
Non-U.S.
Customers**

For European Customers

The Mach1™-T1^R *E. coli* strain is genetically modified to carry the *lacZ*ΔM15 *hsdR lacX74 recA endA tonA* genotype. As a condition of sale, this product must be in accordance with all applicable local legislation and guidelines including EC Directive 90/219/EEC on the contained use of genetically modified organisms.

For All Non-U.S. Customers

The parental strain of Mach1™-T1^R *E. coli* is the non-K-12, wild-type W strain (ATCC #9637, S. A. Waksman). Although the parental strain is generally classified as Biosafety Level 1 (BL-1), we recommend that you consult the safety department of your institution to verify the Biosafety Level.

Methods

Overview

Introduction

The MultiShot™ StripWell Mach1™-T1^R Chemically Competent *E. coli* is chemically competent Mach1™-T1^R *E. coli* packaged in a 96-stripwell plate to simplify medium and high throughput bacterial transformation. The stripwell format allows you to use only the number of wells you need for your particular application.

Features of the Strain

The Mach1™-T1^R *E. coli* strain is modified from the wild-type W strain (ATCC #9637, S. A. Waksman) and has a faster doubling time compared to other standard cloning strains. With Mach1™-T1^R cells, you can visualize colonies 8 hours after plating on ampicillin selective plates. You can also prepare plasmid DNA 4 hours after inoculating a single, overnight-grown colony in the selective media of choice. Note that this feature is not limited to ampicillin selection.

The Mach1™-T1^R *E. coli* strain possesses several additional features that make it an ideal strain to use for most cloning applications. These features include:

- *lacZ*ΔM15 for blue / white color screening of recombinants
 - *hsdR* mutation for efficient transformation of unmethylated DNA from PCR applications
 - *ΔrecA1398* mutation for reduced occurrence of homologous recombination in cloned DNA
 - *endA1* mutation for increased plasmid yield and quality
 - *tonA* mutation to confer resistance to T1 and T5 phage
-

Transforming Competent Cells

Introduction

A procedure is provided below to transform MultiShot™ StripWell Mach1™-T1^R Chemically Competent *E. coli*. Depending on your application of interest, you may modify this procedure to fit your needs.



Important

Do not electroporate chemically competent cells. The salt content of the buffer will cause arcing and kill the cells.

Before Starting

Perform the following before starting the transformation procedure:

- Prepare a container of ice large enough to chill the number of wells you will be using.
- Warm the vial of S.O.C. Medium (supplied with the kit) to room temperature.
- Preheat a water bath to 42°C.
- Warm the selective plates in a 37°C incubator for 30 minutes (use one plate for each transformation). If you are including the pUC19 control, make sure that you have one LB agar plate containing 100 µg/ml ampicillin.

Note: For optimal growth of Mach1™-T1^R *E. coli* cells, it is essential that selective plates are prewarmed to 37°C prior to spreading.

Transformation Procedure

We recommend including the pUC19 control plasmid DNA supplied with the kit in your transformation experiment to verify the efficiency of the competent cells.

1. Remove a MultiShot™ StripWell plate from the freezer and remove the number of wells you need. Return any unused wells to the freezer. Place the wells in the container of ice. Cells should thaw within 1 minute.
 2. Carefully remove the strip of caps from each set of 8 wells and keep them for further use.
-

continued on next page

Transforming Competent Cells, continued

Transformation Procedure, continued

3. Use a multi-channel pipet to add 2-5 μl DNA (2 pg-20 ng) to the wells. Keep the volume consistent between wells for uniform results. If you are transforming the pUC19 control, add 1 μl (10 pg) into one well.
 4. After adding the DNA, cover the wells with the caps and incubate the cells and DNA on ice for 30 minutes.
 5. Transfer the wells to the water bath and heat-shock for 30 seconds at 42°C. **Note:** Be careful not to contaminate the cells.
 6. Transfer the wells back to the ice and allow the wells to cool for 1 minute.
 7. Remove the caps and add 250 μl of room temperature S.O.C. Medium to each well. **Re-cap the wells tightly.**
 8. Incubate the wells at 37°C for 1 hour with shaking (225 rpm). We turn the wells on their side to increase aeration and secure them to the shaker.
 9. For each sample, spread 25-100 μl on a **prewarmed** selective plate. Store the remaining transformation mix at +4°C. Additional cells may be plated out the next day, if desired.
 10. Invert the plates and incubate at 37°C. If you are using ampicillin selection, visible colonies should appear within 8 hours, and blue/white screening can be performed after 12 hours. If you are selecting transformants with an antibiotic other than ampicillin, incubate plates overnight.
 11. Select overnight-grown colonies and analyze by plasmid isolation, PCR, or sequencing. For plasmid isolation, inoculate a single, overnight-grown colony in 2 ml of **prewarmed** selective media (*e.g.* LB + ampicillin, LB + kanamycin, LB + Zeocin™, *etc.*). For optimal results, we recommend inoculating as much of the single colony as possible. Shake at 37°C for 4 hours before isolating the plasmid.
-

Appendix

Technical Service

World Wide Web



Visit the Invitrogen Web Resource using your World Wide Web browser. At the site, you can:

- Get the scoop on our hot new products and special product offers
- View and download vector maps and sequences
- Download manuals in Adobe® Acrobat® (PDF) format
- Explore our catalog with full color graphics
- Obtain citations for Invitrogen products
- Request catalog and product literature

Once connected to the Internet, launch your web browser (Internet Explorer 5.0 or newer or Netscape 4.0 or newer), then enter the following location (or URL):

<http://www.invitrogen.com>

...and the program will connect directly. Click on underlined text or outlined graphics to explore. Don't forget to put a bookmark at our site for easy reference!

Contact Us

For more information or technical assistance, call, write, fax, or email. Additional international offices are listed on our Web page (www.invitrogen.com).

Corporate Headquarters:

Invitrogen Corporation
1600 Faraday Avenue
Carlsbad, CA 92008 USA
Tel: 1 760 603 7200
Tel (Toll Free): 1 800 955 6288
Fax: 1 760 602 6500
E-mail:
tech_service@invitrogen.com

European Headquarters:

Invitrogen Ltd
Inchinnan Business Park
3 Fountain Drive
Paisley PA4 9RF, UK
Tel: +44 (0) 141 814 6100
Tech Fax: +44 (0) 141 814 6117
E-mail:
eurotech@invitrogen.com

continued on next page

Technical Service, continued

MSDS Requests

To request an MSDS, visit our Web site at www.invitrogen.com. On the home page, go to 'Technical Resources', select 'MSDS', and follow instructions on the page.

Limited Warranty

Invitrogen is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, contact our Technical Service Representatives.

Invitrogen warrants that all of its products will perform according to specifications stated on the certificate of analysis. The company will replace, free of charge, any product that does not meet those specifications. This warranty limits Invitrogen Corporation's liability only to the cost of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. Invitrogen reserves the right to select the method(s) used to analyze a product unless Invitrogen agrees to a specified method in writing prior to acceptance of the order.

Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore Invitrogen makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Service Representatives.

Invitrogen assumes no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.

Purchaser Notification

**Limited Use
Label License
No. 144:
Mach1™
Competent
Cells**

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for commercial purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any commercial purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for commercial purposes. Commercial purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of patents owned by Invitrogen and claimed in this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Invitrogen Corporation, 1600 Faraday Avenue, Carlsbad, California 92008. Phone (760) 603-7200. Fax (760) 602-6500.

©2003 -2004 Invitrogen Corporation. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.

Notes

Notes



Corporate Headquarters:

*Invitrogen Corporation
1600 Faraday Avenue
Carlsbad, California 92008
Tel: 1 760 603 7200
Tel (Toll Free): 1 800 955 6288
Fax: 1 760 603 7229
Email: tech_service@invitrogen.com*

European Headquarters:

*Invitrogen Ltd
3 Fountain Drive
Inchinnan Business Park
Paisley PA4 9RF, UK
Tel (Free Phone Orders): 0800 269 210
Tel (General Enquiries): 0800 5345 5345
Fax: +44 (0) 141 814 6287
Email: eurotech@invitrogen.com*

International Offices:

*Argentina 5411 4556 0844
Australia 1 800 331 627
Austria 0800 20 1087
Belgium 0800 14894
Brazil 0800 11 0575
Canada 800 263 6236
China 10 6849 2578
Denmark 80 30 17 40*

*France 0800 23 20 79
Germany 0800 083 0902
Hong Kong 2407 8450
India 11 577 3282
Italy 02 98 22 201
Japan 03 3663 7974
The Netherlands 0800 099 3310
New Zealand 0800 600 200
Norway 00800 5456 5456*

*Spain & Portugal 900 181 461
Sweden 020 26 34 52
Switzerland 0800 848 800
Taiwan 2 2651 6156
UK 0800 838 380
For other countries see our Web site*

www.invitrogen.com

