

Electrocomp™ Kits

Catalog Numbers C664-55, C664-11, C664-24, C665-11, C665-24

Publication Part Number 280078

Publication Number MAN0000621

Revision A.0

Contents

- 5, 10, or 30 × 100 µL electrocompetent *E. coli*
- 10 pg/µL pUC19 supercoiled vector in 10 mM Tris-HCl, 1 mM EDTA
- SOC medium, 15 mL

Genotypes

TOP10F⁺: F⁺{*lacI^qTn10*(Tet^R)} *mcrA* Δ(*mrr-hsdRMS-mcrBC*) Φ80*lacZ*Δ*M15* Δ*lacX74* *recA1* *araD139* Δ(*ara-leu*)7697 *galU* *galK* *rpsL* *endA1* *nupG*

TOP10: F⁻ *mcrA* Δ(*mrr-hsdRMS-mcrBC*) Φ80*lacZ*Δ*M15* Δ*lacX74* *recA1* *araD139* Δ(*ara-leu*)7697 *galU* *galK* *rpsL* *endA1* *nupG*

Prepare for electroporation

You must keep cells and cuvettes cold until after the pulse is applied for high efficiency electroporation. After the pulse is applied, you must immediately add SOC to the cells to ensure good results. Use the following steps to prepare for electroporation.

1. Pre-chill 0.1-cm cuvettes on ice or at -20°C.
2. Set up the electroporation device for bacterial electroporation following the directions provided by the manufacturer. See the following table of settings for selected electroporators. If you are using an electroporator not listed in the table, you may need to vary your settings to achieve optimal transformation efficiency.

Unit	Capacitance	Resistance	Voltage
Bio-Rad® Gene Pulser® II	25 µF	200 Ω	2.5 kV
BTX® ECM™ 630	50 µF	150 Ω	1.5 kV

3. Prepare an ice bucket for the cells and the DNA.
4. Thaw the SOC medium to add to the cells after the pulse. Keep the SOC at room temperature.
5. Prepare sterile 15-mL tubes for use after the pulse has been delivered

Electroporation

6. Place 1–5 µL of DNA in a sterile microcentrifuge tube and chill it on ice. For the pUC19 control DNA, use 1 µL (10 pg).

Caution: Resuspend the DNA in water (rather than TE buffer) to keep the ionic strength to a minimum. Otherwise, arcing might occur during the pulse; arcing can damage the machine and results in no transformation and cell death. Precipitate or desalt DNA in ligation or restriction buffer before electroporation.

7. Gently thaw the cells on ice. Use these cells immediately; do not leave them on ice for an extended period of time. You may re-freeze unused cells for later use, but they will suffer a significant loss of efficiency.
8. Transfer the desired amount of cells to the pre-chilled microcentrifuge tube containing the DNA. We recommend using 20–40 µL (40–80 µL if using 0.2 cm cuvettes) for a small-scale electroporation (< 100 ng of DNA), and 160 µL for a large-scale or library electroporation. The volume of DNA added should not exceed 5% of the total cell/DNA mixture. Mix the cells gently. *Do not mix by pipetting up and down.* Leave the cells on ice for 1 minute.

9. After 1 minute, transfer the cells to a chilled cuvette and gently shake them to the bottom of the cuvette. Check the cuvette by looking at the contents from both sides to make sure the cells make contact all the way across the bottom of the chamber without any air bubbles. Do this as quickly as possible, being careful not to warm up the cuvette and cells. Remove the condensation from the outside surfaces of the cuvette with a tissue.
10. Electroporate your sample using the settings listed or using the manufacturer's recommended settings.
11. Immediately add 480 µL (for small-scale) or 960 µL (for large-scale) of SOC medium to the cells.
12. Transfer the suspension to a 15-mL tube. Incubate the tube at 37°C in a rotary shaking incubator at 225 rpm for 1 hour to allow expression of the antibiotic resistance.
13. Plate the cells on prewarmed selective plates and incubate overnight at 37°C. We recommend plating 2 different volumes to ensure that at least 1 plate has well-spaced colonies.

Note: For the pUC19 control DNA, prepare a 1:10 dilution in SOC medium and plate 25 µL and 100 µL of cells on LB agar with 100 µg/mL of ampicillin.

Safety Data Sheets (SDSs)

Safety Data Sheets (SDSs) are available at www.lifetechnologies.com/support.

Certificate of Analysis

The Certificate of Analysis provides detailed quality control and product qualification information for each product. Certificates of Analysis are available on our website. Go to www.lifetechnologies.com/support and search for the Certificate of Analysis by product lot number, which is printed on the box.

Important licensing information

This product may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.



DISCLAIMER: LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

© 2014 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. Bio-Rad and Gene Pulser are registered trademarks of Bio-Rad Laboratories, Inc.

BTX and ECM are registered trademarks of Harvard Bioscience, Inc.

For support visit lifetechnologies.com/support or email techsupport@lifetech.com

lifetechnologies.com

18 June 2014

