


SuperTaq™ Plus Polymerase

Catalog Number AM2054, AM2056

Pub. No. 4393891 Rev. B

Contents	Quantity	Storage conditions
SuperTaq™ Plus Polymerase, 5 U/μL	AM2054: 50 Units AM2056: 250 Units	Store at -20°C. <i>Do not store in a frost-free freezer.</i>
10X Long PCR Buffer: 100 mM Tris-HCl (pH 9.0), 500 mM KOAc, 15 mM MgSO ₄ , stabilizers	1.25 mL	
10X Long PCR Buffer (-)MgSO ₄ : 100 mM Tris-HCl (pH 9.0), 500 mM KOAc, stabilizers	1.25 mL	
25 mM MgSO ₄	1 mL	
PCR dNTP Mix: 2.5 mM each dATP, dCTP, dGTP, dTTP	200 μL	

 **WARNING!** Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from www.lifetechnologies.com/support.

Product description

SuperTaq™ Plus Polymerase is an extended-range Taq polymerase containing a proof-reading activity that increases the fidelity of DNA replication and reduces error within products. SuperTaq™ Plus Polymerase produces a high yield of PCR products for fragments greater than 1 kb (up to 20 kb), allowing the use of less enzyme per reaction.

SuperTaq™ Plus Polymerase is supplied with dNTPs and two 10X Long PCR Buffers, one with and one without MgSO₄. A MgSO₄ solution is also supplied for titration into the 10X Long PCR Buffer (-)MgSO₄.

Source: A mixture of an *E. coli* strain overexpressing *Thermus aquaticus* DNA polymerase and a proprietary thermostable DNA polymerase with proofreading activity.

Unit (U) definition: One unit of SuperTaq™ Plus DNA Polymerase incorporates 10 nmol of deoxynucleotides into acid insoluble form in 30 minutes at 74°C. Assay Conditions: 25 mM TAPS (pH 9.3 at 25°C), 50 mM KCl, 2 mM MgCl₂, 200 μM each dATP, dCTP, dGTP, dTTP, 0.5 mg/mL activated salmon sperm DNA, 1 mM 2-mercaptoethanol.

Storage Buffer (not included): 50 mM Tris-HCl (pH 8.0), 1 mM EDTA, 1 mM DTT, 50% (v/v) glycerol, stabilizers.

Using SuperTaq™ Plus Polymerase

Note: The elongation step must be carried out at 68°C to maximize production of long PCR products.

Unsatisfactory PCR reactions carried out under ordinary buffer conditions can often be enhanced by increasing the concentration of Mg²⁺ to 1.75 mM in the final reaction volume. Dilution of the enzyme in a smaller volume, thus increasing its concentration, is advisable for difficult PCR reactions such as random amplification of polymorphic DNA (RAPD), PCR of tailed DNA, and single-cell copy PCR.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

For Research Use Only. Not for use in diagnostic procedures.

The information in this guide is subject to change without notice.

DISCLAIMER: LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

NOTICE TO PURCHASER: LIMITED USE LABEL LICENSE: Research Use Only : The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

© 2012 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners. SuperTaq is a trademark of Enzyme Technologies Ltd..

For support visit lifetechnologies.com/support or email techsupport@lifetech.com

lifetechnologies.com